

Test Code	Test Name	Components	Specimen	Method	Comments	Report
Z785	ABORTION PANEL	*CD3 *CD19 *CD16 + CD56*TNF-alpha *HLA DR, DQA1 & DQB1 (Class II) Typing	2 mL (1 mL min.) serum from 1 SST AND 6 mL (4 mL min.) whole blood in 2 Lavender Top (EDTA) tubes AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry, CLIA, PCR		Sample by Mon 4 pm; Report Wed
H022	ABSOLUTE EOSINOPHIL COUNT; AEC		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedance, VCS		Daily
H052	ABSOLUTE LYMPHOCYTE COUNT		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedance, VCS		Daily
H051	ABSOLUTE NEUTROPHIL COUNT		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedance, VCS		Daily
RW002	ACETYLCHOLINE RECEPTOR (ACHR) BLOCKING ANTIBODY		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
S220	ACETYLCHOLINE RECEPTOR (ACHR) BINDING ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Avoid general anaesthetic or muscle relaxant drugs 24 hours prior to sampling.	Radioimmunoassay	This assay should not be requested for patients who have recently received radio-isotopes therapeutically or diagnostically due to potential assay interference. Also see MYASTHENIA GRAVIS PANEL.	Sample by Wed / Sat 9 am; Report Same day

RW003	ACETYLCHOLINE RECEPTOR (ACHR) MODULATING ANTIBODY		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Avoid general anaesthetic or muscle relaxant drugs 24 hours prior to sampling.	Radiobinding Assay		Sample by 7th of the month; Report after 2–3 weeks
RW004	ACETYLCHOLINE STERASE, AMNIOTIC FLUID		3 mL (1 mL min.) amniotic fluid in a sterile screw capped vial. Ship refrigerated or frozen. Indicate gestational age, collection date, clinical indications and the AFP and MoM results if available on the Test Request Form.	Gel Electrophoresis	Sample should be drawn between 14-21 weeks gestation	Sample by 7th of the month; Report after 2–3 weeks
B054	ACID PHOSPHATASE, PROSTATIC	*Acid Phosphatase Total*Acid Phosphatase Prostatic	2 mL (0.5 mL min) serum from 1 SST. Keep whole specimen refrigerated until serum can be separated. Separate serum, add 2 drops (60 uL) 5% Acetic acid and mix. Ship refrigerated or frozen.	Spectrophotometry		Daily
B053	ACID PHOSPHATASE, TOTAL		2 mL (0.5 mL min.) serum from 1 SST. Keep whole specimen refrigerated until serum can be separated. Separate serum, add 2 drops (60 uL) 5% Acetic acid and mix. Ship refrigerated or frozen.	Spectrophotometry		Daily
R040	ACTH		2 mL (0.5 mL min.) plasma from 1 Lavender Top (EDTA) tube. Draw blood between 8–10 a.m. Mix well, separate and transfer plasma IMMEDIATELY into a sterile screw capped vial and FREEZE. Ship frozen. DO NOT THAW.	Chemiluminescent Immunoassay		Sample Daily by 4 pm; Report Same day

E016	ACTIVATED PROTEIN C - RESISTANCE; APCR		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hours. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician	Sample by 2nd/4th Sat by 9 am; Report Next Mon
OS100	ACTN3 (SPORTS GENE), GENOTYPING		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Thu 11 am; Report 10 days

Z774	ACUTE MYELOID LEUKEMIA (AML), CYTOGENETICS PANEL	*t (8;21) or LSI ETO/AML1 * inv (16) (p13q22) or t(16;16) (p13;q22) or CBFB * t (15;17) or LSI PML/RARA * t(variable;11q23) ; MLL gene breakapart *Chromosome analysis, Hematological malignancy	8 mL (6 mL min.) whole blood from 2 Green Top (Sodium Heparin) tubes AND 4 mL (2 mL min.) Bone Marrow from 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE.	FISH & Conventional Karyotyping		Sample Daily by 6 pm; Report 7 Days
Z805	ACUTE KIDNEY INJURY PANEL	*Cystatin C *Neutrophil Gelatinase - Associated Lipocalin (NGAL) *Creatinine, serum *Fractional excretion of Sodium	5 mL (3 mL min.) serum from 1 SST AND 10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	CMIA, Nephelometry, Jaffe's reaction, Indirect ISE		Sample by Mon / Wed / Fri 4 pm; Report Same Day

Z775	ACUTE LYMPHOBLASTIC LEUKEMIA (ALL), CYTOGENETICS PANEL	*BCR/ABL or Philadelphia translocation *t(variable;11q23) ; MLL gene breakapart * t (12;21) (p13;q22); TEL/AML1 *Chromosome analysis, Hematological malignancy	8 mL (6 mL min.) whole blood from 2 Green Top (Sodium Heparin) tubes AND 4 mL (2 mL min.) Bone Marrow from 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE.	FISH & Conventional Karyotyping		Sample Daily by 6 pm; Report 7 Days
R156	ADENOSINE DEAMINASE ; ADA		2 mL (1 mL min.) serum from 1 SST OR 2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube OR 2 mL (1 mL min.) Body Fluid / CSF in a sterile screw capped container. Ship refrigerated or frozen. Hemolysed samples are not acceptable.	Spectrophotometry		Daily
M211	ADENOVIRUS ANTIGEN FOR GASTROENTERI TIS. STOOL		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pmReport Same day
N126	ADENOVIRUS, QUALITATIVE, PCR		3 ml (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
M213	ADENOVIRUS ANTIGEN FOR RESPIRATORY INFECTIONS		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pmReport Same day

OS102	ADRB2 GENOTYPING FOR BETA-2 AGONIST RESPONSIVENES S		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Sat 11 am; Report 10 days
M094	AFB IDENTIFICATION , RAPID; PNB TEST		Submit pure growth of AFB on appropriate media in screw capped container in sealed plastic bag. Ship refrigerated OR when used with other AFB culture tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Mycobacterium TB Complex / MOTT indentification.	Automated Fluorescent	This test can be used in conjunction with other AFB Cultures on specific request as additional test. For Nucleic Acid Probes, use test code N001; for AFB Antibiotic Sensitivity, use appropriate test code. Request as additional tests.	Sample Daily by 3.30 pm; Report 15 days
P030	AFB STAIN, MISCELLANEOU S	Auramine and Ziehl Neelsen staining	Miscellaneous specimens (Urine, Pus, Body Fluids like pleural, peritoneal, pericardial, gastric, bronchial, semen etc.) in a sterile, screw capped container OR Swabs in sterile normal saline OR Bone marrow in 1 Green Top (Sodium Heparin) tube OR smeared onto 4 clean glass slides fixed immediately in 95% ethanol for 30 minutes and air dried. For gastric lavage , patient should be fasting for previous 12 hours. Transport the gastric lavage specimens within 4 hrs of collection. Specify type of specimen on container and test request form. Ship refrigerated.	Fluorescent and conventional microscopic examination	For tissue specimens use test code J050 AFB Stain	Sample Daily by 3.30 pm; Report Same day

R084	AFP (ALPHA FETOPROTEIN), AMNIOTIC FLUID		3 mL (1 min.) Amniotic fluid. Avoid contamination with maternal blood. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Also see MATERNAL SERUM SCREENING PANELS.	Daily
F001	AFP (ALPHA FETOPROTEIN), TUMOR MARKER		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Give brief clinical history.	Chemiluminescent Immunoassay	Recommended for cancer follow up only.	Daily
R030	AFP (ALPHA FETOPROTEIN), MATERNAL		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	The test is recommended only between 14–22 weeks of pregnancy. Also see MATERNAL SERUM SCREENING PANELS.	Daily
P094	ALBERT STAIN		Submit Throat swabs in Aerobic transport medium available from LPL OR 2 air dried smears on sterile glass slides. Ship refrigerated. DO NOT FREEZE.	Gram staining and Metachromatic staining		Daily
B094	ALBUMIN, CSF		2 mL (1 mL min.) Cerebrospinal fluid (CSF) in a sterile screw capped container. Ship refrigerated or frozen.	Nephelometry		Sample Mon through Sat by 4 pm; Report Same day
P086	ALBUMIN, FLUID		2 mL (1 mL min.) Body fluid in a sterile screw capped container. Ship refrigerated or frozen.	BCG		Daily
Z088	ALBUMIN & IgG, CSF	*CSF IgG/ Albumin Ratio	2 mL (1 mL min.) Cerebrospinal fluid (CSF) in a sterile screw capped container. Ship refrigerated or frozen.	Nephelometry	Also see MULTIPLE SCLEROSIS PANELS.	Sample Mon through Sat by 4 pm; Report Same day
B039	ALBUMIN, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	BCG		Daily

G184	ALCOHOL ETHYL GLUCURONIDE (EtG) & ETHYL SULPHATE (EtS), QUANTITATIVE ASSAY, URINE		5 mL (2 mL min.) aliquot of random urine without preservative. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 5 pm; Report Next day
K017	ALCOHOL (ETHANOL) SCREEN, SERUM		3 mL (2 mL min.) serum from 1 SST. Cleanse collection site with NON-ALCOHOLIC agent before drawing blood. Ship refrigerated or frozen.	Alcohol Dehydrogenase UV		Daily
S024	ALDEHYDE CHOPRA TEST FOR KALA AZAR		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Tube clotting		Daily
B052	ALDOLASE		2 mL (0.5 mL min.) serum from 1 SST OR 2 mL (0.5 mL min) plasma from 1 Green Top (Sodium Heparin) / Lavender Top (EDTA) tube. Hemolysed specimens are not acceptable. Ship refrigerated or frozen.	Spectrophotometry		Daily
R049	ALDOSTERONE, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Patient should be ambulatory / upright 2 hours prior to sampling. Ship refrigerated or frozen. Overnight fasting is preferred.	CLIA	A complete investigation of Aldosterone pathophysiology must include assessment of PLASMA RENIN, ACTIVITY.	Sample by Fri 9 am; Report Same day

B175	ALDOSTERONE / PLASMA RENIN DIRECT RATIO		2 mL (1 mL min.) serum from 1 SST AND 3 mL (2.5 mL min.) plasma from 1 Lavender Top (EDTA) tube. Collect blood between 8.00-10.00 am. Patient should be ambulatory / upright 2 hours prior to sampling. Separate plasma in a plastic transport vial and freeze. Ship frozen. DO NOT THAW. Overnight fasting is preferred.	CLIA		Sample by Fri 9 am; Report Same day
U039	ALDOSTERONE, 24-HOUR URINE		5 mL (1 mL min.) aliquot of 24-hour urine. Keep refrigerated during collection. Add boric acid (1g/100 ml urine) to the 24-hour collection to maintain pH below 7.5. Mix well and aliquot. Ship refrigerated or frozen. Record 24-hour volume on urine container and test request form.	CLIA	A complete investigation of Aldosteronopathophysiology must include assessment of PLASMA RENIN, ACTIVITY.	Sample by Thu 4 pm; Report Next day
E011	ALKALINE PHOSPHATASE (ALP) ISOENZYMES	Total alkaline phosphatase & Isoenzymes	2 mL (1.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Hemolysed specimens are not acceptable. Overnight fasting is preferred.	Agarose Gel Electrophoresis	Also see NTx (COLLAGEN CROSS-LINKED N-TELOPEPTIDE);PYRIDINIUM COLLAGEN CROSS -LINKS	Sample by Thu 9 am; Report Next day
R121	ALKALINE PHOSPHATASE, BONE-SPECIFIC; BSAP		2 mL (0.5 mL min.) serum from 1 SST. Hemolysed specimens are not acceptable. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see NTx (COLLAGEN CROSS-LINKED N-TELOPEPTIDE); PYRIDINIUM COLLAGEN CROSS LINKS	Sample by Tue 9 am; Report Same day
B041	ALKALINE PHOSPHATASE; ALP		2 mL (0.5 mL min.) serum from 1 SST. Hemolysed specimens are not acceptable. Ship refrigerated or frozen. Overnight fasting is preferred.	IFCC		Daily

G182	ALKAPTONURIA, URINE, QUANTITATIVE	Homogentisic Acid	15 mL (10 mL min.) aliquot of random urine in a sterile screw capped container without preservative. Ship refrigerated or frozen.	Gas Chromatography/Ma ss Spectrometry		Sample Daily by 5 pm; Report 5 days
G102	ALKAPTONURIA, URINE, QUALITATIVE		10 mL (5 mL min.) aliquot of freshly voided urine in a sterile screw capped container. No preservative required. Wrap container in aluminium foil to protect from light. Ship refrigerated or frozen. Clinical and drug history must accompany sample.	Chemical		Daily
Z684	ALLERGY ASTHMA / RHINITIS SCREENING PANEL ADULT	Common Ragweed,Dermat ophagoides pteronysinus, Dermatophagoid es farinae, Alterneria tenuis,Aspergillus fumigatus, Candida albicans,Penicilliu m notatum, House dust & Cockroach	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z695	ALLERGY NON- VEGETARIAN PANEL 2	Fish, Shrimp, Lobster Mussel, Tuna & Salmon	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
A038	ALLERGY, BUCKWHEAT (Kutu)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A077	ALLERGY, CASEIN		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A030	ALLERGY, CHICKEN (Murgi)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A007	ALLERGY, CLADOSPORIU M HERBARUM		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A071	ALLERGY, COCOA		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A014	ALLERGY, COW DANDER		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A015	ALLERGY, DOG DANDER		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A090	ALLERGY, ELM (Papdi Chibil)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A088	ALLERGY, GOOSEFOOT (Bathua)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A008	ALLERGY, HOUSE DUST		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A066	ALLERGY, LENTIL (Dal)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A034	ALLERGY, MELALEUCA		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A006	ALLERGY, MILK (Doodh)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A096	ALLERGY, MULBERRY (Shahtoot)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A035	ALLERGY, OAT		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A046	ALLERGY, ONION (Pyaz)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A067	ALLERGY, PEA (Matar)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A003	ALLERGY,PENICILLIUM NOTATUM		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
Z531	ALLERGY VEGETABLE COMPREHENSIVE PANEL	Spinach, Cabbage, Paprika Tomato, Mushroom, Garlic, Onion, Celery, Carrot, Potato, Cucumber & Brinjal	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z687	ALLERGY VEGETABLE PANEL 2	Tomato, Spinach, Cabbage & Paprika	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days

Z526	ALLERGY WEED PANEL	Ragweed, Goosefoot, Common pigweed, Cocklebur, Mugwort & Rapeseed	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z161	ALLERGY DUST PANEL	House dust, Mites likeDermatophag oides Pteronyssinus, Dermatophagoid es Farinae & Cockroach	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z162	ALLERGY TREE PANEL	Acacia, Olive, Willow, Melaleuca, White pine & Eucalyptus	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
A093	ALLERGY, ALDER (Bhidur)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11 am; Report Same Day
A049	ALLERGY, APPLE (Seb)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A005	ALLERGY, ASPERGILLUS FUMIGATUS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A050	ALLERGY, BANANA (Kela)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A040	ALLERGY, BEEF		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A069	ALLERGY, CAROB		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A059	ALLERGY, CARROT (Gajar)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A047	ALLERGY, CELERY (Ajmoda)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A073	ALLERGY, CHEDDAR (Cheese)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A070	ALLERGY, CHICKPEA (Kabuli chana)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A010	ALLERGY, COCKROACH		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A089	ALLERGY, COMMON PIGWEED (Kantewali Chauli)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A095	ALLERGY, COTTONWOOD		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A061	ALLERGY, CUCUMBER (Kheera)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A041	ALLERGY, EGG YOLK (Anda)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A022	ALLERGY, BERMUDA GRASS (Cynodon dactylon / Doob)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A097	ALLERGY, COCKLEBUR (Bichha / Kutta)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A016	ALLERGY, EGG WHITE (Anda)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A045	ALLERGY, GARLIC (Lasoon)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A019	ALLERGY, GLUTEN (Lasa)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A013	ALLERGY, HORSE DANDER		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A023	ALLERGY, JOHNSON GRASS (Sorghum halepense / Jowar)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A057	ALLERGY, LEMON (Neembu)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A080	ALLERGY, MANGO (Aam)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A084	ALLERGY, MEADOW GRASS (Lawn grass)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A091	ALLERGY, MESQUITE (Pahari Keekar)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A009	ALLERGY, MITES	Dermatophagoides pteronyssinus, Dermatophagoides farinae	3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A025	ALLERGY, MUG WORT (Artemisia vulgaris / Barna)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A079	ALLERGY, MUTTON (Meat)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A094	ALLERGY, OAK		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A032	ALLERGY, OLIVE		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A048	ALLERGY, ORANGE (Santara)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A054	ALLERGY, PAPRIKA		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A056	ALLERGY, PEAR (Nashpati)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A062	ALLERGY, PECAN		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A063	ALLERGY, CASHEWNUT (Kaju)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A076	ALLERGY, COCONUT (Nariyal)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A027	ALLERGY, EUCALYPTUS (Nilgiri / Safeda)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A078	ALLERGY, LOBSTER		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A109	ALLERGY, MUSHROOM (Kukurmutta)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A051	ALLERGY, PEACH (Aaru)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A058	ALLERGY, PINEAPPLE (Ananas)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A039	ALLERGY, PORK		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A060	ALLERGY, POTATO (Aloo)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A098	ALLERGY, RAPESEED (Sarson)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A074	ALLERGY, RYE		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A082	ALLERGY, RYEGRASS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A037	ALLERGY, SESAME (Til)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A017	ALLERGY, SOYBEAN		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A052	ALLERGY, SPINACH (Palak)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A083	ALLERGY, TIMOTHY GRASS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A121	ALLERGY, TUNA		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A042	ALLERGY, TURKEY		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A018	ALLERGY, WHEAT (Gehu)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A068	ALLERGY, WHITE BEAN		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A033	ALLERGY, WILLOW		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A011	ALLERGY, CAT EPITHELIUM & DANDER		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A072	ALLERGY, COFFEE		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A031	ALLERGY, FISH (Machli)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A036	ALLERGY, MAIZE (Bhutta)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A119	ALLERGY, MUSSEL		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A021	ALLERGY, PEANUT (Moongfali)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A029	ALLERGY, RICE (Chawal)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A043	ALLERGY, TOMATO (Tamatar)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

Z089	ALLERGY COMPREHENSIV E PROFILE	Allergy Screen, 7 Allergy Panels & 32 Individual allergens	6 mL (4 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z092	ALLERGY POST SCREENING PROFILE	7 Allergy Panels & 32 Individual allergens	6 mL (4 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z164	ALLERGY ANIMAL PANEL	Cat epithelium. Cow, Dog & Horse Dander	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z685	ALLERGY ASTHMA / RHINITIS SCREENING PANEL PEDIATRIC	Common Ragweed, Dermatophagoid es pteronysinus, Dermatophagoid es farinae, Aspergillus fumigatus, Candida albicans, Penicillium notatum, Egg white, Milk & Cockroach	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days

Z690	ALLERGY ECZEMA PANEL FOR NON- VEGETARIANS	Wheat, Rice, Peanut, Soybean, Potato, Milk, Dermatophagoid es pteronysinus, Dermatophagoid es farinae, Fish, Shrimp, Egg white & Meat	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z689	ALLERGY ECZEMA PANEL FOR VEGETARIANS	Wheat, Rice, Peanut, Lentil, Soybean, Potato, Milk, Dermatophagoid es pteronyssinus & Dermatophagoid es farinae	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z691	ALLERGY FOOD PANEL2	Egg white, Milk, Wheat, Soybean, Fish & Peanut	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z524	ALLERGY FRUIT COMPREHENSIV E PANEL	Orange, Apple, Banana, Peach, Mango, Strawberry, Pear, Lemon, Pineapple & Coconut	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen	ImmunoCAP		Sample Daily by 11 am; Report 3 days

Z693	ALLERGY, FRUIT PANEL 2	Strawberry, Pear, Lemon & Pineapple	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z525	ALLERGY GRASS PANEL	Bermuda grass, Rye grass, Timothy grass, Meadow grass, Johnson grass & Bahia grass	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z160	ALLERGY MOLD PANEL	Penicillium notatum, Cladosporium herbarum, Aspergillus fumigatus, Candida albicans & Alternaria tenuis	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z523	ALLERGY NUT PANEL	Pecan, Cashew nut, Pistachio, Walnut, Almond & Peanut	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z522	ALLERGY PULSES PANEL	Lentils, White bean & Chick pea	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
A112	ALLERGY PHADIATOP, PEDIATRIC		3 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11am; Report Same Day
A001	ALLERGY SCREEN	Total IgE and Phadia Top	3 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report Same Day

Z527	ALLERGY ECZEMA COMPREHENSIVE PANEL	Egg white, Cow's milk, Fish, Shimp, Wheat, Rice, Peanut, Lentil, Soybean, Chick pea, Lemon, Banana, Mites, Dog dander, Cow dander, Horse dander & Cat epithelium	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z091	ALLERGY FOOD PANEL 1	Egg white, Milk, Wheat, Soybean, Fish , Chicken, Rice, Gluten, Shrimp & Peanut	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z520	ALLERGY NON- VEGETARIAN PANEL 1	Pork, Beef, Egg, Turkey & Chicken	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z686	ALLERGY VEGETABLE PANEL 1	Tomato, Yeast, Garlic, Onion & Celery	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z688	ALLERGY VEGETABLE PANEL 3	Carrot, Spinach, Potato & Cucumber	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days

Z529	ALLERGY WHEEZE / RHINITIS COMPREHENSIV E PANEL ADULT	Elm, Bermuda grass, Rye grass, Timothy grass, Johnson grass, Meadow grass, Bahia grass, Alder, Mugwort, Common pigweed, Gluten, Pea, Casein, Housedust, Cocoa, Coffee, Cheddar, Mesquite, Ragweed, Goosefoot, Cat epithelium, Dog dander, Cockroach, AlternariaÇ	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
------	--	---	---	-----------	--	--

Z528	ALLERGY WHEEZE / RHINITIS COMPREHENSIV E PANEL PEDIATRIC	Elm, Lemon, Egg white, Milk, Rice, Lentil, Pea, White bean, Mesquite, Ragweed, Goosefoot, Cat epithelium, Dog dander, Mites, Cock roach, Alternaria tenuis, Aspergillus fumigatus, Candida albicans, Penicillium notatum & Cladosporium herbarum	6 mL (4 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
Z521	ALLERGY CEREAL PANEL	Oat, Maize, Sesame, Buckwheat, Rye, Wheat, Rice & Carob	6 mL (4 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
A024	ALLERGY, ACACIA (Babool)		3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11 am; Report Same Day
A075	ALLERGY, ALMOND (Badam)		3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11 am; Report Same Day

Z545	ALLERGY, ANY 10 ALLERGENS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
Z544	ALLERGY, ANY 5 ALLERGENS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A086	ALLERGY, BAHIA GRASS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A092	ALLERGY, BRINJAL (Baingan)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A053	ALLERGY, CABBAGE (Bandgobi)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A004	ALLERGY, CANDIDA ALBICANS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A111	ALLERGY PHADIATOP, ADULT		3 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report Same Day
Z692	ALLERGY FRUIT PANEL 1	Orange, Apple, Banana & Peach	6 mL (4 mL min.) serum from 2 SST's.Ship refrigerated or frozen.	ImmunoCAP		Sample Daily by 11 am; Report 3 days
A028	ALLERGY, ALTERNARIA TENUIS		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11 am; Report Same Day
A064	ALLERGY, PISTACHIO (Pista)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A120	ALLERGY, SALMON		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day

A020	ALLERGY, SHRIMP (Prawn / Jhinga)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A065	ALLERGY, WALNUT (Akhrot)		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A026	ALLERGY, WHITE PINE		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A087	ALLERGY, RAGWEED		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A055	ALLERGY, STRAWBERRY		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
A044	ALLERGY, YEAST		3 mL (2 ml min.) serum from 1 SST. Ship refrigerated or frozen.	ImmunoCAP Specific IgE		Sample Daily by 11am; Report Same Day
B096	ALPHA-1- ANTITRYPSIN QUANTITATION; AAT		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum from cells immediately. Ship refrigerated or frozen. Overnight fasting is preferred.	Nephelometry		Sample by Mon/Thu 4 pm; Report Same day
B194	ALPHA-2- MACROGLOBUL IN		2 mL (0.5 mL min.) serum from 1 SST. Avoid hemolysis. Overnight fasting is preferred. Ship refrigerated or frozen.	Nephelometry		Sample by Tue / Fri 4 pm; Report Same Day
RW009	ALPHA SUBUNIT		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Radioimmunoassay		Sample by 7th of the month; Report after 2–3 weeks

OS103	ALPHA THALASSEMIA MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Tue 11 am ; Report 9 days
OS104	ALPHA THALASSEMIA MUTATION DETECTION, PRENATAL		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube each from both the parents AND 10 mL Amniotic fluid / CVS. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Thu 11 am; Report 9 days
RW006	ALPHA-1-ACID GLYCOPROTEIN		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Overnight fasting is preferred. Ship refrigerated or frozen.	Fixed Rate Time Nephelometry		Sample by 7th of the month; Report after 2–3 weeks
RW007	ALPHA-1- ANTITRYPSIN (AAT) PHENOTYPE		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Isoelectrofocusing		Sample by 7th of the month; Report after 2–3 weeks
C001	ALUMINIUM, DIALYSIS FLUID		10 mL (5 mL min) dialysis fluid in metal free container available from LPL. Collect specimen directly into the container. Ship refrigerated or frozen. Use powderless gloves during specimen collection. Patients should refrain from taking antacids and aluminium containing medications prior to specimen collection.	ICPMS	Use only specified containers available from LPL.	Sample Daily by 6 pm; Report Next day

C002	ALUMINIUM, PLASMA		3 mL (1 mL min.) plasma from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection. Patients should refrain from taking antacids and aluminium containing medications prior to specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C003	ALUMINIUM, SERUM		3 mL (1.0 mL min.) serum from 1 metal free White Top (Z-No Additive) tube available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection. Patients should refrain from taking antacids and aluminium containing medications prior to specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
Z094	AMENORRHOEA PROFILE	*Estradiol *FSH *LH *Prolactin *Progesterone* TSH	3 mL (1 mL min.) serum from 1 SST. Collect specimen 4 hrs after the patient has awakened. Overnight fasting is preferred. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	For male patients use test code Z099.	Daily
Z479	AMENORRHOEA / ANDROLOGY SCREEN	*FSH *LH *Prolactin	3 mL (1 mL min.) serum from 1 SST. Collect specimen 4 hrs after the patient has awakened. Overnight fasting is preferred. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily

RW270	AMIKACIN		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
G175	AMINO ACIDS QUANTITATIVE, PLASMA, 45 AMINO ACIDS / FULL PANEL		2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Transfer plasma into a sterile screw capped vial. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	LC-MS/MS		Sample by Mon / Wed 5 pm; Report 3 days
U099	5-AMINOLEVULINIC ACID (5 -ALA), QUANTITATIVE, URINE		10 mL (5 mL min.) aliquot of 24 hour urine collected in a dark coloured container / wrapped completely with black paper or aluminium foil. Ship refrigerated or frozen. Record 24 hour volume on test request form and urine container. Uncovered vials exposed to light are unacceptable.	Column Chromatography	Also see LEAD; PORPHYRINS, TOTAL; PORPHOBILINOGEN; LEAD POISONING PANEL	Sample by Tue / Fri 9 am; Report Nextday
G178	AMINO ACID, NON-KETOTIC HYPERGLYCINEMIA PANEL QUANTITATIVE, CSF & PLASMA		1 mL (0.5 mL min.) CSF in a sterile screw capped vial AND 2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Transfer plasma into a sterile screw capped vial. Ship refrigerated or frozen. Clinical details and drug history must accompany sample. CSF & Plasma specimens should be drawn at the same time.	LC-MS/MS		Sample by Mon / Wed 5 pm; Report 3 days
G176	AMINO ACID QUANTITATIVE, CSF, 45 AMINO ACIDS / FULL PANEL		1 mL (0.5 mL min.) CSF. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	LC-MS/MS		Sample by Mon / Wed 5 pm; Report 3 days

G179	AMINO ACIDS, MAPLE SYRUP URINE DISEASE (MSUD) PANEL, PLASMA		2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Transfer plasma into a sterile screw capped vial. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	LC-MS/MS		Sample by Mon / Wed 5 pm; Report 3days
G106	AMINO ACIDS, QUALITATIVE, CSF		1 mL (0.5 mL min.) CSF. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	Thin Layer Chromatography, One dimensional		Sample by Mon / Wed / Fri 4 pm; Report 2 days
G107	AMINO ACIDS, QUALITATIVE, TWO DIMENSIONAL, URINE		5 mL (2 ml min.) aliquot of first morning urine in a sterile screw capped container. No preservative required. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	Thin Layer Chromatography, Two dimensional		Sample by Mon / Wed / Fri 4 pm; Report 2 days
G109	AMINO ACIDS, QUALITATIVE, URINE & PLASMA		5 mL (2 ml min.) aliquot of first morning urine in a sterile screw capped container. No preservative required AND 2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Transfer plasma into a sterile screw capped vial. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	Thin Layer Chromatography, One dimensional		Sample by Mon / Wed / Fri 4 pm;Report 2 days
G174	AMINO ACIDS, QUANTITATIVE URINE, 45 AMINO ACIDS / FULL PANEL		5 mL (2 ml min.) aliquot of first morning urine in a sterile screw capped container. No preservative required. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	LC-MS/MS		Sample by Mon / Wed 5 pm; Report 3days

U033	5-AMINOLEVULINIC ACID (5 -ALA), QUALITATIVE, URINE		10 mL (5 mL min.) random urine in a dark coloured container / wrapped completely with black paper or aluminium foil. Ship refrigerated or frozen.	Column Chromatography	Also see LEAD; PORPHYRINS, TOTAL; PORPHOBILINOGEN; LEAD POISONING PANEL	Sample by Tue / Fri 9 am; Report Nextday
G104	AMINO ACIDS, QUALITATIVE, ONE DIMENSIONAL, URINE		5 mL (2 ml min.) aliquot of first morning urine in a sterile screw capped container. No preservative required. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	Thin Layer Chromatography, One dimensional		Sample by Mon / Wed / Fri 4 pm; Report 2 days
G165	AMINO ACIDS QUANTITATIVE, BLOOD, 10 AMINO ACIDS	*Alanine *Arginine *Citrulline *Glycine *Leucine*Methionine *Ornithine *Phenylalanine *Tyrosine *Valine	1 drop of heel prick blood each on 3 spots of filter paper available from LPL OR 2 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Clinical details and drug history must accompany sample. Ship refrigerated. DO NOT FREEZE.	Tandem Mass Spectrometry		Sample by Mon / Wed / Fri 4 pm; Report 2 days
G105	AMINO ACIDS, QUALITATIVE, PLASMA		2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Transfer plasma into a sterile screw capped vial. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	Thin Layer Chromatography, One dimensional		Sample by Mon / Wed / Fri 4 pm; Report 2 days

G177	AMINO ACIDS, QUANTITATIVE, URINE & PLASMA , 45 AMINO ACIDS / FULL PANEL		5 mL (2 ml min.) aliquot of first morning urine in a sterile screw capped container. No preservative required. AND 2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Transfer plasma into a sterile screw capped vial. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	LC-MS/MS		Sample by Mon / Wed 5 pm; Report 3days
RW012	AMITRIPTYLINE	Amitriptyline & metabolite Nortriptyline	3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: >12 hours post oral dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
Z535	AML CATEGORIZATION, PCR	*PML RARA t(15;17)(q22;q12)) GENE REARRANGEMENT *AML ETO t(8;21) GENE REARRANGEMENT *Inv 16 (p13 q22) / t(16;16)(p13;q22)) GENE REARRANGEMENT	3 mL (2 mL min.) Whole blood / Bone marrow from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N045	AML ETO t(8;21) GENE REARRANGEMENT, PCR QUALITATIVE		3 mL (2 mL min.) Whole blood / Bone marrow from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat

B061	AMMONIA, BLOOD		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Separate plasma IMMEDIATELY and transfer onto ice. Ship refrigerated or frozen. Avoid hemolysis.	Spectrophotometry	Patients to report to National Reference Lab, Rohini, New Delhi for testing. Outstation samples are not acceptable.	Daily
S014	AMOEBIC SEROLOGY, IgG, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed / / Sat 9 am; Report Same day
K013	AMPHETAMINE SCREEN, URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily
U017	AMYLASE (DIASTASE), URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated. DO NOT FREEZE	Spectrophotometry		Daily
RW013	AMYLASE ISOENZYMES		3 mL (1 mL min.) serum from 1 SST. Avoid hemolysis. Ship refrigerated or frozen.	Enzymatic		Sample by 7th of the month; Report after 2–3 weeks
P053	AMYLASE, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen. Specify type of fluid on test request form and sample container.	Spectrophotometry		Daily
B050	AMYLASE, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Avoid hemolysis. Ship refrigerated or frozen.	Spectrophotometry		Daily

Z015	ANALYZER 26; SMA 26 TEST PANEL	*Glucose (F&PP) *Urea *Creatinine *UricAcid*Choles terol *Triglycerides *CPK *HBDH *LDH *SGOT *SGPT *GGTP *Bilirubin, Total *Bilirubin, Direct *Bilirubin, Indirect *Protein, Total*Albumin *A:G Ratio *Alkaline Phosphatase *Calcium *Phosphorus *Sodium *PotasÇ	3 mL (2 mL min.) serum from 1 SST, AND 2 mL (1.5 mL min) plasma from 1 Grey Top (Sodium Fluoride) tube EACH for fasting and PP glucose. Avoid use of tourniquet and fist clenching during sampling. 12 hours overnight fasting essential. Ship refrigerated or frozen.	Spectrophotometry, Indirect ISE, Hexokinase		Daily
------	--------------------------------------	---	--	---	--	-------

Z018	ANALYZER 18; SMA 18 TEST PANEL	*Glucose (F) *Urea *Creatinine *Uric Acid*Cholesterol *SGOT *SGPT *GGTP *Bilirubin,Total *Protein, Total *Albumin *A:G Ratio*Alkaline Phosphatase *Calcium *Phosphorus*So dium *Potassium *Chloride	3 mL (2 mL min.) serum from 1 SST AND 2 mL (1.5 mL min) plasma from 1 Grey Top (Sodium Fluoride) tube. Avoid use of tourniquet and fist clenching during sampling. 12 hours overnight fasting essential. Ship refrigerated or frozen.	Spectrophotometry, Indirect ISE, Hexokinase		Daily
S171	ANCA (ANTINEUTROPHIL CYTOPLASMIC ANTI- BODIES), IFA	Includes differentiation of perinuclear (p-ANCA) and cytoplasmic (c-ANCA) pattern.	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	For SERIAL DILUTIONS, use test code S172. Also see ANCA-MPO, ANCA PR3.	Sample Daily by 11 am; Report Same Day
S172	ANCA (ANTINEUTROPHIL CYTOPLASMIC ANTI- BODIES), IFA IN DILUTIONS		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	Also see ANCA-MPO, ANCA PR3.	Sample Daily by 11 am; Report Next day

S048	ANCA-MPO (ANTI MYELOPEROXID ASEANTIBODIES) , EIA	Corresponds with p-ANCA	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay	Also see ANCA-PR3; ANTINEUTROPHIL CYTOPLASMIC ANTIBODIES by immunofluorescence .	Sample by Tue / Thu / Sat 9 am;Report Same day
S047	ANCA-PR3 (SERIN PROTEINASE 3ANTIBODIES), EIA	Corresponds with c-ANCA	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay	Also see ANCA- MPO; ANTINEUTROPHIL CYTOPLASMIC ANTIBODIES by immunofluorescence .	Sample by Tue / Thu / Sat 9 am; Report Same day
Z009	ANDROSTENEDI ONE STIMULATION BY ACTH		3 mL (1.5 mL min.) serum from 1 SST for each timed specimen. Draw fasting baseline sample; administer 250 ug of synthetic ACTH (Synacthen) IV, draw additional sample at 90 min. Ship refrigerated or frozen. Both tubes must be clearly marked with time drawn. Submit both tubes with 1 test request form.	Chemiluminescent Immunoassay	Synthetic ACTH (Synacthen) is notsupplied by LPL.Refer to Y Test codes :Y014 Basal Y015 90 min. Y017	Sample by Mon / Thu 4 pm ; ReportSame day
Z553	ANDROGEN INDEX	*Sex Hormone Binding Globulin *Testosterone, Total	3 mL (2 mL min.) serum from 1 SST. Specify age and sex on test request form. Ship refrigerated or frozen.	CLIA		Sample by Mon / Thu 4 pm; Report Same day

R060	ANDROSTENEDIONE		3 mL (1.5 mL min.) serum from 1 SST. An early morning specimen is preferred. Ship refrigerated or frozen. Specify age and sex on test request form.	Chemiluminescent Immunoassay	In females, preferred sampling time is 1week before or after the menstrual period. Also See ANDROSTENEDIONES TIMULATION BY ACTH.	Sample by Mon / Thu 4 pm; ReportSame day
Z099	ANDROLOGY PANEL	*FSH *LH *Testosterone, Total & Free	3 mL (2 mL min.) serum from 1 SST. Specify age and sex on test request form. Ship refrigerated or frozen.	CLIA, RIA	In females, preferred sampling time isDay 2 / Day 3 of the menstrual cycle.	Sample by Mon through Sat 9 am; Report Same day
Z288	ANEMIA PANEL BASIC	*Complete Blood Count *Serum Iron Studies*Ferritin *CRP	3 mL (1.5 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. Overnight fasting is preferred.	Electrical Impedence, VCS, Spectrophotometry, CLIA, Immunoturbidimetry		Daily
Z289	ANEMIA PANEL 1	*Complete Blood Count *Serum Iron Studies*Ferritin *Protein Total *Albumin *A:G Ratio*Transferrin *CRP *Hemoglobin Electrophoresis	3 mL (1.5 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. Overnight fasting is preferred.	Electrical Impedence, VCS, Spectrophotometry, CLIA, Immunoturbidimetry , HPLC		Sample by 9 am; Report Next day

Z513	ANEMIA PANEL 2	*Complete Blood Count *Serum Iron Studies*Ferritin *Protein Total *Albumin *A:G Ratio*Transferrin *CRP *Lead Blood	3 mL (1.5 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube AND 2 mL (1 mL min.) whole blood from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. Overnight fasting is preferred.	Electrical Impedence, VCS, Spectrophotometry, CLIA, Immunoturbidimetry, AAS		Sample Daily by 6 pm; Report Next day
N053	ANEUPLOIDY DETECTION, AMNIOTIC FLUID		5 mL (3 mL min.) Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Analysis Test request form.	Quantitative Fluorescence PCR		Sample by Tue / Fri 11 am; Report Fri / Tue
N054	ANEUPLOIDY DETECTION, BLOOD		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Analysis Test request form.	Quantitative Fluorescence PCR		Sample by Tue / Fri 11 am; Report Fri / Tue
N055	ANEUPLOIDY DETECTION, PRODUCTS OF CONCEPTION		Submit curretted tissue in normal saline. Fetus not accepted. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Analysis Test request form.	Quantitative Fluorescence PCR	Formalin fixed specimen not preferred	Sample by Tue / Fri 11 am; Report Fri/ Tue
OS105	ANGELMAN SYNDROME		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Methylation specific PCR		Sample by Sat 11 am; Report 10 days
RW015	ANGIOTENSIN II		2 mL (1 mL min) plasma from 1 Lavender Top (EDTA) tube. Ship refrigerated or frozen.	Immunoassay	See also PLASMA RENIN, ACTIVITY.	Sample by 7th of the month; Report after 2–3 weeks
B063	ANGIOTENSIN CONVERTING ENZYME; ACE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Spectrophotometry		Daily

Z290	ANTENATAL PANEL BASIC	*Complete Blood Count *Blood Group + Rh Factor *Urine R/E *HIV-1 & 2 *HBsAg *VDRL*Glucose F / PP / R	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 3 mL (1.5 mL min.) serum from 1 SST AND 10 mL (5 mL min.) aliquot of random urine. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence,VCS, EMT, Dipstick, CMIA, Flocculation, Hexokinase		Daily
Z546	ANTENATAL PANEL 1	*Complete Blood Count *Blood Group + Rh Factor *Urine R/E *HIV-1 & 2 *HBsAg *VDRL*Glucose F / PP / R *Hemoglobin electrophoresis	6 mL (3 mL min.) whole blood in 2 Lavender Top (EDTA) tubes AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 3 mL (1.5 mL min.) serum from 1 SST AND 10 mL (5 mL min.) aliquot of random urine. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence, VCS, EMT, Dipstick, CMIA, Flocculation, Hexokinase, HPLC		Sample by 9 am; Report Next day
H046	ANTI A TITRE		2 mL (1 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	Tube Agglutination	Also see ANTI B TITRE & COOMBS TEST.	Sample Daily by 9 am; Report Next day
H047	ANTI B TITRE		2 mL (1 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	Tube Agglutination	Also see ANTI A TITRE & COOMBSTEST.	Sample Daily by 9 am; Report Next day
B149	ANTI CYCLIC CITRULLINATED PEPTIDE; ANTI CCP		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay		Daily

S096	ANTI ds DNA ANTIBODY, CRITHIDIA IFA	Double Stranded, Native	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	For SERIAL DILUTIONS use testcode S119.	Sample Daily by 11 am; Report Same Day
S029	ANTI - ds DNA ANTIBODY, EIA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay		Sample by Mon through Sat 9 am; Report Same day
B075	ANTI DNase B; ANTIDEOXYRIBONUCLEASE B		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Latex Enhanced Nephelometry		Sample by Mon / Thu 4 pm; Report Same day
S119	ANTI - ds DNA ANTIBODY, CRITHIDIA IFA IN DILUTIONS		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 11 am; Report Nextday
S246	ANTI IA-2; INSULINOMA ASSOCIATED ANTIGEN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Fri 9 am; Report Next day.
B171	ANTI MULLERIAN HORMONE; AMH		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample Daily by 4 pm; Report Same Day
S239	ANTI NMDA RECEPTOR / ANTI-GLUTAMATE ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 9 am; Report 5 days

S044	ANTI NUCLEAR ANTIBODY / FACTOR (ANA / ANF), IFA (HEP-2)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay (Hep-2 Substrate)	For SERIAL DILUTIONS, use test code S109.	Sample Daily by 11 am; Report Same Day
S109	ANTI NUCLEAR ANTIBODY / FACTOR (ANA / ANF), IFA (HEP-2) IN DILUTIONS		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay (Hep-2 Substrate)		Sample Daily by 11 am; Report Nextday
S027	ANTI NUCLEAR ANTIBODY / FACTOR (ANA/ANF), EIA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay		Sample Daily by 9 am; Report Same day
S235	ANTI NUCLEAR ANTIBODY / FACTOR (ANA / ANF), IFA (HEP-2) REFLEX TESTING		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay (Hep-2 Substrate)		Sample Daily by 11 am; Report Nextday
S238	Anti-NMO (NEUROMYELITIS OPTICA) ANTIBODY /AQUAPORIN 4		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 9 am; Report 5 days

B154	ANTI SOLUBLE LIVER ANTIGEN; SLA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Sat 9 am; Report Same day
S097	ANTI - ss DNA ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S012	ANTI SPERM ANTIBODY, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed / Sat 9 am; Report Same day
B073	ANTI STREPTOLYSIN - O (ASO) TITRE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry	Also see ANTI DNase B	Daily
S095	ANTI THYROGLOBULIN ANTIBODY; ANTI Tg		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Combined assay with Anti TPO Antibody (Test Code Z201) is recommended	Daily
Z201	ANTI THYROID ANTIBODIES PANEL	*Antithyroglobulin Antibody *Antithyroid Peroxidase Antibody	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	This is the recommended test for Thyroid Antibodies.	Daily

S094	ANTI THYROID PEROXIDASE ANTIBODY; ANTI TPO		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Combined assay with Anti Tg Antibody (Test Code Z201) is recommended. Anti TPO replaces the earlier Anti Thyroid Microsomal Antibody (TMA) because of greater specificity for the thyroid.	Daily
RW016	ANTITHROMBIN III, ANTIGENIC		3 mL (1.5 mL min.) Platelet Poor Plasma (PPP) from 2 Blue Top (Sodium Citrate) tubes. Prepare PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory Patient should avoid Anabolic steroids, Gemfibrozil, Heparin therapy, Asparaginase, Estrogens, Gestodene and Oral contraceptives 3 days prior to specimen collection.	Fixed Rate Time Nephelometry		Sample by 7th of the month; report after 2–3 weeks.

H101	ANTITHROMBIN ACTIVITY, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical clot detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Thu 9 am; Report Next day
OS106	APO E GENOTYPING		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Hybridization		Sample by Mon 11 am; Report 10 days
Z212	APOLIPOPROTEI NS A1 & B	Includes Apo B / Apo A1 Ratio	2 mL (1 mL min.) serum from 1 SST. Minimum 12 hour fasting is mandatory . Ship refrigerated or frozen.	Immunoturbidimetry	Also see LIPID PROFILE; LIPOPROTEIN FRACTIONATION, ULTRACENTRIFUGATION; HDL CHOLESTEROL SUBCLASSES; LIPOPROTEIN (a).	Daily

RW017	ARGININE VASOPRESSIN / ANTI DIURETIC HORMONE (ADH), PLASMA		4 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Draw blood in 2 pre-chilled tubes, transfer in an ice bath to refrigerated centrifuge. Separate plasma immediately and freeze within 1 hour. Ship frozen. DO NOT THAW	Extraction, Radioimmunoassay	Patient must be fasting for 6 hours prior to sampling. No liquids including water allowed.	Sample by 7th of the month; Report after 2–3 weeks
G149	ARGINASE		4 mL (3 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Spectrophotometry		Sample by Mon 4 pm; Report Thu
C006	ARSENIC, BLOOD		4 mL (2 mL min.) whole blood in 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Patient should avoid eating sea food and herbal supplements for 3 days prior to sample collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

C004	ARSENIC, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Patient should avoid eating sea food and herbal supplements for 3 days prior to sample collection. Use powderless gloves during specimen collection.	ICPMS	Recommended test for Arsenic estimation. Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C005	ARSENIC, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Patient should avoid eating sea food and herbal supplements for 3 days prior to sample collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

Z261	ARTHRITIS PANEL 2	*CBC *CRP *RA Factor*Anti Nuclear Antibody/ Factor *Uric Acid *ASO*Compleme nt C3 *Complement C4	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE. Overnight fasting is preferred.	Electrical Impedance, VCS, Spectrophotometry, I mmunoturbidimetry, EIA	Recommended as a follow up test if Anti Nuclear Antibody/Factor is positive in Z260 ARTHRITIS PANEL 1. Also see RHEUMATOID tests, LUPUS ERYTHEMATOSUS PANEL, SJOGREN'S SYNDROM E ANTIBODIES, COLLAGEN DISEASES PANEL, PHOSPHOLIPID SYNDROME PANEL and DEEP VEIN Ç	Sample Daily by 9 am; Report Same Day
------	----------------------	---	--	--	--	---

Z260	ARTHRITIS PANEL 1	*CBC *CRP *RA Factor*Anti Nuclear Antibody/ Factor *ASO *Uric Acid	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE. Overnight fasting is preferred.	Electrical Impedence,VCS, Spectrophotometry,I mmunoturbidimetry, EIA	This panel is recommended as initial work up for arthritis patients. Also see RHEUMATOID tests, LUPUS ERYTHEMATOSUS PANEL, SJOGREN'S SYNDROME ANTIBODIES, COLLAGEN DISEASES PANEL, PHOSPHOLIPID SYNDROME PANEL and DEEP VEIN THROMBOSIS PANEL.	Sample Daily by 9 am; Report Same day
Z316	ARTHRITIS PANEL BASIC	*CBC *CRP *RA Factor *ASO *Uric Acid	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE. Overnight fasting is preferred.	Electrical Impedence, VCS, Spectrophotometry,I mmunoturbidimetry		Daily
Z675	ASCA (ANTI SACCHAROMYC ES CEREVISIAE ANTIBODY), IgA & IgG PANEL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day

S219	ASCA (ANTI SACCHAROMYC ES CEREVISIAE ANTIBODY), IgA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S218	ASCA (ANTI SACCHAROMYC ES CEREVISIAE ANTIBODY), IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S102	ASPERGILLUS ANTIBODY, IgG, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see ALLERGY, INDIVIDUAL,ASPERGI LLUS.	Sample by Wed 9 am; Report Same day
S155	ASPERGILLUS ANTIBODY, IgM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see ALLERGY, INDIVIDUAL,ASPERGI LLUS.	Sample by Wed 9 am; Report Same day
Z215	ASPERGILLUS ANTIBODIES PANEL, IgG & IgM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see ALLERGY, INDIVIDUAL,ASPERGI LLUS.	Sample by Wed 9 am; Report Same day
RW259	AUTOIMMUNE CEREBELLAR ATAXIA PANEL		3 mL(1.5 mL min.) serum from 1 Red Top (No Additive) Tube.Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Immunoassay, Immunofluorescence Assay (IFA), Immunoturbidimetry , Radiobinding Assay		Sample by 7th of the month; Reportafter 2–3 weeks
Z619	AUTOIMMUNE HEPATITIS PANEL	*ANA *AMA *ASMA *LKM Antibody	3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample Daily by 11 am; Report Same Day

RW019	BABESIA ANTIBODIES, IgG & IgM		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Immunofluorescence Assay	Individuals infected with malaria may falsely test positive	Sample by 7th of the month; Report after 2–3 weeks
OS107	BACTERIAL MENINGITIS CONFIRMATORY PANEL	*H.influenzae *N.meningitidis *S.pneumoniae	4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube OR 1 mL CSF OR Throat swab. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed / Sat 11 am; Report Sat / Wed
OS108	BACTERIAL PNEUMONIA PANEL , PCR	*Mycoplasma pneumoniae *S. pneumoniae *Chlamydia genus	Submit Throat Swab / Nasopharyngeal swab / Respiratory secretions / BAL fluid in a sterile screw capped container Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily ; Report 7 days
S103	BACTERIAL MENINGITIS SCREEN	*Group B Streptococcus *H.influenzae Type B *Streptococcus pneumoniae *Neisseria meningitidis Group A, C, Y, W135 & B *E.coli K1	2 mL (1 mL min.) serum from 1 SST OR 2 mL (1 mL min.) Cerebrospinal fluid OR 5 mL (3 mL min.) random urine in a sterile screw capped container. Ship refrigerated or frozen.	Latex Agglutination		Daily
K014	BARBITURATE SCREEN, QUALITATIVE, URINE	Amobarbital, Butabarbital, Butalbital, Pentobarbital, Phenobarbital and Secobarbital	10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily

OS109	B-CELL REARRANGEMENT DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Sat 11 am; Report 10 days
N035	BCR-ABL PCR, QUANTITATIVE		5 mL (3 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Wed / Fri 9 am; Report Tue / Thu / Sat
N037	BCR-ABL PCR, QUALITATIVE		5 mL (3 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Wed / Fri 9 am; Report Tue / Thu / Sat
U103	BENCE JONES PROTEIN, URINE, QUALITATIVE		15 mL (10 mL min.) Urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Chemical		Daily
K002	BENZODIAZEPINE SCREEN, QUALITATIVE, URINE	Chlordiazepoxide , Clonazepam, Demoxepam, Desalkylflurazepam Diazepam, Flunitrazepam, Flurazepam, Lorazepam, Nitrazepam, Nordiazepam, Oxazepam & Triazolam	10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily
R168	BETA CROSSLAPS; BETA CTx		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Overnight fasting is mandatory. Ship refrigerated or frozen.	ECLIA		Sample by Mon 4 pm; Report Same day

B178	BETA 2 GLYCOPROTEIN 1, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Thu 9 am; Report Same day
Z547	BETA 2 GLYCOPROTEIN 1 PANEL, IgG, IgM, IgA		:2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Thu 9 am; Report Same day
B179	BETA 2 GLYCOPROTEIN 1, IgA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Thu 9 am; Report Same day
B177	BETA 2 GLYCOPROTEIN 1, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Thu 9 am; Report Same day
RW020	6 BETA HYDROXYCORTI SOL, 24-HOUR URINE	Includes creatinine	10 mL (5 mL min.) aliquot of 24-hour urine collected with 10g boric acid or keep urine refrigerated if collected without preservative. Record 24-hour volume on urine container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
U051	BETA - 2 - MICROGLOBULI N, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen. Record 24 hour volume on Test Request Form and urine container.	CLIA		Daily
R061	BETA - 2 - MICROGLOBULI N, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Avoid hemolysis. Overnight fasting is preferred.	CLIA		Daily

OS111	BETA THALASSEMIA MUTATION DETECTION , PRENATAL	*5 common mutations	4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube from both the Parents AND 10 mL amniotic fluid / CVS. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Tue 11am; Report 10 days
RW245	BETA -2- TRANSFERRIN, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen.	Electrophoresis, Immunofixation		Sample by 7th of the month; Report after 2–3 weeks
B162	BICARBONATE, SERUM		2 mL (1 mL min) serum from 1 SST Ship refrigerated or frozen.	Ion Selective Electrode	Also see BLOOD GAS ANALYSIS	Daily
B197	BILE ACIDS, TOTAL, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzymatic Colorimetric		Sample by Mon/Wed/Fri 9 am; Report same day
Z106	BILIRUBIN TOTAL, DIRECT & INDIRECT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	DPD, Calculated		Daily
P099	BILIRUBIN, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen.	DPD, Calculated		Daily
B035	BILIRUBIN, TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	DPD		Daily
B036	BILIRUBIN, DIRECT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	DPD		Daily
G129	BIOTINIDASE ACTIVITY, QUANTITATIVE, BLOOD		2 mL (1 mL min.) serum from 1 SST OR 2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Ship refrigerated or frozen.	Enzyme Assay		Sample Daily by 4 pm; Report Same Day
R165	BIOTINIDASE, NEWBORN SCREEN		1 drop of heel prick blood each on 3 spots of filter paper available from LPL. Ship refrigerated or frozen.	Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm; Report 2 days

RW021	BISMUTH, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid washed (metal free) container kit available from LPL. Ship refrigerated or frozen. To avoid potential contamination, do not measure 24-hour volume. Avoid mineral supplements and bismuth preparations like Pepto-Bismal for at least 1 week prior to specimen collection.	Inductively-Coupled Plasma/Mass Spectrometry		Sample by 7th of the month; Report after 2–3 weeks
RW022	BISMUTH, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Avoid mineral supplements and bismuth preparations like Pepto-Bismal for at least 1 week prior to specimen collection.	Inductively-Coupled Plasma/Mass Spectrometry		Sample by 7th of the month; Report after 2–3 weeks
N069	BK VIRUS, PCR, QUALITATIVE		3 mL (2 mL min.) Plasma from 1 Lavender top (EDTA) tube OR 15 mL (5 mL min.) Urine in a sterile screw capped container. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N096	BK VIRUS, PCR, QUANTITATIVE		3 mL (2 mL min.) plasma from 1 Lavender Top (EDTA) tube OR 15 mL (5 mL min.) random urine in a sterile screw capped container Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
RW023	BLASTOMYCES ANTIBODY		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunodiffusion		Sample by 7th of the month; Report after 2–3 weeks
H015	BLEEDING TIME; BT		Conducted on patient.	Ivy's method	Test conducted only in LPL Laboratories	Daily

P001	BLOOD GAS ANALYSIS, ARTERIAL		2 mL arterial blood in a preheparinised (1000 IU/mL) disposable syringe and needle. Avoid air bubbles. Mix gently 10 times. Seal needle with rubber stopper. Sample to be drawn by referring doctor. Ship refrigerated. DO NOT FREEZE. Outstation samples not accepted.	ISE, Amperometry, Reflectance photometry		Daily
P002	BLOOD GAS ANALYSIS, VENOUS		2 mL venous blood in a pre-heparinised (1000 IU/mL) disposable syringe and needle. Avoid air bubbles. Mix gently 10 times. Seal needle with rubber stopper. Ship refrigerated. DO NOT FREEZE. Outstation samples not accepted.	ISE, Amperometry, Reflectance photometry		Daily
H029	BLOOD PICTURE; PERIPHERAL BLOOD SMEAR EXAMINATION		2 Peripheral blood smears. Give brief clinical history. Ship in a slide mailer.	Microscopy		Daily
H206	BLOOD GROUP ABO & Rh FACTOR	*Forward & Reverse grouping	4 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Erythrocyte Magnetized Technology / Tube Agglutination		Daily
B147	BNP; B - TYPE NATRIURETIC PEPTIDE		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily

Z169	BOH (BAD OBSTETRIC HISTORY) PANEL	*Cardiolipin Antibodies Panel *Lupus Anticoagulant Panel	3 mL (1.5 mL min.) serum from 1 SST AND 4 mL (3 mL min.) Platelet Poor Plasma (PPP) from 3 Blue Top (Sodium Citrate) tubes. Prepare PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	EIA, Electromechanical Clot Detection	Also see DEEP VEIN THROMBOSIS PANEL.	Sample Daily by 3 pm; Report Next day
Z476	BONE MARROW EXAMINATION PANEL 1	*Microscopy of Bone Marrow Aspirate & Peripheral Blood *Special Stains *CBC	2 mL (1 mL min.) Bone Marrow Aspirate in 1 Lavender Top (EDTA) tube AND 4 Bone Marrow Aspirate Smears AND 4 mL (2 mL min) whole blood in 1 Lavender Top (EDTA) tube AND 2 peripheral blood smears. Air dry & fix all smears in 100% Methanol for 30 minutes. Ship refrigerated. DO NOT FREEZE. Give clinical history on Bone Marrow Examination Test request form	Microscopy of Leishman / Jenner-Giemsa stained smears, Electrical Impedence, VCS	Special stain includes Prussian Blue Reaction for Iron.	Sample Daily by 1 pm; Report Next day
Z477	BONE MARROW EXAMINATION PANEL 2	*Microscopy of Bone Marrow Aspirate & Peripheral Blood *Special Stains *CBC * Histopathology	2 mL (1 mL min.) Bone Marrow Aspirate in 1 Lavender Top (EDTA) tube AND 4 Bone Marrow Aspirate Smears AND 4 mL (2 mL min) whole blood in 1 Lavender Top (EDTA) tube AND 2 peripheral blood smears. Air dry & fix all smears in 100% Methanol for 30 minutes AND Bone Marrow Trepine biopsy in 10% formal saline. Ship refrigerated. DO NOT FREEZE. Give clinical history on Bone Marrow Examination Test request form	Microscopy of Leishman / Jenner-Giemsa stained smears, Electrical Impedence, VCS, HPE	Special stains include Prussian Blue Reaction for Iron & Reticulin stain.	Sample Daily by 6 pm; Report 5 days

Z499	BONE MARROW, TREPHINE BIOPSY PANEL	*Histopathology *Reticulin stain *Iron stain	Submit specimen in 10% Formal-saline. Ship at room temperature. Give clinical history on Bone Marrow Examination Test request form.	Processing, Staining, Microscopy		Sample Daily by 6 pm; Report 5 days
H227	BONE MARROW, IRON STAIN		Minimum 4 smears of bone marrow aspirate AND 2 peripheral blood smears. Fix smears in 95% ethanol for 30 min and air dry. Give clinical history on test request form. Ship in a slide mailer.	Prussian Blue Reaction		Sample Daily by 12 noon; Report Next day
N113	BRAF MUTATION ANALYSIS		Submit formalin fixed paraffin embedded tissue block. Ship at room temperature. Block should have at least 10% tumor tissue.	PCR, Sequencing		Sample by Tue 11 AM, Report Mon
N081	BRCA1 & BRCA2 MUTATION COMPREHENSIVE PANEL	*All BRCA1 & BRCA2 gene mutations mapped	Submit 3 mL (2 mL min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Consent & Clinical Information forms available from LPL are mandatory.	Sequencing	Genetic counselling will be provided with prior appointment.	Sample by 1st / 16th of every month; Report 15th / 30th of same month
N080	BRCA1 & BRCA2 MUTATION SCREEN	*BRCA1 exon2 (187delAG) * BRCA1 exon2(5385insc) * BRCA2 exon11 (6174delT) * 5 Experimental markers	Submit 3 mL (2 mL min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Consent & Clinical Information forms available from LPL are mandatory.	SNP Genotyping using Microarray	Genetic counselling will be provided with prior appointment.	Sample Daily by 9 am; Report 8 weeks
Z240	BREAST CANCER PANEL 1	*ER *PR *Includes pathologist review for presence of malignant cells.	Tumour biopsy in 10% Formal-Saline OR Formalin fixed tumour paraffin block OR 5 unstained sections / alcohol fixed cytology smears on Poly-L-Lysine coated slides. Ship at room temperature.	Immunohistochemistry	Use test code J113 for Photomicrograph. Also see CA 15.3.	Sample Daily by 6 pm Report 5 days: Blocks / Sections / Smears 7 days: Tissue

Z242	BREAST CANCER PANEL 2	*ER *PR *Cathepsin D* Includes pathologist review for presence of malignant cells.	Tumour biopsy in 10% Formal-Saline OR Formalin fixed tumour paraffin block OR 6 unstained sections / alcohol fixed cytology smears on Poly-L-Lysine coated slides. Ship at room temperature.	Immunohistochemist ry	Use test code J113 forPhotomicrograph.	Sample Daily by 6 pmReport5 days: Blocks / Sections / Smears7 days: Tissue
Z293	BREAST CANCER PANEL 3	*ER *PR *HER-2 * Includes Pathologist review for presence of malignant cells.	Tumour biopsy in 10% Formal-Saline OR Formalin fixed tumour paraffin block OR 6 unstained sections / alcohol fixed cytology smears on Poly-L-Lysine coated slides. Ship at room temperature.	Immunohistochemist ry	Use test code J113 forPhotomicrograph. Also see CA 15.3.	Sample Daily by 6 pmReport 5 days: Blocks / Sections / Smears7 days: Tissue
Z388	BREAST MONITOR PANEL	*CA 15.3 *CEA	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
Z387	BREAST METASTATIC MARKER	*ER *PR *CA 15.3 *Includes pathologist review for presence of malignant cells.	Tumour biopsy in 10 % Formal-Saline OR Formalin fixed tumour paraffin block OR 5 unstained sections / alcohol fixed cytology smears on Poly-L-Lysine coated slides. Ship at room temperature AND 2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunohistochemist ry, Chemiluminescent Immunoassay	Use test code J113 forPhotomicrograph.	Sample Daily by 6 pmReport 5 days: Blocks / Sections / Smears7 days: TissueSame Day : CA 15.3
S222	BRUCELLA ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue / Fri 9 am; ReportSame day

S003	BRUCELLA AGGLUTINATION TEST	Includes antibody to B.abortus and B.melitensis	2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Tube Agglutination	Also see FEBRILE AGGLUTININ PANEL, WEIL FELIX TEST and WIDAL TEST.	Sample Daily by 6 pm; Report Next day
Z676	BRUCELLA ANTIBODIES, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue / Fri 9 am; Report Same day
S221	BRUCELLA ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue / Fri 9 am; Report Same day
J003	BUCCAL SMEAR FOR BARR BODIES (SEX CHROMATIN)		Rinse mouth with water. Prepare 4 slides of scrapings from the buccal mucosa, 2 from left and 2 from right cheek. Fix immediately in 95% ethanol for 30 min. Air dry and ship at room temperature. Provide clinical details on test request form.	Aceto-orcein stain	Recommended test is CHROMOSOME ANALYSIS, BLOOD (KARYOTYPE)	Sample Daily by 6 pm; Report Next day
B022	BUN; BLOOD UREA NITROGEN		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Urease UV		Daily
RW025	BUPROPION		3 mL (1 mL min.) serum in 1 Red Top (No additive) tube. Separate serum from cells immediately. Ship frozen. DO NOT THAW. Ideal sampling time is trough level prior to next dose.	HPLC-MS/MS		Sample by 7th of the month; Report after 2-3 weeks
RW035	CA 27.29; BREAST CANCER MARKER		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoassay	Also see CA 15.3, BREAST CANCER PANELS.	Sample by 7th of the month; Report after 2-3 weeks

F010	CA 15.3; BREAST CANCER MARKER		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide brief clinical history.	Chemiluminescent MicroparticleImmunoassay	Also see CA 27.29, BREAST CANCERPANELS	Daily
Z306	CA 27.29 & CA 15.3 BREAST CANCER MARKERS PANEL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CMIA, Immunoassay		LPL: DailyUSA: Sample by 7th of the month;report after 2–3 weeks
F015	(CA) CYFRA 21- 1; LUNG CANCER MARKER		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide brief clinical history.	Electrochemiluminescence	Also see CEA & NEURONE SPECIFICENOLASE (NSE).	Sample by Mon 4 pm; Report Same Day
C008	CADMIUM, 24- HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid-washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Patient should avoid eating seafood at least 3 days prior to specimen collection. Use powderless gloves during specimen collection. Avoid catheterized urine specimen as rubber contains trace amounts of cadmium.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

C009	CADMIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid-washed (metal free) container available from LPL. Ship refrigerated or frozen. Patient should avoid eating seafood at least 3 days prior to specimen collection. Use powderless gloves during specimen collection. Avoid catheterized urine specimen as rubber contains trace amounts of cadmium.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C007	CADMIUM, BLOOD		4 mL (2 mL min.) whole blood from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Patient should avoid eating seafood at least 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
F009	CA 72.4; GASTRIC CANCER MARKER		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide brief clinical history.	Electrochemiluminescence		Sample by Mon / Thu 4 pm; Report Same day
U011	CALCIUM, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine. Collect urine with 25 mL 50% concentrated HCl to maintain pH below 3. Record 24-hour urine volume & patient body weight on test request form and urine vial. Ship refrigerated or frozen.	Arsenazo III		Daily

U096	CALCIUM / CREATININE RATIO, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine in 2 mL of 6N HCl to maintain pH below 3. Ship refrigerated or frozen.	Spectrophotometry		Daily
B055	CALCIUM, IONIZED		3 mL (1 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube OR 2 mL (0.5 mL min.) serum from 1 SST. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship refrigerated. DO NOT FREEZE.	Ion Selective Electrode		Daily
U073	CALCIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine in 2 mL of 6N HCl to maintain pH below 3. Ship refrigerated or frozen.	Arsenazo III		Daily
B042	CALCIUM, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship refrigerated or frozen.	Arsenazo III		Daily
R062	CALCITONIN		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Ship frozen. DO NOT THAW. Overnight fasting is preferred.	Chemiluminescent Immunoassay	Also see CALCITONIN STIMULATIONBY CALCIUM PENTAGASTRIN	Sample by Mon / Thu 4 pm; ReportSame day
Z046	CALCITONIN STIMULATION BY CALCIUM - PENTAGASTRIN		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube for each timed specimen. Ship frozen. DO NOT THAW. All containers must be clearly marked with time drawn and submitted with 1 test request form. Draw overnight fasting baseline specimen. With patient supine, administer Calcium Gluconate IV 2 mg/ kg body weight over 1 min, followed by Pentagastrin IV bolus 0.5 ug/kg body weight over 5 sec. Draw additional specimens at 1, 2, 5 and 10 min.	Chemiluminescent Immunoassay	Refer to Y CodesY018 BasalY121 1 minY019 2 min.Y020 5 min.Y021 10 min.Y022 Customized	Sample by Mon / Thu 4 pm; ReportSame day

U097	CALCIUM EXCRETION 2 HOUR FASTING URINE		10 mL (5 mL min.) aliquot of 2 hour urine collection in a sterile screw capped container. Follow detailed instructions as below:1) No food or drink after 6 pm a day prior to the test.2) On the day of the test do not void urine after 8 am.3) At 10 am void urine and discard.4) Drink a glass of water at 10 am and 11 am. Drink only deionized / distilled water.5) Collect urine voided at 12 noon in the container with 6N HCl as the preservative.Ship refrigerated or frozen. Record 2 hour volume on the urine vial and test request form.	Arsenazo III		Daily
Z117	CALCIUM, TOTAL & IONIZED		2 mL (0.5 mL min.) serum from 1 SST OR 3 mL (1 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship refrigerated. DO NOT FREEZE.	ISE, Arsenazo III		Daily
N112	CALR MUTATION DETECTION		3 ml (2 ml min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	PCR, Sequencing		Sample by Mon 11 am; Report Sat
M215	CAMPYLOBACT ER ANTIGEN, STOOL		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily 1 pmReport Same day
G172	CANAVAN DISEASE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen. Provide brief clinical history	Gas Chromatography / MassSpectrophotom etry		Sample Daily by 9 am: Report 5 days
N097	CANCER 50 GENE PANEL		Submit formalin fixed paraffin embedded tissue block. Ship at room temperature. Block should have at least 10% tumor tissue.	Next Generation Sequencing		Sample by 1st / 15th of every month; Report on 15th / 30th of every month

OS114	CANDIDA ALBICANS DETECTION, PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube OR 1 mL CSF/Sputum/Fluids OR Tissue in normal saline. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Thu 11 am; Report 9 days
K003	CANNABINOIDS (Delta-9-THC Metabolites) / MARIJUANA, SCREEN, URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily
F007	CA 125; OVARIAN CANCER MARKER		2 mL (0.5 mL min) serum from 1 SST. Ship refrigerated or frozen. Provide brief clinical history.	Chemiluminescent MicroparticleImmunoassay	Also see ROMA (Risk of OvarianMalignancy)	Daily
F008	CA 19.9; PANCREATIC CANCER MARKER		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide brief clinical history.	Chemiluminescent MicroparticleImmunoassay		Daily
D051	CARBAMAZEPINE		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Chemiluminescent Immunoassay		Daily
RW036	CARBAMAZEPINE, FREE		4 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
Z003	CARDIAC / MYOCARDIAL INFARCTION PANEL	*CPK-MB Mass * Troponin I * Troponin T	2 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	CLIA, ECLIA		Daily

S053	CARDIOLIPIN ANTIBODY, IgA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see PHOSPHOLIPID SYNDROME PROFILE & DEEP VEIN THROMBOSIS PANEL.	Sample by Mon through Sat 9 am Report Same day
G181	CARDIOMETABOLIC PANEL		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Ship refrigerated or frozen. Overnight fasting is mandatory.	LC-MS/MS		Sample Daily by 5 pm; Report 5 days
S054	CARDIOLIPIN ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see PHOSPHOLIPID SYNDROME PROFILE .	Sample by Mon through Sat 9 am Report Same day
Z102	CARDIOLIPIN ANTIBODIES PANEL, IgG, IgA & IgM		4 mL (3 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see PHOSPHATIDYLSERINE ANTIBODIES PANEL, PHOSPHOLIPID PANEL, DEEP VEIN THROMBOSIS PANEL & LUPUS ANTICOAGULANT PANEL.	Sample by Mon through Sat 9 am Report Same day
S055	CARDIOLIPIN ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see PHOSPHOLIPID SYNDROME PROFILE & DEEP VEIN THROMBOSIS PANEL.	Sample by Mon through Sat 9 am Report Same day
G170	CARNITINE, SERUM	*Carnitine *Acylcarnitine	2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 9 am; Report 2 days

RW038	CAROTENE		3 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Transfer into sterile screw capped vial. Wrap in aluminium foil to protect from light. Ship refrigerated or frozen. Overnight fasting for 12-14 hours is mandatory. Patient must not consume alcohol or Vitamin A supplement 24 hours prior to sampling.	HPLC	Also see VITAMIN A.	Sample by 7th of the month; Report after 2–3 weeks
U024	CATECHOLAMINES, 24-HOUR URINE	Includes Epinephrine, Norepinephrine and Dopamine.	50 mL (10 mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl, to maintain pH between 1-2. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen. Patient should strictly avoid Alpha-one blockers, Aminophylline, Amphetamines, Ampicillin, Beta-blockers, Ephedrine, Imipramine, Methyldopa, MAO inhibitors, Nicotine, Phenacetin, Phenothiazine, Theophylline, Vasodilators, Aspirin, PAS, Alpha-2 agonists, Bromocriptine, Calcium channel blockers (long term use), Clofibrate, Propranolol, Penicillin, Reserpine, Thyroxine, Alpha-methyldopa, Isoproterenol, Labetalol, Mandelamine, Paracetamol, Cimetidine, Metoclopramide, vigorous exercise, alcoholic beverages, tea/coffee and tobacco at least 72 hours before and during specimen collection.	High Performance Liquid Chromatography	Also see PHEOCHROMOCYTOMAPROFILE	Sample by Wed / Sat 9 am; Report 2 days

Z118	CATECHOLAMINES & VMA, 24-HOUR URINE		50 mL (10 mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl to maintain pH between 1-2. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen. Patient should strictly avoid Alpha-one blockers, Aminophylline, Amphetamines, Ampicillin, Beta-blockers, Ephedrine, Imipramine, Methyldopa, MAO inhibitors, Nicotine, Phenacetin, Phenothiazine, Theophylline, Vasodilators, Aspirin, PAS, Alpha-2 agonists, Bromocriptine, Calcium channel blockers (long term use), Clofibrate, Propranolol, Penicillin, Reserpine, Thyroxine, Alpha-methyldopa, Isoproterenol, Labetalol, Mandelamine, Paracetamol, Cimetidine, Metoclopramide, vigorous exercise, alcoholic beverages, tea/coffee, tobacco, chocolate, vanilla & banana at least 72 hours before and during specimen collection.	High Performance Liquid Chromatography	Also see PHEOCHROMOCYTOMAPROFILE.	Sample by Wed / Sat 9 am; Report 2 days
------	-------------------------------------	--	--	--	-----------------------------------	---

R097	CATECHOLAMINES, PLASMA	Includes Epinephrine, Norepinephrine and Dopamine.	8 mL (4 mL min.) plasma from 4 special Green Top (Sodium Heparin tubes containing 5 mg Sodium metabisulphite) available from LPL on request. Adhere to REQUESTED specimen volume. Draw specimen as follows: Patient should be at rest for 30 minutes after insertion of an IV cannula. Collect specimen with patient supine. Immerse specimen immediately in ice water. Separate plasma within 10 minutes, preferably in cold centrifuge, transfer into sterile screw capped vial and FREEZE immediately. Ship frozen. DO NOT THAW. Patient should strictly avoid drugs like Aldomet (Alpha methyldopa), Labetalol, Isoproterenol, Isoetharine, Alpha-one blockers, Aminophylline, Amphetamines, Ampicillin, Beta-blockers, Ephedrine, Imipramine, Nicotine, Phenacetin, Theophylline, Phenothiazine, Vasodilators (Minoxidil, Nitrates & Hydralazine), Aspirin, PAS, Alpha-2 agonists, Bromocriptine, Calcium channel blockers (long term use), L-dopa, Cimetidine, clofibrate, MAO inhibitors, Propranolol, Penicillin, Reserpine, Thyroxine, vigorous exercise, severe stress and cigarette smoking at least 72 hrs before and during specimen collection.	High Performance LiquidChromatography	The recommended test is CATECHOLAMINES, 24-HOUR URINE. Also seePHEOCHROMOCYTOMA PROFILE.	Sample by Wed 9 am; Report Fri
------	------------------------	--	--	---------------------------------------	--	--------------------------------

U043	CATECHOLAMINES, RANDOM URINE	Includes Epinephrine, Norepinephrine and Dopamine.	50 mL (10 mL min.) aliquot of random urine. Collect urine with appropriate volume of 50% HCl , to maintain pH between 1–2. Do not use concentrated HCl. Ship refrigerated or frozen. Patient should strictly avoid Alpha-one blockers, Aminophylline, Amphetamines, Ampicillin, Beta-blockers, Ephedrine, Imipramine, Methyldopa, MAO inhibitors, Nicotine, Phenacetin, Phenothiazine, Theophylline, Vasodilators, Aspirin, PAS, Alpha-2 agonists, Bromocriptine, Calcium channel blockers (long term use), Clofibrate, Propranolol, Penicillin, Reserpine, Thyroxine, Alpha-methyldopa, Isoproterenol, Labetalol, Mandelamine, Paracetamol, Cimetidine, Metoclopramide, vigorous exercise, alcoholic beverages, tea/coffee and tobacco at least 72 hours before and during specimen collection.	High Performance LiquidChromatography	Also see PHEOCHROMOCYTOMAPROFILE.	Sample by Wed / Sat 9 am; Report 2 days
Z108	C3 & C4 COMPLEMENT PANEL		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunoturbidimetry	Also see HISTOPATHOLOGY, DIRECTIMMUNOFLUORESCENCE TESTS.	Daily
RW031	C6 COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radial Immunodiffusion		Sample by 7th of the month; Reportafter 2–3 weeks
RW034	C9 COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radial Immunodiffusion		Sample by 7th of the month; Reportafter 2–3 weeks

H105	C3 COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunoturbidimetry	Also see HISTOPATHOLOGY, DIRECTIMMUNOFLU ORESCENCE TESTS.	Daily
RW033	C8 COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radial Immunodiffusion		Sample by 7th of the month; Report after 2–3 weeks
H106	C4 COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunoturbidimetry	Also see HISTOPATHOLOGY, DIRECTIMMUNOFLU ORESCENCE TESTS.	Daily
RW032	C7 COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radial Immunodiffusion		Sample by 7th of the month; Report after 2–3 weeks
RW029	C2 COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen.	Radial Immunodiffusion		Sample by 7th of the month; Report after 2–3 weeks
RW030	C5 COMPLEMENT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radial Immunodiffusion		Sample by 7th of the month; Report after 2–3 weeks
T006	C. DIFFICILE TOXIN A & B, STOOL		2 g (0.5 g min.) fresh stool in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. For patients using diapers, first line the diaper with clean plastic to prevent absorption.	Immunochromatogr aphy	Also see CULTURE, CLOSTRIDIUMDIFFICI LE.	Daily

F003	CEA; CARCINO EMBRYONIC ANTIGEN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Give brief clinical history.	Chemiluminescent MicroparticleImmunoassay		Daily
OS115	CEBPA MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed 11 am; Report 8 days
S079	CENTROMERE ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred	Enzyme Immunoassay	Also see EXTRACTABLE NUCLEARANTIGENS.	Sample by Tue / Fri 9 am; Report Same day
B087	CERULOPLASMIN		2 mL (0.5 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	Nephelometry	Also see COPPER.	Sample Mon through Sat by 4 pm;Report Same day Sample Sun by 9 am: Report Same day
Z770	CERVICAL SCREEN	*Cervical PAP Smear, Genital, Female (Liquid based) * Human Papilloma Virus (HPV) Detection & Genotyping, PCR, Qualitative	Collect Endocervical smears with a brush after inserting the speculum. Dip & swirl the brush completely into the LBC vial with liquid fixative. Mix thoroughly. Use Special Collection Kit available from LPL only. Ship refrigerated. DO NOT FREEZE. Brief clinical history in Cytology Requisition Form is mandatory.	Microscopy of Papanicolaou stained smears, RT-PCR	Do not take smears during menstruation.	Sample by Wed / Sat 11 am; Report Fri/ Tue

H102	C1 ESTERASE INHIBITOR, PROTEIN QUANTITATION; C1 INHIBITOR, QUANTITATIVE		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use gel barrier tubes. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Overnight fasting is preferred.	Radial Immunodiffusion		Sample by Wed 9 am; Report next Mon
RW027	C1 ESTERASE INHIBITOR, FUNCTIONAL ACTIVITY; C1 INACTIVATOR, FUNCTIONAL		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use gel barrier tubes. Separate serum within 1 hour of collection and freeze. Ship frozen. DO NOT THAW.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
N131	CHIKUNGUNYA VIRUS, PCR QUALITATIVE		3 mL (2 mL min.) Whole Blood / Plasma in 1 Lavender Top (EDTA) Tube OR 2 mL (1 mL min.) Serum from 1 SST OR 2 mL (1 mL min.) CSF in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE	Real Time PCR		Sample Daily by 11 am; Report Next day
S209	CHIKUNGUNYA VIRUS ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunochromatography		Daily
N110	CHIMERISM SPLIT CELL ANALYSIS		3 ml (2 ml min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR, Fragment Analysis		Sample Daily by 11 am; Report 7 days
N058	CHIMERISM, POST-ENGRAFTMENT		4 ml (2 ml min) whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR, Fragment analysis	Indicate date of transplant on Test request form. Whole blood is the preferred sample.	Sample Daily by 11 am; Report 7 days

N059	CHIMERISM, PRE- ENGRAFTMENT		4 ml (2 ml min) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	PCR, Fragment analysis	Whole blood is the preferred sample. Buccal swab/Nail clippings also acceptable.	Sample Daily by 11 am; Report 7 days
N057	CHIMERISM, PRE- ENGRAFTMENT, DONOR & RECIPIENT		4 ml (2 ml min) whole blood in 1 Lavender Top (EDTA) tube each for Donor and Recipient. Ship refrigerated. DO NOT FREEZE. Indicate Donor and Recipient on tubes appropriately.	PCR, Fragment analysis	Whole blood is the preferred sample. Buccal swab/Nail clippings also acceptable.	Sample Daily by 11 am; Report 7 days
S237	CHLAMYDOPHIL A PNEUMONIAE IgG		3 mL (1.5 mL min.) serum from 1 SST.Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Fri 9 am; Report Same day
RW039	CHLAMYDIA / CHLAMYDOPHIL A SEROLOGY	*C.pneumoniae IgA, IgG & IgM *C.psittaci IgA,IgG & IgM *C.trachomatis IgA, IgG & IgM	3 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks
Z779	CHLAMYDIA TRACHOMATIS AND NEISSERIA GONORRHOEA , PCR QUALITATIVE		15 ml (8 mL min.) first morning urine in a sterile screw capped container OR Vaginal/Endocervical /Conjunctival swab in normal saline. Ship refrigerated or frozen	Real Time PCR		Sample Daily by 11 am; Report Next day

S156	CHLAMYDIA IgG ANTIBODY PANEL	*Chlamydia trachomatis, IgG *Chlamydophila pneumoniae, IgG	3 mL (1.5 mL min.) serum from 1 SST.Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Fri 9 am; Report Same day
N056	CHLAMYDIA TRACHOMATIS, PCR		10 mL (5 mL min.) urine in a sterile screw capped container OR Genital / Conjunctival swabs in normal saline. Ship refrigerated or frozen.	Polymerase Chain Reaction		Sample Daily by 11 am; Report Next day
S236	CHLAMYDIA TRACHOMATIS IgG		3 mL (1.5 mL min.) serum from 1 SST.Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Fri 9 am; Report Same day
U044	CHLORIDE, 24- HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine. Refrigerate during collection. Record 24-hour volume on test request form and urine container. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
U015	CHLORIDE, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
P088	CHLORIDE, FLUID		Submit 5 mL (2 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
B046	CHLORIDE, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
P054	CHOLESTEROL, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen. Specify type of fluid on test request form and sample container.	CHO-POD		Daily
B025	CHOLESTEROL, TOTAL		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CHO-POD	Also see LIPID PROFILE TESTS.	Daily
B065	CHOLINESTERAS E, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Spectrophotometry		Daily

R098	CHROMOGRANIN A; CGA		2 mL (0.5 mL min.) plasma from 1 Lavender Top (EDTA) tube. Ship frozen. DO NOT THAW	Enzyme Immunoassay		Sample by Tue / Fri 9 am; Report Same day
Q015	CHROMOSOME ANALYSIS, PRODUCTS OF CONCEPTION REFLEX TESTING TO FISH FOR ANEUPLOIDY DETECTION		Submit curetted tissue in normal saline. (Fetus not accepted). Ship immediately at 18–22°C. Specimens must reach LPL, New Delhi within 24 hours. Formalin fixed specimen not acceptable. Give clinical history on Chromosome Analysis Test request form	Culture, Microscopy, Karyotype, FISH	If culture is unsuccessful, FISH will be performed for Aneuploidy detection.	Sample Daily by 6 pm; Report 4 weeks
N091	CHROMOSOME ANALYSIS FOR AUTISM, MICROARRAY 180K		4 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Analysis Test request form.	Microarray 180K with focus on Autism specific genes		Sample Daily by 11 am; Report 25 days
Q018	CHROMOSOME ANALYSIS, CHORIONIC VILLUS		40 mg (20 mg min.) Chorionic villus biopsy collected aseptically in 10 mL transport medium available from LPL. Sample to be collected under strict aseptic precautions. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL New Delhi within 48 hours. Test is recommended between 10-13 weeks of gestation. Consent form for Prenatal Genetic Testing is mandatory.	Culture, Microscopy, Karyotype		Sample Daily by 6 pm; Report 21 days

Q019	CHROMOSOME ANALYSIS, CORD BLOOD		4 mL (2 mL min.) Cord blood in 1 Green Top (Sodium Heparin) tube. Avoid clot formation during sampling. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL, New Delhi within 24 hrs. Test is recommended between 18-20 weeks of gestation. Consent form for Prenatal Genetic Testing is mandatory. Provide gestational age of patient.	Culture, Microscopy, Karyotype		Sample Daily by 9 am; Report 10 days
C012	CHROMIUM, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
Q014	CHROMOSOME ANALYSIS, SOLID TUMOR		Submit tumor tissue in a sterile screw capped container in Hank's / Ringer's solution. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL, New Delhi within 24 hrs. Formalin fixed specimen not acceptable. Give clinical history on Chromosome Analysis Test request form.	Culture, Microscopy, Karyotype	Trim tissue free of fat, normal and necrotic tissue.	Sample Daily by 9 am; Report 21 days
Q020	CHROMOSOME ANALYSIS, HIGH RESOLUTION, NEONATAL		4 ml (2 ml min). Whole blood in 1 Green top (Sodium Heparin) tube. Ship refrigerated immediately. DO NOT FREEZE. Give clinical history on Chromosome Analysis Test Request Form.	Culture, Microscopy, High Resolution Karyotype		Sample Daily by 6 pm; Report 21 days

C013	CHROMIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
Q023	CHROMOFIC CHROMOSOME SNP MICROARRAY 750K, HIGH RESOLUTION		4 mL (2 mL min.) whole blood in 1 lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Genomic Microarray Requisition form.	Affymetrix CytoScan™ 750K Microarray		Sample Daily by 4 pm; Report 10 days
Q025	CHROMOTOUCH CHROMOSOME SNP MICROARRAY OPTIMA, PRODUCTS OF CONCEPTION		Submit 5 mg (2 mg min.) curretted tissue in normal saline. Give clinical history on Genomic Microarray Requisition form.	Affymetrix Optima Suite™ Microarray	Test can also be performed on Formalin-fixed, paraffin-embedded tissue or Products of conception (POC) specimens that fail to grow in culture	Sample Daily by 4 pm; Report 10 days

C052	CHROMIUM, BLOOD		4 mL (2 mL min.) whole blood from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
Q024	CHROMOTOUCH CHROMOSOME SNP MICROARRAY OPTIMA, PRENATAL		Amniotic fluid: 15 mL (10 mL min.) in a sterile screw capped container. Chorionic villus: 30 mg (20 mg min.) biopsy collected aseptically in 10 mL transport medium available from LPL. Umbilical Cord blood: 4 mL (2 mL min.) cord blood in 1 Lavender top (EDTA) tube. Avoid clot formation during sampling. Ship refrigerated immediately. DO NOT FREEZE. Give clinical history on Genomic Microarray Requisition form. Consent form for Prenatal genetic testing is mandatory.	Affymetrix Optima Suite™ Microarray		Sample Daily by 4 pm; Report 10 days
Q005	CHROMOSOME ANALYSIS FOR HEMATOLOGIC MALIGNANCY		2–3 mL (1.5 mL min.) heparinized Bone Marrow OR 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL, New Delhi within 24 hrs. Give clinical history on Chromosome Analysis Test request form.	Culture, Robotic Microscopy, Karyotype	Preferred sample is Bone Marrow .	Sample Daily by 9 am; Report 4 days

Q008	CHROMOSOME ANALYSIS (KARYOTYPE), BLOOD		5 mL (3 mL min.) whole blood in 2 Green Top (Sodium Heparin) tubes. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL, New Delhi within 24 hrs. Give clinical history on Chromosome Analysis Test request form.	Culture, Robotic Microscopy, Karyotype		Sample Daily by 9 am; Report 12 days
N092	CHROMOSOME ANALYSIS FOR UNIPATERNAL DISOMY (UPD), MICROARRAY 180K		4 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Analysis Test request form.	Microarray 180K with focus to detect UPD		Sample Daily by 11 am; Report 25 days
N074	CHROMOSOME ANALYSIS, MICROARRAY 60K		4 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Analysis Test request form.	Microarray		Sample Daily by 11 am; Report 25 days
Q016	CHROMOSOME ANALYSIS, PHILADELPHIA		2–3 ml (1.5 mL min.) heparinized bone marrow in 1 Green Top (Sodium Heparin) tube. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL, New Delhi within 24 hrs. Give clinical history on Chromosome Analysis Test request form.	Culture, Microscopy, Karyotype		Sample Daily by 9 am; Report 7 days

Q017	CHROMOSOME ANALYSIS, AMNIOTIC FLUID		20 mL (15 mL min.) Amniotic fluid in a sterile screw-capped container. Sample to be collected under strict aseptic precautions. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL New Delhi within 48 hours. Test is recommended between 15-20 weeks of gestation. Consent form for Prenatal Genetic Testing is mandatory.	Culture, Microscopy, Karyotype	It is recommended to discard the initial 2 - 3 mL of aspirate, since this "first draw" may contain some maternal cells trapped within the needle during insertion.	Sample Daily by 6 pm; Report 21 days
Q022	CHROMOSOME INTERPHASE PROFILING, AMNIOTIC FLUID		20 mL (15 mL min.) Amniotic fluid in a sterile screw-capped container. Sample to be collected under strict aseptic precautions. Ship refrigerated immediately. DO NOT FREEZE. Specimen must reach LPL New Delhi within 48 hours. Test is recommended between 15-20 weeks of gestation. Consent form for Prenatal Genetic Testing is mandatory.	FISH		Sample Daily by 4 pm ; Report 7 days
Q021	CHROMOSOME INTERPHASE PROFILING, PRODUCTS OF CONCEPTION		Submit 5 mg (2 mg. min) Placental tissue OR Fetal tissue (Facia lata / diaphragm / tendon / skin / tissue from internal organs (if fresh) / chest wall cartilage particularly if macerated) in Normal saline / 10% Formalin.. Ship refrigerated. DO NOT FREEZE.	FISH		Sample Daily by 4 pm; Report 7 days
Z799	CHRONIC FATIGUE SYNDROME	*CBC * C-Reactive Protein * LFT *KFT *CPK * Glucose, fasting * Iron Studies * TSH, Ultrasensitive	5 mL (3 mL min.) serum from 1 SST AND 3 mL (1 mL min.) whole blood from 1 Lavender Top (EDTA) tube AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. 12 hours overnight fasting essential. Ship refrigerated or frozen.	Electrical Impedence, Immunoturbidimetry , Spectrophotometry, IFCC, Hexokinase		Daily

P072	CHYLE EXAMINATION		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Staining, Microscopy		Sample Daily by 6 pm; Report Next day
RW247	CIRCULATING TUMOR CELLS (CTC) FOR COLORECTAL CANCER		15 mL (7.5 mL min.) whole blood in 2 Special Cell save tubes available from LPL. Ship ambient immediately.	Cell search system	Patients on Adriamycin must wait a minimum of 7 days post administration before blood is drawn for this test.	Contact National Reference Lab, NewDelhi
RW248	CIRCULATING TUMOR CELLS (CTC) FOR PROSTATE CANCER		15 mL (7.5 mL min.) whole blood in 2 Special Cell save tubes available from LPL. Ship ambient immediately.	Cell search system	Patients on Adriamycin must wait a minimum of 7 days post administration before blood is drawn for this test.	Contact National Reference Lab, NewDelhi
RW246	CIRCULATING TUMOR CELLS (CTC) FOR BREAST CANCER		15 mL (7.5 mL min.) whole blood in 2 Special Cell save tubes available from LPL. Ship ambient immediately.	Cell search system	Patients on Adriamycin must wait a minimum of 7 days post administration before blood is drawn for this test.	Contact National Reference Lab, NewDelhi
G147	CITRATE, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine. Ship frozen. DO NOT THAW.	Spectrophotometry, Enzymatic	Also see KIDNEY STONEFORMATION tests.	Sample Daily by 4 pm; Report Same Day
OS118	C-KIT MUTATION DETECTION , PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed 11 am; Report 9 days

RW044	CLOMIPRAMINE	Includes metabolite N-Desmethyldomipramine	4 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: >12 hours post oral dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
RW045	CLONAZEPAM		3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: Before the administration of the next dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
RW046	CLORAZEPATE	Measured as metabolite Nordiazepam.	3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Draw sample before next dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
OS119	CLOSTRIDIUM DIFFICILE DETECTION, PCR		5g (3g min.) Stool in a leak proof screw capped container. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed / Sat 11 pm; Report Mon / Tue
H037	CLOT RETRACTION TEST		5 mL whole blood in a special graduated tube available from LPL. Insert a wooden stick immediately and permit blood to clot.	Manual		Daily
H016	CLOTTING TIME; CT		Conducted on patient.	Lee & White	Test conducted only in LPL Laboratories	Daily
RW047	CLOZAPINE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Recommended specimen: 0.5 to 1 hour before next oral dose at steady state.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks

Z041	COAGULATION PROFILE 2	*CBC *BT *CT *Prothrombin studies *APTT*Thrombin Time*Factor VIII Functional	3 mL (1.5 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE and 3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electrical Impedence, Photo- opticalClot Detection, Ivy, Lee & White	BT and CT will be performed at LPL Laboratories. It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from tÇ	Sample by Wed / Sat 9 am; Report Fri/ Tue.
Z042	COAGULATION PROFILE 1.	*CBC *BT *CT *Prothrombin studies *APTT*Thrombin TimeDoes not include Factor VIII assay	3 mL (1.5 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE AND 3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electrical Impedence, Photo- opticalClot Detection, Ivy, Lee & White	BT and CT will be performed at LPL Laboratories	Daily

H243	COAGULATION SCREEN	*Platelet Count *Prothrombin studies *APTT	3 mL (1.5 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE AND 3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: entrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection,Electrical Impedence		Daily
C014	COBALT, BLOOD		4 mL (2 mL min.) whole blood in 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Patient should avoid mineral supplements, Vitamin B12 or Vitamin B Complex 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

C016	COBALT, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Patient should avoid mineral supplements, Vitamin B12 or Vitamin B Complex 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C015	COBALT, 24- HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in metal free container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Patient should avoid mineral supplements, Vitamin B12 or Vitamin B Complex 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
K011	COCAINE & METABOLITE SCREEN, URINE	Measured primarily as metabolite Benzoyllecgonine.	10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily

RW048	COCCIDIOIDES ANTIBODY		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunodiffusion	False-negatives may occur in patients with solitary pulmonary lesions. Patients with active Histoplasmosis may generate cross reactions. Elevated serum antibodies may be present in patients with Meningeal coccidiomycosis.	Sample by 7th of the month; Report after 2–3 weeks
Z155	COLLAGEN DISEASE ANTIBODIES PANEL	*ANA *Anti-ds-DNA *Sm *SS-A/Ro *SS-B/La *Scl-70 *Jo-1 *Centromere *U1RNP	4 mL (2.5 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see ARTHRITIS PANELS.	Sample by Tue / Fri 9 am; Report Same day
RW049	COMPLEMENT, TOTAL (CH50)		3 mL (1.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection and freeze immediately Ship frozen. DO NOT THAW	Automated liposome lysis assay	Also see individual Complement Components by name (C1 to C9)	Sample by 7th of the month; Report after 2–3 weeks

Z048	COMPLETE BODY PANEL	*SMA 26 *Hemogram *Lipid Profile, Extended*Urine R/E	2 mL (1 mL min.) plasma each from 2 Grey Top (Sodium Fluoride) tubes for Glucose F & PP AND 3 mL (2 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube AND 10 mL (5 mL min.) aliquot of random urine. Ship refrigerated. DO NOT FREEZE. 12 hours overnight fasting is mandatory.	See individual assays		Sample Daily by 9 am; Report 3days
Z021	COMPLETE BLOOD COUNT; CBC	*Hemoglobin *PCV *TLC *RBC Count *MCV *MCH *MCHC *Platelet Count *Automated DLC *Absolute Differential Counts *RDWDoes not include ESR & Peripheral Smear	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence, VCS	For Peripheral smear register code H029 separately.	Daily
OS121	COMT GENOTYPING		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Mon 11 am ; Report 8 days
OS122	CONNEXIN 26 MUTATION DETECTION		4 mL (2 mL min.)Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed 11 am; Report 9 days
OS123	CONNEXIN 30 MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Wed 11 am; Report 9 days

H027	COOMBS TEST, INDIRECT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Erythrocyte Magnetised Technology		Sample Daily by 9 am; Report Same Day
H028	COOMBS TEST, DIRECT		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Erythrocyte Magnetised Technology		Sample Daily by 9 am; Report Same Day
C018	COPPER, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If barium containing contrast media has been administered, the specimen should not be collected for 96 hours. Also see CERULOPLASMIN.	Sample Daily by 6 pm; Report Next day
C019	COPPER, SERUM		3 mL (1 mL min.) serum from 1 metal free White Top (Z-No Additive) tube available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium, iodine or barium containing contrast media has been administered, the specimen should not be collected for 96 hours. Also see CERULOPLASMIN.	Sample Daily by 6 pm; Report Next day

C017	COPPER, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid-washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If barium containing contrast media has been administered, the specimen should not be collected for 96 hours. Also see CERULOPLASMIN.	Sample Daily by 6 pm; Report Next day
G199	CORTICOSTERONE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Diurnal variation present; take sample between 8–10 a.m OR 4–6 p.m. Specify time and date on specimen container and test request form. Overnight fasting is preferred. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days
U021	CORTISOL, FREE, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine. Collect urine with 10 g of boric acid to maintain sample integrity. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
R043	CORTISOL, MIDNIGHT, SERUM		2 mL (1 mL min.) serum from 1 SST. Take sample between 11 p.m–12 midnight. Specify time and date on sample container and test request form. Ship refrigerated or frozen. Four hours fasting is recommended.	Chemiluminescent Immunoassay		Daily

Z075	CORTISOL, MORNING & EVENING, SERUM		2 mL (1 mL min.) serum from 1 SST for each timed specimen. Recommended timings: Morning: 7-9 a.m; Evening: 3-5 p.m. Overnight fasting is preferred for the morning specimen, four hours fasting is recommended for the evening specimen. Specify time and date of collection on container and test request form. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Also see CORTISOL, STIMULATION BY ACTH; CORTISOL SUPPRESSION BY DEXAMETHASONE	Daily
Z010	CORTISOL STIMULATION BY ACTH		2 mL (0.5 mL min.) serum from 1 SST for each timed specimen . Draw fasting baseline sample; administer 250 ug of synthetic ACTH (Synacthen) IV, draw additional samples at 30 and 60 min. Ship refrigerated or frozen. All tubes must be clearly marked with time drawn. Submit all tubes with 1 test request form.	Chemiluminescent Immunoassay	Synthetic ACTH (Synacthen) is not supplied by LPL. Refer to Y Test codes :Y028 Basal Y029 30 min.Y030 60 min.Y031 Customized	Daily
R016	CORTISOL, SERUM (MORNING SAMPLE)		2 mL (1 mL min.) serum from 1 SST. Diurnal variation present; take specimen between 7-9 a.m. Specify time and date on sample container and test request form. Ship refrigerated or frozen. Overnight fasting is preferred.	Chemiluminescent Immunoassay	Also see CORTISOL, STIMULATION BY ACTH; CORTISOL SUPPRESSION BY DEXAMETHASONE.	Daily
G198	CORTISONE, SERUM		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Diurnal variation present; take sample between 8–10 a.m OR 4–6 p.m. Specify time and date on specimen container and test request form. Overnight fasting is preferred. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days

Z058	CORTISOL SUPPRESSION BY DEXAMETHASONE, OVERNIGHT LOW DOSE		Administer tablet 1 mg Dexamethasone orally at 11 p.m. at night. Draw 2 mL (1 mL min.) serum from 1 SST at 8 am next morning.	Chemiluminescent Immunoassay	Patient should not be on steroid medication.	Daily
RW054	CORTISONE (COMPOUND E), 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine in 10 g boric acid in a screw capped container OR keep urine refrigerated during collection if preservative is not used. Record 24- hour volume on test request form and specimen container. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; report after 2–3 weeks.
G200	CORTISONE : CORTISOL RATIO	*Cortisone *Cortisol *Ratio	2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Diurnal variation present; take sample between 8–10 a.m OR 4–6 p.m. Specify time and date on specimen container and test request form. Overnight fasting is preferred. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days
RW053	CORTISOL FREE, SERUM		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	LC-MS/MS, Equilibrium Dialysis		Sample by 7th of the month; Report after 2–3 weeks
R175	CORTISOL, SALIVA		5ml (2ml min) saliva in a sterile screw capped container. Collect specimen before brushing teeth. Avoid stress, exercise, eating & drinking 30 minutes prior to sample collection. Rinse mouth thoroughly before collecting sample. Ship refrigerated or frozen.	ECLIA	The ideal sampling time is late at night between 11pm to 12 midnight.	Sample by Mon 9 am; Report Same day

Z059	CORTISOL SUPPRESSION BY DEXAMETHASONE, OVERNIGHT HIGH DOSE		Day 1: 2 mL (1 mL min.) serum from 1 SST. Draw baseline Cortisol sample at 8 a.m. Administer tablet 4 mg Dexamethasone orally at 11 p.m. at night. Day 2: Draw 2nd specimen at 8 a.m. next morning. Both containers must be marked with time and date drawn, and submitted with one test request form. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Patient should not be on steroid medication. Refer to Y codes Y038 Basal Y039 Post Dexamethasone	Daily
R017	CORTISOL, SERUM (EVENING SAMPLE)		2 mL (1 mL min.) serum from 1 SST. Diurnal variation present; take sample between 3-5 p.m. Specify time and date on sample container and test request form. Ship refrigerated or frozen. Four hours fasting is recommended.	Chemiluminescent Immunoassay	Also see CORTISOL, STIMULATION BY ACTH; CORTISOL SUPPRESSION BY DEXAMETHASONE	Daily
S110	COXSACKIE ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same Day
S176	COXSACKIE ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same Day
RW056	COXSACKIE A ANTIBODY PANEL	Includes 2, 4, 7, 9, 1 & 16	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Complement Fixation		Sample by 7th of the month; Report after 2-3 weeks
Z217	COXSACKIE ANTIBODIES, IgG & IgM PANEL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same Day

Z276	C-PEPTIDE STIMULATION BY GLUCAGON		2 mL (0.5 mL min.) serum from 1 SST for each timed specimen. Draw baseline fasting specimen. Administer 1 mg Glucagon I.V. Draw additional specimens at 5, 10 and 15 minutes post-Glucagon. Note time drawn on each container label and test request form. 12 hours overnight fasting mandatory . Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Glucagon is not supplied by LPL.Refer to Y Codes:Y041 Basal Y042 5 min. Y043 10 min. Y044 15 min. Y045 Customized	Daily
R055	C-PEPTIDE, FASTING, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is mandatory.	Chemiluminescent Immunoassay	Also see GLUCOSE TOLERANCE TEST WITH C-PEPTIDE (Z037).	Daily
R157	C-PEPTIDE, POST PRANDIAL, SERUM		2 mL (0.5 mL min.) serum from 1 SST 2 hours after starting meal. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Also see GLUCOSE TOLERANCE TEST WITH C-PEPTIDE (Z037).	Daily
Z037	C-PEPTIDE RESPONSE TO GLUCOSE, 5 HOURS		2 mL (1 mL min.) plasma from 1 GreyTop (Sodium Fluoride) tube AND 2 mL (0.5 mL min.) corresponding serum from 1 SST for each timed specimen. Draw baseline fasting specimen. Dissolve 75g Glucose in adults / 1.75 g / Kg body weight in children (not exceeding 75g) in 300 mL water. Administer orally over a period of 5 minutes. Draw additional specimens at 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 4.0 & 5.0 hours. Note time drawn and test required on each container label and test request form. Overnight fasting mandatory. Ship frozen. DO NOT THAW. Urine specimens NOTREQUIRED.	Chemiluminescent Immunoassay,Hexokinase	refer to Y CodesY074 BasalY075 0.5 hrY076 1.0 hrY077 1.5 hrY078 2.0 hrY079 2.5 hrY080 3.0 hrY081 4.0 hrY082 5.0 hrY083 Customized	Daily

U095	C-PEPTIDE, 24-HOUR URINE		5 mL (1 mL min.) aliquot of 24 hour urine. Do not use preservatives. Record 24 hr urine volume on container and Test Request Form. Ship frozen. DO NOT THAW.	Chemiluminescent Immunoassay		Daily
B029	CPK; CREATINE KINASE		2 mL (0.5 mL min.) serum from 1 SST. Hemolysed specimens are not acceptable. Ship refrigerated or frozen.	IFCC	Also see CPK ISOENZYME ELECTROPHORESIS, ALDOLASE & CARNITINE.	Daily
E007	CPK ISOENZYME ELECTROPHORESIS	Includes Total CPK.	2 mL (1.5 mL min.) fresh serum from 1 SST. Hemolysed specimens are not acceptable. Ship refrigerated or frozen.	Agarose gel Electrophoresis	See ALDOLASE & CARNITINE for Muscular Dystrophies.	Sample Wed by 9 am; Report Same day
R169	CPK-MB MASS		2 mL (1 mL min.) serum from 1 SST. Hemolysed specimens are not acceptable. Ship refrigerated or frozen.	ECLIA		Daily
B051	CPK-MB		2 mL (1 mL min.) serum from 1 SST. Hemolysed specimens are not acceptable. Ship refrigerated or frozen.	Enzymatic Immunoinhibition	Also see CPK ISOENZYME ELECTROPHORESIS.	Daily
RW028	C1q COMPLEMENT		2 mL (0.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen.	Radial Immunodiffusion		Sample by 7th of the month; Report after 2–3 weeks
B074	C-REACTIVE PROTEIN; CRP		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunoturbidimetry		Daily
B134	C-REACTIVE PROTEIN, CARDIO; hsCRP		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry		Daily
U045	CREATININE, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container.. Ship refrigerated or frozen.	Compensated Jaffe's Reaction IDMS Traceable		Daily

Z107	CREATININE CLEARANCE TEST	Includes serum and urine creatinine.	10 mL (5 mL min.) aliquot of 24-hour urine AND 2 mL (0.5 mL min.) serum from 1 SST. Record patient height, weight and total 24-hour urinary volume on test request form and urine container. Ship refrigerated or frozen.	Compensated Jaffe's Reaction IDMS Traceable		Daily
B023	CREATININE, SERUM		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Compensated Jaffe's Reaction IDMS Traceable		Daily
P082	CREATININE, FLUID		10 mL (5 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen.	Compensated Jaffe's Reaction IDMS Traceable		Daily
U007	CREATININE, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine. Record 24-hour volume and patient weight on test request form and urine container. Ship refrigerated or frozen.	Compensated Jaffe's Reaction IDMS Traceable		Daily
H115	CRYOGLOBULIN S QUALITATIVE TEST		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship at 18-22° C. DO NOT FREEZE.	Precipitation method		Sample Daily by 9 am; Report 5 days
H117	CRYOGLOBULIN S PANEL	Immunofixation Electrophoresis for classification of Cryoglobulins	3 mL (1.5 mL min.) serum. Collect blood in a PREWARMED Red Top (No Additive) tube using a PRE-WARMED syringe and needle at 37°C. Allow to clot and centrifuge preferably at 37°C. Separate serum and transport at room temperature. DO NOT FREEZE.	Immunofixation Electrophoresis	Strictly adhere to temperature requirements for sampling.	Sample Daily by 9 am; Report 5 days
S112	CRYPTOCOCCUS ANTIGEN, SERUM		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed 9 am; Report Same day
T011	CRYPTOSPORIDIUM / ISOSPORA IDENTIFICATION, STOOL	Includes Modified Ziehl Neelsen stain.	5 g (1 g min.) of fresh stool in special transport medium available from LPL. Ship at 18–22°C.	Conventional microscopic examination		Sample Daily by 3.30 pm; Report Same Day

S111	CRYPTOCOCCUS ANTIGEN, CSF		2 mL (0.5 mL min.) Cerebrospinal fluid in a sterile screw capped container. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed 9 am; Report Same day
M202	CRYPTOSPORIDIUM ANTIGEN, STOOL		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pm; Report Same day
RW057	CRYPTOCOCCUS ANTIBODY		2 mL (0.5 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum immediately. Ship refrigerated or frozen.	Tube Agglutination	Also see CULTURE, FUNGUS	Sample by 7th of the month; Report after 2–3 weeks
OS124	CRYPTOCOCCUS NEOFORMANS DETECTION, PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube OR 1 mL CSF/Sputum/Fluids OR Tissue in normal saline. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am; Report 18 days
E004	CSF PROTEIN ELECTROPHORESIS / OLIGOCLONAL BANDS		2 mL (1.5 mL min.) CSF in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Agarose gel Electrophoresis	Also see MULTIPLE SCLEROSIS PANELS; MYELIN-ASSOCIATED GLYCOPROTEIN; MYELIN ASSOCIATED GLYCOPROTEIN SULFATED GLUCURONIC PARAGLOBOSIDE AUTOANTIBODY, IgM.	Sample Daily by 9 am; Report Next day

E018	CSF & SERUM PROTEIN ELECTROPHORE SIS / OLIGOCLONAL BANDS		2 mL (1.5 mL min.) CSF in a sterile screw capped container AND 2 mL (1.5 mL min.) serum from 1 SST. Serum should be separated within 1 hour of collection. Ship refrigerated. DO NOT FREEZE.	Agarose gel / Capillary Electrophoresis	Also see MULTIPLE SCLEROSIS PANELS; MYELIN- ASSOCIATED GLYCOPROTEIN; MYELIN ASSOCIATED GLYCOPROTEIN SULFATED GLUCURONICPARAG LOBOSIDE AUTOANTIBODY, IgM.	Sample Daily by 9 am; Report Next day
M010	CULTURE, AEROBIC, BLOOD, RAPID	Includes ANTIBIOTIC SENSITIVITY if indicated.	Collect 8-10 mL blood aseptically in special Plus Aerobic bottle available from LPL. For children less than 12 yrs collect 2 mL blood in Special Paeds Plus bottle available from LPL. Mix by gentle swirling. Ship at room temperature.	Rapid Automated Aerobic Culture, Identification and Sensitivity	Also see CULTURE, ANAEROBIC,BLOOD, RAPID.	Sample Daily by 3.30 pmReport Interim: 24 & 48 hrs Final: 6 days
M015	CULTURE, AEROBIC, EAR SWAB	Includes ANTIBIOTIC SENSITIVITY if indicated.	Ear swab inoculated in special Aerobic Transport Medium available from LPL. Specify left or right ear swab on transport medium vial and test request form. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3 days

M133	CULTURE, AEROBIC, BODY FLUIDS, RAPID	Includes ANTIBIOTIC SENSITIVITY if indicated.	Collect 8-10 mL body fluids aseptically in special Plus Aerobic bottle available from LPL. For children less than 12 yrs OR low fluid volumes, collect 3 mL (0.5 mL min.) in Special Paeds Plus bottle available from LPL. Mix by gentle swirling. Ship at room temperature.	Rapid Automated Aerobic Culture, Identification and Sensitivity	This test code includes Bone Marrow & body fluids from sterile sites (Pleural/Pericardial/Synovial/Ascitic/CSF/ / Ocular). This test code does not include Urine / Pus specimens.	Sample Daily by 3.30 pm Report Interim: 24 & 48 hrs Final: 6 days
M014	CULTURE, AEROBIC, EYE SWAB	Includes ANTIBIOTIC SENSITIVITY if indicated.	Eye swab inoculated in special Aerobic Transport Medium available from LPL. Specify left or right eye swab on transport medium vial & test request form. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3days
P058	CULTURE, AEROBIC, MISCELLANEOUS WITH GRAM STAIN	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit catheter tip in sterile screw capped container / 5 mL (1 mL min.) body effusion, aspirates, etc., in sterile screw capped container OR swab in Aerobic Transport medium. Specify collection site on container and test request form.	Conventional culture and Automated Sensitivity, Gram stain		Sample Daily by 3.30 pm; Report 3days
M007	CULTURE, URINE	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit first morning mid stream urine in a special urine collection vacutainer kit available from LPL as per instructions enclosed OR sterile screw capped container. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3days
M145	CULTURE, AEROBIC, CSF	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit 2-5 mL (1 mL min.) CSF in a sterile screw capped container. Ship at 18–22°C.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3 days

M078	CULTURE, DIPHTHERIA	Does not include ANTIBIOTIC SENSITIVITY.	Submit nasal / throat swab (mucus membrane of pharynx / tonsils) obtained by rubbing over affected / localised lesion, in Amie's transport medium available from LPL. Ship refrigerated.	Conventional culture		Sample Daily by 3.30 pm; Report 3days
M018	CULTURE, GONOCOCCUS	Does not include ANTIBIOTIC SENSITIVITY.	Submit swab from urethra, anus, eye, endocervix or pharynx in Amie's charcoal medium available from LPL, Specify collection site on transport container and test request form. Ship at room temperature.	Conventional culture and Automated identification		Sample Daily by 3.30 pm; Report 3days
M072	CULTURE, ACTINOMYCES	Does not include ANTIBIOTIC SENSITIVITY.	Submit pus, intrauterine devices, exudate, CSF, pleural or other fluids, tissues in an Anaerobic transport medium available from LPL. Ship at 18–22°C.	Anaerobic culture on specific media		Sample Daily by 3.30 pm; Report Interim : 1 week Final : 3 weeks
M017	CULTURE, AEROBIC, NASAL SWAB	Includes ANTIBIOTIC SENSITIVITY if indicated.	Nasal swab inoculated in special Aerobic Transport Medium available from LPL. Specify left or right nasal swab on transport medium vial and test request form. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3days
M012	CULTURE, AEROBIC, PUS	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit pus in sterile screw capped container OR pus swab inoculated in special Aerobic Transport Medium available from LPL. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3days

M011	CULTURE, AEROBIC, SPUTUM	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit 5–10 mL (1 mL min.) sputum in a sterile screw capped container. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity	If AFB is suspected, use test code M050 CULTURE, AFB, RAPID. If fungus is suspected, use test code M104 CULTURE, FUNGUS, and/or M096 CULTURE, FUNGUS, NOCARDIA.	Sample Daily by 3.30 pm; Report 3 days
M016	CULTURE, AEROBIC, VAGINAL SWAB	Includes ANTIBIOTIC SENSITIVITY if indicated.	Vaginal swab inoculated in Amie's charcoal transport medium available from LPL. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity	If N. gonorrhoeae is suspected, use test code M018 CULTURE, GONOCOCCUS.	Sample Daily by 3.30 pm; Report 3 days

M139	CULTURE, AFB (MYCOBACTERIUM), RAPID WITH PNB TEST	Does not include ANTI-TUBERCULAR DRUG SENSITIVITY.	Blood/Bone Marrow: Collect 8 mL (5 mL min.) Blood/Bone marrow in Myco F / Lytic bottle available from LPL. Mix by gentle inversion 8–10 times. Ship at 18–22°C. CSF: Collect 2 mL (1 mL min.) CSF in a Red Top (No Additive) tube or sterile screw capped container. Do not use SST gel barrier tubes. Ship refrigerated. Pus/Body fluids (Pleural/ Pericardial /Ascitic/Synovial/Ocular) Aspirates/ Semen /BAL / Bronchial washings: Submit as much as possible (1 mL min.) in sterile screw capped container. Ship refrigerated. Endometrial curettings/Tissue: Submit in sterile normal saline in a sterile screw capped container. Ship refrigerated. Swabs: Submit swabs in 1 mL sterile normal saline in sterile screw capped container. Ship refrigerated. Sputum/Urine: Submit 2 Spot (random) morning samples, 5–10 mL (1 mL min.) sputum / 10 mL (5 mL min.) urine in sterile screw capped container. Do not pool the samples. Ship refrigerated. Gastric Lavage: Submit 5–10 mL (2 mL min.) gastric lavage in a sterile screw capped container. Ship refrigerated . Overnight fasting is mandatory.	Rapid Automated Fluorescent AFB Culture and Automated AFB Identification by PNB test.	For AFB Antitubercular Sensitivity, use appropriate test code. Request as additional tests.	Sample Daily by 3.30 pm Report Interim: 14 days & 4 weeks Final: 6 weeks PNB test: 10 days after AFB culture is positive
M068	CULTURE, ANAEROBIC	Does not include ANTIBIOTIC SENSITIVITY.	Submit as much as possible (1 mL min.) CSF OR Pus/Body Fluids (Pleural/Pericardial/Ascitic/Synovial/Ocular)/ Aspirates OR Endometrial curettings/ Tissues OR Bone Marrow in Anaerobic transport medium available from LPL . Ship at room temperature.	Conventional anaerobic culture and Automated Identification		Sample Daily by 3.30 pm; Report 8 days

M134	CULTURE, ANAEROBIC, BODY FLUID, RAPID	Does not include ANTIBIOTIC SENSITIVITY.	Collect 8-10 mL body fluids aseptically in special Plus Anerobic bottle available from LPL. Mix by gentle swirling. Ship at room temperature.	Rapid Automated Anaerobic Culture and Identification	This test includes body fluid from sterile sites. This test code does not include Urine /Pus specimens.	Sample Daily by 3.30 pm; Report 6days
M076	CULTURE, CHOLERA	Includes ANTIBIOTIC SENSITIVITY if indicated.	5 mL (1 mL min.) stool in a sterile screw capped container within 6 hours. If there is a delay of more than 6 hours collect specimen in special Stool Culture Transport Medium available from LPL. Ship refrigerated.DO NOT FREEZE	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3days
M020	CULTURE, AEROBIC, MISCELLANEOUS	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit tissue in normal saline /catheter tip in sterile screw capped container / 5 mL (1 mL min.) body effusion, aspirates, etc., in sterile screw capped container OR swab in Aerobic Transport medium. Specify collection site on container and test request form.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3days
M021	CULTURE, ANAEROBIC, BLOOD, RAPID	Does not include ANTIBIOTIC SENSITIVITY	Collect 8-10 mL blood aseptically in special Plus Anerobic bottle available from LPL. Mix by gentle swirling. Ship at room temperature.	Rapid Automated Anaerobic Culture and Identification	Also see CULTURE, AEROBIC, BLOOD,RAPID.	Sample Daily by 3.30 pm; Report 6days
M195	CULTURE ESBL (EXTENDED SPECTRUM BETA- LACTAMASE), SCREENING		Rectal / Wound SWABS / Pure GROWTH in Amies transport medium. Ship refrigerated only. DO NOT FREEZE.	Chromogenic		Sample Daily by 11 am;Report Next day
M024	CULTURE, FUNGUS, BLOOD / BODY FLUIDS RAPID	Includes ANTIBIOTIC SENSITIVITY for Yeast only if indicated.	Collect 8-10 mL whole blood / body fluids aseptically in a Myco / F-Lytic bottle available from LPL. Mix by gentle swirling. Ship at room temperature.	Rapid Automated Fungus Culture		Sample Daily by 3.30 pmReport Interim: 2 days Final: 15 days

M095	CULTURE, HAEMOPHILUS	Does not include ANTIBIOTIC SENSITIVITY.	Submit Throat / Eye / Ear swab OR Sputum in Amie's Transport Medium available from LPL OR CSF in sterile screw capped container. Ship at room temperature.	Conventional culture and Automated identification		Sample Daily by 3.30 pm; Report 3days
M198	CULTURE KPC (CARBAPENEMASE), SCREENING		Rectal / Wound SWABS / Pure GROWTH in Amies transport medium. Ship refrigerated only. DO NOT FREEZE.	Chromogenic		Sample Daily by 11 am Report Next day
M197	CULTURE MRSA (METHICILLIN RESISTANT STAPHYLOCOCCUS AUREUS), SCREENING		Nasal / Axilla / Groin / Wound SWABS / Pure GROWTH in Amies transport medium. Ship refrigerated only. DO NOT FREEZE.	Chromogenic		Sample Daily by 11 am Report Next day
M096	CULTURE, NOCARDIA	Does not include ANTIBIOTIC SENSITIVITY.	Submit pus from discharging sinus / sputum in sterile screw capped container without fixatives OR tissue/ biopsy material / swabs in 1 mL sterile normal saline. Ship refrigerated.	Conventional culture		Sample Daily by 3.30 pm; Report Interim : 1 week Final: 3 weeks
M035	CULTURE, WATER; PRESUMPTIVE COLIFORM COUNT		180 mL potable water in (3 × 60 mL) sterile, screw capped containers. Ship refrigerated. Specify water source on containers and test request form. Strictly adhere to requested volumes.	Multiple tube test		Sample Daily by 3.30 pm; Report 3days
M009	CULTURE, AEROBIC, THROAT SWAB	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit throat swab in Amie's transport medium available from LPL. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity	If C. diphtheriae is suspected, use test code M078 CULTURE, DIPHTHERIA.	Sample Daily by 3.30 pm; Report 3days

M141	CULTURED AEROBIC BACTERIA IDENTIFICATION AND SENSITIVITY		Submit pure isolate on appropriate media. Ship refrigerated. DO NOT FREEZE. Specify collection site and provide brief clinical history.	Sub-culture, Automated identification and Sensitivity	Stability of the isolate will be determined by the laboratory.	Sample Daily by 3.30 pm; Report 3days
M140	CULTURED AEROBIC BACTERIA SENSITIVITY		Submit pure known isolate on appropriate media. Ship refrigerated. DO NOT FREEZE. Specify collection site and name of the isolate grown.	Sub-culture and Automated Sensitivity	Mixed bacteria culture sensitivity or fungus sensitivity or AFB sensitivity are not included in this test.	Sample Daily by 3.30 pm; Report Nextday
M182	CULTURED AFB ANTITUBERCULAR DST, RAPID 10 SECOND LINE DRUGS PANEL	*Amikacin *Kanamycin *Capreomycin *Clofazimine *Ethionamide *Ofloxacin *Levofloxacin *Moxifloxacin *PAS *Cycloserine	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M183	CULTURED AFB ANTITUBERCULAR DST, RAPID 5 SECOND LINE DRUGS PANEL	*Ofloxacin *Levofloxacin *Moxifloxacin *PAS *Cycloserine	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M181	CULTURED AFB OTHER THAN TUBERCULOSIS (MOTT), DST FOR RAPID GROWERS		Submit pure growth of MOTT on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests the pure growth (MOTT) obtained in the Laboratory specimens already submitted will be processed for MOTT Drug Sensitivity.	Broth Microdilution (MIC)	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 7days
M013	CULTURE, AEROBIC, SEMEN	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit semen in a sterile screw capped container. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity		Sample Daily by 3.30 pm; Report 3days

M050	CULTURE, AFB (MYCOBACTERIUM), RAPID	Does not include ANTI-TUBERCULAR DRUG SENSITIVITY.	Blood/Bone Marrow: Collect 8 mL (5 mL min.) Blood/Bone marrow in Myco F / Lytic bottle available from LPL. Mix by gentle inversion 8–10 times. Ship at 18–22°C. CSF: Collect 2 mL (1 mL min.) CSF in a Red Top (No Additive) tube or sterile screw capped container. Do not use SST gel barrier tubes. Ship refrigerated. Pus/Body fluids (Pleural/Pericardial/Ascitic/ Synovial/Ocular) Aspirates /Semen/ BAL Bronchial washings: Submit as much as possible (1 mL min.) in sterile screw capped container. Ship refrigerated. Endometrial curettings/Tissue: Submit in sterile normal saline in a sterile screw capped container. Ship refrigerated. Swabs: Submit swabs in 1 mL sterile normal saline in sterile screw capped container. Ship refrigerated. Sputum/Urine: Submit 2 Spot (random) morning samples, 5-10 mL (1 mL min.) sputum/10 mL (5 mL min.) urine in sterile screw capped container. Do not pool the samples. Ship refrigerated. Gastric Lavage: Submit 5-10 mL (2 mL min.) gastric lavage in a sterile screw capped container. Ship refrigerated . Overnight fasting is mandatory.	Automated Fluorescent, ICT, MALDI-TOF	For AFB Antitubercular Sensitivity, use appropriate test code. Request as additional tests.	Sample Daily by 3.30 pm Report Interim: 14 days & 4 weeks Final: 6 weeks Positive results between 2-42 days shall be informed telephonically.
M070	CULTURE, CLOSTRIDIUM DIFFICILE	Does not include ANTIBIOTIC SENSITIVITY.	5 g (1 g min.) of fresh stool in Anaerobic transport medium available from LPL. Ship refrigerated.	Conventional anaerobic culture and Automated Identification	Also see C. DIFFICILE TOXIN A and B, STOOL.	Sample Daily by 3.30 pm; Report 8 days

M104	CULTURE, FUNGUS	Does not include ANTIFUNGAL SENSITIVITY.	CSF: 2 mL (1 mL min.) in a Red Top (No Additive) tube or sterile screw capped vial. Do not use SST gel barrier tubes. Ship refrigerated. Bone Marrow: 2 mL (0.5 mL min.) Bonemarrow in 1 Green Top (Sodium Heparin) tube. Ship refrigerated. Pus/Body fluids (Pleural/ Pericardial/ Ascitic/ Synovial/ Ocular/ Bronchial washings/BAL). Aspirates /Semen/ Stool: Submit as much as possible (1 mL/1g min.) in sterile screw capped container. Ship refrigerated. Endometrial curettings/Tissue: Submit in sterile normal saline in a sterile screw capped container. Ship refrigerated. Swabs: Submit swabs in 1 mL sterile normal saline in sterile screw capped container. Ship refrigerated. Sputum/Urine: Submit 2 Spot (random) morning samples, 5–10 mL (1 mL min.) sputum/10 mL (5 mL min.) urine in sterile screw capped container. Ship refrigerated. Skin: Submit 2–3 mm scrapings from periphery of active skin lesion in sterile container or between 2 sterile glass slides. Wrap slides in black paper. Ship refrigerated. Hair: Submit plucked 10–12 infected, lustreless hair including root in a sterile container. Ship refrigerated. Nail: Submit clippings and powdery material from friable and discolored region of affected nail. Ship refrigerated.	Conventional Fungus Culture	Contact lab within 10 days of release of result if Antifungal drug sensitivity is required.	Sample Daily by 3.30 pm; Report Interim 1: 1 week Interim 2: 2 weeks Interim 3: 3 weeks Final: 4 weeks
M008	CULTURE, STOOL	Includes ANTIBIOTIC SENSITIVITY if indicated.	Submit 2–5 g (1 g min.) stool in a Special Stool Culture Transport Medium available from LPL. Ship refrigerated.	Conventional Culture, Automated Identification & Sensitivity	Also see C. DIFFICILE TOXIN A and B and CULTURE, CLOSTRIDIUM DIFFICILE.	Sample Daily by 3.30 pm; Report 3 days

M069	CULTURED AEROBIC BACTERIA IDENTIFICATION	Does not include ANTIBIOTIC SENSITIVITY.	Submit pure isolate on appropriate media. Ship refrigerated. DO NOT FREEZE. Specify collection site and provide brief clinical history.	Sub-culture and Automated identification	Mixed bacteria culture identification or fungus identification or AFB identification are not included in this test. For ANTIBIOTIC SENSITIVITY, use test code M140 CULTURED, AEROBIC, BACTERIA SENSITIVITY. Request as additional test.	Sample Daily by 3.30 pm; Report Nextday
M180	CULTURED AFB OTHER THAN TUBERCULOSIS (MOTT), DST FOR SLOW GROWERS		Submit pure growth of MOTT on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests the pure growth (MOTT) obtained in the Laboratory specimens already submitted will be processed for MOTT Drug Sensitivity.	Broth Microdilution (MIC)	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 14days

M093	CULTURED FUNGUS IDENTIFICATION	Does not include ANTIFUNGAL SENSITIVITY	Submit pure isolate on appropriate medium. Ship refrigerated. DO NOT FREEZE.	Conventional fungus culture and identification		Sample Daily by 3.30 pm; Report 3–28days
M176	CULTURED AFB ANTITUBERCULAR DST, 15 DRUGS PANEL	Streptomycin, Isoniazid, Rifampicin, Ethambutol Pyrazinamide Clofazimine Ethionamide Capreomycin Kanamycin Amikacin Ofloxacin +PAS Cycloserine Levofloxacin Moxifloxacin	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M224	CULTURED AFB ANTITUBERCULAR DST, 16 DRUGS PANEL	Streptomycin, Isoniazid, Rifampicin, Ethambutol Pyrazinamide Clofazimine Ethionamide Capreomycin Kanamycin Amikacin Ofloxacin +PAS Cycloserine Levofloxacin Moxifloxacin * Linezolid	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M109	CULTURED AFB ANTITUBERCUL AR DST, STREPTOMYCIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M110	CULTURED AFB ANTITUBERCUL AR DST, ISONIAZID		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M177	CULTURED AFB ANTITUBERCUL AR DST, CYCLOSERINE		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M115	CULTURED AFB ANTITUBERCUL AR DST, ETHIONAMIDE		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M092	CULTURED ANAEROBIC BACTERIAIDENT IFICATION	Does not include ANTIBIOTIC SENSITIVITY.	Submit pure isolate on appropriate medium. Ship refrigerated. DO NOT FREEZE. Specify collection site and provide brief clinical history.	Sub-Culture and Automated Identification		Sample Daily by 3.30 pm; Report 8days

M107	CULTURED AFB ANTITUBERCUL AR DST, 10 DRUGS PANEL	Streptomycin,Iso niazid,Rifampicin, Ethambutol Pyrazinamide +Clofazimine Ethionamide Capreomycin Kanamycin Amikacin	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M123	CULTURED AFB ANTITUBERCUL AR DST, 11 DRUGS PANEL	Streptomycin, Isoniazid, Rifampicin, Ethambutol Pyrazinamide Clofazimine Ethionamide Capreomycin Kanamycin Amikacin +Ofloxacin	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M105	CULTURED AFB ANTITUBERCUL AR DST, 4 DRUGS PANEL	Streptomycin,Iso niazid,Rifampicin, Ethambutol	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M106	CULTURED AFB ANTITUBERCUL AR DST, 5 DRUGS PANEL	Streptomycin,Iso niazid,Rifampicin, Ethambutol + Pyrazinamide	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M225	CULTURED AFB ANTITUBERCUL AR DST, LINEZOLID		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M117	CULTURED AFB ANTITUBERCUL AR DST, CAPREOMYCIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M114	CULTURED AFB ANTITUBERCUL AR DST, AMIKACIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M120	CULTURED AFB ANTITUBERCUL AR DST, PAS		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M111	CULTURED AFB ANTITUBERCUL AR DST, RIFAMPICIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	u	Sample Daily by 3.30 pm; Report 21days
M112	CULTURED AFB ANTITUBERCUL AR DST, ETHAMBUTOL		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M113	CULTURED AFB ANTITUBERCUL AR DST, KANAMYCIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M119	CULTURED AFB ANTITUBERCUL AR DST, OFLOXACIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M124	CULTURED AFB ANTITUBERCUL AR DST, 12 DRUGS PANEL	Streptomycin, Isoniazid, Rifampicin, Ethambutol Pyrazinamide Clofazimine Ethionamide Capreomycin Kanamycin Amikacin +Ofloxacin + PAS	Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M108	CULTURED AFB ANTITUBERCUL AR DST, PYRAZINAMIDE		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M179	CULTURED AFB ANTITUBERCUL AR DST, MOXIFLOXACIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M178	CULTURED AFB ANTITUBERCUL AR DST, LEVOFLOXACIN		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days

M118	CULTURED AFB ANTITUBERCUL AR DST, CLOFAZIMINE		Submit pure growth of Mycobacterium tuberculosis complex on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the Laboratory from specimens already submitted will be processed for Antitubercular Sensitivity.	Automated Fluorescent	Specimens submitted in the laboratory for AFB culture are preserved for a period of 15 days if the culture is positive. If Antitubercular sensitivity is required on this sample, kindly order the test within 15 days	Sample Daily by 3.30 pm; Report 21days
M223	CULTURED MOLD ANTIFUNGAL DRUG SENSITIVITY (DST) TEST		Submit Pure GROWTH of Mold on appropriate media in sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Micro-broth		Sample Daily by 9 am; Report 7days
G159	CYCLOSPORINE A		3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. Ideal Sampling Time: Just before the next dose (Trough levels). Provide date and time of sampling, amount & time of dose.	LC-MS/MS		Sample Daily by 5 pm; Report Same day
OS126	CYP3A4*22 GENOTYPING		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample by Sat 11 am; Report 10 days
OS125	CYP2C19 GENOTYPING		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample by Wed 11 am; Report 11 days

B173	CYSTATIN C		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Nephelometry		Sample by Mon / Wed / Fri 4 pm; Report Same day
S033	CYSTICERCOSIS (TAENIA SOLIUM) ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
OS127	CYSTIC FIBROSIS MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube OR 10 mL Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 11 am; Report 13 days
R162	CYSTIC FIBROSIS, NEWBORN, SCREEN		1 drop of heel prick blood each on 3 spots of filter paper available from LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen	Fluoroimmunoassay		Sample Tue / Thu / Sat by 5 pm; Report 2 days
J006	CYTOLOGY (PAP SMEAR), GENITAL, FEMALE, CONVENTIONAL		Collect smears onto 2 glass slides and fix immediately in Absolute alcohol for 30 min and air dry. Do not allow smears to air dry before fixation. Brief clinical history in Cytology Requisition Form is mandatory.	Microscopy of Papanicolaou stained smears	Do not take smears during menstruation.	Sample Daily by 6 pm; Report Next day
J152	CYTOLOGY (PAP SMEAR), GENITAL, FEMALE, LIQUID BASED		Collect Endocervical smears with a brush after inserting the speculum. Dip & swirl the brush completely into the container with liquid fixative. Mix thoroughly. Use Special Collection Kit available from LPL only. Ship at room temperature. Brief clinical history in Cytology Requisition Form is mandatory.	Microscopy of Papanicolaou stained smears	Do not take smears during menstruation.	Sample Daily by 6 pm; Report Next day

J011	CYTOLOGY EXAMINATION, BODY FLUIDS		Submit Body fluids (Ascitic / Pleural / Pericardial / CSF/ Synovial / Prostatic / Urethral) OR Bronchial washings / Lavage OR Sputum OR second morning Urine specimen OR Nipple discharge in a sterile screw capped container. Mix with equal quantities of 50% alcohol. Alternatively centrifuge, make smears from the deposit and fix in Absolute alcohol for 30 minutes and air dry. Do not allow smears to air dry before fixation. Specimens in syringes are not acceptable. Brief clinical history in Cytology Requisition Form is mandatory.	Microscopy of Papanicolaou / MGG stained smears		Sample Daily by 6 pm; Report Next day
J074	CYTOLOGY, CELL BLOCK EXAMINATION		Submit cytopathological material in a screw capped container and fix immediately in equal quantity of 50% alcohol for 30 min. Give brief clinical history on Cytopathology Test Requisition form.	Cell block		Sample Daily by 6 pm; Report 3 days
J004	CYTOLOGY, FINE NEEDLE ASPIRATION; FNAC		Aspirate material from Breast / Superficial lumps / CT or USG guided deep lumps. Smear on 4 glass slides and fix immediately in Absolute alcohol for 30 minutes. Air dry completely. Do not allow smears to air dry before fixation. Specimens in syringes are not acceptable. Brief clinical history in Cytology Requisition Form is mandatory.	Microscopy of Papanicolaou / MGG stained smears		Sample Daily by 6 pm; Report Next day
J008	CYTOLOGY CONSULTATION , SECOND OPINION		Submit cytopathological slides previously reported upon. Previous report(s), clinical history to be submitted with the Cytopathology Test Requisition form.	Microscopy		Sample Daily by 6 pm; Report Next day

Z122	CYTOMEGALOVIRUS (CMV) ANTIBODIES PANEL, IgG & IgM		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
S041	CYTOMEGALOVIRUS (CMV) ANTIBODY, IgG		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
S022	CYTOMEGALOVIRUS (CMV) ANTIBODY, IgM		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
B156	CYTOMEGALOVIRUS (CMV) AVIDITY, IgG		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample Daily by 9 am; Report Same Day
N017	CYTOMEGALOVIRUS (CMV), QUALITATIVE PCR		5 mL (3 mL min.) Amniotic fluid in a sterile screw capped container collected after 21 weeks of pregnancy for Prenatal testing OR 10 mL(5 mL min) Urine/Saliva in a sterile screw capped container OR 2 mL(1 mL min) Plasma from 1 Lavendar top (EDTA) tube OR 2 mL(1 mL min) Serum from 1 Red top (No additive) tube OR 2 mL(1 mL min) CSF/BAL fluid in a sterile screw capped container OR 5 g (min 3 g) Fresh Tissue in normal saline.Ship refrigerated. Separate serum / plasma aseptically within 2 hours of collection. 12 hours overnight fasting is recommended. Brief clinical history and nature of specimen on Test Request Form is mandatory.	Real time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat

N031	CYTOMEGALOVIRUS (CMV), DNA, QUANTITATIVE, PCR		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma aseptically within 6 hours of collection. Ship refrigerated or frozen. 12 hours overnight fasting is recommended.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
H134	D – DIMER, QUANTITATIVE		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to one labelled clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Immunoturbidimetry	Also see DEEP VEIN THROMBOSIS PANEL & FIBRINOGEN DEGRADATION PRODUCTS.	Daily
Z056	DEEP VEIN THROMBOSIS PANEL; DVT PANEL	*Antithrombin Activity Functional *Cardiolipin Antibodies IgG and IgM *dRVVT *Protein C Functional *Protein S Functional	2 mL (1 mL min.) serum from 1 SST AND 6 mL Whole blood in 2 Blue Top (Sodium Citrate) tubes. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical & Photo-optical Clot Detection, EIA	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician	Sample by Wed 7 am; Report Fri

Z674	DENGUE FEVER ANTIBODIES IgG & IgM, EIA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Recommended test is NS1 Antigen by ELISA in the first 5 days of fever. After 7-10 days of fever, the recommended test is Dengue feverantibodies IgG & IgM by ELISA.	Sample Daily by 12 noon; Report same day
S151	DENGUE FEVER ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	ICT	Recommended test is NS1 Antigen by ELISA in the first 5 days of fever. After 7-10 days of fever, the recommended test is Dengue fever antibodies IgG & IgM by ELISA.	Daily
Z510	DENGUE FEVER COMBINED PANEL	*Dengue fever antibodies IgG & IgM *Dengue fever NS1 antigen	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	ICT, Lateral Flow Immuno- chromatography	Recommended test is NS1 Antigen by ELISA in the first 5 days of fever. After 7-10 days of fever, the recommended test is Dengue feverantibodies IgG & IgM by ELISA.	Daily

S217	DENGUE FEVER IgG ANTIBODY, EIA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Recommended test is NS1 Antigen by ELISA in the first 5 days of fever. After 7-10 days of fever, the recommended test is Dengue feverantibodies IgG & IgM by ELISA.	Sample Daily by 12 noon; Report same day
H063	DENGUE FEVER NS1 ANTIGEN, RAPID		2 mL (1 mL min.) serum from 1 SST OR 2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Ship refrigerated or frozen.	Lateral Flow Immunochromatography	Recommended test is NS1 Antigen by ELISA in the first 5 days of fever. After 7-10 days of fever, the recommended test is Dengue feverantibodies IgG & IgM by ELISA.	Daily
S225	DENGUE FEVER NS1 ANTIGEN, EIA		:2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample Daily by 12 noon; Report same day

S216	DENGUE FEVER IgM ANTIBODY, EIA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Recommended test is NS1 Antigen by ELISA in the first 5 days of fever. After 7-10 days of fever, the recommended test is Dengue fever antibodies IgG & IgM by ELISA.	Sample Daily by 12 noon; Report same day
Z853	DENGUE & CHIKUNGUNYA VIRUSES, PCR QUALITATIVE		3 mL (2 mL min.) Whole Blood / Plasma in 1 Lavender Top (EDTA) Tube OR 2 mL (1 mL min.) Serum from 1 SST OR 2 mL (1 mL min.) CSF in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE	Real Time PCR		Sample Daily by 11 am; Report Next day
N132	DENGUE VIRUS, PCR QUALITATIVE		3 mL (2 mL min.) Whole Blood / Plasma in 1 Lavender Top (EDTA) Tube OR 2 mL (1 mL min.) Serum from 1 SST OR 2 mL (1 mL min.) CSF in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE	Real Time PCR	NA	Sample Daily by 11 am; Report Next day
RW058	DEOXYCORTICO STERONE		3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum from cells immediately. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
RW059	DESIPRAMINE		4 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum and transfer into screw capped vial. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: >12 hours post oral dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
S099	DESMOGLEIN (DSG) 1 AND 3 ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	IFA		Sample by Mon 9 am; Report Wed 6 pm

RW060	DEXAMETHASONE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Draw specimen between 8–10 A.M. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
R100	DHEA; DEHYDROEPIANDROSTERONE		2 mL (0.5 mL min.) serum from 1 SST. Overnight fasting is preferred. Specify age and sex on test request form. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see DHEA STIMULATION BY ACTH.	Sample by Tue 9 am; Report Same day
R044	DHEA SULPHATE; DEHYDROEPIANDROSTERONE SULPHATE		2 mL (0.5 mL min.) serum from 1 SST. Overnight fasting is preferred. Specify age and sex on test request form. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay		Daily
Z022	DHEA STIMULATION BY ACTH		2 mL (0.5 mL min.) serum from 1 SST for each timed specimen. Draw fasting baseline sample; administer 250 ug of synthetic ACTH (Synacthen) IV, draw additional sample at 90 min. Ship refrigerated or frozen. Both tubes must be clearly marked with time drawn. Submit both tubes with 1 test request form.	Enzyme Immunoassay	Synthetic ACTH (Synacthen) is not supplied by LPL. Refer to Y CODES Y046 Basal Y047 90 min. Y049 Customized	Sample by Tue 9 am; Report Same day
Z312	DIABETES PANEL BASIC	*Glucose F / PP / R *HbA1c	2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 3 mL (1.5 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Hexokinase, HPLC		Daily
Z359	DIABETES SCREEN	*Glucose F + PP *Urine Glucose	2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube each AND 10 mL (5 mL min.) aliquot of first morning urine. Ship refrigerated. DO NOT FREEZE. Overnight fasting is mandatory.	Hexokinase, Dipstick		Daily

Z285	DIABETES PANEL 1	*Glucose F + PP *Cholesterol * Triglycerides*Hb A1c *Urea *Creatinine	2 mL (1 mL min.) plasma each from 2 Grey Top (Sodium Fluoride) tubes for Glucose F + PP AND 3 mL (1.5 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Minimum 12 hours fasting is mandatory.	Spectrophotometry, HPLC		Daily
Z286	DIABETES PANEL 2	*Glucose F + PP *HbA1c *Uric Acid *Creatinine *Cholesterol *HDL *LDL Calculated *Triglyceride *Microalbumin *Urine R/E	2 mL (1 mL min.) plasma each from 2 Grey Top (Sodium Fluoride) tubes for Glucose F + PP AND 3 mL (2 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube AND 10 mL (5 mL min.) aliquot of first morning urine Ship refrigerated. DO NOT FREEZE. 12 hours strict overnight fasting is mandatory.	Spectrophotometry, HPLC, Dipstick		Daily
Z833	DIABETES TYPE 1, AUTOIMMUNE PROFILE	*GAD-65* Anti- IA2* Insulin Autoantibody* Islet Cell Antibody	5 mL (3 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay, IFA		Sample by Fri 9 am; Report Next day
RW061	DIAZEPAM	Includes metabolite Nordiazepam.	3 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
RW062	DIGITOXIN		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Separate serum from cells immediately. Ship refrigerated or frozen. Ideal sampling time: just before the next dose is due at steady state.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks

D008	DIGOXIN		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Separate serum within 2 hours of collection. Ship refrigerated or frozen. Trough: just before the next dose is due at steady state. Peak : 6–8 hours after dose.	Chemiluminescent Immunoassay		Daily
R103	5α-DIHYDROTESTOSTERONE; 5α-DHT		2 mL (1 mL min.) serum from 1 SST. Specify age and sex on test request form. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same Day
RW063	DIHYDROTESTOSTERONE (DHT), FREE		5 mL (2.5 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	LC-MS / MS		Sample by 7th of the month; report after 2–3 weeks.
RW064	DISOPYRAMIDE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum from cells immediately. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
H245	DISSEMINATED INTRAVASCULAR COAGULATION (DIC) PROFILE	*CBC *Prothrombin studies *APTT *Fibrinogen* D-Dimer	3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. AND 6 mL Whole blood in 2 Blue Top (Sodium Citrate) tubes. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection, Electrical Impedance	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Daily

H023	DLC (DIFFERENTIAL LEUCOCYTE COUNT), AUTOMATED		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electronic Cell Cytometry, VCS & Microscopy	Also see ABSOLUTE COUNTS.	Daily
P084	DLC (DIFFERENTIAL LEUCOCYTE COUNT), FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Manual count		Daily
N107	DNA ANALYSIS FOR KINSHIP TESTING		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube for each Parent and Child. Ship refrigerated. DO NOT FREEZE. Paternity Consent form with photograph and identification document for each individual is mandatory.	PCR, STR analysis	For 2 individuals price is Rs.10,000/-	Sample Daily by 11 am; Report 10 days
N106	DNA ANALYSIS FOR PATERNITY TESTING		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube for each Parent and Child. Ship refrigerated. DO NOT FREEZE. Paternity Consent form with photograph and identification document for each individual is mandatory.	PCR, STR analysis	For 2 individuals price is Rs.10,000/-	Sample Daily by 11 am; Report 10 days
N115	DNA ANALYSIS FOR PATERNAL Y CHROMOSOME		3 mL (2 mL min.) Whole blood from 1 Lavender Top (EDTA) tube for each individual. Ship refrigerated. DO NOT FREEZE. Paternity Consent form with photograph and identification document for each individual is mandatory.	PCR, STR analysis	For 2 individuals price is Rs.20,000/-	Sample Daily by 11 am; Report 12 days
G140	DNPH, URINE		10 mL (5 mL min.) aliquot of first morning urine OR post feed urine. Ship refrigerated or frozen.	Biochemical		Sample Daily by 5 pm; Report Same Day

RW065	DOXEPIN	Includes metabolite Nordoxepin.	3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: >12 hours post oral dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
OS130	DPD Gene Mutations (5-FU toxicity) detection		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Sat 9 am; Report 10 days
R122	DPYD (DEOXYPYRIDINOLINE), URINE		20 mL (10 mL min.) aliquot of first morning urine. Ship refrigerated or frozen. Wrap the specimen in aluminium foil to protect from light.	Chemiluminescent Immunoassay	Also see ALKALINE PHOSPHATASE ISOENZYMES; ALKALINE PHOSPHATASE, BONE SPECIFIC; COLLAGEN CROSS-LINKED N-TELOPEPTIDE; OSTEOPOROSIS PANEL. Deoxypyridinoline, when used for monitoring therapy, should be collected at the same time on each occasion as for the Ç	Sample by Mon 4 pm; Report Same Day
OS132	DRPLA GENE ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days

V617	DRUGS OF ABUSE, 7 DRUGS URINE SCREEN	* Opiates *Benzodiazepine *Cannabinoids* Cocaine & metabolites *Amphetamine* Barbiturate * Phencyclidine (PCP)	10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily
N116	DUCHENNE / BECKER MUSCULAR DYSTROPHY (DMD / BMD) GENE MUTATION		3 mL (2 mL min.) Whole blood in 1 lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	MLPA		Sample by Mon 11 AM; Report Sat
S036	ECHINOCOCCUS (HYDATID SEROLOGY) IgG, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see ECHINOCOCCUS DETECTION FOR SCOLICES.	Sample by Tue / Sat 9 am; Report Same day.
P096	ECHINOCOCCUS DETECTION FOR SCOLICES IN HYDATID CYST FLUID		1 mL (0.5 mL min.) Cyst fluid in a sterile screw capped container. Ship refrigerated or frozen.	Microscopy		Sample Daily by 6 pm; Report Next day
Z218	ECHOVIRUS ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
S177	ECHOVIRUS ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day

S125	ECHOVIRUS ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
N072	EGFR, MUTATIONAL ANALYSIS		Submit formalin fixed paraffin embedded block containing tumor tissue. Ship at room temperature. Provide a copy of the Histopathology report. Indicate site of biopsy and provide clinical history	PCR, Fragment Analysis, Sequencing	Paraffin block must contain minimum 10% tumor cells	Sample by Tue 11 am; Report Mon
Z380	ELECTROLYTES, RANDOM URINE	*Sodium *Potassium *Chloride	10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
J154	ELECTRON MICROSCOPY (TEM)		Submit Kidney/ Muscle/ Nerve /Other tissue biopsies in 3% Buffered Glutaraldehyde solution available from LPL. Ship refrigerated.	Transmission Electron microscopy		Sample Daily by 6 pm; Report 7 days
Z383	ELECTROLYTES, 24 HOUR URINE	*Sodium *Potassium *Chloride	10 mL (5 mL min.) aliquot of 24-hour urine in a sterile screw capped container. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
Z020	ELECTROLYTES, SERUM	*Sodium *Potassium *Chloride	3 mL (1.5 mL min.) serum from 1 SST. Avoid hemolysis. Ship refrigerated or frozen. Separate serum within 1 hour of collection.	Indirect Ion Selective Electrode		Daily
S078	ENDOMYSIAL ANTIBODY, IgA, IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	For SERIAL DILUTIONS use test code S105.	Sample Daily by 11 am; Report Same Day
S105	ENDOMYSIAL ANTIBODY, IgA, IFA IN DILUTIONS		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 11 am; Report Next day

Z819	ENHANCED LIVER FIBROSIS (ELF) PANEL	*Hyaluronic acid (HA) * Amino terminal Propeptide type III Procollagen (PIIINP) *Tissue inhibitor Metalloproteinase 1 (TIMP-1) * ELF score	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is mandatory.	Chemiluminescent Microparticle Immunoassay		Sample by Mon / Thu 4 pm; Report Next day
M210	ENTEROVIRUS ANTIGEN, STOOL		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pmR Same day
M204	ENTEROHEMOR RHAGIC E.COLI, VEROTOXIN 1 & 2 ANTIGEN		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pmR Same day
N125	ENTEROVIRUSE S, QUALITATIVE, PCR		3ml (2mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 3ml (2mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
S195	EPIDERMAL (PEMPHIGUS) ANTIBODY, IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	For SERIAL DILUTIONS use test code S198.	Sample Daily by 11 am; Report Same day
S198	EPIDERMAL (PEMPHIGUS) ANTIBODY, IFA IN DILUTIONS		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 11 am; Report Nextday

N127	EPSTEIN BARR VIRUS (EBV), QUALITATIVE PCR		3 ml (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
RW067	EPSTEIN-BARR VIRUS, ANTIBODY TO EARLY ANTIGEN (D), IgG		2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
S070	EPSTEIN-BARR VIRUS ANTIBODY TO VIRAL CAPSID ANTIGEN (VCA), IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed / Sat 9 am; Report Same day
Z172	EPSTEIN BARR VIRUS ANTIBODIES TO VIRAL CAPSID ANTIGEN (VCA), IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed / Sat 9 am; Report Same day
S069	EPSTEIN-BARR VIRUS ANTIBODY TO VIRAL CAPSID ANTIGEN (VCA), IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed / Sat 9 am; Report Same day

RW068	EPSTEIN-BARR VIRUS, ANTIBODY TO NUCLEAR ANTIGEN (EBNA), IgG		2 mL (0.5 mL min.) serum from 1 SST. Separate serum immediately. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
R018	ERYTHROPOIETIN; EPO		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample by Mon / Thu 4 pm; Report Same day
H009	ESR; ERYTHROCYTES SEDIMENTATION RATE		3 mL (1.5 mL min.) fasting whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Capillary Photometry / Westergren		Daily
R019	ESTRADIOL (E2), SERUM		3 mL (1 mL min.) serum from 1 SST. Specify age and sex on test request form. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
R101	ESTRIOL UNCONJUGATED / FREE (μ E3)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Also see MATERNAL SERUM SCREEN3 and 4.	Daily
J009	ESTROGEN RECEPTOR; ER	Includes pathologist review for the presence of malignant cells.	Submit tumour tissue in 10% Formal saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate / Cytology Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report.	Immunohistochemistry	Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report Block / Section / Smears : 5 days Tissues : 7 days
RW069	ESTROGEN, TOTAL, SERUM		3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Extraction, Radioimmunoassay		Sample by 7th of the month; Report after 2–3 weeks

RW070	ETHOSUXIMIDE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Ideal sampling time: Trough level, just prior to next dose.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
G161	EVEROLIMUS		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. Ideal sampling time: 1 hour before next dose (12 Hour Trough)	LC-MS/MS		Sample Daily by 5 pm. Report Same Day
Z011	EXECUTIVE OFFICERS PANEL	*Glucose (F) *Creatinine *Uric Acid *Cholesterol *Triglycerides *SGOT *SGPT *GGTP *Bilirubin, Total *Protein, Total *Albumin *A:G Ratio *Alkaline Phosphatase *CBC *Blood Group & Rh factor *TSH Ultrasensitive	2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube for Glucose, F AND 4 mL (2 mL min.) serum from 1 SST AND 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. 12 hours overnight fasting is mandatory.	Spectrophotometry, Electrical Impedance, VCS, EMT		Daily
S118	EXTRACTABLE NUCLEAR ANTIGENS (ENA), QUALITATIVE PROFILE	Combined group screening test for 15 antigens	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Line Immunoassay	Also see EXTRACTABLE NUCLEAR ANTIGENS (ENA), QUANTITATIVE PROFILE	Sample Mon through Sat by 9 am; Report Same day

Z101	EXTRACTABLE NUCLEAR ANTIGENS (ENA), QUANTITATIVE PROFILE	Individual antibodies to *Sm *U1RNP *SS-A/Ro *SS- B/La *Jo-1 *Scl- 70 *Centromere	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue / Fri 9 am; Report Same day
G127	FABRY DISEASE, QUANTITATIVE, BLOOD		10 mL (5–7 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium Heparin) tubes. Ship refrigerated. DO NOT FREEZE. Clinical details must accompany sample.	Enzyme Assay	This tests deficiency of the enzyme Alpha Galactosidase A.	Sample Daily by 4 pm; Report 4 Days
H173	FACTOR II, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Thu 9 am; Report Next day

H174	FACTOR V, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Thu 9 am; Report Next day
H240	FACTOR VIII FUNCTIONAL		3 mL Whole blood in 1 Blue Top(Sodium Citrate) tube. Mix thoroughlyby inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-Optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Also see COAGULATION PC	Sample by Wed / Sat 9 am; Report Fri/ Tue

H176	FACTOR X, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Thu 9 am; Report Next day
H177	FACTOR XI, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician	Sample by Thu 9 am; Report Next day

H179	FACTOR XII, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Thu 9 am; Report Next day
H041	FACTOR XIII (CLOT SOLUBILITY), FUNCTIONAL, QUALITATIVE		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Clot Dissolution with 5M Urea	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample Daily by 1 pm; Report Next day

H175	FACTOR IX, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Also see HEMOPHILIA PANEL.	Sample by Wed / Sat 9 am; Report Fri/ Tue
H209	FACTOR V LEIDEN MUTATION ANALYSIS	*Factor V Leiden Mutation *Prothrombin Gene Mutation *MTHFR Gene Mutation	3 mL (2 mL min.) Whole blood in 1Lavender Top (EDTA) tube. Shiprefrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Thu 11 am; ReportWed / Sat

H178	FACTOR VII FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Thu 9 am; Report Next day
H247	FACTOR VIII INHIBITOR		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection, EIA	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Wed / Sat 9 am; Report Fri/ Tue

XX025	FANCONI'S ANEMIA, STRESS CYTOGENETICS		5 mL (3 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE	Stress Cytogenetics		Sample Daily by 6 pm; Report 7 Days
Z837	FATTY LIVER INDEX	*GGTP *Triglycerides	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is mandatory. Record patient height, weight, waist circumference on test request form.	IFCC,GPO-POD		Sample by Mon through Sat 9 am; Report same day
Z124	FEBRILE AGGLUTININS PANEL Combine d Brucella, Widal and Weil Felix tests panel	*Brucella abortus *Brucella melitensis *Salmonella paratyphi Type AH and Type BH *Salmonella typhi H & O *OX2 *OXK *OX19	3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Latex & Tube Agglutination	Also see BRUCELLA ANTIBODY, IgG& IgM	Sample Daily by 6 pm; Report Next day
S247	FECAL ELASTASE		Submit 5g (2g min.) Stool in a clean dry container. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day.
S231	FECAL CALPROTECTIN		Submit 5g (2g min.) Stool in a clean dry container. Ship refrigerated. DO NOT FREEZE.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
RW071	FELBAMATE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Ideal sampling time: 1 hour before next dose. Do not use SST gel barrier tubes.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks

G141	FERRIC CHLORIDE TEST, URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Chemical	Also see AMINO ACID PANELS, QUALITATIVE AND QUANTITATIVE.	Sample Daily by 5 pm; Report Same day
H060	FERRITIN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
Z390	FEVER PANEL 1	*Hemogram *MP *Widal *Urine R/E	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 2 mL (1 mL min.) serum from 1 SST AND 10 mL (5 mL min.) aliquot of random urine. Ship refrigerated. Overnight fasting is preferred	Electrical impedance, VCS, Microscopy, Agglutination, Dipstick		Daily
Z391	FEVER PANEL 2	*Hemogram *MP *Widal *Urine R/E *SGPT * Aerobic Blood Culture	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 2 mL (1 mL min.) serum from 1 SST AND 10 mL (5 mL min.) aliquot of random urine. Ship refrigerated AND 10 mL Blood collected aseptically in Special Plus Aerobic bottle available from LPL. Ship at 18–22°C. Overnight fasting is preferred.	Electrical impedance, VCS, Microscopy, Agglutination, Dipstick, IFCC, Conventional culture & Automated sensitivity		Sample Daily by 6 pm; Report Interim : 24 hrs & 48 hrs Final : 6 Days
B088	FIBRINOGEN, CLOTTING ACTIVITY		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection		Daily

H026	FIBRINOGEN DEGRADATION PRODUCTS (FDP), BLOOD		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Latex Agglutination	Also see D-DIMER & DEEP VEIN THROMBOSIS PANEL	Daily
RW206	FIBROBLAST GROWTH FACTOR; FGF2		2 mL (1 mL min.) Plasma from 1 Lavender Top (EDTA) tube. Separate plasma within 6 hours of collection in a cold centrifuge. Freeze immediately. Ship frozen. DO NOT THAW.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
Z800	FIBROMYALGIA PANEL	*CMV, IgG *Selenium, Serum * RA factor * Anti Nuclear Antibody (ANA) * Creatine Phosphokinase (CPK), TSH, Ultrasensitive *C-Reactive Protein	4 mL (3 mL min.) serum from 1 SST AND 2 mL (1 mL min.) serum from 1 White top (Z-No Additive) tube available from LPL. Ship refrigerated or frozen.	ICPMS, CLIA, Immunoturbidimetry , IFCC, EIA	If Gadolinium or Iodine containing contrast media has been administered, specimen should not be collected for 96 hours Use only specified containers available from LPL.	Sample Daily by 6 pm; Report Next day
S170	FILARIA ANTIGEN		2 mL (1 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE.	Immunochromatography		Daily

OS135	FIP1L1-PGDFRA GENE REARRANGEMENT DETECTION		4 ml (2 ml min.) Whole blood / Bone marrow from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	PCR		Sample by Tue 11 am; Report 10 days
XX031	FISH - ALK-1 Breakapart Rearrangement		Formalin fixed paraffin embedded tissue block. Ship at room temperature	FISH		Sample Daily by 4 pm; Report 7 Days.
XX015	FISH - Amnio one probe (Trisomy 21) / Down syndrome		5 mL Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Clinical history mandatory. Consent form for Prenatal Genetic Testing is mandatory.	FISH	Consent form for Prenatal Genetic Testing is mandatory.	Sample Daily by 4 pm; Report 7 Days.
XX017	FISH - Amnio three probes (Trisomy 18, X, Y)		5 mL Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Clinical history mandatory. Consent form for Prenatal Genetic Testing is mandatory.	FISH	Consent form for Prenatal Genetic Testing is mandatory.	Sample Daily by 4 pm; Report 7 Days.
XX016	FISH - Amnio two probes (Trisomy 13; 21)		5 mL Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Clinical history mandatory. Consent form for Prenatal Genetic Testing is mandatory.	FISH	Consent form for Prenatal Genetic Testing is mandatory.	Sample Daily by 4 pm; Report 7 Days.
XX036	FISH - Aneuploidy detection, Products of Conception (POC) using chromosomes 13,18, 21, X and Y		Submit Placenta (villi) or Fascia lata, diaphragm, tendon, skin, tissue from internal organs (if fresh), chest wall cartilage (particularly if macerated) in normal saline. Ship immediately at 18-22°C. Formalin fixed specimen not acceptable. Give clinical history on Chromosome Analysis Test Request Form.	FISH		Sample Daily by 4 pm; Report 7 Days.

XX001	FISH - BCR / ABL or Philadelphia translocation		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX024	FISH - Bladder cancer recurrence: Chromosomes 3,7, 17 & p16 gene on chromosome 9		Submit 10 mL Urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX023	FISH - CLL Panel	*del11q23(ATM) *del17p13.1(TP53) *del13q14.3 *trisomy 12 *IGH Gene breakapart	5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX029	FISH - C-MYC		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX005	FISH - del 13q or LSI D 13S319		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX012	FISH - EGFR Gene amplification		Formalin fixed paraffin embedded tissue block. Ship at room temperature	FISH		Sample Daily by 4 pm; Report 7 Days.

XX051	FISH- EQUIVOCAL HER2 (ERBB2)		Formalin fixed paraffin embedded tissue block / three 4µ sections on Poly-L-Lysine coated slides. Ship at room temperature. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 5 days
XX048	FISH-EWING SARCOMA (EWS), 22q12 (EWSR1) REARRANGEMENT		Formalin fixed paraffin embedded tissue block / three 4µ sections on Poly-L-Lysine coated slides. Ship at room temperature. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX045	FISH- Eosinophilic Leukemia Panel (PDGFRA, PDGFRB, FGFR1 and CFBF)		5 mL (3 mL min.) Bone marrow OR 5 mL (3 mL min.) Whole blood in 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX044	FISH-FGFR1 (8p12) GENE REARRANGEMENT		5 mL (3 mL min.) Bone marrow OR 5 mL (3 mL min.) Whole blood in 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX011	FISH - HER2 (ERBB2) amplification		Formalin fixed paraffin embedded tissue block. Ship at room temperature	FISH		Sample Daily by 4 pm; Report 7 Days.
XX004	FISH - inv (16) or LSI CFBF		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX022	FISH - MDS Panel - Chromosomes 5q, 7q, 8q & 20q.		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.

XX050	FISH-MET (7q31) AMPLIFICATION		Formalin fixed paraffin embedded tissue block / three 4µ sections on Poly-L-Lysine coated slides. Ship at room temperature. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX027	FISH - Multiple Myeloma (5 Probes)	Del 13q14.3,del 17p13.1, t(11;14) CCND1/IGH, t(4;14) FGFR3/IGH, and t(14;16) IGH/MAF	5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX046	FISH-MYCN AMPLIFICATION		Submit formalin fixed paraffin embedded tissue block / three 4 µ sections on Poly-L-Lysine coated slides. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX038	FISH - Microdeletion detection for Williams syndrome		5 mL (3 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX032	FISH - Multiple Myeloma with CD138 positive selection cells (5 Probes)	Del 13q14.3,del 17p13.1, t(11;14) CCND1/IGH, t(4;14) FGFR3/IGH, and t(14;16) IGH/MAF	5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
Z787	FISH - NSCLC (Non-small cell lung cancer) ALK1 & ROS1		Formalin fixed paraffin embedded tissue block / three 4u sections on Poly-L-Lysine coated slides. Ship at room temperature.	FISH		Sample Daily by 4 pm; Report 10 Days

XX010	FISH - Opposite Sex BMT (XX / XY)		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX018	FISH - Prader Willi Syndrome SNRPN		5 mL (3 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX030	FISH - Prenatal screening panel (13,18,21,X and Y)		6 mL Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Clinical history mandatory. Consent form for Prenatal Genetic Testing is mandatory.	FISH	Consent form for Prenatal Genetic Testing is mandatory.	Sample Daily by 4 pm; Report 7 Days.
XX042	FISH-PDGFR alpha (4q12) GENE REARRANGEMENT		4 mL (3 mL min.) Bone marrow OR 4 mL (3 mL min.) Whole blood in 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX043	FISH-PDGFR beta (5q33) GENE REARRANGEMENT		4 mL (3 mL min.) Bone marrow OR 4 mL (3 mL min.) Whole blood in 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX013	FISH - 1p and 19q Codeletion		Formalin fixed paraffin embedded tissue block. Ship at room temperature	FISH		Sample Daily by 4 pm; Report 7 Days.
XX026	FISH - 17 p		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.

XX014	FISH - 22q deletion or LSI Di George / VCFS		5 mL (3 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX007	FISH - 11q 23 or LSI MLL Gene breakapart		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX009	FISH - 14q 32.3 or LSI IGH Gene breakapart		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX049	FISH- ROS1 (6Q22) REARRANGEMENT		Formalin fixed paraffin embedded tissue block / three 4µ sections on Poly-L-Lysine coated slides. Ship at room temperature. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX047	FISH - SYNOVIAL SARCOMA (SS18 / SYT) 18q11.2 GENE REARRANGEMENT		Submit formalin fixed paraffin embedded tissue block / three 4 µ sections on Poly-L-Lysine coated slides. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 days
XX028	FISH - SRYGENE		5 mL (3 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX003	FISH - t(8;21) or LSI ETO/AML1		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.

XX039	FISH - t(8;14) or MYC/IGH for Burkitt's & Non-Hodgkin lymphoma		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX006	FISH - t(11;14) or LSI IGH/CCND1		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX008	FISH - t (12;21) or LSI TEL/AML1		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX002	FISH - t (15;17) or LSI PML/RARA		5 mL (3 mL min.) whole blood OR 4 mL (2 mL min.) Bone Marrow from 2 Green Top (Sodium Heparin) tubes. Ship at 18-22°C. DO NOT FREEZE. Clinical history mandatory.	FISH		Sample Daily by 4 pm; Report 7 Days.
XX037	FISH - Trisomy 21 / Down syndrome		5 mL (3 mL min.) Whole blood / Cord blood from 2 Green Top (Sodium Heparin) tubes OR 5 mL (3 mL min.) Amniotic fluid in a sterile screw capped container OR 2g Chorionic villus biopsy in normal saline. Ship at specified temperatures.Consent form for Prenatal Genetic Testing is mandatory.	FISH	Consent form for Prenatal Genetic Testing is mandatory.	Sample Daily by 4 pm; Report 7 Days.
RW072	FLECAINIDE		3 mL (1.5 mL min) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Ideal sampling time: 1 hour before next dose.	LC-MS/MS		Sample by 7th of the month; Reportafter 2–3 weeks

H151	FLOW CYTOMETRY, CD 25		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H232	FLOW CYTOMETRY, CD 123		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H149	FLOW CYTOMETRY, CD 22		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H222	FLOW CYTOMETRY, CD 79a		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H142	FLOW CYTOMETRY, CD 11c		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H143	FLOW CYTOMETRY, CD 13		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H229	FLOW CYTOMETRY, CD 138		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H144	FLOW CYTOMETRY, CD 14		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H145	FLOW CYTOMETRY, CD 15		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H156	FLOW CYTOMETRY, CD 41a		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H231	FLOW CYTOMETRY, CD 64		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H139	FLOW CYTOMETRY, CD 7		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H223	FLOW CYTOMETRY, CD 117		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H138	FLOW CYTOMETRY, CD 5		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H162	FLOW CYTOMETRY, HLA DR		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H237	FLOW CYTOMETRY, CD 34 (STEM CELL) ENUMERATION		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) Bone marrow in 1 Green Top (Sodium heparin) tube OR 1 mL(0.5mL min) Cord blood / Apheresis sample in a sterile screw capped container. Ship immediately at 2 –8°C. DO NOT FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H167	FLOW CYTOMETRY, IgG Heavy Chain (Cytoplasmic)		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H141	FLOW CYTOMETRY, CD 10 (CALLA)		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H181	FLOW CYTOMETRY, CD 103		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
Z166A	FLOW CYTOMETRY, CD 16 + 56		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H147	FLOW CYTOMETRY, CD 19		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
Z219A	FLOW CYTOMETRY, CD 4		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H236	FLOW CYTOMETRY, CD 79b		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H203	FLOW CYTOMETRY, CD MPO		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H160	FLOW CYTOMETRY, FMC-7		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H136	FLOW CYTOMETRY, CD 3		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H154	FLOW CYTOMETRY, CD 34		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H230	FLOW CYTOMETRY, CD 36		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H135	FLOW CYTOMETRY, CD 2		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H157	FLOW CYTOMETRY, CD 45		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H161	FLOW CYTOMETRY, GLYCOPHORIN		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H153	FLOW CYTOMETRY, CD 33		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H155	FLOW CYTOMETRY, CD 38		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
Z147A	FLOW CYTOMETRY, - CD 8		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H172	FLOW CYTOMETRY, TdT		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H228	FLOW CYTOMETRY, ZAP-70		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H158	FLOW CYTOMETRY, CD 61		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H168	FLOW CYTOMETRY, IgG Heavy Chain (Surface)		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H169	FLOW CYTOMETRY, IgM Heavy Chain (Cytoplasmic)		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H170	FLOW CYTOMETRY, IgM Heavy Chain (Surface)		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day

H171	FLOW CYTOMETRY, Kappa & Lambda Light Chains		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H148	FLOW CYTOMETRY, CD 20		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H234	FLOW CYTOMETRY, CD 1a		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
H150	FLOW CYTOMETRY, CD 23		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone marrow. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am. Report Same day
N052	FLT3 GENE MUTATION		5 mL (3 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory.	Real Time PCR, Fragment Analysis		Sample by Mon / Thu 11 am ; Report Wed / Sat

P012	FLUID EXAMINATION, ROUTINE, CSF/ NASAL FLUID		2 mL (1mL min.) Cerebrospinal fluid / Nasal fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Biochemistry and microscopy of Giemsa stained smears		Daily
P013	FLUID EXAMINATION, ROUTINE, BODY FLUIDS		5 ml (2 mL min.) fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Specify type of fluid on Test Request Form and container.	Biochemistry and microscopy of Giemsa stained smears		Daily
RW073	FLUOXETINE	Includes metabolite Norfluoxetine.	3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Ideal sampling time: >12 hours post oral dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
Z781	FMF FIRST TRIMESTER SCREEN		2ml(1 ml min.) serum form 1SST. Ship refrigerated or frozen. Provide maternal Date of Birth (dd/mm/yy), LMP, USG report between 11-13 weeks gestation including CRL, NT & Nasal Bone, number of fetuses; Diabetic status, body weight; IVF, Smoking & Previous h/o Trisomy 21. Enclose Patient demographic details in Maternal Serum Screen Request Form. Valid between 9-13 weeks gestation (Ideal 10-13 weeks).	ECLIA	The test is done using FMF approved kits. Enclose demographic details in maternal serum screen request form	Sample by Tues / Thu / Sat 9 am; Report same day
R048	FOLATE (FOLIC ACID), RBC	Includes Serum Folate.	2 mL (1 mL min.) whole blood from 1 Lavender Top (EDTA) tube AND 2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Wrap tubes in aluminium foil to protect from light. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample Daily by 2 pm; Report Same Day
R046	FOLATE (FOLIC ACID), SERUM		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Wrap tube in aluminium foil to protect from light. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily

A110	FOOD INTOLERANCE TEST		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Protein Microarray		Sample Daily by 1 pm; Report 2 days
B191	FRACTIONAL EXCRETION OF SODIUM	*Creatinine, serum * Creatinine, urine *Sodium, serum *Sodium, urine	3 mL (2 mL min.) serum from 1 SST AND 10 mL (5 mL min) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	Jaffe's reaction, Indirect ISE		Daily
OS137	FRAGILE X (FMR1) MUTATION SCREEN		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Tue 11 am; Report 10 days
B151	FREE BETA HCG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Also see MATERNAL SERUM SCREEN 2.	Daily
OS138	FREIDREICHS ATAXIA MUTATION ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days
B071	FRUCTOSAMINE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Spectrophotometry	Also see GLYCOSYLATED/GLYCATED HEMOGLOBIN, LIPID PROFILE & MICROALBUMIN URINE.	Daily
R080	FSH; FOLLICLE STIMULATING HORMONE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age, sex and day of menstrual cycle on test request form.	Chemiluminescent Immunoassay	In females, preferred sampling time is Day 2 / Day 3 of the menstrual cycle.	Daily

Z079	FSH & LH STIMULATION BY GONADOTROPI NRELEASING HORMONE (GnRH)		2 mL (1 mL min.) serum from 1 SST for each timed specimen . Draw fasting baseline specimen; administer GnRH 2.5 ug/kg body weight (100 ug maximum) as IV bolus; draw additional specimens at 30, 45 & 60 min for FSH and LH. Ship refrigerated or frozen. All containers must be clearly marked with time drawn and submitted with one test request form.	Chemiluminescent Immunoassay	refer to Y Codes:Y050 BasalY051 30 min.Y052 45 min.Y053 60 min.Y054 Customized	Daily
Z126	FSH & LH		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age, sex & day of menstrual cycle on test request form.	Chemiluminescent Immunoassay	In females, preferred sampling time is Day 2 / Day 3 of the menstrual cycle.	Daily
P021	FUNGUS EXAMINATION, ROUTINE, KOH CALCOFLUOR WHITE PREPARATION		Submit Skin, Nail scrapings, Hair (including root), Cerebrospinal fluid, in a sterile screw capped container OR 2 mL whole blood / bone marrow in 1 Green Top (Sodium Heparin) tube. Ship at 18–22°C. Submit Body fluid, Sputum, Bronchial washings / Tissue in a sterile screw capped container OR Oral/Genital/wound swabs in 0.5 mL sterile normal saline. Ship refrigerated.	Fluorescent microscopy		Sample Daily by 3.30 pm; Report Same day
RW074	GABAPENTIN		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum immediately. Ship refrigerated or frozen. Ideal sampling time: 2 hours after last dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
B150	GAD-65 (GLUTAMIC ACID DECARBOXYLAS E- 65), IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Sat 9 am: Report Same day

G119	GALACTOSEMIA PANEL 1	*Galactosemia Screening *Galactosemia Classical (Transferase)	4 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube AND 4 mL (2 mL min.) CONTROL blood from 1 Green Top (Sodium Heparin) tube. Ship refrigerated. DO NOT FREEZE. Avoid sample collection for 60 days post transfusion. Clinical and drug history must accompany sample.	Spot Test, Enzyme assay	CONTROL blood should not be taken from blood relatives.	Sample Daily by 4 pm; Report 5 Days
G120	GALACTOSEMIA PANEL 2	*Galactosemia Screening *Galactosemia Classical (Transferase) *Galactose Plasma	4 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube AND 4 mL (2 mL min.) CONTROL blood from 1 Green Top (Sodium Heparin) tube AND 2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Ship refrigerated. DO NOT FREEZE. Avoid sample collection for 60 days post transfusion. Clinical and drug history must accompany sample.	Spot Test, Enzyme assay	CONTROL blood should not be taken from blood relatives.	Sample Daily by 4 pm; Report 5 days
G121	GALACTOSEMIA PANEL 3	*Galactosemia Classical (Transferase)*Galactosemia (Epimerase)	4 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship refrigerated. DO NOT FREEZE. Avoid sample collection for 60 days post transfusion. Clinical and drug history must accompany sample.	Spot Test, Enzyme assay		Sample Daily by 4 pm; Report 5 days
S242	GALACTOMANNAN (ASPERGILLUS ANTIGEN), SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	EIA		Sample Mon through Sat by 12 Noon; Report Same day
R164	GALACTOSEMIA , NEWBORN SCREEN		1 drop of heel prick blood each on 3 spots of filter paper available from LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen.	Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm; Report 2 days

G118	GALACTOSEMIA (EPIMERASE), QUANTITATIVE, BLOOD		4 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship refrigerated. DO NOT FREEZE. Avoid sample collection for 60 days post transfusion. Clinical and drug history must accompany sample.	Enzyme assay		Sample Daily by 4 pm; Report 5 Days
OS140	GALACTOSEMIA (GALT) GENE MUTATION DETECTION		4 mL (2 mL min.) Whole Blood from 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Tue 11am; Report 11 days
G117	GALACTOSEMIA CLASSICAL(TRANSFERASE), QUANTITATIVE, BLOOD		4 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship refrigerated. DO NOT FREEZE. Avoid sample collection for 60 days post transfusion. Clinical and drug history must accompany sample.	Enzyme assay		Sample Daily by 4 pm; Report 5 Days
G115	GALACTOSE -1-PHOSPHATE, QUANTITATIVE, BLOOD		4 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Immerse specimen immediately in ice water. Ship refrigerated. DO NOT FREEZE. Clinical and drug history must accompany sample.	Fluorometry	Sample must be kept cold as Galactose -1-phosphate levels fall upto 50% within 4 hours if kept at room temperature.	Sample Daily by 4 pm; Report 5 Days
G116	GALACTOSE, QUANTITATIVE, PLASMA		2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Ship refrigerated or frozen. Clinical / drug history & diet must accompany sample.	Fluorometry		Sample Daily by 4 pm; Report 3 days

G114	GALACTOSEMIA SCREENING, BLOOD		4 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube AND 4 mL (2 mL min.) CONTROL blood from 1 Green Top (Sodium Heparin) tube. Ship refrigerated. DO NOT FREEZE. Avoid sample collection for 60 days post transfusion. Clinical and drug history must accompany sample.	Beutler Fluorescence Spot Test	CONTROL blood should not be taken from blood relatives.	Sample Daily by 4 pm; Report Next day
OS141	GANCICLOVIR RESISTANCE DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Tue 11 am; Report 9 days
RW075	GANGLIOSIDE GD1a ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting preferred.	Enzyme Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
RW076	GANGLIOSIDE GD1a ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting preferred.	Enzyme Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
RW080	GANGLIOSIDE GQ1B ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
S234	GANGLIOSIDE ANTIBODY PROFILE, IgM	*GM1 *GM2 *GM3 *GD1a *GD1b *GT1b *GQ1b	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoblot		Sample by Tue 9 am; Report Same day
RW077	GANGLIOSIDE GD1b ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting preferred.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
RW078	GANGLIOSIDE GD1b ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting preferred.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks

RW079	GANGLIOSIDE ANTIBODIES, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting preferred.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks.
S233	GANGLIOSIDE ANTIBODY PROFILE, IgG	*GM1 *GM2 *GM3 *GD1a *GD1b *GT1b *GQ1b	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoblot		Sample by Tue 9 am; Report Same day
R020	GASTRIN		2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW. It is recommended that patient discontinues proton pump inhibitor drugs like Omeprazole, Esomeprazole, Pantoprazole, Lansoprazole & Rabeprazole for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Overnight fasting is mandatory.	Chemiluminescent Immunoassay		Sample by Mon / Thu 4 pm; Report Same day
Z061	GASTRIN STIMULATION BY SECRETIN		2 mL (1 mL min.) serum from 1 SST for each timed specimen. Draw baseline specimen; administer Secretin 2 U/kg body weight IV as a bolus injection, and draw additional post secretin specimens at 5, 10, 15, 20, 25 and 30 minutes for Gastrin. Ship Frozen. DO NOT THAW. All containers must be clearly marked with time drawn and submitted with one test request form. Overnight fasting is mandatory	Chemiluminescent Immunoassay	Refer to Y codes: Y058 Basal Y059 5 min. Y060 10 min. Y061 15 min. Y062 20 min. Y122 25 min. Y123 30 min. Y063 Customized	Sample by Mon / Thu 4 pm; Report Same day
G125	GAUCHER DISEASE, QUANTITATIVE, BLOOD		10 mL (7.5 mL min.) whole blood in 3 Lavender Top (EDTA) OR Green Top (Sodium heparin) tubes. Ship refrigerated. DO NOT FREEZE. Clinical details must accompany sample.	Enzyme assay	This tests the deficiency of enzyme Beta Glucosidase.	Sample Daily by 4 pm; Report 4 Days

N076	GENETIC MAPPING FOR WHOLE GENOME		Submit 3 mL (2 mL min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Consent & Clinical Information forms available from LPL are mandatory.	SNP Genotyping using Microarray	Genetic counselling will be provided with prior appointment.	Sample Daily by 9 am; Report 8 weeks
N077	GENETIC MAPPING FOR ONCOLOGY / CANCER		Submit 3 mL (2 mL min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Consent & Clinical Information forms available from LPL are mandatory.	SNP Genotyping using Microarray	Genetic counselling will be provided with prior appointment.	Sample Daily by 9 am; Report 8 weeks
G173	GENETIC COUNSELLING		Patient & Relative to report to National Reference Laboratory, Rohini, New Delhi OR Telephonically with prior appointment.	Counselling		Daily
N078	GENETIC MAPPING FOR CARDIOVASCULAR DISORDERS		Submit 3 mL (2 mL min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Consent & Clinical Information forms available from LPL are mandatory.	SNP Genotyping using Microarray	Genetic counselling will be provided with prior appointment.	Sample Daily by 9 am; Report 8 weeks
N079	GENETIC MAPPING FOR GYNECOLOGICAL DISORDERS		Submit 3 mL (2 mL min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Consent & Clinical Information forms available from LPL are mandatory.	SNP Genotyping using Microarray	Genetic counselling will be provided with prior appointment.	Sample Daily by 9 am; Report 8 weeks
Z518	GENEXPERT.COM	*Genexpert with Rifampicin resistance * Culture AFB, Rapid	6 mL (3 mL min.) first morning Sputum OR 6 mL (3 mL min) concentrated sediments prepared from induced or expectorated Sputum OR 6 mL (3 mL min) BAL fluid OR 2 mL (1 mL min.) CSF OR 6 mL (3 mL min) Body fluids / Aspirates OR 2 mL (1 mL min.) Pus in a sterile screw capped container OR Fresh tissue biopsies including endometrial aspirate / curettage in sterile normal saline Ship refrigerated. DO NOT FREEZE.	Real Time PCR, Automated fluorescent, ICT	Blood & Swabs are not acceptable	Sample Daily by 3.30 pm; Report: Interim: 1 day, 14 days & 4 weeks Final: 6 weeks

N061	GENEXPERT MTB WITH RIFAMPICIN RESISTANCE, QUALITATIVE, PCR		3 mL (1 mL min.) Sputum OR 3 mL (1 mL min) concentrated sediments prepared from induced or expectorated Sputum OR 3 mL (1 mL min) BAL fluid OR 2 mL (1 mL min.) CSF OR 5 mL (3 mL min) Body fluids/ Aspirates OR Pus in a sterile screw capped container OR Fresh tissue biopsies including endometrial aspirate / curettage in sterile normal saline. Ship refrigerated. DO NOT FREEZE.	Real Time PCR	Blood & Swabs are not acceptable	Sample Daily by 11 am; Report Same Day
Z813	GERM CELL TUMOR PANEL	*AFP * HCG, Total *LDH	3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CMIA, Spectrophotometry		Daily
B168	GFR (GLOMERULAR FILTRATION RATE)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Patient's age & sex is mandatory.	Spectrophotometry, Calculation	GFR calculation by this formula is not applicable for patient's less than 18 years. Race also plays an important role.	Daily
B034	GGTP (GAMMA GLUTAMYL TRANSPEPTIDASE)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	IFCC		Daily
H211	GIEMSA STAIN, PERIPHERAL BLOOD / BONE MARROW		Submit Peripheral Blood / Bone Marrow OR Smears fixed in absolute alcohol for 30 minutes and air dried. Ship refrigerated.	Microscopy		Sample Daily by 6 pm; Report NextDay
P211	GIEMSA STAIN, BODY FLUIDS		Submit Body fluids/Pus/Exudates/ FNA aspirates OR Smears fixed in absolute alcohol for 30 minutes and air dried. Ship refrigerated.	Microscopy		Sample Daily by 6 pm; Report NextDay

S158	GLIADIN DEAMIDATED ANTIBODY, IgA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
Z204	GLIADIN DEAMIDATED ANTIBODIES PANEL, IgA & IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S159	GLIADIN DEAMIDATED ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S098	GLOMERULAR BASEMENT MEMBRANE (GBM) ANTIBODY, IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 11 am; Report Same Day
RW081	GLUCAGON		3 mL (1.5 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma in a plastic transport vial & freeze immediately. Ship frozen. DO NOT THAW. Overnight fasting is mandatory.	Extraction, RIA		Sample by 7th of the month; Report after 2–3 weeks
B001	GLUCOSE, FASTING (F)		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Overnight fasting is mandatory. Ship refrigerated or frozen.	Hexokinase		Daily

G001	GLUCOSE ESTIMATION, 2.5 HOURS		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 10 mL (5 mL min.) aliquot of corresponding urine for each timed specimen. Collect fasting baseline plasma and urine specimen (0 hr). Dissolve 75g Glucose in adults / 1.75 g / Kg body weight in children (not exceeding 75g) in 300mL water. Administer orally over a period of 5 minutes. Collect additional plasma and urine samples at 0.5, 1.0, 1.5, 2.0 & 2.5 hours. Note time drawn on each tube, urine vial and test request form. Specify amount of glucose administered on test request form. Overnight fasting mandatory. Ship refrigerated or frozen.	Hexokinase, Strip test	As per American Diabetes Association (ADA) & WHO recommended test for diagnosis of Diabetes Mellitus is G010. For pregnant patient use test code G009. Also see INSULIN & C-PEPTIDE TESTS. Refer TO Y CODES: Y064 BasalY065 0.5 hrY066 1.0 hrÇ	Daily
G002	GLUCOSE ESTIMATION EXTENDED, 5 HOURS		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 10 mL (5 mL min.) aliquot of corresponding urine for each timed specimen. Collect fasting baseline plasma and urine specimens (0 hr). Dissolve 75g Glucose in adults / 1.75 g / Kg body weight in children (not exceeding 75g) in 300mL water. Administer orally over a period of 5 minutes. Collect additional plasma and urine samples at 0.5, 1.0, 1.5, 2.0, 2.5. 3.0, 4.0 and 5.0 hours. Note time drawn on each tube, urine vial and test request form. Specify amount of glucose administered on test request form. Overnight fasting mandatory . Ship refrigerated or frozen.	Hexokinase, Strip test	As per American Diabetes Association (ADA) & WHO recommended test for diagnosis of Diabetes Mellitus is G010. Also see INSULIN & C-PEPTIDE TESTS. Refer to Y CODES: Y064 BasalY065 0.5 hrY066 1.0 hrY067 1.5 hrY068 2.0 hrY069 Ç	Daily

B196	GLUCOSE TOLERANCE TEST (GTT), PREGNANCY, 100g GLUCOSE		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube for each timed specimen. Draw baseline fasting specimen. Dissolve 100 g of glucose in 300 ml water. Administer orally over a period of 5 minutes. Draw additional specimens at 1.0 , 2.0 & 3.0 hrs. Note time drawn on each container label and test request form. Overnight fasting is mandatory. Ship refrigerated or frozen.	Hexokinase	Two step screening procedure to Diagnose Gestational Diabetes as Recommended by American Diabetes Association Refer to Y CODES: Y124 Basal Y125 1 hr Y126 2 hrs	Daily
B195	GLUCOSE CHALLENGE TEST (GCT), PREGNANCY, 75G GLUCOSE		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube drawn 2 hours after oral administration of 75 g glucose. Overnight fasting is NOT mandatory. Ship refrigerated or frozen.	Hexokinase	Procedure to Diagnose Gestational Diabetes as Recommended by Diabetes in Pregnancy Study Group India (DIPSI)	Daily

G009	GLUCOSE TOLERANCE TEST (GTT), PREGNANCY, 75g GLUCOSE		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube for each timed specimen. Draw baseline fasting specimen. Dissolve 75g Glucose in 300 mL water. Administer orally over a period of 5 minutes. Draw additional specimens at 1.0 & 2.0 hours. Note time drawn on each container label and test request form. Overnight fasting mandatory. Ship refrigerated or frozen.	Hexokinase	This is the recommended test for pregnant patients. If results are normal in a clinically suspected case of Gestational Diabetes mellitus, a repeat test in third trimester is recommended. Refer to Y CODES: Y124 Basal Y125 1 hr Y126 2 hrs	Daily
Z131	GLUCOSE, FASTING (F) AND POST MEAL (PP), 2 HOURS		2 mL (1 mL min.) plasma each for Fasting and Post Prandial specimens from 1 Grey Top (Sodium Fluoride) tube. Overnight fasting mandatory. Collect fasting specimen. Ask patient to have a normal to heavy meal. Note the time patient starts the meal. Collect PP sample exactly after 2 hours. Ship refrigerated or frozen	Hexokinase		Daily
P085	GLUCOSE, FLUID		2 mL (1 mL min.) Body fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Hexokinase		Daily

B002	GLUCOSE, POST PRANDIAL (PP), 2 HOURS		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Collect sample exactly 2 hours after starting meal. Ship refrigerated or frozen.	Hexokinase	Also see GLUCOSE TOLERANCE TEST. Use the following test codes if glucose is administered OR PP time is altered. Code Meal/GlucoseB187 0.5 Hr after 75 gB143 1.0 Hr after 75 gB006 1.5 Hrs after Ç	Daily
G004	GLUCOSE CHALLENGE TEST (GCT); GLUCOSE 1 HOUR OBSTETRICAL		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube drawn 1 hour after oral administration of 50 g glucose. Does not include a fasting specimen . Ship refrigerated or frozen.	Hexokinase		Daily
G010	GLUCOSE TOLERANCE TEST, ORAL (OGTT)		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube each for Fasting and Post glucose specimens. Collect fasting specimen. Dissolve 75g Glucose in adults / 1.75 g / Kg body weight in children (not exceeding 75g) in 300mL water. Administer orally over a period of 5 minutes. Collect post glucose sample exactly after 2 hours. Overnight fasting is mandatory. Ship refrigerated or frozen. Urine Specimens NOT REQUIRED.	Hexokinase	As per American Diabetes Association (ADA) & WHO this is the recommended test for diagnosis of Diabetes Mellitus.	Daily
B003	GLUCOSE, RANDOM (R)		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Ship refrigerated or frozen.	Hexokinase		Daily

OS142	GLUTATHIONE-S-TRANSFERASE MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 9 am; Report 29 days
B198	GLYCEMARK (1,5 - ANHYDROGLUCITOL) SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzymatic Colorimetric	The test can help to identify patients who may need therapy changes despite a "good" HbA1c	Sample by Mon/Thu 9 am; Report same day
G123	GM1 GANGLIOSIDOSIS, QUANTITATIVE, BLOOD		10 ml (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium heparin) tubes. Ship refrigerated. DO NOT FREEZE. Clinical details must accompany sample.	Enzyme Assay	This tests deficiency of enzyme Beta Galactosidase.	Sample Daily by 4 pm; Report 4 Days
G124	GM2 GANGLIOSIDOSIS, QUANTITATIVE, BLOOD; TAY SACHS AND SANDHOFFDISEASE		10 ml (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium heparin) tubes. Ship refrigerated. DO NOT FREEZE. Clinical details must accompany sample.	Enzyme Assay	This tests deficiency of enzyme Total Hexosaminidase A & B.	Sample Daily by 4 pm; Report 4 days
Z200	GOODPASTURE'S SYNDROME PANEL	*Glomerular Basement Membrane Antibody*ANCA-PR3 / c-ANCA	3 ml (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Wed 9 am; ReportNext day
OS139	G6PD GENE MUTATION DETECTION		4 mL (2 mL min.)Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Tue 11 am; Report 8 days

R158	G-6-PD, NEWBORN SCREEN		1 drop of heel prick blood each on 3 spots of filter paper available from LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen.	Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm;Report 2 days
H018	G-6-PD, QUANTITATIVE		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Spectrophotometry		Sample Mon through Sat by 4 pm;Report Same day
P025	GRAM STAIN		Submit specimens in a sterile screw capped container OR smeared onto 2 sterile glass slides in a slide mailer.Ship refrigerated.	Microscopy of Gram stained smears		Daily
M214	GROUP A STREPTOCOCCA L ANTIGEN		Throat Swab in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily 1 pmReport Same day
Z069	GROWTH HORMONE SUPPRESSION BY GLUCOSE		2 mL (1 mL min.) serum from 1 SST for each timed specimen. Draw baseline specimen after an overnight fast. Administer Glucose orally 1.75 g/kg body weight (maximum 75g). Draw additional specimens at 60 & 120 min for GH. Ship refrigerated or frozen. All containers must be clearly marked with time drawn and submitted with 1 test request form.	Chemiluminescent Immunoassay	Patient should be at rest for 30 minutes prior to basal sample collectionrefer to Y Codes:Y100 Basal Y102 60 min.Y104 120 min. Y105 Customized	Daily

R031	GROWTH HORMONE; GH		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is mandatory.	Chemiluminescent Immunoassay	Patient should be at rest for 30 minutes prior to sample collection. Also see GROWTH HORMONE STIMULATION TEST, GROWTH HORMONE SUPPRESSION TEST & GROWTH DISORDER PANEL.	Daily
RW082	GROWTH HORMONE ANTIBODY		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Radiobinding Assay		Sample by 7th of the month; Report after 2–3 weeks
Z208	GROWTH DISORDER PANEL	*Growth Hormone *IGFBP-3 *IGF-I	3 mL (1.5 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW. Overnight fasting is mandatory.	Chemiluminescent Immunoassay	Patient should be at rest for 30 minutes prior to sample collection.	Sample by Mon / Thu 4 pm; Report Same day
RW083	GROWTH HORMONE RELEASING HORMONE; GHRH		3 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum into a sterile screw capped vial and freeze immediately. Ship frozen DO NOT THAW.	Extraction, Radioimmunoassay	Patient to avoid all medications that may influence pituitary secretion.	Sample by 7th of the month; Report after 2–3 weeks

Z300	GROWTH HORMONE STIMULATION TEST		2 mL (1 mL min.) serum from 1 SST for each timed specimen. Draw baseline specimen after an overnight fast. Administer stimulating agent Clonidine (Arkamine) 300 ug OR 5ug / Kg body weight orally. Draw additional specimens at 30, 60 & 90 mins. for Growth Hormone. Ship refrigerated or frozen. All containers must be clearly marked with time drawn and submitted with one test request form. Other stimulating agents that can be used are: Insulin (0.10 - 0.15 U / Kg body weight IV followed by additional samples at 30, 60 & 90 min). Exercise (20 min. vigorous exercise followed by immediate sampling)	Chemiluminescent Immunoassay	Patient should be at rest for 30 minutes prior to basal sample collection. Refer to Y Codes: Y094 Basal Y095 30 min. Y096 60 min. Y097 90 min. Y098 Customized	Daily
Z835	HAINS MYCOBACTERIU M TUBERCULOSIS FIRST & SECOND LINE DRUG RESISTANCE PROFILE		Submit 3mL (2mL min.) Sputum / bronchoalveolar lavages / bronchial aspirates / Culture growth in a sterile screw capped container.	PCR, Line probe assay		Sample by Sat 11 am, Report Fri
N068	HAINS MYCOBACTERIU M TUBERCULOSIS FIRST LINE DRUG RESISTANCE PROFILE		Submit 2 mL (1 mL min) Sputum / Bronchoalveolar Lavage / Pleural aspirates / CSF / Body fluid OR Cultured growth on appropriate media in a sterile screw capped container in a sealed plastic bag. Ship refrigerated. DO NOT FREEZE	PCR, Line probe assay		Sample by Tue / Fri 11 am, Report Thu / Mon

N117	HAINS MYCOBACTERIU M TUBERCULOSISS ECOND LINE DRUG RESISTANCE PROFILE		Submit 3mL (2mL min.) Sputum / bronchoalveolar lavages / bronchial aspirates / Culture growth in a sterile screw capped container.	PCR, Line probe assay		Sample by Sat 11 am, Report Fri
RW084	HALOPERIDOL		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Collect specimen 11–17 hours after the last dose. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
H119	HAM TEST; PAROXYSMAL NOCTURNALHE MOGLOBINURIA (PNH), SCREENING TEST		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Gel Technique	Also see CD 55 & CD 59 - Paroxysmal Nocturnal Hemoglobinuria confirmatory test by Flow Cytometry.	Sample Daily by 9 am; Report Same Day
OS144	HAPTOGLOBIN GENOTYPING		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed/Sat 11 am; Report 6 days
B100	HAPTOGLOBIN		2 mL (1 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	Immunoturbidimetry		Daily
B080	HbA1c; GLYCOSYLATED HEMOGLOBIN		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	High Performance Liquid Chromatography	Also see FRUCTOSAMINE, LIPID PROFILE & MICROALBUMINURIA.	Daily

B030	HBDH; ALPHA - HYDROXYBUTYR ATEDEHYDROG ENASE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	GSCC		Daily
F002	HCG, BETA, TOTAL, QUANTITATIVE, TUMOR MARKER		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
R021	HCG, BETA, TOTAL, QUANTITATIVE, MATERNAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay	Also see MATERNAL SERUM SCREEN 3 & 4.	Daily
U075	HCG, BETA, TOTAL, QUANTITATIVE, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
U022	HCG, BETA, TOTAL, QUANTITATIVE, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine in a sterile screw capped container. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
RW085	HDL CHOLESTEROL SUBCLASSES	*HDL-2 *HDL-3	4 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Minimum 12 hours overnight fasting is mandatory. Ship refrigerated or frozen	Precipitation, Colorimetric, Enzymatic		Sample by 7th of the month; Reportafter 2–3 weeks
B128	HDL CHOLESTEROL, DIRECT		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzymatic Immunoinhibition	Also see LIPID PROFILE TESTS.	Daily

Z615	HEART SCREEN	*hsCRP *Uric acid *Glucose Fasting *Lipid profile basic *Homocysteine	3 mL (2 mL min.) serum from 1 SST AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Ship refrigerated or frozen. Minimum 12 hours fasting is mandatory.	Immunoturbidimetry , Hexokinase,Spectrophotometry, CMIA		Daily
F024	HE4; HUMAN EPIDIDYMIS PROTEIN 4		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CMIA	Also see ROMA (Risk of Ovarian Malignancy)	Sample by Mon / Thu 4 pm; Report Same day
H204	HEINZ BODIES, PERIPHERAL BLOOD		2 air dried peripheral blood smears. Ship in a slide mailer at room temperature. Give brief clinical history.	Microscopy		Daily
Z121	HELICOBACTER PYLORI ANTIBODIES PANEL, IgA & IgG		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Giemsa staining for H. pylori also available.	Sample by Thu 9 am; Report Same day
S082	HELICOBACTER PYLORI, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Giemsa staining for H. pylori also available.	Sample by Thu 9 am; Report Same day
S083	HELICOBACTER PYLORI, IgA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Giemsa staining for H. pylori also available.	Sample by Thu 9 am; Report Same day
S227	HELICOBACTER PYLORI ANTIGEN, STOOL		5g (3g min.) stool in a leakproof screw capped container. Do not use matchbox or absorbent material. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
OS143	HEMOCHROMA TOSIS 2 MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 9 am; Report 32 days

OS146	HEMOGLOBIN D IRAN MUTATION ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample by Tue/ Thu 9 am; Report 8 days
OS148	HEMOGLOBIN E (GLU 26LYS) MUTATION ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample by Tue/ Thu 9 am; Report 8 days
H001	HEMOGLOBIN; Hb		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Photometry		Daily
OS147	HEMOGLOBIN D PUNJAB MUTATION ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample by Tue/ Thu 9 am; Report 8 days
E002	HEMOGLOBIN HPLC/ELECTROPHORESIS	Includes Hematological Indices.	4 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Must Record date of birth, sex, race and religion on the test request form.	High Performance Liquid Chromatography/ Electrophoresis	HPLC is the method of choice. Hemoglobin electrophoresis (Alkaline pH) on cellulose acetate will be carried out, if necessary, in selected cases, at no additional charge and reported after 3 days.	Sample Daily by 10 am; Report Same Day
U047	HEMOGLOBIN, FREE, URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated.	Strip Test		Daily

Z019	HEMOGRAM	*Hemoglobin *PCV *TLC *RBC Count *MCV*MCH *MCHC *Platelet Count *Automated DLC*Absolute Differential counts *RDW *ESR	4 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Overnight fasting is preferred.	Electrical Impedence, VCS, Capillary Photometry / Westergren	Peripheral smear reported in abnormal cases only	Daily
H241	HEMOPHILIA A CARRIER DETECTION	*Factor VIII *vWF Antigen *Factor VIII / vWF Ratio	3 mL Whole blood in 1 Blue Top(Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection, EIA		Sample by 1st & 3rd Mon 11 am; Report Next day

Z374	HEMOPHILIA PANEL	*APTT *Factor VIII *Factor IX	3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical & Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Also see COAGULATION PROFİÇ	Sample by Wed / Sat 9 am; Report Fri / Tue
U046	HEMOSIDERIN, URINE		10 mL (5 mL min.) aliquot of random urine. First morning specimen is preferred. Ship refrigerated. DO NOT FREEZE OR Centrifuge urine at 3200 rpm for 15 min. Smear deposit on 4 glass slides, air dry, fix in 95% ethanol for 30 min. Ship at room temperature.	Prussian Blue Reaction		Sample Daily by 9 am; Report Same day

RW086	HEPARIN, ANTI Xa ASSAY		2 mL (1 mL min.) platelet poor plasma from 2 Blue Top (Sodium Citrate) tubes. Make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred.	Chromogenic		Sample by 7th of the month; Report after 2–3 weeks
Z372	HEPATITIS ACUTE VIRUS CONFIRMATION	*Anti-HAV IgM *HBsAg *HBe Ag * Anti-HBe*Anti-HBc IgM *Anti-HCV *Anti-HEV IgM	4 mL (2 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay, EIA		Sample by Mon through Sat 9 am; Report Same day
S057	HEPATITIS B SURFACE ANTIBODY; ANTI HBs		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay		Daily
S116	HEPATITIS B SURFACE ANTIGEN (HBsAg) CONFIRMATION		2 mL (1 mL min.) serum from 1 SST. Indicate on test request form that specimen has been found positive and that a confirmation procedure is needed. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay		Daily
S215	HEPATITIS B SURFACE ANTIGEN (HBsAg), QUANTITATIVE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay		Daily

S056	HEPATITIS B SURFACE ANTIGEN; HBsAg; AUSTRALIA ANTIGEN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay	For Hepatitis B surface antigen confirmation by neutralization, use test code S116 HEPATITIS B SURFACE ANTIGEN CONFIRMATION.	Daily
N041	HEPATITIS B VIRAL DRUG RESISTANCE & GENOTYPING		3 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen.	Sequencing	The minimum viral load required is 2000 IU/mL	Sample Daily by 11 am; Report 7 days
N011	HEPATITIS B VIRUS, PCR, QUALITATIVE;H BV DNA PCR, QUALITATIVE		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. Overnight fasting is recommended.	Polymerase Chain Reaction		Sample by Tue / Thu / Sat 11 am; Report Wed / Fri / Mon
S059	HEPATITIS Be ANTIGEN; HBeAg		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
Z137	HEPATITIS Be PANEL	* Anti-HBe * HBeAg	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
S062	HEPATITIS C ANTIBODY (Anti-HCV)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily

N028	HEPATITIS C VIRAL RNA (HCV RNA) QUANTITATIVE, REAL TIME PCR		3 mL (2 mL min.) serum from 2 SST's OR 3 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Separate serum/plasma aseptically within 2 hours of collection. 12 hours fasting is recommended. Ship refrigerated or frozen.	Real Time PCR	Avoid freeze thaw cycles.	Sample by Tue / Fri 11 am; Report Thu / Mon
N032	HEPATITIS C VIRAL RNA, GENOTYPE		2 mL (1 mL min.) serum from 1 SST OR 2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate serum/plasma aseptically within 1 hour of collection. 12 hours fasting is recommended. Ship refrigerated or frozen.	Real Time PCR		Sample by Tue / Fri 11 am; Report Thu / Mon
S051	HEPATITIS E ANTIBODY (Anti-HEV), IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon through Sat 9 am; Report Same day
N109	HEPATITIS E VIRUS DETECTION, PCR		2 ml (1 ml min.) serum from 1 SST OR 2 ml (1 ml min.) plasma from 1 Lavender top (EDTA) tube. Separate plasma / serum within 4 hours of collection. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
OS150	HEPATITIS G VIRUS RNA DETECTION , PCR		4 mL (2 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended.	PCR		Sample by Thu 11 am; Report Wed
Z365	HEPATITIS A & B VIRUS EVALUATION	*Anti-HAV, IgG *Anti-HAV IgM *Anti-HBc Total*Anti-HBc IgM *Anti-HBs *HBsAg *Anti- HBe*HBeAg	4 mL (2 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immun oassay		Daily
Z133	HEPATITIS ACUTE DIAGNOSTIC PANEL	*Anti-HAV IgM *HBsAg *Anti-HB Core IgM	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immun oassay		Daily

Z371	HEPATITIS ACUTE VIRUS SCREEN	*Anti-HAV IgM *HBsAg *Anti-HB Core IgM*Anti-HCV *Anti-HEV IgM	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CMIA, EIA		Sample by Mon through Sat 9 am; Report Same day
Z373	HEPATITIS ATYPICAL VIRUS CONFIRMATION	*CMV IgM *HSV IgM *Varicella Zoster IgM*EBV (VCA), IgM	4 mL (2 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	CMIA, EIA		Sample by Mon / Thu 9 am; Report Thu / Mon
S058	HEPATITIS B CORE ANTIBODY (Anti-HBc), IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmunoassay		Daily
OS145	HEPATITIS A VIRUS RNA DETECTION, PCR		4 mL (2 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended.	PCR		Sample by Fri 11 am; Report Wed
Z017	HEPATITIS B PROFILE	*HBsAg *Anti-HBs *HBeAg *Anti-HBe *Anti-HBc IgM *Anti-HBc Total	4 mL (2 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmunoassay		Daily
Z602	HEPATITIS B VIRUS GENOTYPING, DRUGRESISTANCE & QUANTITATIVE PCR		3 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended	Sequencing, Real Time PCR	Avoid freeze thaw cycles	Sample Daily by 11 am; Report 7 days

Z511	HEPATITIS C VIRAL COMBO	*HCV RNA Quantitative Real Time PCR * HCV Genotype	6 mL (4 mL min.) serum from 2 SSTs OR 6 mL (4 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Separate serum/plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended.	Real Time PCR	This test is not intended for use in the initial diagnosis or confirmation of HCV infection.	Sample by Tue / Fri 11 am; Report Thu / Mon
S202	HEPATITIS E ANTIBODY (Anti-HEV), IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Wed / Fri 9 am; Report Same day
Z138	HEPATITIS, VIRAL, COMPREHENSIV E PANEL	*Anti-HAV, IgG *Anti-HAV IgM *Anti-HBc Total*Anti-HBc IgM *Anti-HBs *HBsAg *Anti- HBe*HBe Ag *Anti-HCV *Anti- HEV IgG *Anti- HEV IgM	6 mL (4 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immun oassay, EIA		Sample by Mon / Wed / Fri 9 am; Report Same day
S090	HEPATITIS A ANTIBODY (Anti-HAV), IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immun oassay		Daily
S061	HEPATITIS A ANTIBODY (Anti-HAV), IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immun oassay		Daily
Z132	HEPATITIS A DIAGNOSTIC PANEL	*Anti-HAV, IgG *Anti-HAV, IgM	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CMIA		Daily
Z134	HEPATITIS B CHRONIC PANEL	*HBsAg * HBeAg * Anti-HBe	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay		Daily

S115	HEPATITIS B CORE ANTIBODY (Anti-HBc), TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Microparticle Immunoassay		Daily
Z135	HEPATITIS B IMMUNITY SCREEN	*Anti-HBc Total *Anti-HBs *HBsAg	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
Z136	HEPATITIS B INFANT FOLLOW-UP PANEL	*HBsAg *Anti- HBs	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
N027	HEPATITIS B VIRAL DNA (HBV DNA) QUANTITATIVE, REAL TIME PCR		3 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes OR 2 mL (1 mL min.) serum from 1 SST. Separate plasma / serum aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended.	Real Time PCR	Avoid freeze thaw cycles	Sample by Tue / Thu / Sat 11 am; Report Wed / Fri / Mon
S060	HEPATITIS Be ANTIBODY; Anti-HBe		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay		Daily
N004	HEPATITIS C VIRAL, PCR, QUALITATIVE; HCV RNA, PCR, QUALITATIVE		2 mL (1 mL min.) serum from 1 SST OR 2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate serum/plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. Overnight fasting is recommended.	Polymerase Chain Reaction	Avoid freeze thaw cycles.	Sample by Tue / Fri 11 am;Repot Thu/ Mon
OS149	HEPATITIS D VIRUS DETECTION, PCR		4 mL (2 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended.	PCR		Sample by Thu 11 am; Report Wed

S117	HEPATITIS DELTA ANTIBODY (Anti-HDV), IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Wed 11 am; Report Same Day
S021	HERPES SIMPLEX VIRUS (HSV) 1+2, IgG, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
S071	HERPES SIMPLEX VIRUS (HSV) 1+2, IgM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
S214	HERPES SIMPLEX VIRUS (HSV) 2 IgM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	For detection of Acute HSV infections, HSV PCR is the recommended test.	Sample by Tue / Fri 9 am; Report Same Day
Z556	HERPES SIMPLEX VIRUS (HSV) 1 ANTIBODIES PANEL, IgG & IgM, SERUM		3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay Chemiluminescent Immunoassay	For detection of Acute HSV infections, HSV PCR is the recommended test.	Sample by Tue / Fri 9 am; Report Same Day
S211	HERPES SIMPLEX VIRUS (HSV) 1 IgG, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	For detection of Acute HSV infections, HSV PCR is the recommended test.	Sample by Tue / Fri 9 am; Report Same Day

S213	HERPES SIMPLEX VIRUS (HSV) 1 IgM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	For detection of Acute HSV infections, HSV PCR is the recommended test.	Sample by Tue / Fri 9 am; Report Same Day
Z557	HERPES SIMPLEX VIRUS (HSV) 2 ANTIBODIES PANEL, IgG & IgM, SERUM		3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	For detection of Acute HSV infections, HSV PCR is the recommended test.	Sample by Tue / Fri 9 am; Report Same Day
S212	HERPES SIMPLEX VIRUS (HSV) 2 IgG, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	For detection of Acute HSV infections, HSV PCR is the recommended test.	Sample by Tue / Fri 9 am; Report Same Day
Z303	HERPES SIMPLEX VIRUS (HSV) TYPE 1 & 2, PCR, QUALITATIVE		2 mL (1 mL min.) CSF/EDTA Plasma in a sterile screw capped container OR Blister fluid/Swab from base of ulcer/ Scrapings from base of ulcer in 1 mL sterile normal saline or phosphate buffered saline (PBS) OR Fresh tissue in sterile normal saline OR Paraffin embedded tissue. Ship refrigerated or frozen.	Polymerase Chain Reaction, Amplification		Sample by Mon / Thu 11 am; ReportWed / Sat
N013	HERPES SIMPLEX VIRUS (HSV) TYPE 2, PCR, QUALITATIVE		2 mL (1 mL min.) CSF/EDTA Plasma in a sterile screw capped container OR Blister fluid/Swab from base of ulcer/ Scrapings from base of ulcer in 1 mL sterile normal saline or phosphate buffered saline (PBS) OR Fresh tissue in sterile normal saline OR Paraffin embedded tissue. Ship refrigerated or frozen.	Polymerase Chain Reaction, Amplification		Sample by Mon / Thu 11 am; ReportWed / Sat

N012	HERPES SIMPLEX VIRUS (HSV) TYPE I, PCR, QUALITATIVE		2 mL (1 mL min.) CSF/EDTA Plasma in a sterile screw capped container OR Blister fluid / Swab from base of ulcer/ Scrapings from base of ulcer in 1 mL sterile normal saline or phosphate buffered saline (PBS) OR Fresh tissue in sterile normal saline OR Paraffin embedded tissue. Ship refrigerated or frozen.	Polymerase Chain Reaction, Amplification		Sample by Mon / Thu 11 am; Report Wed / Sat
Z205	HERPES SIMPLEX VIRUS (HSV) 1+2 ANTIBODIES PANEL, IgG & IgM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
RW088	HERPES VIRUS 6 (HHV-6) ANTIBODY PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks
G185	HHH SYNDROME; HYPERORNITHI NEMIA - HYPERAMMON EMIA - HOMOCITRULLI NURIA SYNDROME		15 mL (10 mL min.) aliquot of random urine in a sterile screw capped container. No preservative required. Ship refrigerated or frozen.	LC-MS/MS, GC-MS		Sample Daily by 5 pm; Report 5 days

U034	5-HIAA (5-HYDROXY INDOLE ACETIC ACID), 24-HOUR URINE		50 mL (10 mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl to maintain pH between 1–2. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen. Patient should strictly avoid banana, kiwi, walnut, avocado, brinjal, pineapple, plum, tomato & drugs like Fluorouracil, Melphalan, Paracetamol & Acetaminophen, Heparin, L-dopa, Reserpine, Salicylates, Chlorpromazine, Imipramine, Isoniazid, MAO inhibitors, Phenothiazines, Promethazine, alcoholic beverages, caffeine, tea/coffee, tobacco and strenuous exercise at least 72 hrs before and during specimen collection.	High Performance LiquidChromatography		Sample by Wed / Sat 9 am; Report Fri / Tue
U038	5 - HIAA (5-HYDROXY INDOLE ACETIC ACID), RANDOM URINE		50 mL (10 mL min.) aliquot of random urine. Add 50% HCl drop by drop to bring the pH between 1-2. Do not use concentrated HCl. Ship refrigerated or frozen. Patient should strictly avoid banana, kiwi, walnut, avocado, brinjal, pineapple, plum, tomato & drugs like Fluorouracil, Melphalan, Paracetamol & Acetaminophen, Heparin, L-dopa, Reserpine, Salicylates, Chlorpromazine, Imipramine, Isoniazid, MAO inhibitors, Phenothiazines, Promethazine, alcoholic beverages, caffeine, tea/coffee, tobacco and strenuous exercise at least 72 hrs before and during specimen collection.	High Performance Liquid Chromatography		Sample by Wed / Sat 9 am; Report Fri / Tue
Z040	HIRSUTISM PANEL	*Androstenedione *DHEA Sulphate *Free Testosterone	5 mL (3 mL min.) serum from 2 SST's. An early morning specimen is preferred. Ship refrigerated or frozen. Specify age and sex on test request form.	CLIA, CMIA, RIA		Sample by Mon / Thu 9 am; Report Same Day

RW089	HISTAMINE, 24-HOUR URINE		5 mL (3 mL min.) aliquot of 24-hour urine. Collect urine with 10 mL of 50% 6N HCl to maintain pH between 1–3. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen. Avoid direct sunlight, antihistamines, oral corticosteroids and H2 receptor blockers, 24 hours prior to specimen collection.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
J113	HISTO / CYTOPATHOLOGY / IMMUNOHISTOCHEMISTRY (IHC), PHOTOMICROGRAPH		Submit stained Histopathology sections/ Cytology smears / IHC stained sections. Ship at room temperature in a slide mailer. Brief clinical history in Test Requisition Form is mandatory.	Digital Photography	2 Magnifications provided on Photographic paper or CD ROM.	Sample Daily by 6 pm; Report Next day
J119	HISTO/CYTOPATHOLOGY, SLIDE(S) / BLOCK(S), ISSUE		Submit written request for issue of slide(s) / block(s) to department of Histo / Cytopathology.	NA	Only representative slides / blocks will be issued on advise of the referring consultant.	Slides / Blocks shall be issued from the Histopathology department 72 hours after receiving the request.
S019	HISTONE ANTIBODIES		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay	Also see LUPUS ERYTHEMATOSUS PANEL & EXTRACTABLE NUCLEAR ANTIGENS.	Sample by Thu 9 am; Report Same day

Z158	HISTOPATHOLOGY, DIRECT IMMUNOFLOUORESCENCE (DIF),SKIN / CONJUNCTIVAL BIOPSY,	* C3 *C1q *IgA *IgG *IgM	Submit Skin / Conjunctival Biopsy specimen in buffered normal saline. Biopsy should measure minimum 3 mm. Ship refrigerated.	Cryoprocessing, Fluorescence Microscopy	Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report 4 days
J040	HISTOPATHOLOGY, DIRECT IMMUNO-FLUORESCENCE (DIF),SKIN / CONJUNCTIVAL BIOPSY, COMPLEMENT C1q		Submit Skin / Conjunctival Biopsy specimen in buffered normal saline. Biopsy should measure minimum 3 mm.Ship refrigerated.	Cryoprocessing, Fluorescence Microscopy	. Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report 4 days
Z829	HISTOPATHOLOGY, DIRECT IMMUNOFLOUORESCENCE (DIF), KIDNEY BIOPSY	*Complement C3 *Complement C1q *IgA * IgG *IgM *Kappa Light Chains *Lambda Light Chains	Submit Kidney Biopsy specimen in buffered normal saline. Kidney biopsy should measure minimum 1 cm in length. Ship refrigerated. Brief clinical history in Renal Biopsy form for Native / Transplant kidneys is mandatory.	Cryoprocessing, Fluorescence Microscopy	Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report 3 days
J024	HISTOPATHOLOGY, DIRECT IMMUNO-FLUORESCENCE (DIF),SKIN / CONJUNCTIVAL BIOPSY, COMPLEMENT C3		Submit Skin / Conjunctival Biopsy specimen in buffered normal saline. Biopsy should measure minimum 3 mm. Ship refrigerated.	Cryoprocessing, Fluorescence Microscopy	Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report 4 days

J026	HISTOPATHOLOGY, DIRECT IMMUNO-FLUORESCENCE (DIF),SKIN / CONJUNCTIVAL BIOPSY, IgG		Submit Skin / Conjunctival Biopsy specimen in buffered normal saline. Biopsy should measure minimum 3 mm. Ship refrigerated.	Cryoprocessing, Fluorescence Microscopy	Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report 4 days
J027	HISTOPATHOLOGY, DIRECT IMMUNO-FLUORESCENCE (DIF),SKIN / CONJUNCTIVAL BIOPSY, IgM		Submit Skin / Conjunctival Biopsy specimen in buffered normal saline. Biopsy should measure minimum 3 mm. Ship refrigerated.	Cryoprocessing, Fluorescence Microscopy	Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report 4 days
J007	HISTOPATHOLOGY, BIOPSY - SECOND OPINION		Histopathology sections or paraffin embedded tissue blocks previously reported upon. Ship at room temperature. Previous report(s), clinical history and other details to be submitted along with the test request form.	Microscopy	Special stains performed only on request, as additional tests.	Sample Daily by 6 pm; Report 3 days
J118	HISTOPATHOLOGY, BIOPSY, BONE		Submit specimen in 10% Formal-saline. Formal-saline should cover the specimen completely. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory. For patients with multiple specimens kindly register all specimens under single Lab number.	Processing, Staining & Microscopy	Use Test Code J116 for Large Bone Resections. Bone specimens require decalcification & reporting time can vary. Special stains performed only on request, as additional tests.	Sample Daily by 6 pm; Report maximum 12 days

J116	HISTOPATHOLOGY, BIOPSY, LARGE COMPLEX / CANCER SPECIMEN	*Esophagectomy *Gastrectomy *Mastectomy *Hemi / Total Colectomy *Large Bone Resection*Ovarian Tumor Resection *Radical Nephrectomy for Cancer *Radical Neck Dissection *Radical Hysterectomy *Radical Orchidectomy *Soft tissue tumor resections *Head & Neck C	Submit specimen in 10% Formal-saline. Formal-saline should cover the specimen completely. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory. For patients with multiple specimens kindly register all specimens under single Lab number.	Processing, Staining & Microscopy	Bone specimens require decalcification and reporting will take 14-16 days. Special stains performed only on request, as additional tests.	Sample Daily by 6 pm; Report 5 days. Bone specimens 14 - 16 days
------	---	---	---	-----------------------------------	---	--

J001	HISTOPATHOLOGY, BIOPSY, SMALL SPECIMEN	*Endometrium *Cervical biopsy *Endoscopic biopsies *Trucut biopsies *Appendix *FallopianTubes *Conjunctival biopsy *Small diagnostic / Incision biopsies	Submit specimen in 10% Formal-saline. Formal-saline should cover the specimen completely. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory. For patients with multiple specimens kindly register all specimens under single Lab number.	Processing, Staining & Microscopy	Special stains performed only on request, as additional tests.	Sample Daily by 11 am; Report Nextday Sample Daily 11 am to 6 pm; Report 2days
Z353	HISTOPATHOLOGY, KIDNEY BIOPSY PANEL 1	*Histopathology *Special Stains *DIF for IgG, IgA, IgM, C3 & C1q *IHC for C4d & SV40 (Transplant biopsy) *Photomicrograph	Submit 2 Kidney Biopsy specimens; 1 specimen in buffered normal saline for DIF studies AND 1 specimen in 10% formal saline for Histopathology. Ship refrigerated. Brief clinical history in Renal Biopsy form for Native / Transplant kidneys is mandatory.	Cryoprocessing, Fluorescence and Light Microscopy, IHC, Digital Photography	Each Biopsy should measure minimum 1 cm in length. Special stains include PAS, Masson Trichrome, Silver & Congo Red.	Sample Daily by 6 pm; Report 5 days
Z355	HISTOPATHOLOGY, KIDNEY BIOPSY PANEL 2	*Histopathology *Special Stains *IHC for C4d & SV40 (Transplant biopsy) *Photomicrograph	Submit Kidney Biopsy specimen in 10% formal saline. Ship at room temperature. Brief clinical history in Renal Biopsy form for Native / Transplant kidneys is mandatory.	Light Microscopy, IHC, DigitalPhotography	Each Biopsy should measure minimum 1 cm in length. Special stains include PAS, Masson Trichrome, Silver & Congo Red.	Sample Daily by 6 pm; Report 5 days

Z830	HISTOPATHOLOGY, KIDNEY BIOPSY PANEL 3	*Histopathology *Special Stains *DIF for IgG,IgA, IgM, C3 & C1q *IHC for C4d & SV40 (Transplant biopsy) *Photomicrograph * Electron Microscopy	Submit 3 Kidney Biopsy specimens; 1 specimen in buffered normal saline for DIF studies AND 1 specimen in 10% formal saline for Histopathology AND 1 specimen in 3% Buffered Glutaraldehyde solution available from LPL. Ship refrigerated. Brief clinical history in Renal Biopsy form for Native / Transplant kidneys is mandatory.	Cryoprocessing, Fluorescence and Light Microscopy, IHC, Digital Photography, Transmission Electron Microscopy	Each Biopsy should measure minimum 1 cm in length. Special stains include PAS, Masson Trichrome, Silver & Congo Red.	Sample Daily by 6 pm; Report 7 days
Z850	HISTOPATHOLOGY, KIDNEY BIOPSY SECOND OPINION		Histopathology sections or paraffin embedded tissue blocks previously reported upon. Ship at room temperature. Previous report(s), clinical history and other details to be submitted along with the test request form.	Light Microscopy	Special Stains Performed On Request Depending On Availability Of Unstained Slides.	Sample Daily by 6 pm; Report 4 days
J025	HISTOPATHOLOGY, DIRECT IMMUNOFLOUORESCENCE (DIF), SKIN / CONJUNCTIVAL BIOPSY, IgA		Submit Skin / Conjunctival Biopsy specimen in buffered normal saline. Biopsy should measure minimum 3 mm. Ship refrigerated.	Cryoprocessing, Fluorescence Microscopy	Use test code J113 for Photomicrograph.	Sample Daily by 6 pm; Report 4 days

J092	HISTOPATHOLOGY, BIOPSY, MEDIUM SPECIMEN	*Breast lump *Pilonidal sinus *Fistula / Sinus *Lymph node *Ovarian cyst *Gall bladder *Prostate (TURP / Enucleation) *Superficial lumps *Brain & Spinal cord tumors * Small excision biopsies *Fibroids (Enucleated) *Products of conception	Submit specimen in 10% Formal-saline. Formal-saline should cover the specimen completely. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory. For patients with multiple specimens kindly register all specimens under single Lab number.	Processing, Staining & Microscopy	Special stains performed only on request, as additional tests.	Sample Daily by 6 pm; Report 3 days
Z386	HISTOPATHOLOGY, LIVER BIOPSY PANEL	*Histopathology *Special Stains *Photomicrograph	Submit Liver Biopsy specimen in 10% formal saline. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Light Microscopy, Digital Photography	Biopsy should measure minimum 1.5 cm in length. Special stains include PAS, PAS-D, Masson Trichrome, Iron & Reticulin. Congo Red Staining is done if required.	Sample Daily by 6 pm; Report 5 days

J114	HISTOPATHOLOGY, POLARISINGMICROSCOPY		Submit Histopathology Sections OR Paraffin embedded tissue blocks. Ship at room temperature. Previous report(s), clinical history and other details to be submitted along with the test request form.	Polarising Microscopy		Sample Daily by 6 pm; Report Next day
J050	HISTOPATHOLOGY, SPECIAL STAIN, ACID FAST (Z.N STAIN)		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J084	HISTOPATHOLOGY, SPECIAL STAIN, ALCIAN BLUE		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J052	HISTOPATHOLOGY, SPECIAL STAIN, AMYLOID (CONGO RED)		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J095	HISTOPATHOLOGY, SPECIAL STAIN, COPPER (RHODANINE)		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days

J078	HISTOPATHOLOGY, SPECIAL STAIN, FONTANA MASSON		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J079	HISTOPATHOLOGY, SPECIAL STAIN, GOMORI'S METHENAMINE SILVER		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J053	HISTOPATHOLOGY, SPECIAL STAIN, IRON (PERL'S STAIN)		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J081	HISTOPATHOLOGY, SPECIAL STAIN, MASSON TRICHROME		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days

J080	HISTOPATHOLOGY, SPECIAL STAIN, MUCICARMINE		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J054	HISTOPATHOLOGY, SPECIAL STAIN, PERIODIC ACID SCHIFF (PAS)		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J082	HISTOPATHOLOGY, SPECIAL STAIN, PHOSPHOTUNGSTIC ACID HEMATOXYLIN (PTAH)		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
J059	HISTOPATHOLOGY, SPECIAL STAIN, RETICULIN		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days

J083	HISTOPATHOLOGY, SPECIAL STAIN, VAN GIESONS		Submit tissue biopsy in 10% formal saline OR formalin fixed paraffin embedded tissue block OR FNAC specimen, alcohol fixed, smeared onto 4 glass slides and air dried. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory.	Special stain, Microscopy		Sample Daily by 6 pm; Report 3 days
Z701	HISTOPATHOLOGY, SKIN BIOPSY PANEL	Histopathology examination *Special stains *Photomicrograph	Submit specimen in 10% Formal-saline. Formal-saline should cover the specimen completely. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory. For patients with multiple specimens kindly register all specimens under single Lab number.	Processing, Staining & Microscopy		Sample Daily by 6 pm; Report 5 days
RW269	HISTOPLASMA ANTIBODY, IMMUNODIFFUSION		1 mL (0.5 mL min.) Serum from 1 SST. Ship refrigerated or frozen.	Immunodiffusion		Sample by 7th of the month; Report after 2–3 weeks
S010	HIV 1 & 2 ANTIBODIES, SCREENING TEST	Includes p24 Antigen	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Patient consent on HIV Consent form is mandatory.	Chemiluminescent Microparticle Immunoassay		Daily
S063	HIV 1 & 2 ANTIBODIES, WESTERN BLOT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Western Blot		Sample by Tue / Thu / Sat 9 am; Report Same day
OS151	HIV-2 DETECTION, PCR		4 mL (2 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended.	PCR		Sample by Mon 11 am; Report Fri

N060	HIV 1 GENOTYPING & DRUG RESISTANCE		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma within 6 hours of collection in a plastic transport vial and freeze immediately. Ship frozen. DO NOT THAW.	Reverse Transcriptase PCR and DNA Sequencing		Sample by 1st / 15th; Report on 15th / 30th of every month;
H216	HIV 1 RNA QUANTITATIVE, REAL TIME PCR WITH CD3, CD4 / CD8 (IMMUNE DEFICIENCY PANEL 2) MARKERS		PCR: 3 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended. FLOW: 4 mL (3 mL min.) whole blood from 1 Lavender Top (EDTA) AND 1 Green top (Sodium Heparin) tube each. Ship at 18-22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on Test Request Form.	Real Time PCR, Flow Cytometry		Sample by Wed / Sat 11 am; Report Fri / Tue
N029	HIV 1 RNA QUANTITATIVE, REAL TIME PCR		3 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended.	Real Time PCR	Avoid freeze thaw cycles.	Sample by Wed / Sat 11 am; Report Fri/ Tue
H215	HIV 1 RNA QUANTITATIVE, REAL TIME PCR WITH CD3 / CD4 (IMMUNE DEFICIENCY PANEL 4) MARKERS		PCR: 3 mL (2 mL min.) plasma from 2 Lavender Top (EDTA) tubes. Separate plasma aseptically within 2 hours of collection. Ship refrigerated or frozen. 12 hours fasting is recommended. FLOW: 4 mL (3 mL min.) whole blood from 1 Lavender Top (EDTA) AND 1 Green top (Sodium Heparin) tube each. Ship at 18-22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on Test Request Form.	Real Time PCR, Flow Cytometry		Sample by Wed / Sat 11 am; Report Fri / Tue
L021	HLA - A2 (A*02)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday

L025	HLA - ABC (CLASS 1) TYPING		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex	Doctor's prescription on the letter head with stamp AND 1 photograph of recipient / donor is mandatory.	Sample Daily by 9 a.m; Report Nextday
L005	HLA - AB (CLASS 1) TYPING		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex	Doctor's prescription on the letter head with stamp AND 1 photograph of recipient / donor is mandatory.	Sample Daily by 9 a.m; Report Nextday
L019	HLA - B22 (B*22)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L026	HLA - B7 (B*07)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L009	HLA - B38 (B*38)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L008	HLA - B5 (B*51/52)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L007	HLA - B8 (B*08)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday

L010	HLA - B15 (B*15)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L004	HLA-B27, FLOW CYTOMETRY		4 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) / Green Top (Sodium Heparin) tube. Ship at 18–22°C. DO NOT REFRIGERATE OR FREEZE.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
L028	HLA B LOCUS TYPING - ANY ANTIGEN		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
N044	HLA - B27, QUALITATIVE PCR		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample Daily by 9 am; Report Same Day
L038	HLA CADAVERIC CROSS MATCH		RECIPIENT: 6 mL (3 mL min.) whole blood in 1 Yellow Top (ACD) tube AND 3 mL (2 mL min.) serum from 1 SST. Collect specimen preferably 48 hours after dialysis.CADAVERIC DONOR: 12 mL (10 mL min.) whole blood in 2 Yellow Top (ACD) tubes. Ship refrigerated. DO NOT FREEZE. Brief history in Organ Transplant Clinical information & Consent form is mandatory.	Complement Dependent Cytotoxicity with DTT treatment	Doctor's prescription on the letter head with stamp AND photograph of recipient is mandatory. Above specimen requirements are for a single donor being cross matched with a single recipient. For additional cross matches, send additional specimens of donor	Sample Daily; Report 8 hrs

L017	HLA - CELIAC DISEASE (DR-DQB1-DQA1 HAPLOTYPE ASSOCIATION INCLUDING DQ2,DQ8)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	SSP - PCR	For Pediatric cases use Microtainers available from LPL	Sample Daily by 9 a.m; Report Nextday
L032	HLA C LOCUS TYPING, ANY ANTIGEN		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L020	HLA - Cw6 (C*06)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
N024	HLA DNA TYPING FOR BONE MARROW TRANSPLANT, 10 ANTIGENS, PCR	*A, B, C, DR & DQ Low Resolution Typing	4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	SSP - PCR	Also see OPPOSITE SEX BMT MISMATCH AND CHIMERISM tests.	Sample Daily by 9 a.m; Report 2 days
N086	HLA DNA HIGH RESOLUTION TYPING - A LOCUS		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.
N070	HLA DNA HIGH RESOLUTION TYPING A, B, C, DR & DQ		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 8 days.

N023	HLA DNA TYPING FOR KIDNEY / LIVER TRANSPLANT, PCR	*A, B & DR Low Resolution Typing	4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex	Doctor's prescription on the letter head with stamp AND 1 photograph each of recipient and donor is mandatory.	Sample Daily by 9 a.m; Report 2 days
N087	HLA DNA HIGH RESOLUTION TYPING - B LOCUS		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.
N083	HLA DNA HIGH RESOLUTION TYPING - C LOCUS		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.
N085	HLA DNA HIGH RESOLUTION TYPING - DQB1 LOCUS		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.
N084	HLA DNA HIGH RESOLUTION TYPING - DRB1 LOCUS		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.

N025	HLA DNA TYPING FOR BONE MARROW TRANSPLANT, 6 ANTIGENS, LUMINEX	*A, B & DR Low Resolution Typing	4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report 2 days
L037	HLA DONOR SPECIFIC IgG ANTIBODIES & PANEL REACTIVE ANTIBODIES SCREEN PROFILE		RECIPIENT: 3 mL (2 mL min.) serum from 1 SST. Collect specimen preferably 48 hours after dialysis. DONOR: 12 mL (10 mL min.) whole blood in 2 Yellow Top (ACD) tubes. Ship refrigerated. DO NOT FREEZE. Brief history in Organ Transplant Clinical information & Consent form form is mandatory.	Luminex Fluoro - Beads X-Map technology		Sample by Tue / Thu / Sat 9 amReport Wed / Fri / Mon 6 pm
L036	HLA - DQB1 TYPING		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex	Doctor's prescription on the letter head with stamp AND 1 photograph of recipient / donor is mandatory.	Sample Daily by 9 a.m; Report Nextday
L016	HLA - DR7 (DRB1*07)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L018	HLA - DR & DQB1 (CLASS II) TYPING		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex	Doctor's prescription on the letter head with stamp AND 1 photograph of recipient / donor is mandatory.	Sample Daily by 9 a.m; Report Nextday

L006	HLA - DR (CLASS II) TYPING		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex	Doctor's prescription on the letter head with stamp AND 1 photograph of recipient / donor is mandatory.	Sample Daily by 9 a.m; Report Nextday
L014	HLA - DR5 (DRB1*11/12)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L012	HLA - DR3 (DRB1*03)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L013	HLA - DR4 (DRB1*04)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday
L011	HLA - DR2 (DRB1*02)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Reverse SSOP - Luminex		Sample Daily by 9 a.m; Report Nextday

L022	HLA DSA (DONOR SPECIFIC IgG ANTIBODIES), CLASS I & II		RECIPIENT: 3 mL (2 mL min.) serum from 1 SST. Collect specimen preferably 48 hours after dialysis. DONOR: 12 mL (10 mL min.) whole blood in 2 Yellow Top (ACD) tubes. Ship refrigerated. DO NOT FREEZE. Brief history in Organ Transplant Clinical information & Consent form form is mandatory.	Luminex X-Map technology	Doctor's prescription on the letter head with stamp AND 1 photograph each of recipient and donor is mandatory. Above specimen requirements are for a single donor being cross matched with a single recipient. For additional cross matches, send additional sç	Sample by Tue / Thu / Sat 9 am; Report Wed / Fri / Mon
L029	HLA HYPERSENSITIVITY TO ABACAVIR (B*57;01)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.
L030	HLA HYPERSENSITIVITY TO CARBAMAZEPINE & PHENYTOIN (B*15;02)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.
L031	HLA HYPERSENSITIVITY TO ALLOPURINOL (B*58;01)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube OR Buccal swab collected in a sterile container. Ship refrigerated. DO NOT FREEZE.	PCR, Sequencing		Sample Daily by 9 a.m; Report 7 days.

L027	HLA - NARCOLEPSY (DRB1*15, DQB1*06:02, DQA1*01:02)		4 mL (3 mL min.) Whole blood in 1 Lavender Top (EDTA) tube OR 6 mL (3 mL min.) Whole blood in 1 Yellow Top (ACD) tube. Ship refrigerated. DO NOT FREEZE.	Sequencing Based Typing (SBT)		Sample Daily by 9 a.m; Report 7 days.
L024	HLA - PANEL REACTIVE IgG ANTIBODIES (PRA), CLASS II, QUANTITATIVE		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Collect specimen preferably 48 hours after dialysis. Brief history in Organ Transplant Clinical information & Consent form is mandatory.	Luminex X-Map technology		Sample by Thu; Report Sat
L023	HLA - PANEL REACTIVE IgG ANTIBODIES (PRA), CLASS I, QUANTITATIVE		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Collect specimen preferably 48 hours after dialysis. Brief history in Organ Transplant Clinical information & Consent form is mandatory.	Luminex X-Map technology		Sample by Thu; Report Sat
L003	HLA - PANEL REACTIVE IgG ANTIBODIES (PRA), SCREEN	*Class I *Class II	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Collect specimen preferably 48 hours after dialysis. Brief history in Organ Transplant Clinical information form is mandatory.	Luminex X-Map technology		Sample by Mon / Wed / Fri 9 am Report Tue / Thu / Sat
L035	HLA SINGLE ANTIGEN BEAD ASSAY FOR CLASS II IgG ANTIBODIES		2 mL (1 mL min) serum from 1 SST .Ship refrigerated or frozen.	Luminex X-Map technology		Sample by Thu 9 am; Report Sat

L034	HLA SINGLE ANTIGEN BEAD ASSAY FOR CLASS I IgG ANTIBODIES		2 mL (1 mL min) serum from 1 SST . Ship refrigerated or frozen.	Luminex X-Map technology		Sample by Thu 9 am; Report Sat
L002	HLA T & B CELL LYMPHOCYTE CROSS MATCH		RECIPIENT: 6 mL (3 mL min.) whole blood in 1 Yellow Top (ACD) tube AND 3 mL (2 mL min.) serum from 1 SST. Collect specimen preferably 48 hours after dialysis.DONOR: 12 mL (10 mL min.) whole blood in 2 Yellow Top (ACD) tubes. Ship refrigerated. DO NOT FREEZE. Brief history in Organ Transplant Clinical information & Consent form is mandatory.	Complement Dependent Cytotoxicity with DTT treatment	Doctor's prescription on the letter head with stamp AND 1 photograph each of recipient and donor is mandatory. Above specimen requirements are for a single donor being cross matched with a single receipient. For additional cross matches, send additional sÇ	Sample Daily by 1 pm; Report Next day

L033	HLA T & B CELL FLOW CYTOMETRY, CROSSMATCH		RECIPIENT: 2 mL (1 mL min) serum from 1 SST AND 4 mL (2 mL min) whole blood from 1 Green Top (Sodium Heparin) tube. DONOR: 3 mL (2 mL min) whole blood in 3 Green Top (Sodium Heparin) tubes. Ship immediately at 18- 22° C. DO NOT FREEZE. Brief history in Organ Transplant Clinical information & Consent form is mandatory. Doctor's prescription on the letter head with stamp AND 1 photograph each of recipient and donor is mandatory. Above specimen requirements are for a single donor being cross matched with a single recipient. For additional cross matches, send additional specimens of donors. Request as additional tests. Collect specimen preferably 48 hours after dialysis.	Flow cytometry		Sample Daily by 9 am; Report Next day
------	--	--	--	----------------	--	---------------------------------------

J012	HISTOPATHOLOGY, BIOPSY, LARGE SPECIMEN	*Uterus with cervix *Uterus with cervix and adnexae / fibroids *Ovarian tumors / Large cysts *Eyeball (tumorous) *Non tumorous conditions of Thyroid gland / Testes / Kidney - nephrectomy / Intestinal resection *Lymph node (Block dissection) *Placenta *EyeÇ	Submit specimen in 10% Formal-saline. Formal-saline should cover the specimen completely. Ship at room temperature. Brief clinical history in Histopathology Requisition Form is mandatory. For patients with multiple specimens kindly register all specimens under single Lab number.	Processing, Staining & Microscopy	Special stains performed only on request, as additional tests.	Sample Daily by 6 pm; Report 4 days
N034	H1N1 (SWINE FLU) QUALITATIVE, REAL TIME PCR		Submit Nasal / Throat / Nasopharyngeal Swab inoculated in special viral transport medium available from LPL. Ship refrigerated in a special transport box available from LPL. Clinical History is mandatory	Real Time PCR		Sample Daily by 11 am: Report Same Day (except Sundays)

Z828	HOMA-IR; INSULIN RESISTANCE INDEX, C- PEPTIDE MODEL	* Glucose, Fasting *C- Peptide, Fasting *Beta cell function *Insulin sensitivity *HOMA-IR index	2 mL (1 mL min.) plasma from 1 Grey Top (Sodium fluoride) tube AND 2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is mandatory.	Hexokinase, CLIA		Sample Daily by 4 pm; Report Same Day
Z827	HOMA-IR ; INSULIN RESISTANCE INDEX	*Glucose-F * Insulin-F*Beta cell function* Insulin sensitivity*HOM A-IR index	2 mL (1 mL min.) plasma from 1 Grey Top (Sodium fluoride) tube AND 2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is mandatory.	Hexokinase, Chemiluminescent MicroparticleImmun oassay		Sample Daily by 4 pm; Report Same Day
R143	HOMOCYSTEINE , QUANTITATIVE, SERUM		2 mL (1 mL min.) serum from 1 SST. Place the sample immediately on ice & separate serum within 1 hour of collection. Ship refrigerated or frozen. Overnight fasting is preferred.	Chemiluminescent MicroparticleImmun oassay		Daily
U036	HOMOVANILLIC ACID (HVA), 24- HOUR URINE		50 mL (10 mL min.) aliquot of 24 Hour urine. Add 50% HCl drop by drop to bring the pH between 1–2. Do not use concentrated HCl. Ship refrigerated or frozen. Patient should strictly avoid L-Dopa at least 72 hours before and during collection.	High Performance LiquidChromatograp hy	Also see NEURONE SPECIFIC ENOLASE and VMA.	Sample by Wed / Sat 9 am; Report Fri / Tue
U049	HOMOVANILLIC ACID (HVA), RANDOM URINE		50 mL (10 mL min.) aliquot of random urine. Add 50% HCl drop by drop to bring the pH between 1–2. Do not use concentrated HCl. Ship refrigerated or frozen. Patient should strictly avoid L-Dopa at least 72 hours before and during collection.	High Performance LiquidChromatograp hy	Also see NEURONE SPECIFIC ENOLASE and VMA.	Sample by Wed / Sat 9 am; Report onFri / Tue

Z794	HPV, CHLAMYDIA TRACHOMATIS & NEISSERIA GONORRHEA PCR PANEL, QUALITATIVE		Submit cervical swab in a special container available from LPL. Ship refrigerated. DO NOT FREEZE	Real Time PCR		Sample by Wed / Sat 11 am. Report Fri /Tue
N040	HPV (HUMAN PAPILLOMA VIRUS) DETECTION & GENOTYPING, QUALITATIVE, PCR		Submit Cervical swab OR Oral swab in a special container available from LPL OR Paraffin embedded tissue block OR Fresh tissue in normal saline. Ship refrigerated. DO NOT FREEZE	Real Time PCR		Sample by Wed / Sat 11 am; Report Fri / Tue
RW097	HTLV - I / II ANTIBODY CONFIRMATION		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen in plastic screw capped vial. DO NOT USE GLASS TUBES.	Line Immunoassay		Sample by 7th of the month; Report after 2–3 weeks.
RW096	HTLV - I / II ANTIBODY, SCREEN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen in plastic screw capped vial. DO NOT USE GLASS TUBES.	Chemiluminescence, Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
N130	HUMAN HERPESVIRUS 6 (HHV6), QUALITATIVE PCR		3 ml (2 mL min.) Whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
N129	HUMAN PARECHOVIRUS QUALITATIVE, PCR		3 ml (2 mL min.) Whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat

N123	HUMAN HERPES VIRUS 7 (HHV7), QUALITATIVE PCR		3 ml (2 mL min.) Whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
N118	HUNTINGTON DISEASE MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days
R163	17 - HYDROXYPROG ESTERONE, NEWBORN SCREEN; CAH SCREEN		1 drop of heel prick blood each on 3 spots of filter paper available from LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen	Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm; Report 2 days
U032	17 - HYDROXY- CORTICOSTEROI DS (17-OHCS), 24-HOUR URINE		20 mL (10 mL min.) well mixed aliquot of 24-hour urine. Add appropriate volume of 50% HCl to maintain pH below 5. Record 24-hour urine volume on test request form and urine container. Ship refrigerated or frozen.	Column Chromatography		Sample Daily by 4 pm; Report 3 days
Z086	17 - HYDROXY- CORTICOSTEROI DS & 17 - KETOSTEROIDS, 24-HOUR URINE		20 mL (10 mL min.) well mixed aliquot of 24-hour urine. Add appropriate volume of 50% HCl to maintain pH below 5. Record 24 hour urine volume on test request form and urine container. Ship refrigerated or frozen.	Column Chromatography		Sample Daily by 4 pm; Report 3 days
RW100	17 - HYDROXYPREG NENOLONE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks

R053	17 - HYDROXYPROG ESTERONE (17- OHP)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. An early morning specimen is preferred.	Radioimmunoassay	Also see 17- HYDROXY- PROGESTERONESTIM ULATION BY ACTH.	Sample Tue / Sat by 9 am; Report Same Day
Z024	17- HYDROXYPROG ESTERONE (17-OHP) STIMULATION BY ACTH		2 mL (0.5 mL min.) serum from 1 SST for each timed specimen. Draw fasting early morning baseline sample; administer 250 ug of synthetic ACTH (Synacthen) IV, draw additional sample at 60 min. Ship refrigerated or frozen. Both tubes must be clearly marked with time drawn. Submit both tubes with 1 test request form.	Radioimmunoassay	Synthetic ACTH (Synacthen) is not supplied by LPL. Refer to Y Test codes :Y106 Basal Y108 60 min. Y109 Customized	Sample Tue / Sat by 9 am; Report Same Day
G186	HYDROXYPROLI NE, PLASMA		2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Separate plasma into sterile screw capped vial. Ship refrigerated or frozen. Overnight fasting is mandatory.	LC/MS	Non fasting samples acceptable for pediatric patients	Sample by Mon / Wed 5 pm; Report 3 days
Z165	HYPERCALCEMI A PANEL 2	*PTH (Intact) *Calcium, Serum *Phosphorus, Serum * Vitamin D 25 Hydroxy	6 mL (4 mL min.) serum from 2 SST's collected between 8-10 a.m. Avoid use of tourniquet and fist clenching during sampling. Separate serum and freeze immediately. Ship frozen. DO NOT THAW. Overnight fasting is preferred.	CMIA, Arsenazo III, CLIA, Molybdate UV	Also see PTH - RELATED PEPTIDE.	Daily
Z311	HYPERTENSION PANEL	*Urea *Creatinine *Glucose, Fasting *Electrolytes, serum *Lipid Profile Basic	2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 2 mL (1 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE. Minimum 12 hours fasting is mandatory.	Hexokinase, Spectrophotometry, I ndirect ISE		Daily

Z798	HYPERTHYROID / GRAVE'S DISEASE PANEL	*Free T4 * TSH * Total T3 *TSH Receptor Antibody	4 mL (3 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA, ECLIA		Sample Daily by 4 pm; Report Same day
Z175	HYPERCALCEMIA PANEL 1	*PTH (Intact) *Calcium, Serum *Calcium 24-hour Urine	3 mL (1 mLmin.) serum from 1 SST collected between 8-10 a.m. Separate serum and freeze immediately. Avoid use of tourniquet and fist clenching during sampling AND 10 mL (5 mL min.) aliquot of 24-hour urine collected with 10 mL of 50% HCl to maintain pH below 3. Overnight fasting is preferred. Ship frozen. DO NOT THAW. Record 24 hour urine volume & body weight on test request form and urine container.	CMIA, Arsenazo III		Daily
RW102	HYPERSENSITIVITY PNEUMONITIS SCREEN	*Aspergillus fumigatus *Micropolyspora faeni*Pigeon serum *S.viridis *Thermoactinomyces candidus *Thermoactinomyces vulgaris	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunodiffusion		Sample by 7th of the month; Reportafter 2–3 weeks
G201	HYPOGONADISM PANEL, MALES	*FSH *LH *Testosterone - LC-MS/MS	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Morning sample preferred.	CLIA, LC-MS/MS		Sample Daily by 4 pm; Report 2 days
G202	HYPOGONADISM PANEL, FEMALES	*FSH *LH *Estradiol- LC-MS/MS	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	CLIA, LC-MS/MS	Prefered sampling time is Day 2 or Day 3 of menstrual cycle in the morning	Sample Daily by 4 pm; Report 2 days

Z797	HYPOTHYROIDISM / HASHIMOTO'S DISEASE PANEL	*Free T4 * TSH *Anti TPO	4 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
RW103	IgA SUBCLASSES		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Nephelometry		Sample by 7th of the month; Report after 2–3 weeks
RW104	IGF BINDING PROTEIN–1; IGFBP–1		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Overnight fasting is mandatory.	Radioimmunoassay		Sample by 7th of the month; Report after 2–3 weeks
R105	IGF BINDING PROTEIN - 3; IGFBP - 3		2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW.	Chemiluminescent Immunoassay	Also see GROWTH DISORDER PANEL.	Sample by Mon / Thu 4 pm; Report Same day
RW105	IGF BINDING PROTEIN - 2; IGFBP - 2		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Overnight fasting is preferred.	Radioimmunoassay		Sample by 7th of the month; Report after 2–3 weeks
Z176	IGF - I & GROWTH HORMONE PANEL		3 mL (1.5 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW. Overnight fasting is mandatory.	Chemiluminescent Immunoassay	Also see GROWTH DISORDER PANEL.	Daily
R106	IGF - I; SOMATOMEDIN - C		2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW.	Chemiluminescent Immunoassay	Also see GROWTH DISORDER PANEL.	Daily
OS154	IgVH GENE MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Wed 11 am; Report 8 days

N073	IMATINIB RESISTANCE MUTATION ANALYSIS; IRMA		4 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Submit latest BCR-ABL Quantitative report along with sample.	PCR, Sequencing		Sample by Mon 11 am; Report Sat
G135	IMD PANEL, QUANTITATIVE, BLOOD	*Amino Acids *Organic Acids *Fatty Acids	1 drop of heel / finger prick blood each on 3 spots of filter paper available from LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen.	Tandem Mass Spectrometry		Sample by Mon / Wed / Fri 4 pm; Report 2 days
RW107	IMIPRAMINE	Includes metabolite Desipramine.	3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: >12 hours post oral dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
H250	IMMATURE PLATELET FRACTION (IPF)		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Fluorescent Flow Cytometry		Daily
RW108	IMMUNE COMPLEX DETECTION BY C1q BINDING		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum and freeze within 1 hour of collection. Ship frozen. DO NOT THAW. Do not transport in glass tubes. Do not use SST gel barrier tubes.	Enzyme Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
J132	IMMUNOHISTOCHEMISTRY, C4d		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J096	IMMUNOHISTO CHEMISTRY, CD 5		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J174	IMMUNOHISTO CHEMISTRY, CD 21		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
Z147	IMMUNE DEFICIENCY PANEL 2;	*CD3 *CD4 *CD8 Includes % and Absolute Counts	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry	This is the recommended test for CD4 / CD8 Counts and Ratio.	Sample Daily by 9 am; Report Same Day
Z148	IMMUNE DEFICIENCY PANEL 3; CD 3 COUNT	Includes % and Absolute Counts.	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am; Report Same Day

Z219	IMMUNE DEFICIENCY PANEL 4; CD 4 COUNTS	*CD3 *CD4Includes % and Absolute Counts	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18-22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
Z220	IMMUNE DEFICIENCY PANEL 5	*CD3 *CD19Includes % and Absolute Counts.	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18-22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
E005	IMMUNOFIXATION ELECTROPHORESIS (IFE), SERUM		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Grossly hemolysed specimens are not acceptable.	Agarose Gel Electrophoresis and Immunofixation	Includes qualitative identification of IgA/IgG/IgM Heavy Chains and Kappa/Lambda Light Chains. Does not include identification of IgD/IgE Heavy Chains. For IgD/IgE Heavy Chain typing, use test code E012.	Sample Daily by 9 am; Report Next day
E012	IMMUNOGLOBULIN IgD & IgE TYPING, SERUM		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Grossly hemolysed specimens are not acceptable.	Agarose Gel Electrophoresis and Immunofixation	For qualitative identification of IgD or IgE Heavy Chains. To be used as a follow-up to test codes E005 or E006, when indicated.	Sample Daily by 9 am; Report Next day

E013	IMMUNOGLOBULIN IgD & IgE TYPING, URINE		20 mL (10 mL min.) aliquot of 24-hour urine without preservative. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen.	Agarose Gel Electrophoresis and Immunofixation	For qualitative identification of IgD or IgE Heavy Chains. To be used as a follow-up to test codes E005 or E006, when indicated	Sample Daily by 9 am; Report 4 days
RW244	IMMUNOGLOBULIN IgD, SERUM		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Nephelometry		Sample by 7th of the month; Report after 2–3 weeks
S028	IMMUNOGLOBULIN IgE, SERUM		2 mL (1 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	ImmunoCAP	Also see ALLERGY tests & HISTAMINE tests.	Daily
RW109	IMMUNOGLOBULIN IgG SUBCLASSES		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Nephelometry		Sample by 7th of the month; Report after 2–3 weeks
B091	IMMUNOGLOBULIN IgG, CSF		2 mL (1 mL min.) Cerebrospinal fluid in a sterile screw capped container. Ship refrigerated or frozen.	Nephelometry		Sample Mon through Sat by 4 pm; Report Same day
J153	IMMUNOHISTOCHEMISTRY, CD 1a		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J097	IMMUNOHISTO CHEMISTRY, CD 10		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
Z146	IMMUNE DEFICIENCY PANEL 1	*CD3 *CD4 *CD8 *CD19 *CD16+/CD56 Includes % and Absolute Counts.	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
E006	IMMUNOFIXATION ELECTROPHORESIS (IFE), 24-HOUR URINE		20 mL (10 mL min.) aliquot of 24-hour urine without preservative. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen.	Agarose Gel Electrophoresis and Immunofixation	Includes qualitative identification of IgA/IgG/IgM Heavy Chains and Kappa/Lambda Light Chains. Does not include identification of IgD/IgE Heavy Chains. For IgD/IgE Heavy Chain typing, use test code E013.	Sample Daily by 9 am; Report 4 days
B078	IMMUNOGLOBULIN IgA, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry		Daily

B190	IMMUNOGLOBULIN IgG SUBCLASS 4		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Nephelometry		Sample by Mon/Thu 4 pm; Report Same day
Z180	IMMUNOGLOBULIN IgG SYNTHESIS INDEX & RATE		2 mL (1 mL min.) Cerebrospinal fluid in a sterile screw capped container AND 2 mL (1 mL min.) serum from 1 SST. The collection date and time must be the same for both specimens. Ship refrigerated or frozen.	BCG, Nephelometry, Immunoturbidimetry	Also see MULTIPLE SCLEROSIS PANELS; MYELIN-ASSOCIATED GLYCOPROTEIN; MYELIN-ASSOCIATED GLYCOPROTEIN-SULFATED GLUCURONIC PARAGLOBOSIDE AUTO-ANTIBODY, IgM.	Sample Daily by 4 pm; Report Next day
RW111	IMMUNOGLOBULIN PROFILE, CSF	*IgG *IgA *IgM	2 mL (1 mL min.) Cerebrospinal fluid in a sterile screw capped container. Ship refrigerated or frozen. Overnight fasting is preferred.	Fixed Time Nephelometry		Sample by 7th of the month; Report after 2–3 weeks
Z178	IMMUNOGLOBULIN PROFILE, SERUM	*IgG *IgA *IgM	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry		Daily

J150	IMMUNOHISTO CHEMISTRY, ANY 3 MARKERS		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J076	IMMUNOHISTO CHEMISTRY, ACTIN (Smooth Muscle Actin)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J137	IMMUNOHISTO CHEMISTRY, AFP (Alpha Feto Protein)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days

J180	IMMUNOHISTO CHEMISTRY, AMACR		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 days Tissues : 7 days
J043	IMMUNOHISTO CHEMISTRY, HER-2 / ERBB2		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J029	IMMUNOHISTO CHEMISTRY, CD 45 (LCA)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J157	IMMUNOHISTO CHEMISTRY, CK 5 / 6		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J127	IMMUNOHISTO CHEMISTRY, CYTOKERATIN 20		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J107	IMMUNOHISTO CHEMISTRY, FSH		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J044	IMMUNOHISTO CHEMISTRY, p53		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J045	IMMUNOHISTO CHEMISTRY, pS2		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
RW110	IMMUNOGLOBULIN IgM, CSF		2 ml (1 mL min.) Cerebrospinal fluid in a sterile screw capped container. Ship refrigerated or frozen.	Rate Nephelometry		Sample by 7th of the month; Report after 2–3 weeks
B079	IMMUNOGLOBULIN IgM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry		Daily
J151	IMMUNOHISTO CHEMISTRY, ANY 5 MARKERS		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

Z700	IMMUNOHISTO CHEMISTRY, CUSTOM PANEL		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology/ Cytopathology report. Indicate site of biopsy and provide Clinical history.	Imunohistochemistry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J196	IMMUNOHISTO CHEMISTRY, Ber-Ep4 (EpCAM)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J143	IMMUNOHISTO CHEMISTRY, 34-BETA E12		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J186	IMMUNOHISTO CHEMISTRY, CALDESMON		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J141	IMMUNOHISTO CHEMISTRY, CALRETININ		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J163	IMMUNOHISTO CHEMISTRY, CD 4		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J100	IMMUNOHISTO CHEMISTRY, CD 34		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J037	IMMUNOHISTO CHEMISTRY, CD 45 Ro (PAN T)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J065	IMMUNOHISTO CHEMISTRY, CD 79a		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J102	IMMUNOHISTO CHEMISTRY, CD 99 (MIC-2)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J146	IMMUNOHISTO CHEMISTRY, CD 138		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J111	IMMUNOHISTO CHEMISTRY, CHROMOGRANI N A		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J187	IMMUNOHISTO CHEMISTRY, CDX-2		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J172	IMMUNOHISTO CHEMISTRY, IgG		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J046	IMMUNOHISTO CHEMISTRY, Ki67		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J128	IMMUNOHISTO CHEMISTRY, TdT		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J178	IMMUNOHISTO CHEMISTRY, BCL-6		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J158	IMMUNOHISTO CHEMISTRY, CALCITONIN		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J041	IMMUNOHISTO CHEMISTRY, CATHEPSIN D		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J066	IMMUNOHISTO CHEMISTRY, CD 3		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J062	IMMUNOHISTO CHEMISTRY, CD 20		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J175	IMMUNOHISTO CHEMISTRY, CD 35		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J101	IMMUNOHISTO CHEMISTRY, CD 117		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J130	IMMUNOHISTO CHEMISTRY, CYCLIN-D1		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J017	IMMUNOHISTOCHEMISTRY, CYTOKERATIN (PANCYTOKERATIN)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J188	IMMUNOHISTOCHEMISTRY, CYTOKERATIN 19		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J167	IMMUNOHISTOCHEMISTRY, CYTOMEGALOVIRUS (CMV)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J030	IMMUNOHISTO CHEMISTRY, EPITHELIAL MEMBRANE ANTIGEN (EMA)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J034	IMMUNOHISTO CHEMISTRY, GLIAL FIBRILLARY ACIDIC PROTEIN (GFAP)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J199	IMMUNOHISTO CHEMISTRY, GLYPICAN - 3		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J105	IMMUNOHISTO CHEMISTRY, GROWTH HORMONE		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J191	IMMUNOHISTO CHEMISTRY, HUMAN CHORIONIC GONADOTROPI N (HCG)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J138	IMMUNOHISTO CHEMISTRY, LAMBDA LIGHT CHAINS		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
G169	IMMUNOSUPPR ESSANT DRUG PROFILE 1	*Tacrolimus *Sirolimus	3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated . Ideal sampling time: Just before the next dose (Trough levels). Provide date & time of sampling, dose taken & time of dosing.	LC-MS/MS		Sample Daily by 5 pm; Report Same Day

G168	IMMUNOSUPPRESSANT DRUG PROFILE 4	*Tacrolimus *Everolimus	3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated .Ideal sampling time: Just before the next dose (Trough levels). Provide date & time of sampling, dose taken & time of dosing.	LC-MS/MS		Sample Daily by 5 pm; Report Same Day
J140	IMMUNOHISTOCHEMISTRY, INHIBIN		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J159	IMMUNOHISTOCHEMISTRY, p63		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J093	IMMUNOHISTOCHEMISTRY, PROSTATE SPECIFIC ANTIGEN (PSA)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J166	IMMUNOHISTO CHEMISTRY, SERUM ASSOCIATED AMYLOID (SAA)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J067	IMMUNOHISTO CHEMISTRY, THYROGLOBULIN		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
G167	IMMUNOSUPPRESSANT DRUG PROFILE 3	*Cyclosporine *Tacrolimus	3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated . Ideal sampling time: Just before the next dose (Trough levels). Provide date & time of sampling, dose taken & time of dosing.	LC-MS/MS		Sample Daily by 5 pm; Report Same Day
J126	IMMUNOHISTO CHEMISTRY, CYTOKERATIN 7		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J160	IMMUNOHISTO CHEMISTRY, HEPATOCYTE SPECIFIC ANTIGEN (HSA)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J192	IMMUNOHISTO CHEMISTRY, INI- 1		Submit tumour tissue in 10% formal-saline OR Formalin Fixed paraffin embedded block. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of Biopsy and provide clinical history	Immunohistochemist ry		Sample Daily by 6 pmReport Block:5 daysReport Tissue: 7 Days
J189	IMMUNOHISTO CHEMISTRY, MELAN-A		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J142	IMMUNOHISTO CHEMISTRY, MYELOPEROXID ASE		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J190	IMMUNOHISTO CHEMISTRY, MYOGENIN		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J036	IMMUNOHISTO CHEMISTRY, NEURON SPECIFIC ENOLASE (NSE)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J179	IMMUNOHISTO CHEMISTRY, PAX -5		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J010	IMMUNOHISTO CHEMISTRY, PROGESTERONE RECEPTOR (PR)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J031	IMMUNOHISTO CHEMISTRY, S100		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J193	IMMUNOHISTO CHEMISTRY, UROPLAKIN - II		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J144	IMMUNOHISTO CHEMISTRY, TTF-1		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J136	IMMUNOHISTO CHEMISTRY, WT-1		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
G166	IMMUNOSUPPRESSANT DRUG PROFILE 2	*Cyclosporine *Sirolimus	3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated . Ideal sampling time: Just before the next dose (Trough levels). Provide date & time of sampling, dose taken & time of dosing.	LC-MS/MS		Sample Daily by 5 pm; Report Same Day
J077	IMMUNOHISTO CHEMISTRY, EPIDERMAL GROWTH FACTOR RECEPTOR (EGFR)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J103	IMMUNOHISTO CHEMISTRY, HMB 45		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J139	IMMUNOHISTO CHEMISTRY, KAPPA LIGHT CHAINS		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J125	IMMUNOHISTO CHEMISTRY, SYNAPTOPHYSI N		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J032	IMMUNOHISTO CHEMISTRY, VIMENTIN		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J018	IMMUNOHISTO CHEMISTRY, DESMIN		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J171	IMMUNOHISTO CHEMISTRY, IgA		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J061	IMMUNOHISTO CHEMISTRY, CD 15		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J099	IMMUNOHISTO CHEMISTRY, CD 31		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J064	IMMUNOHISTO CHEMISTRY, CD 68		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J108	IMMUNOHISTO CHEMISTRY, ACTH		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J094	IMMUNOHISTO CHEMISTRY, BCL-2		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J195	IMMUNOHISTO CHEMISTRT, C- MYC		Submit tumour tissue in 10% formal-saline OR Formalin Fixed paraffin embedded block. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of Biopsy and provide clinical history	Immunohistochemist ry		Sample Daily by 6 pmReport Block:5 daysReport Tissue: 7 Days
J098	IMMUNOHISTO CHEMISTRY, CD 23		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days

J063	IMMUNOHISTO CHEMISTRY, CD 30		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J161	IMMUNOHISTO CHEMISTRY, CD 56		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J145	IMMUNOHISTO CHEMISTRY, CEA		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S meas : 5 daysTissues : 7 days
J194	IMMUNOHISTO CHEMISTRY, BETA-CATENIN		Specimen: Submit tumour tissue in 10% formal-saline OR Formalin Fixed paraffin embedded block. Ship at room temperature.Provide a copy of the Histopathology / Cytopathology report. Indicate site of Biopsy and provide clinical history.	Immunohistochemist ry		Sample Daily by 6 pmReport Blocks :5 DaysTissues : 7 Days

J173	IMMUNOHISTO CHEMISTRY, DOG-1		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
J197	IMMUNOHISTO CHEMISTRY, GATA-3		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days
B077	IMMUNOGLOBULIN IgG, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry		Daily
J147	IMMUNOHISTO CHEMISTRY, ALK-1		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemistry	Use test code J113 for Photomicrograph	Sample Daily by 6 pm Report Block/Sections/Smears : 5 days Tissues : 7 days

J162	IMMUNOHISTO CHEMISTRY, CD 8		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J033	IMMUNOHISTO CHEMISTRY, NEUROFILAMEN T POLYPEPTIDE		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J181	IMMUNOHISTO CHEMISTRY, Oct 2		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days

J148	IMMUNOHISTO CHEMISTRY, POLYOMA BK VIRUS (SV 40)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J104	IMMUNOHISTO CHEMISTRY, PROLACTIN		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 5 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
Z680	IMUNOHISTOCH EMISTRY, MICROSATELLIT E INSTABILITY PANEL	*MLH-1 *MSH-2 *MSH-6 *PMS-2	Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block OR Alcohol fixed Fine Needle Aspirate Smears OR 10 unstained sections on Poly-L-Lysine Coated Slides. Ship at room temperature. Provide a copy of the Histopathology/Cytopathology report. Indicate site of biopsy and provide Clinical history.	Immunohistochemist ry	Use test code J113 forPhotomicrograph	Sample Daily by 6 pmReport Block/Sections/S mears : 5 daysTissues : 7 days
J198	IMUNOHISTOCH EMISTRY, NAPSIN-A		Submit tumour tissue in 10% formal-saline OR Formalin Fixed paraffin embedded block. Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of Biopsy and provide clinical history	Immunohistochemist ry		Sample Daily by 6 pmReport Block:5 daysReport Tissue: 7 Days

P068	INDIA INK PREPARATION, FLUID		1 mL (0.5 mL min.) Fluid. Ship at 18–22°C.	Light Microscopy	Also see CULTURE, FUNGUS and FUNGUS EXAMINATION, ROUTINE	Daily
P051	INDIA INK PREPARATION, CSF		1 mL (0.5 mL) Cerebrospinal fluid. Ship at 18-22°C.	Light Microscopy	Also see CULTURE, FUNGUS and FUNGUS EXAMINATION, ROUTINE	Daily
Z812	INFERTILITY COMPREHENSIVE PANEL, FEMALES	*FSH *LH *Prolactin *Testosterone, Total * TSH Ultrasensitive * Estradiol * Anti Mullerian Hormone (AMH)	4 mL (2 mL min.) serum from 1 SST. Collect specimen 4 hrs after the patient has awakened. Overnight fasting is preferred .Ship refrigerated or frozen.	Chemiluminescent Immunoassay, EIA	In females, preferred sampling time is Day 2 / Day 3 of the menstrual cycle.	Sample Daily by 9 am; Report Same day
Z605	INFERTILITY PANEL, MALES	*FSH *LH *Testosterone, Total *TSH Ultrasensitive	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
Z179	INFERTILITY SCREENING PANEL, FEMALES	*FSH *LH *Prolactin *Testosterone, Total * TSH Ultrasensitive	3 mL (1.5 mL min.) serum from 1 SST. Collect specimen 4 hrs after the patient has awakened. Overnight fasting is preferred. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	In females, preferred sampling time is Day 2 / Day 3 of the menstrual cycle.	Daily
F025	INHIBIN B		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Day 2 / 3 sample of the menstrual cycle is preferred to assess Ovarian reserve.	Sample by Mon 9 am; Report Same day

F022	INHIBIN A, REPRODUCTIVE MARKER		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample Daily by 9 am; Report Same Day
F021	INHIBIN A, PREGNANCY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample Daily by 9 am; Report Same Day
J185	IN-SITU- HYBRIDIZATION - EPSTEIN BARR VIRUS EARLY RNA (EBER-ISH)		Submit tumour tissue in 10% Formal-saline OR Formalin fixed paraffin embedded block Ship at room temperature. Provide a copy of the Histopathology / Cytopathology report. Indicate site of biopsy and provide Clinical history.	In-Situ-Hybridization		Sample Daily by 6 pm Report Block : 5 days Tissues : 7 days
S133	INSULIN ANTIBODIES		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Sat 9 am; Report Same day
R042	INSULIN, PP		2 mL (1 mL min.) serum from 1 SST. Samples to be collected 2 hours post meal. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay	Also see INSULIN RESPONSE TO GLUCOSE Test	Daily
Z385	INSULIN, FASTING (F) & POST PRANDIAL (PP), PANEL		2 mL (1 mL min.) serum from 1 SST each for F & PP samples. Overnight fasting is mandatory. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay	Also see INSULIN RESPONSE TO GLUCOSE Test	Daily
R041	INSULIN, RANDOM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay	Also see NSULIN RESPONSE TO GLUCOSE Test	Daily
RW112	INSULIN ANTIBODIES, HIGHLY SENSITIVE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum within 1 hour of collection. Ship refrigerated or frozen.	Radio binding Assay		Sample by 7th of the month; Reportafter 2–3 weeks

Z084	INSULIN RESPONSE TO GLUCOSE, 5 HOURS		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 2 mL (1 mL min.) corresponding serum from 1 SST for each timed specimen set. Draw baseline fasting specimens. Dissolve 75g Glucose in adults / 1.75 g / Kg body weight in children (not exceeding 75g) in 300mL water. Administer orally over a period of 5 minutes. Draw additional samples at 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 4.0 and 5.0 hours. Note time drawn and tests required on each container label and test request form. Overnight fasting mandatory. Ship refrigerated or frozen. Urine specimens NOT REQUIRED.	Chemiluminescent MicroparticleImmun oassay, Hexokinase	Refer to Y Codes:Y084 BasalY085 0.5 hrY086 1.0 hrY087 1.5 hrY088 2.0 hrY089 2.5 hrY090 3.0 hrY091 4.0 hrY092 5.0 hrY093 Customized	Daily
RW113	INSULIN, FREE (BIOACTIVE)		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Overnight fasting is mandatory. Ship refrigerated or frozen.	Immunoassay	The free (bioactive) insulin test is recommended for patients who have detectable levels of insulin antibodies.	Sample by 7th of the month; Reportafter 2–3 weeks
R039	INSULIN, FASTING		2 mL (1 mL min.) serum from 1 SST. Overnight fasting is mandatory. Ship refrigerated or frozen.	Chemiluminescent MicroparticleImmun oassay	Also see INSULIN RESPONSE TO GLUCOSE Test	Daily

RW115	INTERLEUKIN - 6; IL-6		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube OR 2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate serum / plasma within 2 hours of collection in a plastic transport vial and freeze immediately. Ship frozen. DO NOT THAW.	Enzyme Immunoassay	Cytokine levels may demonstrate diurnal variation. Recommend cytokine levels be determined at the same time of day for improved longitudinal comparison.	Sample by 7th of the month; Report after 2–3 weeks
RW256	INTERFERON - BETA IgG, MAID		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Collect sample at least 8 hours after Interferon injection. Verify medication and provide dose & duration of therapy. Ship refrigerated or frozen.	Multianalyte Immunodetection (MAID)		Sample by 7th of the month; Report after 2–3 weeks
RW114	INTERLEUKIN - 2 RECEPTOR; IL-2R		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum within 1 hour of collection in a plastic transport vial and freeze immediately. Ship frozen. DO NOT THAW.	Enzyme Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
N050	INTERLEUKIN 28B rs12979860 & rs 8099917 GENOTYPING, QUALITATIVE, PCR		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Thu / Mon 9 am; Report Sat / Wed
S244	INTRINSIC FACTOR ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	EIA		Sample by Mon 9 am; report Same day

N046	Inv 16 (p13q22) / t(16;16)(p13;q22) GENE REARRANGEMENT QUALITATIVE, PCR		3 mL (2 mL min.) Whole blood / Bone Marrow from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical History is mandatory.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
H020	IRON STUDIES	*SERUM IRON *TOTAL IRON BINDING CAPACITY * TRANSFERRIN SATURATION)	3 mL (1.5 mL min) serum from 1 SST. Avoid hemolysis. Ship refrigerated or frozen. Overnight fasting is preferred. Avoid iron supplements for minimum 7 days prior to specimen collection.	Spectrophotometry	Also see FERRITIN.	Daily
Z783	IRON STUDIES MONITORING PANEL	*Serum Iron *Total Iron Binding Capacity *Transferrin Saturation *Ferritin	4 mL (1.5 mL min) serum from 1 SST. Avoid hemolysis. Ship refrigerated or frozen. Overnight fasting is preferred. Avoid iron supplements for minimum 7 days prior to specimen collection.	Spectrophotometry, CLIA		Daily
S134	ISLET CELL ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Indirect Immunofluorescence Assay		Sample Daily by 11 am; Report Same day
S186	ISLET CELL ANTIBODY, IN DILUTIONS		2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Indirect Immunofluorescence Assay		Sample Daily by 11 am; Report Next day
N122	JAK 2 EXON 12, MUTATION DETECTION		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give brief clinical history.	PCR, Fragment Analysis		Sample by Mon / Thu 11 am; Report Wed / Sat

N039	JAK 2 MUTATION DETECTION, QUALITATIVE PCR		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory.	Real Time PCR		Sample by Mon / Thu 11 am; ReportWed / Sat
Z839	JAK 2 V617F REFLEX TO JAK 2 EXON 12, MUTATION DETECTION		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give brief clinical history	Real Time PCR, Fragment Analysis	0	Sample by Mon 11 am; ReportSat
Z840	JAK 2 V617F, CALR & MPL, MUTATION DETECTION PROFILE		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give brief clinical history	Real Time PCR, Sequencing, Fragment Analysis	0	Sample by Mon 11 am; ReportSat
OS157	JAPANESE ENCEPHALITIS VIRUS DETECTION , PCR		3 mL (2 mL min.) CSF in a sterile screw capped container. Ship refrigerated or frozen.	Real Time PCR		Sample by Wed 11 am; Report 7 days
OS156	JC/BK VIRUS DETECTION , PCR		4 mL (2 mL min.) Serum from 1 SST OR 4 mL (2 mL min.) plasma from 1 Lavendar top (EDTA) tube OR 3 mL (2 mL min.) CSF / Urine in a sterile screw capped container OR Tissue in normal saline. Ship refrigerated or frozen.	PCR		Sample by Wed/Sat 11 am; Report 6 days
S080	Jo-1 ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see EXTRACTABLE NUCLEARANTIGENS.	Sample by Tue / Fri 9 am; Report Same day

B148	KAPPA / LAMBDA LIGHT CHAINS, FREE, SERUM	*Kappa Light Chains * Lambda Light chains* Kappa Lambda Ratio	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Nephelometry		Sample Mon through Sat by 4 pm; Report Same day
U079	KAPPA / LAMBDA LIGHT CHAINS, FREE, URINE	*Kappa Light Chains* Lambda Light chains	10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	Nephelometry		Sample Mon through Sat by 4 pm; Report Same day
U005	17 – KETOSTEROIDS, 24-HOUR URINE		20 mL (10 mL min.) well mixed aliquot of 24-hour urine. Add appropriate volume of 50% HCl to maintain pH below 5. Record 24 hour urine volume on test request form and container. Ship refrigerated or frozen.	Column Chromatography		Sample Daily by 4 pm; Report 3 days
Z379	KIDNEY DIALYSIS PANEL	*Glucose R *Creatinine *Uric Acid *Cholesterol*Pro tein Total *Albumin *A:G Ratio *Alkaline Phosphatase *Calcium *Phosphorus *Electrolytes, serum *BUN *CBC *Urine R/E	3 mL (1.5 mL min.) serum from 1 SST AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred AND 10 mL (5 mL min.) aliquot of first morning urine. Ship refrigerated. DO NOT FREEZE.	Spectrophotometry, Indirect ISE, Hexokinase		Daily

Z784	KIDNEY PANEL EXTENDED	*Urea *Creatinine *Uric Acid *Protein Total *Albumin *A:G Ratio *Alkaline Phosphatase *Calcium *Phosphorus *Electrolytes, serum *Bicarbonate *Glucose Random	4 mL (1.5 mL min) serum from 1 SST AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Avoid use of tourniquet and fist clenching during sampling. Avoid hemolysis. Overnight fasting is preferred. Ship refrigerated or frozen.	Spectrophotometry, Hexokinase, Indirect ISE		Daily
Z007	KIDNEY PANEL; KFT	*Urea *Creatinine *Uric Acid *Protein Total*Albumin *A:G Ratio *Alkaline Phosphatase *Calcium *Phosphorus *Electrolytes, serum	3 mL (1.5 mL min.) serum from 1 SST. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship refrigerated or frozen.	Spectrophotometry, Indirect ISE		Daily

Z230	KIDNEY STONE FORMATION, DIAGNOSTIC PANEL	*pH *Creatinine *Calcium *Magnesium *Citrate*Oxalate *Osmolality *Electrolytes, urine *Uric Acid *Phosphorus	Day 1: 50 mL (30 mL min.) aliquot of 24-hour urine. Do not use preservatives. Mention total volume on Test Request Form & label container as Non-acidified urine. Day 2: 50 mL (30 mL min.) aliquot of 24 hour urine. Collect urine with 25 mL 50% concentrated HCl to maintain pH below 3. Mention total volume on Test Request Form & label container as Acidified urine. Ship both specimens frozen. DO NOT THAW. Record patient weight on Test Request Form.	Spectrophotometry, Enzymatic, Indirect ISE		Sample by Wed / Sat 9 am; Report Same day
Z231	KIDNEY STONE FORMATION, THERAPEUTIC MONITORING PANEL	*pH *Creatinine *Calcium *Magnesium *Citrate*Oxalate *Osmolality *Sodium *Uric Acid *Phosphorus	Day 1: 50 mL (30 mL min.) aliquot of 24-hour urine. Do not use preservatives. Mention total volume on Test Request Form & label container as Non-acidified urine. Day 2: 50 mL (30 mL min.) aliquot of 24 hour urine. Collect urine with 25 mL 50% concentrated HCl to maintain pH below 3. Mention total volume on Test Request Form & label container as Acidified urine. Ship both specimens frozen. DO NOT THAW. Record patient weight on Test Request Form.	Spectrophotometry, Enzymatic, Indirect ISE		Sample by Wed / Sat 9 am; Report Same day
OS158	KRAS CODON 61 MUTATION DETECTION		Submit fresh tissue in normal saline OR formalin fixed paraffin embedded block	PCR		Sample by Mon 11 am; Report 10 days
N067	KRAS MUTATION, CODON 12 & 13 FRAGMENT ANALYSIS		Submit formalin fixed paraffin embedded tissue block. Ship at room temperature. Clinical history & Histopathology report is mandatory.	PCR, Fragment analysis		Sample by Tue 11 am; Report Mon

Z816	KRAS & NRAS MUTATION ANALYSIS PANEL		Submit formalin fixed paraffin embedded tissue block. Ship at room temperature. Block should have at least 10% tumor tissue.	PCR, Fragment analysis and Sequencing		Sample by Tue 11 am; Report Mon
B113	LACTATE, PLASMA, POST EXERCISE		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube without applying tourniquet. Separate plasma immediately. Ship refrigerated or frozen.	Spectrophotometry	Separation of plasma at source is mandatory.	Daily
B099	LACTATE, PLASMA		2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube without applying tourniquet. Separate plasma immediately. Ship refrigerated or frozen.	Spectrophotometry	Separation of plasma at source is mandatory. Use Test Code B113 for Lactate Stimulation by Exercise.	Daily
B098	LACTATE, CSF		1 mL (0.5 mL min) Cerebrospinal fluid. Ship refrigerated or frozen.	Spectrophotometry		Daily
RW118	LAMOTRIGINE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Recommended specimen: 0.5–1 hour before next dose at steady state. Ship refrigerated or frozen. Do not use SST gel barrier tubes.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
P052	LDH (LACTATE DEHYDROGENASE), FLUID		5 mL (2 mL min.) fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	IFCC		Daily
B031	LDH; LACTATE DEHYDROGENASE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE. Hemolysed specimens are not acceptable.	IFCC		Daily
B129	LDL CHOLESTEROL, DIRECT		2 mL (0.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzymatic selective protection	Also see LIPID PROFILE TESTS.	Daily

E014	LDL SUBFRACTIONS		2 mL (1.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated. DO NOT FREEZE. Minimum 12 hours fasting is mandatory. Grossly hemolysed samples are not acceptable.	Polyacrylamide Gel Electrophoresis(PAGE)		Sample by Fri 9 am; Report Same day
C020	LEAD, BLOOD		4 mL (2 mL min.) whole blood from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Use powderlessgloves during specimen collection.	Atomic Absorption Spectrometry withZeeman correction	Use only specified container available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; ReportNext day
C021	LEAD, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in kit. Measure 24 hour total volume and record on vial and Test request form. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	Atomic Absorption Spectrometry withZeeman correction	Use only specified container available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; ReportNext day

Z038	LEAD POISONING PANEL	*5-Amino Levulinic Acid (5- ALA) *Lead Blood	10 mL (5 mL min.) aliquot of 24 hour urine collected in a dark coloured container without preservative Ship frozen. DO NOT THAW. Record 24 hour volume on test request form and urine container. Refrigerate during collection. AND 4 mL (2 mL min.) whole blood from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Use powderless gloves during specimen collection.	Column chromatography, AAS with Zeeman correction	Use only specified container available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample by Tue / Fri 9 am; Report Next day
C022	LEAD, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container kit available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	Atomic Absorption Spectrometry with Zeeman correction	Use only specified container available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
H036	LE CELL PHENOMENON		10 mL (5 mL min.) whole blood in a plastic container with 4 glass beads. Rotate for 30 min to defibrinate. Ship refrigerated. DO NOT FREEZE.	Romanowsky staining	Also see ANTINUCLEAR ANTI BODY/FACTOR.	Sample Daily by 6 pm; Report Next day
M217	LEGIONELLA ANTIGEN, URINE		5 mL aliquot of random urine in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pm Report Same day

RW121	LEGIONELLA ANTIGEN, URINE		5 mL (3 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
RW119	LEGIONELLA PNEUMOPHILA TOTAL ANTIBODY, SERUM	Serogroups 1, 2-6 & 8	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks
S025	LEISHMANIA (KALA AZAR) ANTIBODY IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Sat 9 am; Report Same day
B164	LEPTIN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Sat 9 am; Report Same day
Z302	LEPTOSPIRA ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min) serum from 1 SST. Ship refrigerated or frozen	Enzyme Immunoassay		Sample by Tue / Sat 9 am; Report Same day
S157	LEPTOSPIRA ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue / Sat 9 am; Report Same day
S197	LEPTOSPIRA ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue / Sat 9 am; Report Same day
OS159	LEPTOSPIRA DETECTION, PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube OR 10 mL Urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Mon / Thu 11 am; Report Wed / Sat

H050	LEUCOCYTE ALKALINE PHOSPHATASE SCORE; LAP SCORE		5 fresh blood smears from finger prick, fixed in 95% ethanol for 30 min. Air dry and ship at room temperature immediately.	Cytochemistry	Whole blood in EDTA tubes not accepted. Only fresh blood smears to be sent.	Sample Daily by 9 am; Report Next day
Z151	LEUKEMIA / LYMPHOMA DIAGNOSTIC PANEL: CHRONIC LYMPHOPROLIF ERATIVE DISORDERS, T & B CELL		3 mL (2 mL min.) whole blood each in 1 Lavender Top (EDTA) tube AND 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone Marrow AND Aspirate Smear. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. Clinical History Mandatory.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
Z601	LEUKEMIA DIAGNOSTIC COMPREHENSIV E PROFILE 4	*Chromosome analysis for hematological malignancy *AML categorization by PCR / BCR ABL gene rearrangement, PCR Qualitative	2 mL (1.5 mL min.) Bone Marrow Aspirate / Whole Blood in 1 Green Top (Sodium Heparin) tube AND 3 mL (2 mL min.) Whole blood / Bone Marrow Aspirate in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Test Request Form .	Culture, Microscopy, Karyotype, Real Time PCR		Sample Daily by 11 am; Report 7 days
Z665	LEUKEMIA DIAGNOSTIC COMPREHENSIV E PROFILE 6	*Chromosome analysis for hematological malignancy *BCR ABL gene rearrangement, PCR Qualitative	2 mL (1.5 mL min.) Bone Marrow Aspirate / Whole Blood in 1 Green Top (Sodium Heparin) tube AND 3 mL (2 mL min.) Whole blood / Bone Marrow Aspirate in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Test Request Form.	Culture, Microscopy, Karyotype, Real Time PCR		Sample Daily by 11 am; Report 7 days

Z512	LEUKEMIA DIAGNOSTIC PANEL, CUSTOMIZED (ANY 10 MARKERS)		3 mL (2 mL min.) whole blood each in 1 Lavender Top (EDTA) tube AND 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone Marrow AND Aspirate Smear. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. Clinical History Mandatory.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
Z282	LEUKEMIA DIAGNOSTIC PANEL: CLL/HCL/SLL (BASIC)		3 mL (2 mL min.) whole blood each in 1 Lavender Top (EDTA) tube AND 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone Marrow AND Aspirate Smear. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. Clinical History Mandatory.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
Z538	LEUKEMIA DIAGNOSTIC COMPREHENSIV E PROFILE 3	*Leukemia diagnostic panel, Acute Leukemia - T, B or Myeloid, *Chromosome analysis for hematological malignancy, *AML categorization by PCR / BCR ABL gene rearrangement, PCR Qualitative	2 mL (1.5 mL min.) Bone Marrow Aspirate in 1 Green Top (Sodium Heparin) tube AND 3 mL (2 mL min.) Whole blood / Bone Marrow in 1 Lavender Top (EDTA) tube. Ship at specified temperatures. Give clinical history on Chromosome Test Request Form	Flow Cytometry, Culture, Microscopy, Karyotype, Real Time PCR		Sample Daily by 11 am; Report 7 days

Z780	LEUKEMIA DIAGNOSTIC PANEL: PLASMA CELL DYSCRASIAS		3 mL (2 mL min.) whole blood each in 1 Lavender Top (EDTA) tube AND 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone Marrow AND Aspirate Smear. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. Clinical History Mandatory.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
Z662	LEUKEMIA GENETIC PROFILE - ANY 6 MARKERS, PCR, QUALITATIVE	Any 6 of the following markers can be selected - BCR-ABL, PML-RARA, AML-ETO, INV16, NPM1, FLT3, MLLAF9, MLL-AF4, MLL-ENL, t(12;21), t(1;19).	3 mL (2 mL min.) Whole blood / Bone Marrow Aspirate in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Specify the 6 markers selected on Test Request Form.	Real Time PCR	Any 6 of the following markers can be selected - BCR-ABL, PML-RARA, AML-ETO, INV16, NPM1, FLT3, MLL-AF9, MLL-AF4, MLL-ENL, t(12;21), t(1;19).	Sample by Tue / Fri 11 am; Report Thu/ Mon

Z536	LEUKEMIA DIAGNOSTIC COMPREHENSIV E PROFILE 1	*Microscopy of Bone Marrow Aspirate & Peripheral Blood *Special Stains *CBC * Histopathology *Leukemia diagnostic panel, Acute Leukemia - T, B or Myeloid *Chromosome analysis for hematological malignancy *AML categorization by PCR / BCR ABL gene rearrangeÇ	4 mL (2.5 mL min.) Bone Marrow Aspirate in 2 Green Top (Sodium Heparin) tubes AND 4 Bone Marrow Aspirate Smears. Air dry & fix smears in 100% Methanol for 30 minutes AND 3 mL (2 mL min.) Whole blood / Bone Marrow in 1 Lavender Top (EDTA) tube AND 2 Peripheral blood smears AND Bone Marrow Trephine biopsy in 10% formal salin. Ship at specified temperatures. Give clinical history on Bone Marrow & ChromosomeTest Request Forms.	Flow Cytometry, Culture, Microscopy, Karyotype, Real Time PCR, HPE, Electrical Impedence		Sample Daily by 11 am; Report 7 days
------	---	---	---	---	--	--

Z537	LEUKEMIA DIAGNOSTIC COMPREHENSIVE PROFILE 2	*Microscopy of Bone Marrow Aspirate & Peripheral Blood *Special Stains *CBC *Leukemiadiagnostic panel, Acute Leukemia - T, B or Myeloid, *Chromosome analysis for hematological malignancy, *AML categorization by PCR / BCR ABL gene rearrangement, PCR QualitaÇ	4 mL (2.5 mL min.) Bone Marrow Aspirate in 2 Green Top (Sodium Heparin) tubes AND 4 Bone Marrow Aspirate Smears. Air dry & fix smears in 100% Methanol for 30 minutes AND 3 mL (2 mL min.) Whole blood / Bone Marrow in 1 Lavender Top (EDTA) tube AND 2 Peripheral blood smears AND Bone Marrow Trephine biopsy in 10% formal saline. Ship at specified temperatures. Give clinical history on Bone Marrow & ChromosomeTest Request Forms.	Flow Cytometry, Culture, Microscopy, Karyotype, Real Time PCR, Electrical Impedence		Sample Daily by 11 am; Report 7 days
Z664	LEUKEMIA DIAGNOSTIC COMPREHENSIVE PROFILE 5	*Chromosome analysis for hematologicalma lignancy *AML categorization by PCR	2 mL (1.5 mL min.) Bone Marrow Aspirate / Whole Blood in 1 Green Top (Sodium Heparin) tube AND 3 mL (2 mL min.) Whole blood / Bone Marrow Aspirate in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Test Request Form.	Culture, Microscopy, Karyotype, Real Time PCR		Sample Daily by 11 am; Report 7 days

Z663	LEUKEMIA DIAGNOSTIC COMPREHENSIV E PROFILE 7	*Leukemia diagnostic panel, Acute Leukemia - T, B or Myeloid *Chromosome analysis for hematological malignancy * Leukemia Genetic Profile - any 6 markers, PCR	4 mL (3 mL min.) Bone Marrow Aspirate / Whole Blood in 1 Green Top (Sodium Heparin) tube AND 3 mL (2 mL min.) Whole blood / Bone Marrow Aspirate in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Give clinical history on Chromosome Test Request Form.	Flow Cytometry, Culture, Microscopy, Karyotype, Real Time PCR		Sample Daily by 11 am; Report 7 days
Z268	LEUKEMIA DIAGNOSTIC PANEL:ACUTE LEUKEMIA - T, B or MYELOID		3 mL (2 mL min.) whole blood each in 1 Lavender Top (EDTA) tube AND 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone Marrow AND Aspirate Smear. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. Clinical History Mandatory	Flow Cytometry		Sample Daily by 9 am; Report Same Day
Z269	LEUKEMIA DIAGNOSTIC PANEL:AML CHARACTERIZAT ION		3 mL (2 mL min.) whole blood each in 1 Lavender Top (EDTA) tube AND 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone Marrow AND Aspirate Smear. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. Clinical History Mandatory.	Flow Cytometry		Sample Daily by 9 am; Report Same Day

Z154	LEUKEMIA MONITORING PANEL: HAIRY CELL		3 mL (2 mL min.) whole blood each in 1 Lavender Top (EDTA) tube AND 1 Green Top (Sodium Heparin) tube OR 2 mL (1 mL min.) heparinized Bone Marrow AND Aspirate Smear. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form. Clinical History Mandatory.	Flow Cytometry		Sample Daily by 9 am; Report Same Day
RW205	LEVETIRACETAM		3 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Recommended specimen: Just before next scheduled dose. Do not use SST gel barrier tubes.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
R081	LH; LUTEINISING HORMONE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age and sex on test request form.	Chemiluminescent Immunoassay	In females, preferred sampling time is Day 2 / Day 3 of the menstrual cycle.	Daily
Z207	LH & TESTOSTERONE, TOTAL		3 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age and sex on test request form.	Chemiluminescent Immunoassay	In females, preferred sampling time is Day 2 / Day 3 of the menstrual cycle.	Daily
B059	LIPASE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Spectrophotometry		Daily
RW122	LIPID - ASSOCIATED SIALIC ACID; LASA		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Colorimetry		Sample by 7th of the month; Report after 2–3 weeks

B131	LIPID PROFILE, BASIC	*Cholesterol *Triglycerides *HDL, Direct *LDL, Direct *VLDL, Calculated*Non HDL Cholesterol	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Minimum 12 hours fasting is mandatory. Grossly hemolysed samples are not acceptable.	Spectrophotometry,		Daily
Z060	LIPID PROFILE, EXTENDED	*Apolipoproteins A1, B & Ratio *Lipoprotein (a) *Cholesterol *Triglycerides *HDL/VLDL/LDL Electrophoretic *Cholesterol / HDL Ratio*Chylomicr ons *Non HDL Cholesterol * Cardio CRP	3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen. Minimum 12 hours fasting is mandatory. Grossly hemolysed samples are not acceptable.	Spectrophotometry, Electrophoresis, Imm unoturbidimetry		Sample Daily by 9 am; Report 3 days

E017	LIPID PROFILE, COMPLETE	*Cholesterol *Triglycerides *HDL/VLDL/LDL - Electrophoretic *Cholesterol / HDL Ratio*Chylomicrons*Non HDL Cholesterol	3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen. Minimum 12 hours fasting is mandatory. Grossly hemolysed samples are not acceptable.	Spectrophotometry, Electrophoresis,		Sample Daily by 9 am; Report 3 days
Z279	LIPID PROFILE, COMPREHENSIVE	*Apolipoproteins A1, B & Ratio *Lipoprotein (a) *Cholesterol *Triglycerides *LDL Cholesterol *HDL Cholesterol *Cholesterol / HDL Ratio *Non HDL Cholesterol *Homocysteine *Fibrinogen *Cardio CRP	4 mL (2 mL min.) serum from 2 SST's. Ship refrigerated. DO NOT FREEZE AND 2 mL (1.5 mL min.) Platelet Poor Plasma (PPP) from 1 Blue Top (Sodium Citrate) tube. Prepare PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Duly filled Coagulation Requisition Form is mandatory. Grossly hemolysed samples are not acceptable. Minimum 12 hours fasting is mandatory.	Spectrophotometry, CMIA , Imunoturbidimetry & Photo optical clot detection		Sample Mon through Sat by 4 pm; Report Same day
RW123	LIPOPROTEIN FRACTIONATION, ULTRACENTRIFUGATION		4 mL (2 mL min.) serum from 2 Red Top (No Additive) tubes. Ship refrigerated or frozen. Minimum 12hours fasting is mandatory.	Ultracentrifugation, Colorimetric,Enzymatic	Also see LIPID PROFILE TESTS.	Sample by 7th of the month; Reportafter 2–3 weeks

B101	LIPOPROTEIN (a); Lp(a)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry		Daily
M203	LISTERIA ANTIGEN, STOOL		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily 1 pm Report Same day
B069	LITHIUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Ion Selective Electrode		Daily
Z005	LIVER PANEL 1; LFT	*SGOT *SGPT *GGTP *Bilirubin *Protein, Total *Alkaline Phosphatase	2 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Spectrophotometry		Daily
Z618	LIVER FIBROSIS PANEL	*SGOT *SGPT *GGTP *Bilirubin *Protein, Total * Alkaline Phosphatase * ProthrombinTime * Protein Electrophoresis, Serum	3 mL (2 mL min.) serum from 1 SST AND 3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hours.. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Spectrophotometry, ElectromechanicalClot Detection, Capillary Electrophoresis	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample Daily by 9 am; Report Same Day

Z004	LIVER PANEL 2	*SGOT *SGPT *GGTP *Bilirubin *Protein, Total *Alkaline Phosphatase *Prothrombin Time *HBsAg	3 mL (2 mL min) serum from 1 SST AND 3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Spectrophotometry, CMIA, Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Daily
S165	LIVER KIDNEY MICROSOMAL (LKM) ANTIBODY IN DILUTIONS, IFA		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample Daily by 11 am; Report Nextday
S164	LIVER KIDNEY MICROSOMAL (LKM) ANTIBODY, IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay	For SERIAL DILUTIONS, use testcode S165.	Sample Daily by 11 am; Report Same Day
Z025	LIVER & KIDNEY PANEL	*LFT * KFT	3 mL (1.5 mL min.) serum from 1 SST. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship refrigerated or frozen.	Spectrophotometry, Indirect ISE		Daily

B192	Lp-PLA2 (LIPOPROTEIN ASSOCIATED PHOSPHOLIPAS E A2), ACTIVITY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzymatic		Sample by Tue / Fri 9 am ; Report Same day
RW124	LSD, URINE, QUANTITATIVE		20 mL (10 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen. Wrap container in aluminium foil to protect from light.	GC / MS		Sample by 7th of the month; Report after 2–3 weeks
Z789	LUNG CANCER MUTATION PANEL	*ALK-1 *EGFR	Submit formalin fixed paraffin embedded tissue block OR three 4 micron sections on Poly-L-Lysine coated slides. Ship at room temperature.	FISH, PCR		Sample by Tue 11 am; Report Mon
Z168	LUPUS ANTICOAGULAN T PANEL; LAC PANEL	*APTT *PTT - LA *dRVVT	3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Also see PHOSPHOLIPID SYNDÇ	Sample Daily by 3 pm; Report Next day

H123	LUPUS ANTICOAGULAN T BY dRVVT		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory	Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample Daily by 3 pm; Report Next day
H124	LUPUS ANTICOAGULAN T BY APTT; PTT - LA		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample Daily by 3 pm; Report Next day

S092	LYME DISEASE (BORRELIA BURGDORFERI) ANTIBODIES, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Sat 9 am; Report Same day
RW024	LYME DISEASE (BORRELIA BURGDORFERI) DNA, QUALITATIVE REAL TIME PCR		2 mL (1 mL min.) CSF / Synovial Fluid in a sterile screw capped container. Ship refrigerated or frozen.	Real Time PCR		Sample by 7th of the month; Report after 2–3 weeks
RW125	LYME DISEASE (BORRELIA BURGDORFERI), IgG & IgM, IMMUNOBLOT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoblot		Sample by 7th of the month; Report after 2–3 weeks
RW126	LYME DISEASE (BORRELIA BURGDORFERI), IgG, IMMUNOBLOT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoblot		Sample by 7th of the month; Report after 2–3 weeks
Z149	LYMPHOMA DIAGNOSTIC PANEL		Blood: 3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Bone Marrow: 3 mL (1.5 mL min.) in 1 Green Top (Sodium Heparin) tube. Body Fluids: 50 mL (10 mL min.) fluids in a sterile screw capped container. Lymph Node Aspirate: 1 mL (0.5 mL min.) in normal saline. Ship at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Clinical history is mandatory.	Flow Cytometry		Sample Daily by 9 am; Report Same Day

RW127	LYSOZYME; MURAMIDASE		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of draw in a plastic transport vial. Ship refrigerated or frozen.	Enzymatic		Sample by 7th of the month; Report after 2–3 weeks
Z125	MACROCYTIC ANEMIA PANEL	*Folate, serum * Vitamin B12	3 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Wrap tube in aluminium foil to protect from light. Ship refrigerated or frozen. Overnight fasting is mandatory.	Chemiluminescent Immunoassay		Daily
RW249	MACROPROLACTIN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	PEG, ECLIA		Sample by 7th of the month; Report after 2–3 weeks
U019	MAGNESIUM, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine in a sterile screw capped container. Do not use preservatives. Record 24-hour volume on test request form and urine container. Ship refrigerated or frozen.	Xylidyl blue		Daily
B066	MAGNESIUM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Xylidyl blue		Daily
U060	MAGNESIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen.	Xylidyl blue		Daily
H030	MALARIA PARASITE / BLOOD PARASITE IDENTIFICATION		Submit two thin AND two thick peripheral blood smears.	Microscopic Examination		Daily
S153	MALARIA VIVAX AND FALCIPARUM ANTIGEN		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Immunochromatography		Daily

S189	MALARIA ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunochromatography		Sample by Tue / Fri 9 am; Report Same day
OS160	MALARIA DETECTION , PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Tue / Sat 11 am; Report Sat/Wed
H249	MALARIA PARASITE CONFIRMATION		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Digital Cytometry		Daily
M184	MALDITOF - MS MYCOBACTERIUM SPECIATION		Submit pure growth of Mycobacterium species on a solid medium - LJ or Middlebrook agar. Ship refrigerated. DO NOT FREEZE	MALDITOF-MS (Matrix assisted laser desorption ionisation Time of flight- Mass spectrometry)		Sample by 9 am; Report Same day
M189	MALDITOF-MS BACTERIA IDENTIFICATION		Submit pure growth of bacteria on a solid agar plate. Ship refrigerated. DO NOT FREEZE.	MALDITOF-MS (Matrix assisted laser desorption ionisation Time of flight- Mass spectrometry)		Sample by 9 am; Report Same day
M190	MALDITOF-MS MOLD IDENTIFICATION		Submit pure growth of mold on a solid agar plate. Ship refrigerated. DO NOT FREEZE.	MALDITOF-MS (Matrix assisted laser desorption ionisation Time of flight- Mass spectrometry)		Sample by 9 am; Report Same day

M188	MALDITOF-MS YEAST IDENTIFICATION		Submit pure growth of yeast on a solid agar plate . Ship refrigerated. DO NOT FREEZE.	MALDITOF-MS (Matrix assisted laser desorption ionisation Time of flight- Mass spectrometry)		Sample by 9 am; Report Same day
C058	MANGANESE, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in the kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Patients should refrain from vitamins, mineral or herbal supplements one week prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C064	MANGANESE, BLOOD		4 mL (2 mL min.) whole blood from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Patients should refrain from vitamins, mineral or herbal supplements one week prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next Day

C059	MANGANESE, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Patients should refrain from vitamins, mineral or herbal supplements one week prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report next day
P007	MANTOUX TEST; TUBERCULIN SKIN TEST		Patient to report to NRL / Satellite Labs on Day 1 and Day 3 of the test.	Intradermal skin test		48 hrs post-injection
N075	MATERNAL BLOOD FOR FETAL DNA	*Trisomy detection for 13, 18 & 21 Chromosomes	10 mL Whole blood in a special tube available from LPL. Ship refrigerated. DO NOT FREEZE. Valid between 10-24 weeks of gestation. Give clinical history on Maternal Serum Screen request form available from LPL.	DNA Sequencing		Sample by Tue 9 am; Report next Tue
Z266	MATERNAL SERUM SCREEN 3; TRIPLE TEST	*AFP *BETA HCG *Free Estriol *Risk Evaluation	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide maternal Date of birth (dd/mm/yy); LMP or Ultrasound; IVF, Number of Fetuses (Single/ Twins); Diabetic status and Body Weight in Kg, Smoking & Previous history of Trisomy 21 pregnancy in Maternal Serum Screen Form. Valid between 14-22 weeks gestation (Ideal 15-20 weeks).	Chemiluminescent Immunoassay		Sample Mon through Sat by 9 am; Report Same day

Z514	MATERNAL SERUM SCREEN 4; QUADRUPL TEST	*AFP *BETA HCG *Free Estriol *Inhibin A *Risk Evaluation.	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide maternal Date of birth (dd/mm/yy); LMP or Ultrasound; IVF, Number of Fetuses (Single/ Twins); Diabetic status and Body Weight in Kg, Smoking & Previous history of Trisomy 21 pregnancy in Maternal Serum Screen Form. Valid between 14-22 weeks gestation (Ideal 15-20 weeks).	Chemiluminescent Immunoassay		Sample Mon through Sat by 9 am; Report Same day
Z397	MATERNAL SERUM SCREEN 2; DUAL TEST	*PAPP-A *Free BETA HCG *Risk Evaluation.	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide maternal Date of birth (dd/mm/yy); LMP or Ultrasound; Number of Fetuses (Single/Twins); Diabetic status and Body Weight in Kg, IVF, Smoking & Previous history of Trisomy 21 pregnancy in Maternal Serum Screen Form. Valid between 9-13 weeks gestation (Ideal 10- 13 weeks).	Chemiluminescent Immunoassay	For Combined Risk Assessment (Biochemical risk + Nuchal translucency) kindly enclose Ultrasound report between 11-13 weeks gestation including CRL, NT & Nasal Bone.	Sample Mon through Sat by 9 am; Report Same day
H005	MEAN CORPUSCULAR VOLUME; MCV		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Calculated		Daily
H007	MEAN CORPUSCULAR HEMOGLOBIN CONCENTRATIO N; MCHC		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Calculated		Daily
H006	MEAN CORPUSCULAR HEMOGLOBIN; MCH		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electronic Cell Cytometry, Calculated		Daily

H033	MEAN PLATELET VOLUME; MPV		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence		Daily
Z191	MEASLES (RUBEOLA) ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S072	MEASLES (RUBEOLA) ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S073	MEASLES (RUBEOLA) ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
OS161	MEDULLARY THYROID CARCINOMA, MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am; Report 10 days
OS162	MELAS (MITOCHONDRI AL ENCEPHALOPAT HY & LACTIC ACIDOSIS WITH STROKE LIKE SYNDROME), MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am; Report 10 days

H244	MENNORHAGIA SCREEN	*Prothrombin Time *APTT *Factor VIII Functional *Factor IX	6 mL Whole blood in 2 Blue Top (Sodium Citrate) tubes. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled coagulation Requisition form is mandatory.	Electromechanical & Photo-optical ClotDetection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Thu 11 am; Report Next day.
Z808	MENOPAUSAL SCREENING PANEL	*Estradiol (E2) *TSH, Ultrasensitive * FSH *Anti Mullerian hormone (AMH)	5 mL (3 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA, EIA		Sample Daily by 9 am; Report Same day

C028	MERCURY, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24 hr urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in the kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Patient should refrain from eating predatory fish like Sword Fish, Tuna and Shark at least 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C029	MERCURY, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Patient should refrain from eating predatory fish like Sword fish, Tuna and shark atleast 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified container available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report next day
OS163	MERRF (MYOCLONIC EPILEPSY ASSOCIATED WITH RAGGED RED FIBRES), MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am; Report 10 days

G108	METABOLIC SCREEN, COMPREHENSIV E, URINE	*Reducing Substances *Ferric Chloride Test *DNPH Test *Nitrosonaphthol Test *Nitroprusside Test *Silver nitroprusside Test *Amino Acids, Qualitative	10 ml (5 ml min) aliquot of first morning urine in a sterile screw capped container. No preservative required. Ship refrigerated or frozen. For infants on breast milk, feed the child and then collect urine. Provide a brief clinical and drug history over the last 3 days prior to specimen collection.	Chemical, Thin Layer Chromatography		Sample by Mon / Wed / Fri 4 pm; Report 2 days
G101	METABOLIC SCREEN, BASIC, URINE	*Reducing Substances *Ferric Chloride Test *DNPH Test *Nitrosonaphthol Test *Nitroprusside Test * Silver nitroprusside Test	10 ml (5 ml min) aliquot of first morning urine in a sterile screw capped container. Do not use any preservatives. Ship refrigerated or frozen. For infants on breast milk, feed the child and then collect urine. Provide a brief clinical and drug history over the last 3 days prior to specimen collection.	Chemical		Sample Daily by 4 pm; Report Same day
G122	METACHROMA TIC LEUCODYSTROP HY, QUANTITATIVE, BLOOD		10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) tubes OR 10 mL (7.5 mL min.) whole blood from 3 Green Top (Sodium Heparin) tubes. Ship refrigerated. DO NOT FREEZE. Provide brief clinical history.	Spectrophotometry	This tests deficiency of enzyme Aryl Sulphatase A.	Sample Daily by 4 pm; Report 4 days

Z820	METAL SCREEN 2, WHOLE BLOOD	*Arsenic * Cadmium *Chromium *Cobalt *Lead *Manganese *Thalium	4 mL (2 mL min.) whole blood from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Patient should avoid eating seafood, herbal supplements, mineral supplements, vitamin B12 or vitamin B complex for 3 days prior to specimen collection. Ship refrigerated. DO NOT FREEZE. Use powderless gloves during specimen collection.	ICPMS, AAS with Zeeman correction	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
Z821	METAL SCREEN 3, URINE	*Arsenic * Cadmium *Chromium *Cobalt *Lead *Manganese *Thalium * Copper * Zinc *Selenium *Nickel *Mercury	10 mL (5 mL min.) aliquot of 24 hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in kit. Measure 24 hour total volume and record on vial and Test request form. Patient should avoid eating seafood, herbal supplements, mineral supplements, vitamin B12 or vitamin B complex, sea/fresh water fish for atleast 3 days prior to specimen collection. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	AAS withZeeman correction, ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

Z822	METAL SCREEN 1, SERUM	*Aluminum * Copper * Zinc *Selenium	3 mL (1 mL min.) serum from 1 metal free White Top (Z-No Additive) tube available from LPL. Patients should refrain from taking antacids, Aluminium containing medications, sea/fresh water fish at least 3 days prior to specimen collection. Ship refrigerated. DO NOT FREEZE. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
U035	METANEPHRINE S, 24-HOUR URINE	*Metanephrine *Normetanephri ne	50 mL (10 mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl to maintain pH between 1–2. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen. Patient should strictly avoid Alpha- Methyldopa, Buspirone, Codeine, Isoproterenol metabolite, Mandelamine, L-dopa, Paracetamol, Metoclopramide, pepper, alcoholic beverages, tea/coffee, tobacco and strenuous exercise at least 72 hours before and during specimen collection.	High Performance LiquidChromatograp hy	Also see PHEOCHROMOCYTO MA PROFILE	Sample by Fri 9 am; Report Tue

U052	METANEPHRINE S, RANDOM URINE	*Metanephrine *Normetanephri ne	50 mL (10 mL min.) aliquot of random urine. Add 50% HCl drop by drop to bring the pH between 1–2. Do not use concentrated HCl. Ship refrigerated or frozen. Patient should strictly avoid Alpha-Methyldopa, Buspirone, Codeine, Isoproterenol metabolite, Mandelamine, L-dopa, Paracetamol, Metoclopramide, pepper, alcoholic beverages, tea/coffee, tobacco and strenuous exercise at least 72 hours before and during specimen collection.	High Performance LiquidChromatograp hy	Also see PHEOCHROMOCYTO MA PROFILE	Sample by Fri 9 am; Report on Tue
K019	METHAMPHETA MINE SCREEN, URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen	Lateral Flow Immunoassay		Sample Daily by 4 pm; Report Same day
K021	METHAQUALON E SCREEN, URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Sample Daily by 4 pm; Report Same day
K022	METHADONE SCREEN, URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Sample Daily by 4 pm; Report Same day
D052	METHOTREXAT E		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use gel barrier tubes. Wrap in Aluminium foil to protect from light. Ship refrigerated or frozen.	EMIT	Blood concentrations of Methotrexate are monitored at 24, 48 & 72 hours after a single high dose (>50 mg/m ²) therapy.	Sample by Mon / Wed / Fri 2 pm; ReportSame day

OS164	METHYLGUANINE METHYLTRANSFERASE (MGMT)		Submit fresh tissue in normal saline in a sterile screw capped container. Ship refrigerated . DO NOT FREEZE	PCR		Sample Daily by 11 am; Report 22 days
G189	METHYLMALONIC ACID, QUANTITATIVE, SERUM		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report Next day
G103	METHYLMALONIC ACID, QUALITATIVE, URINE		10 ml (5 ml min.) aliquot of first morning urine in a sterile screw capped container. Do not use any preservatives. Ship refrigerated or frozen. Provide a brief clinical and drug history.	Chemical		Sample Daily by 4 pm; Report Same day
G188	METHYLMALONIC ACID, QUANTITATIVE, URINE		10 ml (5 ml min.) aliquot of first morning urine in a sterile screw capped container. Do not use any preservatives. Ship refrigerated or frozen. Provide a brief clinical and drug history.	LC-MS/MS		Sample Daily by 4 pm; Report Next day
OS165	MICM (MATERNALLY INHERITED CARDIOMYOPATHY), MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am ; Report 10 days
B142	MICROALBUMIN CREATININE RATIO, URINE		15 mL (10 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Immunoturbidimetry , Spectrophotometry		Daily

U027	MICROALBUMIN, 1ST MORNING / RANDOM URINE		10 ml (5 ml min.) aliquot of first morning / random urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Immunoturbidimetry	The recommended test is U031MICROALBUMIN, 24-HOUR URINE.	Daily
U031	MICROALBUMIN, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine in a sterile screw capped container. Do not use preservatives. Record 24 hour volume on test request form and urine container. Ship refrigerated. DO NOT FREEZE.	Immunoturbidimetry		Daily
H024	MICROFILARIA DETECTION		2 air dried peripheral blood smears (1 thick and 1 thin) OR 2 mL (1 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 2 mL (1 mL min.) Hydrocele fluid. Draw specimen between 12 midnight to 2 a.m. Ship refrigerated. DO NOT FREEZE.	Microscopy	Also see FILARIA ANTIGEN	Daily
Z831	MINIMAL RESIDUAL DISEASE (MRD) ANALYSIS FOR B- ALL	*CD10 *CD13 *CD15 *CD19 *CD20 *CD33 *CD34 *CD38 *CD45 *CD58 *CD73 *CD86 *CD123 *CD200 *SYTO-16	3 mL (1 mL min.) First pull Bone Marrow aspirate in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18-22C. DO NOT FREEZE. Following requirements are mandatory:• Date & Time of sampling• Previous immunophenotype report• Clinical history	Flow Cytometry	Following requirements are mandatory:• Time of sampling• Previous immunophenotype report• Clinical history• Date & time of sampling	Sample Daily by 9 am; Report Next day
S046	MITOCHONDRIAL ANTIBODY (AMA), IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	For SERIAL DILUTIONS, use testcode S085.	Sample Daily by 11 am; Report Same Day

S085	MITOCHONDRIAL ANTIBODY (AMA), IFA IN DILUTIONS		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 11 am; Report Nextday
OS166	MITOCHONDRIAL MUTATION DETECTION PANEL 1	*MALAS *MERRF *MM *MMC *ME *PEM *EI	4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am ; Report 10 days
OS167	MITOCHONDRIAL MUTATION DETECTION PANEL 2	*LHON *NARP *CPEO *SNHL *DEAF *AMDF	4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am ; Report 10 days
OS168	MITOCHONDRIAL MUTATION DETECTION COMPREHENSIVE PANEL		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am ; Report 10 days

H242	MIXING STUDIES - APTT WITH PROTHROMBIN TIME STUDIES	*Prothrombin Time *APTT *Thrombin Time *Fibrinogen * Mixing studies	6 mL Whole blood in 2 Blue Top (Sodium Citrate) tubes. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Wed / Sat 11 am; Report Thu / Mon
H184	MIXING STUDIES - APTT	*Fibrinogen *Thrombin Time *APTT *Mixing studies	6 mL Whole blood in 2 Blue Top (Sodium Citrate) tubes. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Wed / Sat 11 am; Report Next day

OS169	MMC (MATERNAL MYOPATHY WITH CARDIOMYOPA THY), MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am ; Report 10 days
OS170	MM (MITOCHONDRI AL MYOPATHY), MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am ; Report 10 days
M187	MODIFIED ACID FAST STAIN FOR LEPROSY		Submit Slit -skin smear, Nasal Smear, Ear lobe Tissue smear and other secretions as applicable. Ship at 18–22°C.	Light Microscopy		Sample Daily by 9 am; Report Same day
S004	MONO TEST (INFECTIOUS MONONUCLEOS IS)		2 mL (1 mL min.) serum from 1 SST. Ship frozen. DO NOT THAW	Latex Agglutination	Also see EPSTEIN-BARR VIRUSTESTS.	Daily
N121	MPL (MYELOPROLIFE RATIVE LEUKEMIA), GENE MUTATION		3 ml (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR, Fragment Analysis		Sample by Mon / Thu 11 am; Report Wed / Sat

N049	MTHFR (METHYLENE TETRAHYDROFO LATE REDUCTASE) GENE MUTATION, QUALITATIVE PCR		3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
G112	MUCOPOLYSAC CHARIDOSIS (MPS) TYPE 1 (HURLER), QUANTITATIVE, BLOOD		10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) tubes OR Green Top (Sodium Heparin) tubes. Ship refrigerated. DO NOT FREEZE. Provide brief clinical history.	Enzyme assay	This tests deficiency of enzyme Alpha Iduronidase.	Sample Daily by 4 pm; Report 4 days
G113	MUCOPOLYSAC CHARIDOSIS (MPS) TYPE VI (MAROTEAUX LAMY), QUANTITATIVE, BLOOD		10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) OR Green Top (Sodium Heparin) tubes. Ship refrigerated. DO NOT FREEZE. Provide brief clinical history.	Enzyme assay	This tests deficiency of enzyme Aryl-Sulphatase B.	Sample Daily by 4 pm; Report 4 days
G110	MUCOPOLYSAC CHARIDOSIS (MPS) SCREEN, URINE		10 ml (5 ml min.) aliquot of first morning urine in a sterile screw capped container. Do not use any preservatives. Ship refrigerated or frozen. Provide brief clinical and drug history.	Toluidine blue spot test		Sample Daily by 4 pm; Report Next day
G111	MUCOPOLYSAC CHARIDOSIS (MPS) TYPING, 24 HOUR URINE		20 mL (10 mL min.) aliquot of 24 hour urine in a sterile screw capped container. Do not use any preservatives. Ship refrigerated or frozen. Provide a brief clinical and drug history.	Two dimensional electrophoresis		Sample Daily by 4 pm; Report 15 days

Z182	MULTIPLE SCLEROSIS PANEL 2	*IgG Synthesis Index & Rate *Oligoclonal Bands, CSF *Myelin Basic Protein, CSF (Code RW135)	5 ml (3 ml min.) Cerebrospinal fluid AND 2 mL (1 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE	BCG, Nephelometry, Immunoturbidimetry, Agarose Gel Electrophoresis, RIA	Strictly adhere to specimen volumes. Also see MYELIN ASSOCIATED GLYCOPROTEIN; MYELIN ASSOCIATED GLYCOPROTEIN - SULFATED - GLUCURONIC PARAGLOBOSIDE AUTOANTIBODY, IgM	LPL: Sample Daily by 9 am; Report Next day USA: Sample by 7th of the Month; Report after 2–3 weeks.
Z181	MULTIPLE SCLEROSIS PANEL 1	*IgG Synthesis Index & Rate *Oligoclonal Bands, CSF	4 ml (3 ml min.) Cerebrospinal fluid AND 2 mL (1 mL min.) serum from 1 SST. Ship refrigerated. DO NOT FREEZE. The collection date and time must be the same for both specimens.	BCG, Nephelometry, Immunoturbidimetry, Agarose Gel Electrophoresis	Strictly adhere to specimen volumes for CSF	Sample Daily by 9 am; Report Next day
Z782	MULTIPLE MYELOMA SCREENING PANEL	*Protein electrophoresis, serum *Immunofixation electrophoresis, serum *Kappa / Lambda Light Chains, Free, serum *Immunoglobulin Profile, serum	5 mL (3 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Capillary electrophoresis, Agarose Gel Electrophoresis, Immunoturbidimetry		Sample Daily by 9 am; Report Next day

Z382	MUMPS VIRUS ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
S037	MUMPS VIRUS ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
S075	MUMPS VIRUS ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon 9 am; Report Same day
S243	MuSK (MUSCLE SPECIFIC KINASE) ANTIBODY		4 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radioimmunoassay		Sample by Fri 9 am; report Next day
Z174	MYASTHENIA GRAVIS PANEL	*Acetylcholine Receptor Binding Antibody *Striated / Skeletal Muscle Antibody in Dilutions	5 mL (3 mL min.) serum from 2 SST's. Ship refrigerated or frozen. Patient should not have general anaesthetic or muscle relaxant drugs, 24 prior to sampling.	Radioimmunoassay, IFA		Sample by Wed / Sat 9 am; Report Same day

N001	MYCOBACTERIUM TUBERCULOSIS, PCR; MYCOSURE; TB PCR	This test includes 4 Drug Sensitivity for MTB Complex if positive. Antitubercular Drug Sensitivity includes Rifampicin, Isoniazid, Ethambutol & Pyrazinamide	Submit first morning Sputum on 3 consecutive days OR 5–10 mL Body fluids / Aspirates OR 2 mL CSF, OR 100 ml (50 ml min) first morning Urine collected on 3 consecutive days OR Fresh tissue biopsies including endometrial aspirate / curettage in sterile normal saline OR formalin fixed paraffin embedded Tissue block OR 3 mL (1 mL min.) Whole blood/Bone Marrow in 1 Lavender Top (EDTA)/ Yellow Top (ACD) tube OR Menstrual blood in a sterile screw capped container OR culture isolates. All specimens to be collected aseptically and shipped refrigerated. DO NOT FREEZE.	Real Time PCR	Antitubercular Drug Sensitivity includes Rifampicin, Isoniazid, Ethambutol & Pyrazinamide	Sample Daily by 11 am ; Report Next day if negative & 4 days if Positive
N094	MYCOBACTERIUM TUBERCULOSIS, rRNA DETECTION, TMA		Submit first morning Sputum (3 consecutive days preferred) OR 5-10 ml Body fluid/Aspirate OR 2 ml CSF OR Fresh tissue biopsies including endometrial aspirate/curettage in sterile normal saline OR 10 ml (5 mL min.) Whole Blood in 3 lavender top (EDTA) tubes OR Menstrual blood in a sterile screw capped container OR culture isolates. Ship refrigerated or frozen.	Transcription Mediated Amplification (TMA) and Hybridization Protection Assay (HPA)		Sample by Mon /Wed/ Fri 11 am, Report Tue/Thu / Sat

Z277	MYCOBACTERIU M COMBINED PANEL;TB. COM	*Mycobacterium Tuberculosis, PCR *Culture, AFB Rapid *Ziehl Neelson Stain	First morning Sputum on 3 consecutive days OR 5–10 mL Body fluids/ aspirate OR 2 mL CSF OR 100 mL first morning Urine collected on 3 consecutive days OR Fresh tissue biopsies including endometrial aspirate / curettage in sterile normal saline OR 1-5 mL Whole blood / Bone Marrow in 1 Lavender Top (EDTA) / Yellow Top (ACD) tube OR 1- 2 mL Menstrual blood in a sterile screw capped container AND 3 mL (1 mL min.) blood in Myco F / Lytic bottle available from LPL.All specimens to be collected aseptically and shipped refrigerated.	Real Time PCR, Automated fluorescent, ICT	Specimens to be collected in 2 separate vials for culture & PCR. Sputum and urine samples should ideally be tested on 3 consecutive days.	Sample Daily by 3.30 pmReport Interim: 14 days & 4 weeksFinal: 6 weeks
Z677	MYCOPLASMA PNEUMONIAE IgG & IgM, ANTIBODIES PANEL, SERUM		2 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample by Mon 9 am; Report Same day
S223	MYCOPLASMA PNEUMONIAE IgG, SERUM		2 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample by Mon 9 am; Report Same day
S224	MYCOPLASMA PNEUMONIAE IgM, SERUM		2 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample by Mon 9 am; Report Same day
G187	MYCOPHENOLIC ACID		3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Separate plasma within 2 hours. Ship refrigerated. Ideal sampling time 1 hour before next dose (Trough level) at steady state. Provide date and time of sampling, amount and time of dose.	LC-MS/MS		Sample Daily by 4 pm; Report Same day

RW135	MYELIN BASIC PROTEIN; MBP		2 mL (1 mL min.) Cerebrospinal fluid in a sterile screw capped container. Ship refrigerated or frozen. Blood contamination and hemolysis may interfere with the results. Ideally submit 4th lumbar puncture collection tube to minimize contamination.	Immunoassay	Also see MULTIPLE SCLEROSIS PANELS; MYELIN ASSOCIATED GLYCOPROTEIN and MAG-SGPG AUTOANTIBODY.	Sample by 7th of the month; Report after 2–3 weeks
RW134	MYELIN ASSOCIATED GLYCOPROTEIN-SULFATED GLUCURONIC PARAGLOBOSIDE AUTOANTIBODY (MAG-SGPG), IgM		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunoassay, Endpoint Titer		Sample by 7th of the month; Report after 2–3 weeks
RW133	MYELIN ASSOCIATED GLYCOPROTEIN (MAG), IgM		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
H225	MYELOPEROXIDASE (MPO) STAIN, BLOOD / BONE MARROW		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 2 mL (1 mL min.) Bone Marrow from 1 Green Top (Sodium heparin) tube. Ship refrigerated. DO NOT FREEZE. Clinical history mandatory.	Cytochemical staining		Sample Daily by 6 pm; Report 3 days
B102	MYOGLOBIN, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Electrochemiluminescence		Daily

S245	MYOSITIS PROFILE	*Mi-2 *Ku* PM- SCL 100 * PM- SCL 75* Jo-1 * SRP* PL-7 * PL- 12* EJ* OJ* Ro- 52	2 mL (1 mL min.) serum from 1 SST.Ship refrigerated or frozen.	Immunoblot		Sample by Thu 9 am; Report Same day
Z836	NAFLD (NON ALCOHOLIC FATTY LIVER DISEASE) FIBROSIS SCORE	*AST *ALT *Albumin* Glucose F & PP *Platelet Count	2 mL (1 mL min.) serum from 1 SST AND 3 mL (1mL min.) whole blood in 1 lavender top (EDTA) tube AND 2 mL (1 mL min.) plasma each for Fasting & Post prandial from 1 Grey Top (Sodium Fluoride) tube. Overnight Fasting is mandatory. Collect fasting Glucose sample. Ask patient to have normal to heavy meal. Note the time when patient starts meal, collect PP sample exactly after 2 hrs. Ship refrigerated. DO NOT FREEZE. Patient's age, height & weight are mandatory.	IFCC, Hexokinase, BCG, Electrical Impedence		Sample Mon through Sat 9 am;Report Same day
Z166	NATURAL KILLER CELLS	*CD3 *CD16 / CD56	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube OR 2 mL(1 mL min.) heparinized Bone Marrow. Ship immediately at 18-22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry	Also see IMMUNE DEFICIENCY PANELS	Sample Daily by 9 am; Report Same day
N088	NEISSERIA GONORRHEA, PCR QUALITATIVE		15 ml (5 mL min.) First morning Urine in a sterile screw capped container OR Vagina / Endocervical Swab. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample Daily by 11 am; Report Next day

F019	NEURON SPECIFIC ENOLASE (NSE), SERUM		2 mL (1 mL min.) serum from 1 SST. Hemolysed specimens not acceptable. Ship refrigerated or frozen.	Electrochemiluminescence	NSE is high in platelets and RBC, therefore plasma and hemolysed specimens are unacceptable.	Sample by Mon 4 pm; Report Same day
Z851	NEUROVIRUSES PANEL, QUALITATIVE PCR	*Adenovirus*Cytomegalovirus*Enteroviruses*Epstein Barr virus*Herpes simplex virus 1 & 2*Varicella Zoster virus*Parvovirus*HHV6*HHV7*Parvovirus B19	3 ml (2 mL min.) Whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
Z184	NEUROBLASTOMA PROFILE, 24-HOUR URINE	*Vanil Mandelic Acid (VMA) *Homovanillic Acid (HVA)	50 mL (10 mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl to maintain pH between 1–2. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen. Patient should strictly avoid L-Dopa, Theophylline, chocolate, vanilla, banana, alcoholic beverages, tea/coffee, tobacco and strenuous exercise at least 72 hrs before and during specimen collection.	High Performance Liquid Chromatography	Also see NEURON SPECIFIC ENOLASE.	Sample by Wed / Sat 9 am; Report Fri / Tue

S232	NEURONAL (PARANEOPLAS TIC) AUTOANTIBODI ES PROFILE	*Amphiphysin *CV2 *PNMA2 (Ma2 / Ta) *ANNA-1 / Hu *ANNA-2 / Ri *PCA-1 / Yo	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoblot		Sample by Mon through Sat 9 am; Report Same day
Z335	NEWBORN SCREENING	*G6PD *TSH *Phenylalanine *Cystic Fibrosis *17- OHP *Galactosemia *Biotinidase	1 drop of heel prick blood each on 3 spots of filter paper available from LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen.	Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm; Report 2 days
G146	NEWBORN SCREENING, EXTENDED	*Biotinidase *Cystic Fibrosis *G6PD *Galactosemia *17-Hydroxy Progesterone *TSH *Tandem Mass Spectrometry - 30 Disorders	1 drop of heel prick blood each on 3 spots of filter paper available from LPL.Clinical details and drug history must accompany sample. Ship refrigerated.	TMS, Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm; Report 2 days
R167	NGAL (NEUTROPHIL GELATINASE - ASSOCIATED LIPOCALIN), URINE		10 mL(5 mL min) aliquot of random urine in a sterile screw capped container. Ship refrigerated or frozen.	CMIA		Sample by Mon / Wed / Fri 4 pm; Report Same Day

C051	NICKEL, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in the kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Patients should refrain from mineral supplements 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C056	NICKEL, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Patients should refrain from mineral supplements 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report next day
RW138	NICOTINE AND COTININE, URINE		10 mL (5 mL min.) random urine in a sterile screw capped container. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks

G126	NIEMANN PICK DISEASE, QUANTITATIVE, BLOOD		10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium Heparin) tubes. Ship refrigerated. DO NOT FREEZE. Provide brief clinical history	Enzyme assay		Sample Daily by 4 pm; Report 4 days
G143	NITROPRUSSIDE TEST, URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Chemical	Also see AMINO ACID PANELS, QUALITATIVE AND QUANTITATIVE.	Sample Daily by 4 pm; Report Same day
G142	NITROSONAPHTHOL, URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Chemical	Also see AMINO ACID PANELS, QUALITATIVE AND QUANTITATIVE.	Sample Daily by 4 pm; Report Same day
M185	NOCARDIA, STAIN	*KOH Calcifluor * Gram stain *Modified ZN stain	Submit Sputum/ BAL /PUS/ Swab /Endotracheal tube secretions / Swabs / Culture Growth of Nocardia. Ship at 18-22°C	Light microscopy		Sample Daily by 9 am; Report Same day
RW140	NORTRIPTYLINE		3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended Specimen: >12 hours post oral dose.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
N051	NPM1 GENE MUTATION		5 mL (3 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N108	NRAS MUTATION CODON 12, 13 & 61		Submit formalin fixed paraffin embedded tissue block. Ship at room temperature. Block should have at least 10% tumor tissue.	PCR, Sequencing		Sample by Tue 11 am; Report Mon

B185	NT-PRO BNP (N-TERMINAL PRO B TYPE NATRIURETIC PEPTIDE)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	ECLIA		Daily
B161	NTx (COLLAGEN CROSS-LINKED N-TELOPEPTIDE), URINE		10 mL (5 mL min.) aliquot of 24-hour / Second void urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen. Record 24 hour volume on vial and Test Request Form.	Enzyme Immunoassay	Also see ALKALINE PHOSPHATASE, BONE SPECIFIC; OSTEOPOROSIS PANEL; PYRIDINIUM COLLAGEN CROSS LINKS.	Sample by Sat 9 am; Report Same day
RW141	NUCLEAR MATRIX PROTEINS; NMP		Submit 10 mL random urine collected between midnight and noon as a single void. Stabilize immediately after collection using the special NMP urine collection kit available from LPL. Follow instructions provided with the kit carefully. Stabilized urine will turn blue / green in color. Keep away from sunlight. Ship refrigerated or frozen.	Immunoassay	Test should not be performed on persons with Total cystectomy. Urine samples should be collected at least 5 days after invasive procedures like Cystoscopy or Catheterization.	Sample by 7th of the month; Report after 2–3 weeks
Z321	OBESEITY PANEL	* Sugar (F) *TSH *Lipid profile Basic *Leptin	3 mL (1.5 mL min.) serum from 1 SST AND 2 mL (1 mL min.) Plasma from 1 Grey Top (Sodium Fluoride) tube. Hemolysed specimens not acceptable. Ship refrigerated or frozen. Minimum 12 hours fasting is mandatory	Spectrophotometry, Chemiluminescent Immunoassay, EIA		Sample by Sat 9 am; Report Same day

Z292	OBSTETRIC FLOW PANEL	*Natural Killer (NK) cells *CD 19	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am; Report Same day
Z818	ONCOGENIC E6/E7 mRNA		Cervical swab / cells in special container with medium available from LPL. Ship refrigerated. DO NOT FREEZE. Store specimen at 2-8 degrees celsius immediately after collection. Exposure of the specimen to high temperatures has adverse effects on sample quality.	Flow Cytometry		Sample by Tue /Thu / Sat 9 am; Report Wed / Fri / Mon
K001	OPIATES SCREEN, RANDOM URINE		15 mL (10 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily
G164	ORGANIC ACIDS, QUANTITATIVE, RANDOM URINE, FULL PANEL		15 mL (10 mL min.) aliquot of random urine in a sterile screw capped container without preservative. Ship refrigerated or frozen. Clinical details and drug history must accompany the sample.	Gas Chromatography / Mass Spectrometry		Sample Daily by 5 pm; Report 5 days
G150	OROTIC ACID		15 mL (10 mL min.) aliquot of random urine. Ship refrigerated or frozen. Collect specimen fasting or 2-3 hours post feed	GC-MS		Sample Daily by 4 pm; Report 2 days
B184	OSMOLAL GAP, SERUM	*Sodium *BUN *Glucose R *Osmolality	3 mL (1.5 mL min.) serum from 1 SST AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Separate serum / plasma within 1 hour of collection. Ship refrigerated or frozen. Grossly hemolysed samples not accepted.	Freezing Point Depression, Spectrophotometry, Calculated		Sample Daily by 4 pm; Report Same day

P064	OSMOLALITY, SERUM		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Grossly hemolysed samples not accepted.	Freezing Point Depression	Also see ARGININE VASOPRESSIN.	Sample Daily by 4 pm; Report Same day
U040	OSMOLALITY, URINE		5 mL (3 mL min.) aliquot of random urine. Do not use preservative. Ship refrigerated. DO NOT FREEZE.	Freezing point depression	Also see ARGININE VASOPRESSIN.	Sample Daily by 4 pm; Report Same day
U100	OSMOLALITY, URINE : SERUM RATIO		2 mL (1.5 mL min.) serum from 1 SST AND 5 mL (3 mL min.) aliquot of random urine. Do not use preservative. Separate serum within 1 hour of collection. Ship refrigerated. DO NOT FREEZE.	Freezing point depression		Sample Daily by 4 pm; Report Same day
H012	OSMOTIC FRAGILITY TEST; RBC FRAGILITY TEST		3 mL (2 mL min.) whole blood from 1 Green Top (Sodium Heparin) tube. Ship refrigerated. DO NOT FREEZE.	Parparts	Oxalate and citrate samples not accepted.	Sample Daily by 12 noon; Report Next day
Z159	OSTEOPOROSIS PANEL	*Calcium *NTX *Estradiol *Osteocalcin *PTH-Intact *Vitamin D, 25-Hydroxy *Phosphorus	5 mL (3 mL min.) serum from 2 SST's. Collect sample between 8-10 a.m, separate serum and freeze within 30 minutes AND 20 mL (10 mL min) aliquot of first morning urine without preservatives. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship frozen. DO NOT THAW. Hemolysed samples are not acceptable.	Arsenazo III, CLIA, ECLIA, Molybdate UV, EIA		Sample by Sat 9 am; Report in two days
Z612	OSTEOSCREEN PANEL	*Calcium *Phosphorus *PTH-Intact *Vitamin D, 25-Hydroxy	4 mL (2 mL min.) serum from 1 SST. Collect sample between 8-10 a.m, in a fasting condition, separate serum and freeze within 30 minutes. Avoid use of tourniquet and fist clenching during sampling. Ship frozen. DO NOT THAW.	CLIA, Spectrophotometry		Daily

R108	OSTEOCALCIN; BONE Gla PROTEIN (BGP)		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour and freeze. Ship frozen. DO NOT THAW. Overnight fasting is preferred.	Electrochemiluminescent Immunoassay		Sample by Mon 4 pm; Report Same day
Z606	OVARIAN RESERVE PANEL	*FSH *Estradiol *Anti Mullerian Hormone	3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age, sex and day of menstrual cycle on test request form.	EIA, CLIA		Sample Daily by 9 am; Report Same day
Z803	OVARIAN CANCER PANEL	*CA 125 *Inhibin B	4 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CMIA, EIA		Sample by Mon 9 am; Report Same day
U055	OXALATE, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine. Collect urine with 25 mL 50% concentrated HCl to maintain pH below 3. For children adjust the volume of acid according to the expected urinary output. Keep urine refrigerated during collection. Mix well before aliquoting. Record body weight & 24 hour volume on container and test request form. Ship refrigerated or frozen. Avoid vitamin C and citrus fruits 48 hours prior to specimen collection.	Spectrophotometry		Sample by Wed / Sat 9 am; Report Same day
U101	OXALATE, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine in 2 mL of 6N HCl to maintain pH below 3. Mix well before aliquoting. Ship refrigerated or frozen. Avoid vitamin C and citrus fruits 48 hrs prior to specimen collection.	Spectrophotometry		Sample by Wed / Sat 9 am; Report Same day
P093	PACKED CELL VOLUME (PCV), BODY FLUIDS		Submit 3 mL (1.5 mL min.) body fluids in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE	Calculated		Daily

H002	PACKED CELL VOLUME (PCV); HEMATOCRIT		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	Calculated		Daily
RW144	PANCREATIC POLYPEPTIDE		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Ship refrigerated or frozen. Overnight fasting is preferred.	Extraction, Radioimmunoassay		Sample by 7th of the month; Report after 2–3 weeks
OS175	PAN FUNGAL DETECTION AND IDENTIFICATION, PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube OR 2 mL CSF/Sputum/ Fluids OR Tissue in Normal saline. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Thu 11 am; Report 8 days
OS174	PAN FUNGAL DETECTION , PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube OR 2 mL CSF/Sputum/ Fluids OR Tissue in Normal saline. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Thu 11 am; Report 8 days
B152	PAPP-A (PREGNANCY ASSOCIATED PLASMA PROTEIN-A)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA	Also see MATERNAL SERUM SCREEN 2; DUAL TEST	Sample Daily by 9 am; Report Same day
S086	PARIETAL CELL ANTIBODY, IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	For SERIAL DILUTIONS, use test code S087.	Sample Daily by 11 am; Report Same day
S087	PARIETAL CELL ANTIBODY, IFA IN DILUTIONS		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 11 am; Report Next day

Z496	PAROXYSMAL NOCTURNAL HEMOGLOBINU RIA (PNH), CONFIRMATOR Y TEST	*CD55 * FLAER on Monocytes and Granulocytes *CD59 on RBC	3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 3 mL (2 mL min.) whole blood in 1 Green Top (Sodium Heparin) tube. Ship immediately at 18–22°C. DO NOT REFRIGERATE OR FREEZE. Specify time, date and clinical details on test request form.	Flow Cytometry		Sample Daily by 9 am; Report Same day
RW145	PAROXETINE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum within 2 hours of collection. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
H014	PARTIAL THROMBOPLAS TIN TIME, ACTIVATED; PTT; APTT		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection	Also see LUPUS ANTICOAGULANT BY APTT and MIXING STUDIES.	Daily
Z185	PARVOVIRUS B19 ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Fri 9 am; Report Same day
S142	PARVOVIRUS B19 ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Fri 9 am; Report Same day

N124	PARVOVIRUS B19, QUALITATIVE PCR		3 ml (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
S143	PARVOVIRUS B19 ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Fri 9 am; Report Same day
H224	PAS (PERIODIC ACID SCHIFF) STAIN, BLOOD / BONE MARROW		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 2 mL (1 mL min.) Bone Marrow from 1 Green Top (Sodium heparin) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory.	Cytochemical staining		Sample Daily by 6 pm; Report 3 days
RW146	PCNA (PROLIFERATING CELL NUCLEAR ANTIGEN), ANTIBODY		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks
B170	PCOD (POLYCYSTIC OVARIAN DISEASE) PANEL	*FSH *LH *Testosterone, Total *Prolactin *Insulin, fasting *Glucose, fasting *TSH Ultrasensitive	4 mL (2 mL min.) serum from 1 SST AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Collect specimen 4 hrs after the patient has awakened. Ship refrigerated or frozen. Overnight fasting is mandatory.	Hexokinase, CLIA, CMIA		Daily

Z838	PCOD (POLYCYSTIC OVARIAN DISEASE) COMPREHENSIV E PANEL	*FSH *LH *Testosterone, Total *Prolactin *Insulin, fasting *Glucose, fasting * TSH Ultrasensitive *DHEAS *AMH *Lipid Profile, Screen *HOMA- IR	2 mL (1 mL min.) plasma from Gray Top (Sodium Fluoride) tubes for Glucose AND 4 mL (2 mL min.) serum from 1 SST; .Ship refrigerated. Overnight fasting is mandatory. Collect specimen 4 hrs after the patient has awakened	Hexokinase, CMIA, CLIA,Spectrophotom etry		Sample Daily by 4 pm; Report Same day
Z517	PERITONEAL EQUILIBRIUM & ADEQUACY TEST		Submit 10 mL (5 mL min.) aliquot of 24 hour Peritoneal dialysate fluid in a sterile screw capped container at 0 hour, 2 hour, 4 hour, overnight & 24 hour intervals AND 10 mL (5 mL min.) aliquot of 24 hour Urine in a sterile screw capped container AND 2 mL (1 mL min.) serum from 1 SST & 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium fluoride) tube	Spectrophotometry		Daily
Z515	PERITONEAL EQUILIBRIUM TEST (PET)		Submit 10 mL (5 mL min.) aliquot of Peritoneal dialysate fluid in a sterile screw capped container at 0 hour, 2 hour, 4 hour & overnight intervals AND 2 mL (1 mL min.) serum from 1 SST & 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium fluoride) tube both along with 2 hour samples only.	Spectrophotometry		Daily

Z826	PERNICIOUS ANEMIA PANEL	*Vitamin B12 *Gastrin *Intrinsic Factor Antibody *Parietal Cell Antibody	5 mL (3 mL min.) serum from 1 SST. Wrap sample in aluminium foil to protect from light. Ship refrigerated or frozen. Overnight fasting is mandatory.	EIA, CLIA, IFA		Sample by Mon 9 am; Report Same Day
K016	PHENCYCLIDINE (PCP), SCREEN, URINE		15 mL (10 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily
D050	PHENOBARBITONE		2 mL (1 mL min.) serum from 1 RedTop (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Chemiluminescent Immunoassay		Daily
V603	PHENOL		10 ml (5 ml min.) urine sample collected at the end of the shift on the last working day of the week. Ship refrigerated or frozen.	HPLC		Sample by Tue 7 am; Report 2 weeks
R160	PHENYLALANINE, NEWBORN SCREEN		1 drop of heel prick blood each on 3 spots of filter paper available from LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen.	Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm; Report 2 days
D049	PHENYTOIN		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Chemiluminescent Immunoassay		Daily

Z144	PHEOCHROMO CYTOMA PROFILE	*VMA *Metanephrines *Catecholamines * Chromogranin A	2 mL (1 mL min.) plasma from 1 Lavender Top EDTA tube AND 50 mL (10 mL min.) aliquot of 24-hour urine. Collect urine with appropriate volume of 50% HCl to maintain pH between 1–2. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container. Ship frozen. DO NOT THAW. Patient should strictly avoid Alpha-one blockers, Aminophylline, Amphetamines, Ampicillin, Beta-blockers, Ephedrine, Imipramine, Methyldopa, Buspiron, Codeine, L-Dopa, MAO inhibitors, Nicotine, Phenacetin, Phenothiazine, Theophylline, Vasodilators, Aspirin,PAS, Alpha-two agonists, Bromocriptine, Calcium channel blockers (long term use), Clofibrate, Propranolol, Penicillin, Reserpine, Thyroxine, Alpha-methyldopa, Isoproterenol, Labetalol, Mandelamine, Paracetamol, Cimetidine, Metoclopramide, vigorous exercise, alcoholic beverages, tea/coffee, tobacco, chocolate, pepper, vanilla & banana at least 72 hours before and during specimen collection.	High Performance Liquid Chromatography, EIA	Also see CHROMOGRANIN A.	Sample by Fri 9 am; Report Tue
Z167	PHOSPHATIDYL SERINE ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue 11 am; Report Fri
Z287	PHOSPHOLIPID ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon through Sat 9 am; Report Same day

Z377	PHOSPHOLIPID SYNDROME PANEL	*Cardiolipin Antibodies IgG & IgM *Phospholipid Antibodies Panel *Lupus Anticoagulant by dRVVT	4 mL (2.5 mL min) serum from 1 SST AND 3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory	EIA, Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Also see Ç	Sample Daily by 3 pm; Report Next day
S120	PHOSPHATIDYL SERINE ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue 11 am; Report Fri
S240	PHOSPHOLIPASE A2 RECEPTOR ANTIBODY (PLA2R)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	IFA		Sample by Mon 9 am; Report Wed
S191	PHOSPHOLIPID ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon through Sat 9 am; Report Same day

U074	PHOSPHORUS, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine. Do not use preservatives. Ship refrigerated or frozen.	Molybdate UV		Daily
B043	PHOSPHORUS, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Molybdate UV		Daily
S121	PHOSPHATIDYL SERINE ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue 11 am; Report Fri
U012	PHOSPHORUS, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine. Do not use preservatives. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen.	Molybdate UV		Daily
S192	PHOSPHOLIPID ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon through Sat 9 am; Report Same day
RW262	PIGEON FEATHERS IgE		3 mL (1.5 mL) serum from 1 Red Top (No Additive) Tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
RW261	PIGEON DROPPINGS, IgE		1 mL (0.5 mL) serum from 1 Red Top (No Additive) Tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
N114	PIK3CA MUTATION ANALYSIS		Submit formalin fixed paraffin embedded tissue block. Ship at room temperature. Block should have at least 10% tumor tissue.	PCR, Sequencing		Sample by Tuey 11 am; Report Mon

OS178	PLASMINOGEN ACTIVATOR INHIBITOR- 1/SERPINE-1, 4G/5G GENOTYPING		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed 11 am; Report 8 days
R072	PLASMA RENIN, DIRECT		3 mL (2.5 mL min.) plasma from 1 Lavender Top (EDTA) tube. Collect blood between 8 am to 10 am. Patient should be ambulatory / upright 2 hours prior to the test. Separate plasma and FREEZE. Ship frozen. DO NOT THAW. Overnight fasting is preferred. When submitting catheterization studies, it is recommended that the referring laboratory retain a portion of each sample. No medication 3 weeks prior to draw is recommended.	CLIA	Also see ANGIOTENSIN II.	Sample by Fri 9 am; Report Same day
RW147	PLASMINOGEN ACTIVATOR INHIBITOR-1 (PAI-1) ACTIVITY		2 mL (1 mL min.) plasma from 1 Blue Top (Sodium Citrate) tube. Separate plasma & freeze immediately. Ship frozen. DO NOT THAW. Overnight fasting is mandatory.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
H251	PLATELET COUNT, FLUORESCENT (PLT-F)		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Fluorescent Flow Cytometry		Daily
H010	PLATELET COUNT		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence		Daily
H210	PLATELET DISTRIBUTION WIDTH; PDW		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence		Daily

N095	PML RARA t(15;17)(q22;q1 2) GENE REARRANGEME NT QUANTITATIVE, PCR		5 mL (3 mL min.) whole blood from 2 Lavender Top (EDTA) tubes. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Wed / Fri 9 am; Report Tue / Thu / Sat
N038	PML RARA t(15;17)(q22;q1 2) GENE REARRANGEME NT QUALITATIVE, PCR		5 mL (3 mL min.) whole blood from 2 Lavender Top (EDTA) tubes. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Wed / Fri 9 am; Report Tue / Thu / Sat
M221	PNEUMOCYSTIS , IFA		Submit Sputum / Induced Sputum / BAL FLUID in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Direct IFA		Sample Daily by 9 am;Report Same day
B165A	P1NP (PROCOLLAGEN TYPE 1 AMINO TERMINAL PROPEPTIDE), TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Electrochemilumes cence		Sample by Mon / Thu 4 pm; Report same day
P097	POLARIZING MICROSCOPY, BODY FLUIDS		5 mL (2 mL min.) body fluids in a sterile screw capped container. Ship refrigerated.	Polarizing microscopy		Sample Daily by 2 pm; Report Same day
G128	POMPE DISEASE, QUANTITATIVE, BLOOD		10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium Heparin) tubes. Ship refrigerated within 48 hrs. Clinical history must accompany sample.	Enzyme assay	This tests deficiency of enzyme Alpha Glucosidase	Sample Daily by 4 pm; Report 4 days

U053	PORPHYRINS TOTAL, QUALITATIVE, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine preferably DURING ACUTE CRISIS in a dark coloured vial / container wrapped completely with black paper or aluminium foil. Do not use preservatives. Ship refrigerated or frozen. Uncovered vials exposed to light are unacceptable.	Column chromatography	The recommended screening test is PORPHYRINS, TOTAL, 24 HOURS URINE. Also see 5-AMINOLEVULINIC ACID, PORPHOBILINOGEN.	Sample Daily by 1 pm.; Report Same day
U098	PORPHOBILINO GEN (PBG), QUANTITATIVE, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine collected in a dark coloured container / wrapped completely with black paper or aluminium foil. Ship refrigerated or frozen. Record 24 hour volume on test request form and urine container. Uncovered vials exposed to light are unacceptable.	Column Chromatography	Also see PORPHYRINS, TOTAL; 5-AMINOLEVULINIC ACID.	Sample by Tue / Fri 9 am; Report Same day
U029	PORPHOBILINO GEN (PBG), QUALITATIVE, URINE		10 mL (5 mL min.) aliquot of random urine collected in a dark coloured vial / container wrapped completely with black paper or aluminium foil. Ship refrigerated or frozen. Uncovered vials exposed to light are unacceptable	Column Chromatography	Also see PORPHYRINS, TOTAL; 5-AMINO-LEVULINIC ACID.	Sample by Tue / Fri 9 am; Report Same day
U003	PORPHYRINS TOTAL, QUANTITATIVE, 24- HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine in a dark coloured vial / container wrapped completely with black paper or aluminium foil. Do not use preservatives. Ship refrigerated or frozen. Record 24 hour volume on test request form and urine container. Uncovered vials exposed to light are unacceptable.	Column chromatography	Also see 5-AMINOLEVULINIC ACID, PORPHOBILINOGEN.	Sample Daily by 1 pm.; Report Same day

U014	POTASSIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine. Do not use preservative. Blood contaminated and highly turbid samples not acceptable. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
U058	POTASSIUM, 24-HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine. Do not use preservative. Record 24 hour volume on container and test request form. Blood contaminated and highly turbid samples not acceptable. Ship refrigerated or frozen	Indirect Ion Selective Electrode		Daily
B045	POTASSIUM, SERUM		2 mL (1 mL min.) serum from 1 SST. Separate serum from cells within 1 hour of collection. Hemolysed samples not acceptable. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
P098	POTASSIUM, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Blood contaminated and highly turbid samples not acceptable.	Indirect Ion Selective Electrode		Daily
OS179	PRADER-WILLI SYNDROME MUTATION DETECTION, METHYLATION SPECIFIC PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Methylation specific PCR		Sample by Sat 11 am; Report 11 days
RW151	PREALBUMIN		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Nephelometry		Sample by 7th of the month; Report after 2–3 weeks

Z841	PRE-ECLAMPSIA PREDICTOR PANEL	*PLGF * sFLT *sFLT-1 : PLGF RATIO	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	ECLIA	NICE (2016) recommends usage of sFlt-1/PLGF ratio, with standard clinical assessment & subsequent clinical follow-up, to help rule-out pre-eclampsia in suspected women between 20 - 34 weeks+ 6 days of gestation.	Sample by Tue / Thu / Sat 9 am; Report Same Day
RW152	PREGNENOLONE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
S006	PREGNANCY TEST, URINE		10 mL (5 mL min.) aliquot of first void specimen of urine. Ship refrigerated.	Immunochromatography		Daily

G133	PRENATAL DIAGNOSIS PANEL 1, CHORIONIC VILLUS BIOPSY	This Panel tests the following diseases: *Metachromatic Leucodystrophy (Code Y556) *GM 1 Gangliosidosis (Code Y557)*GM 2 Gangliosidosis (Code Y558)*Gaucher Disease (Code Y559) *Niemann Pick Disease (Code Y560)*Fabry Disease (Code Y561)*PoÇ	Submit 20-30 mg (15 mg min.) Chorionic villus biopsy. Wash in minimal Essential media and put in Culture medium. Dissect immediately to remove decidua; Wash twice in cold 0.9% normal saline and add 0.2 mL of cold distilled water. Ship frozen within 24 hours. Sample should be taken at 10-13 weeks gestation. It is mandatory to send a Signed Prenatal Consent Form with the specimen.	Enzyme Assay	Kindly specify the disease to be tested. SAMPLE TO BE DISPATCHED WITH PRIOR APPOINTMENT.	Sample Daily by 4 pm; Report 5 days
G134	PRENATAL DIAGNOSIS PANEL 2, AMNIOTIC FLUID	*Mucopolysaccharidoses (MPS)	Submit 10 mL (7 mL min.) Amniotic fluid in a sterile screw capped container. Ship at 18–22°C within 24hrs. DO NOT FREEZE. Sample should be taken after 16 weeks gestation. It is mandatory to send a Signed Prenatal Consent Form with the specimen.	Two Dimensional electrophoresis	SAMPLE TO BE DISPATCHED WITH PRIOR APPOINTMENT.	Sample Daily by 4 pm; Report 10 days

Z614	PRE-OPERATIVE PANEL	*Blood Group & Rh Factor *Urea *Creatinine *Glucose Random *HIV Rapid *HBsAg Rapid *HCV Rapid *Prothrombin Time *Hemoglobin *Urine R/E	3 mL (2 mL min.) serum from 1 SST AND 2 mL (1 mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube AND 3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube AND 10 mL (5 mL min) aliquot of random urine in a sterile screw capped container Ship refrigerated within 48 hrs AND 3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.	See individual assays		Daily
RW153	PROCAINAMIDE	Includes metabolite N- acetylprocainami de	2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
B176	PROCALCITONI N		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen	Electrochemiluminescent Immunoassay		Daily
R011	PROGESTERONE , SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age, sex and menopausal status on test request form.	Chemiluminescent Microparticle Immunoassay		Daily
RW154	PROINSULIN		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Overnight fasting is mandatory.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks

R082	PROLACTIN, SERUM		2 mL (1 mL min.) serum from 1 SST. Collect specimen 4 hours after the patient has awakened. Overnight fasting is preferred. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Prolactin is secreted episodically. Pooling equal volume of serum from 3 samples drawn at 15 minute intervals is recommended.	Daily
RW155	PROPERDIN FACTOR B; C3 PROACTIVATOR		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum within 2 hours of collection. Ship refrigerated or frozen.	Nephelometry		Sample by 7th of the month; Report after 2–3 weeks
K020	PROPOXYPHENE SCREEN, URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen.	Lateral Flow Immunoassay		Daily
Z016	PROTEIN & IMMUNOGLOB ULIN PANEL	*Protein Electrophoresis, Serum *Immunoglobulin profile serum	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated. or frozen. Grossly hemolysed samples not accepted.	Immunoturbidimetry , Capillary electrophoresis		Sample Daily by 9 am; Report Same day
U087	PROTEIN - CREATININE RATIO, URINE		20 mL (10 mL min.) aliquot of first morning / random urine in a sterile screw capped container. .Ship refrigerated. DO NOT FREEZE.	Spectrophotometry		Daily
E001	PROTEIN ELECTROPHORE SIS, SERUM	Includes total protein, electrophoretic fractions with densitometric tracing. Does not include IFE.	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Grossly hemolysed samples not accepted.	Capillary Electrophoresis		Sample Mon through Sat by 9 am; Report Same day

RW158	PROTEIN S, ANTIGEN		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Immunoturbidimetry		Sample by 7th of the month; Report after 2–3 weeks
H061	PROTEIN C, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Chromogenic Assay	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Also see DEEC	Sample by Thu 9 am; Report Next day

Z211	PROTEIN, TOTAL, SERUM	*Total Protein *Albumin *A : G Ratio	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Biuret, BCG		Daily
Z186	PROTEIN ELECTROPHORE SIS & IFE, SERUM	Includes total protein, electrophoretic fractions, densitometric tracing of immunofixation electrophoresis.	2 mL (1.5 mL min) serum from 1 SST. Ship refrigerated or frozen. Grossly hemolysed samples not accepted.	Capillary Electrophoresis, Agarose Gel Immunofixation Electrophoresis		Sample Daily by 9 am; Report Next day
E003	PROTEIN ELECTROPHORE SIS, 24-HOUR URINE	Does not include IFE.	20 mL (10 mL min.) aliquot of 24 hour urine. Do not use preservatives. Record total volume on test request form and urine container. Ship refrigerated or frozen.	Agarose Gel Electrophoresis	For IFE, use test code E006 IMMUNO- FIXATION ELECTROPHORESIS, 24-HOUR URINE.	Sample Daily by 9 am; Report 4 days

H180	PROTEIN S, FUNCTIONAL		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician. Also see DEEC	Sample by Thu 9 am; Report Next day
P089	PROTEIN, TOTAL, BODY FLUIDS		Submit 2 mL (1 mL min.) Synovial fluid in 1 Red Top (No Additive) tube. For other body fluids, submit 10 mL (5 mL min.) in a screw-capped container. Specify source of fluid on the test request form and container. Ship refrigerated. DO NOT FREEZE.	Biuret		Daily
U004	PROTEIN, TOTAL, 24- HOUR URINE		20 mL (10 mL min.) aliquot of 24-hour urine. Do not use preservatives. Record 24 hour volume on test request form and urine container. Ship refrigerated. DO NOT FREEZE.	Spectrophotometry		Daily
U020	PROTEIN, TOTAL, RANDOM URINE		15 mL (10 mL min.) aliquot of random urine without preservative. Ship refrigerated. DO NOT FREEZE.	Spectrophotometry		Daily

RW157	PROTEIN C, ANTIGEN		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Duly filled Coagulation Requisition Form is mandatory.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
B090	PROTEIN, TOTAL, CSF		1 mL (0.5 mL min.) Cerebrospinal fluid. Avoid hemolysis. Ship refrigerated.	Spectrophotometry		Daily
N048	PROTHROMBIN GENE MUTATION ANALYSIS		3 mL (1.5 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
H131	PROTHROMBIN TIME STUDIES	*Prothrombin Time (PT) *Prothrombin Ratio (PR) *International Normalised Ratio (INR)	3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Photo-optical Clot Detection		Daily

F017	PSA (PROSTATE - SPECIFIC ANTIGEN), FREE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Do not draw sample within 7 days of Digital Rectal Examination (DRE) or Rectal Prostatic Ultrasonography.	Chemiluminescent Microparticle Immunoassay		Daily
Z262	PSA PROFILE	*PSA, Total *PSA, Free *Free/Total Percent	2.5 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Do not draw sample within 7 days of Digital Rectal Examination (DRE) or Rectal Prostatic Ultrasonography.	Chemiluminescent Microparticle Immunoassay		Daily
F014	PSA (PROSTATE - SPECIFIC ANTIGEN), TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Do not draw sample within 7 days of Digital Rectal Examination (DRE) or Rectal Prostatic Ultrasonography.	Chemiluminescent Microparticle Immunoassay		Daily
R038	PTH (PARATHYROID HORMONE), INTACT		3 mL (1.5 mL min.) serum from 1 SST drawn between 8-10 a.m in a fasting condition. Allow blood to clot. Separate serum and freeze within 30 minutes. Ship frozen. DO NOT THAW. Hemolysed samples not acceptable.	CLIA	Simultaneous determination of Serum Total Calcium and Phosphorus is recommended.	Daily
RW160	PTHrP; PARATHYROID HORMONE RELATED PEPTIDE		2 mL (1 mL min.) plasma in 1 Green Top (Sodium heparin) tube. Separate plasma immediately and transfer to a screw capped vial. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
RW161	PURKINJE CELL (Yo) ANTIBODY SCREEN, WITH REFLEX TO TITER AND WESTERN BLOT, CSF		3 mL (2 mL min.) CSF in a sterile screw capped container. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	Reflex testing if performed shall be charged additionally.	Sample by 7th of the month; Report after 2–3 weeks

RW162	PURKINJE CELL (Yo) ANTIBODY SCREEN, WITH REFLEX TO TITER AND WESTERN BLOT, SERUM		3 mL (2 mL min.) Serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	Reflex testing if performed shall be charged additionally.	Sample by 7th of the month; Report after 2–3 weeks
RW163	PYRIDINIUM COLLAGEN CROSS-LINKS (PYD), 2-HOUR URINE		10 mL (5 mL min.) aliquot of a 2-hour urine collection. Discard the first morning void. Collect urine for the next two hours, with 2 mL 6N HCl or 1g boric acid. Mix well before aliquoting. Ship refrigerated or frozen	High Performance Liquid Chromatography		Sample by 7th of the month; Report after 2–3 weeks
G148	PYRUVATE; PYRUVIC ACID		3 mL (2 mL min.) whole blood from 1 Grey Top (Sodium Fluoride) tube. Mix whole blood immediately with equal volume of 7% Perchloric Acid. Mix well and centrifuge immediately. Separate supernatant and submit minimum 2 mL. Ship refrigerated or frozen.	Spectrophotometry		Sample Daily by 4 pm; Report Same day
RW165	Q FEVER (COXIELLA BURNETTI) ANTIBODIES, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks

B160	QUANTIFERON-TB GOLD ; GAMMA INTERFERON		Collect 1 mL whole blood each in a set of special QuantiFERON – TB Gold tubes available from LPL. Shake the tubes vigorously at least 10 times to ensure thorough mixing. DO NOT CENTRIFUGE. Ship at room temperature to reach lab within 16 hours OR Incubate the tubes immediately at 37°C for 16-24 hours and centrifuge. Ship refrigerated. DO NOT FREEZE.	Enzyme Immunoassay		Sample by Mon through Sat 4 pm; Report Next Day
M175	RAPID TB ANTIGEN MPT64		Submit pure growth of AFB on appropriate media in sterile screw capped container in a sealed plastic bag. Ship refrigerated OR when used with other AFB culture or Mycobacteria identification tests, the pure growth obtained in the laboratory from specimens already submitted will be processed.	Immunochromatography		Daily
H004	RED BLOOD CELLS; RBC COUNT		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedance		Daily
H019	RED CELL DISTRIBUTION WIDTH; RDW		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedance		Daily
N103	RESPIRATORY PANEL 6	*Staphylococcus aureus *Streptococcus pneumoniae *Chlamydia pneumoniae *Haemophilus influenzae Type B	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat

N100	RESPIRATORY PANEL 3	*Coronaviruses NL63, 229E, OC43 & HKU1	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N102	RESPIRATORY PANEL 5	*Respiratory syncytial virus A/B *Cytomegalovirus *Enterovirus / Parechovirus *Adenovirus	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
Z815	RESPIRATORY COMPREHENSIVE PANEL	*Respiratory Panels 1-8 (33 Pathogens)	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3 mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N101	RESPIRATORY PANEL 4	*Parainfluenza Virus 1 *Human Metapneumovirus A/B *Bocavirus, *Mycoplasma pneumoniae	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
M205	RESPIRATORY SYNCYTIAL VIRUS (RSV), ANTIGEN		Nasal / Nasopharyngeal / BAL / Respiratory secretions in a sterile screw capped container. Ship refrigerated only.DO NOT FREEZE.	ICT		Sample Daily by 1 pmReport Same day
N098	RESPIRATORY PANEL 1	*Influenza A & B *Rhino Virus	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat

N105	RESPIRATORY PANEL 8	*Moraxella catarrhalis *Bordetella pertussis *Influenza C *Haemophilus influenzae	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3 mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N104	RESPIRATORY PANEL 7	*Klebsiella pneumoniae*Leg ionella *Pneumocystis jiroveci *Salmonella species	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N099	RESPIRATORY PANEL 2	*Parainfluenza viruses 2, 3 & 4	Submit Throat/nasal swabs in special collection medium available from LPL OR 5mL (3mL min.) Sputum / Bronchoalveolar lavage fluid. Ship refrigerated or frozen.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
RW166	RETICULIN IgA ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
H011	RETICULOCYTE COUNT		3 mL (2 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence		Daily
OS180	RETT SYNDROME DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Mon 11 am; Report 11 days
B072	RHEUMATOID FACTOR (RA), SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoturbidimetry	Also see RHEUMATOID FACTOR, IgA.	Daily

P070	RHEUMATOID FACTOR (RA), SYNOVIAL FLUID		2 mL (1 mL min.) Synovial Fluid. Ship refrigerated or frozen.	Immunoturbidimetry		Daily
Z097	RHEUMATOID AUTOIMMUNE COMPREHENSIVE PANEL	*ANA *Anti-ds DNA *Sm *U1RNP *Rheumatoid Factor *C3 *C4 *SS-A (Ro) *SS-B (La) *Centromere Antibody	6 mL (4 mL min.) serum from 2 SST's. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay, Immunoturbidimetry	Also see RHEUMATOID FACTOR, IgA.	Sample by Tue / Fri 9 am; Report Same day
Z096	RHEUMATOID AUTOIMMUNE PANEL	*ANA *Anti-ds DNA *Sm *U1RNP *Rheumatoid Factor *C3 *C4	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay, Immunoturbidimetry	Also see RHEUMATOID FACTOR, IgA.	Sample by Tue / Fri 9 am; Report Same day
Z143	RHEUMATOID ARTHRITIS PANEL	*ANA *Rheumatoid Factor *Anti Cyclic Citrullinated Peptide	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay, Immunoturbidimetry, CMIA	Also see RHEUMATOID FACTOR, IgA; Arthritis Panels	Sample Mon through Sat 9 am; Report Same day
S154	RHEUMATOID FACTOR, IgA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Sat 9 am; Report Same day
RW167	RIBOSOME P ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks

RW168	RICKETTSIA ANTIBODY PANEL	*Rocky Mountain Spotted Fever IgG & IgM *Typhus Fever, IgG & IgM	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks
RW169	RISTOCETIN COFACTOR		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.DO NOT THAW. Overnight fasting is mandatory. Duly filled Coagulation Requisition Form is mandatory.	Platelet Agglutination		Sample by 7th of the month; Report after 2–3 weeks
RW170	ROCKY MOUNTAIN SPOTTED FEVER ANTIBODIES	*Rickettsia rickettsii IgG & IgM (Both Screen and Titer)	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks
Z673	ROMA; RISK OF OVARIAN MALIGNANCY ALGORITHM	*CA125 *HE4 *Menopausal status risk estimation	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Provide menopausal status of the patient.	Chemiluminescent Microparticle Immunoassay		Sample by Mon / Thu 4 pm; Report Same day
M212	ROTAVIRUS ANTIGEN, STOOL		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pm Report Same day

S040	RUBELLA (GERMAN MEASLES) ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
B158	RUBELLA (GERMAN MEASLES) AVIDITY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample Daily 9 am. Report Same day.
Z190	RUBELLA (GERMAN MEASLES) ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
S016	RUBELLA (GERMAN MEASLES) ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
OS181	RUBELLA DETECTION , PCR		4 mL (2 mL min.) Plasma from 1 Lavendar top (EDTA) tube OR 3 mL (2 mL min.) CSF/Amniotic fluid in a sterile screw capped container OR Nasopharyngeal swab. Ship refrigerated or frozen.	PCR		Sample by Wed/Sat 11 am; Report Sat/Tue
OS185	SCA-6 (SPINO CEREBELLAR ATAXIA) DETECTION , PCR FRAGMENT ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days

OS184	SCA-3 (SPINO CEREBELLAR ATAXIA) DETECTION , PCR FRAGMENT ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days
OS182	SCA-1 (SPINO CEREBELLAR ATAXIA) MUTATION DETECTION , PCR FRAGMENT ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days
OS183	SCA-2 (SPINO CEREBELLAR ATAXIA) DETECTION , PCR FRAGMENT ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days
OS186	SCA-7 (SPINO CEREBELLAR ATAXIA) DETECTION , PCR FRAGMENT ANALYSIS		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Fragment analysis		Sample by Tue 11 am; Report 10 days
Z192	SCLERODERMA PANEL	*Scl - 70 *Centromere Antibody	3 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Tue / Fri 9 am; Report Same day

S065	ScI-70 (SCLERODERMA) ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see EXTRACTABLE NUCLEAR ANTIGENS.	Sample by Tue / Fri 9 am; Report Same day
V086	SCRUB TYPHUS (TSUTSUGAMUS HI) ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunochromatogr aphy		Sample Daily by 1 pm; Report Same day
C048	SELENIUM, 24- HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in the kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium, iodine or barium containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

C047	SELENIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container kit available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL.If gadolinium, iodine or barium containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C032	SELENIUM, SERUM		2 mL (1.0 mL min.) serum from 1 metal free White Top (Z-No Additive) tube available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium, iodine or barium containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report next day

P006	SEMEN ANALYSIS; SEMINOGRAM		Submit Semen in a sterile screw capped container. Specimen must reach the Lab within 1 hour of collection. Condoms are not acceptable. Print exact time of collection on specimen container and test request form. Patient should refrain from ejaculation for 2-7 days (ideally 3 days) prior to collection. Specimen must be obtained by masturbation only. Immediately prior to collection, patient should urinate, wash genital area with warm water and dry the part. Collect the entire ejaculate. Specimen must be kept at body temperature (shirt pocket- 37°C) after collection.	Physical Examination, Microscopy	Outstation samples not accepted. Use test code P022 for SEMEN FRUCTOSE, QUANTITATIVE available as a separate test.	Sample Daily by 1 pm; Report Same day
P022	SEMEN FRUCTOSE, QUALITATIVE		Submit semen in a sterile screw capped container. Freeze immediately after collection. Ship frozen. DO NOT THAW.	Biochemical analysis		Daily

Z804	SEPSIS PANEL	*Procalcitonin * Lactate * C- Reactive Protein *Blood culture *CBC	5 mL (3 mL min.) serum from 1 SST AND 3 mL (1 mL min.) whole blood from 1 Lavender Top (EDTA) tube AND 2 mL (1mL min.) plasma from 1 Grey Top (Sodium Fluoride) tube. Separate plasma immediately AND Collect 8-10 mL blood aseptically in special Plus Aerobic bottle available from LPL. For children less than 12 yrs collect 2 mL blood in Special Paeds Plus bottle available from LPL. Mix by gentle swirling. Ship refrigerated only. Do not use tourniquet while drawing samples.	ECLIA, Spectrophotometry, Immunoturbidimetry , Rapid Automated Aerobic culture, Electrical Impedence, VCS	Systemic inflammatory response syndrome (SIRS) may have an infectious or non-infectious etiology. If infection is suspected or proven, patients with SIRS are said to have sepsis. Microbial invasion of the blood stream is not essential since local inflammation	Sample Daily by 3.30 pm; Report Interim: 24 & 48 hrsFinal: 6 days
RW173	SEROTONIN (5 - HT), SERUM		3 mL (1.5 mL min) serum from 1 SST. Separate serum within 2 hours of collection in a plastic transport vial and freeze immediately. Ship frozen. DO NOT THAW. Avoid food rich in indoles like avocado, banana, tomato, plum, walnut, pineapple, egg plant (brinjal) and also avoid tobacco, tea and coffee 3 days prior to specimen collection.	HPLC, Fluorescence detection		Sample by 7th of the month; Report after 2–3 weeks
B140	SERUM ASCITES ALBUMIN GRADIENT (SAAG)		2 mL (0.5 mL min.) serum from 1 SST. AND 2 mL (1 mL min.) Ascitic fluid in a sterile screw capped container. Ship refrigerated or frozen.	BCG		Daily

R070	SEX HORMONE BINDING GLOBULIN; SHBG		2 mL (1 mL min.) serum from 1 SST. Specify age and sex on test request form. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample by Mon / Thu 4 pm; Report same day
B032	SGOT; ASPARTATE AMINOTRANSFE RASE (AST)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Hemolysed specimens are not acceptable.	IFCC		Daily
B033	SGPT; ALANINE AMINOTRANSFE RASE (ALT)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Hemolysed specimens are not acceptable.	IFCC		Daily
OS188	SICKLE CELL MUTATION DETECTION , PRENATAL		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube from both parents AND 10 mL (5 mL min.) Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.It is mandatory to send a Signed Prenatal Consent Form with the specimen.	PCR		Sample by Tue 11 am; Report 10 days
H042	SICKLING, SCREENING TEST		2 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Deoxygenation	The recommended test is HEMOGLOBIN HPLC / ELECTROPHORESIS.	Sample Daily by 9 am; Report Next day
OS187	SICKLE CELL MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Tue 11 am; Report 10 days
G144	SILVER NITROPRUSSIDE TEST, URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Chemical	Also see AMINO ACID PANELS, QUALITATIVE AND QUANTITATIVE	Sample Daily by 4 pm; Report Same day

C066	SILVER, SERUM		3 ml (1 ml min) serum from 1 metal free White Top (Z tube- No Additive) tube available from LPL. Ship refrigerated or frozen. Patient should avoid silver containing medication or food-supplements 3 days prior to sample collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Tue / Fri by 6 pm; Report 2 days
G162	SIROLIMUS (RAPAMYCIN)		3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated . Ideal Sampling time: 1 hour before next dose (12 hour trough)	LC-MS/MS		Sample Daily by 5 pm; Report Same day
Z142	SJOGREN'S SYNDROME ANTIBODIES	*SS-A/Ro *SS-B/La	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see EXTRACTABLE NUCLEAR ANTIGENS.	Sample by Tue / Fri 9 am; Report Same Day
Z814	SLE (SYSTEMIC LUPUS ERYTHEMATOSUS) PANEL	*ANA * ds-DNA Antibody * Sm antibody * Phospholipid antibody IgG * Phospholipid antibody IgM	4 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen.	EIA		Sample by Tue / Fri 9 am; Report Same Day
P028	SMEAR EXAMINATION, GONOCOCCUS		Submit in a slide mailer, 2 air dried smears fixed in 95% ethanol for 1 minute. Ship at room temperature.	Gram Stain	Also see CULTURE, GONOCOCCUS.	Daily

P018	SMEAR EXAMINATION, ROUTINE		Submit 2 air dried Nasal/Cervical/ Vaginal/Throat/Pus/Other Smears in a slide mailer. Ship at room temperature. Specify nature of specimen on Test Request Form and container.	Light microscopy, Giemsa stain	Does not include Papanicolaou Stain, ZN Stain and Auramine Stain. For AFB stain use test code P030, AFB STAIN, MISCELLANEOUS.	Daily
S045	SMOOTH MUSCLE ANTIBODY (ASMA), IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay	For SERIAL DILUTIONS use test code S084.	Sample Daily by 11 am; Report Same day
S084	SMOOTH MUSCLE ANTIBODY (ASMA), IFA IN DILUTIONS		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample Daily by 11 am; Report Next day
S038	Sm (SMITH) ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay	Also see EXTRACTABLE NUCLEAR ANTIGENS.	Sample by Tue / Fri 9 am; Report Same day
U025	SODIUM, 24- HOUR URINE		10 mL (5 mL min.) aliquot of 24-hour urine. Do not use preservatives. Record 24 hour volume on container and test request form. Ship refrigerated. DO NOT FREEZE.	Indirect Ion Selective Electrode		Daily
B044	SODIUM, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Indirect Ion Selective Electrode		Daily
B183	SODIUM : OSMOLALITY RATIO, SERUM		2 mL (1.5 mL min.) serum from 1 SST. Separate serum within 1 hour of collection. Ship refrigerated or frozen.	Indirect ISE, Freezing Point Depression, Calculated		Daily

P079	SODIUM, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE. Blood contaminated and highly turbid samples not acceptable	Indirect Ion Selective Electrode		Daily
U013	SODIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine. Do not use preservatives. Ship refrigerated. DO NOT FREEZE.	Indirect Ion Selective Electrode		Daily
RW174	SOMATOSTATIN, PLASMA		2 ml (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Draw blood in prechilled tube, mix by inversion and separate plasma. Freeze immediately. Ship frozen. DO NOT THAW.	Extraction, Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
U018	SPECIFIC GRAVITY, URINE		10 mL (5 mL min.) aliquot of random OR 24 hour urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Automated Strip Test	Outstation specimens not acceptable	Daily
G130	SPHINGOLIPIDOSIS PANEL 1	*Metachromatic Leucodystrophy *GM 1 Gangliosidosis *GM 2 Gangliosidosis *Gaucher Disease *Niemann Pick Disease	10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium Heparin) tubes. Ship refrigerated within 48 hrs. DO NOT FREEZE. Clinical history must accompany sample.	Enzyme assay		Sample Daily by 4 pm; Report 4 days
G132	SPHINGOLIPIDOSIS PANEL 3	*GM 1 Gangliosidosis *Gaucher Disease *Niemann Pick Disease	10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium Heparin) tubes. Ship refrigerated within 48 hrs. DO NOT FREEZE. Clinical history must accompany sample.	Enzyme assay		Sample Daily by 4 pm; Report 4 days

G131	SPHINGOLIPIDOSIS PANEL 2	*Gaucher Disease *Niemann Pick Disease	10 mL (7.5 mL min.) whole blood from 3 Lavender Top (EDTA) / Green Top (Sodium Heparin) tubes. Ship refrigerated within 48 hrs. DO NOT FREEZE. Clinical history must accompany sample.	Enzyme assay		Sample Daily by 4 pm; Report 4 days
N120	SPINAL MUSCULAR ATROPHY (SMA), CARRIER DETECTION		4mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	MLPA		Sample Mon by 11 am; Report Sat
N119	SPINAL MUSCULAR ATROPHY (SMA), MUTATION DETECTION		4mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	MLPA		Sample Mon by 11 am; Report Sat
P027	SPUTUM EXAMINATION, AFB	Includes Auramine and Ziehl Neelsen stain.	Submit Sputum in a sterile screw capped container. Ship refrigerated.	Fluorescent and conventional microscopy		Sample Daily by 3.30 pm; Report Same day
P067	SPUTUM EXAMINATION, OVA & PARASITES		Submit Sputum in a sterile screw capped container. Ship refrigerated.	Light microscopy		Daily
P066	SPUTUM FOR EOSINOPHILS		Submit Sputum in a sterile screw capped container OR smeared onto 4 sterile glass slides, air dried and shipped in a slide mailer. Ship refrigerated.	Light microscopy		Daily

P016	SPUTUM EXAMINATION, ROUTINE		Submit Sputum in a sterile screw capped container OR 2 air dried smears in a slide mailer. Ship refrigerated.	Giemsa stain, Light microscopy	Does not include Gram, Ziehl Neelsen or Auramine Stains. For AFB stain use test code P027 SPUTUM EXAMINATION, AFB.	Daily
S066	SS-A/Ro ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay	Also see EXTRACTABLE NUCLEAR ANTIGENS & SJOGREN'S SYNDROME ANTIBODIES.	Sample by Tue / Fri 9 am; Report Same day
S067	SS-B/La ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Enzyme Immunoassay	Also see EXTRACTABLE NUCLEAR ANTIGENS & SJOGREN'S SYNDROME ANTIBODIES.	Sample by Tue / Fri 9 am; Report Same day
Z314	STD PANEL	*VDRL *HIV 1 & 2, Antibodies Screening test *HSV 1 & 2, IgG & IgM *Chlamydia trachomatis IgG Antibody *TPHA	5 mL (3 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Slide flocculation, CMIA, EIA, Hemagglutination		Sample by Fri 9 am; Report Next day

G191	STEROID PANEL FOR CONGENITAL ADRENAL HYPERPLASIA (CAH)	*Androstenedione *Cortisol *11-Deoxycortisol *DHEA *17- α -Hydroxyprogesterone *Progesterone *Testosterone, Total	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days
G194	STEROID PANEL FOR PREMATURE ADRENARCHE	*Androstenedione *DHEA *DHEAS *17- α -Hydroxyprogesterone *Testosterone, Total *Estradiol	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days
G190	STEROID PANEL, 21-HYDROXYLASE DEFICIENCY/STRESS DIFFERENTIATION	*17-Hydroxyprogesterone *Androstenedione *Cortisol *17-Hydroxyprogesterone + Androstenedione : Cortisol Ratio	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS / MS		Sample Daily by 4 pm; Report 2 days

G192	STEROID PANEL FOR POLYCYSTIC OVARY SYNDROME (PCOS)	*Androstenedione *DHEA *Testosterone, Total	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days
G196	STEROID PANEL 2, 8 STEROIDS	* Androstenedione *DHEA *DHEAS *Dihydrotestosterone (DHT) *Estradiol *17- α -Hydroxyprogesterone * Progesterone *Testosterone, Total	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days

G197	STEROID PANEL 3, 13 STEROIDS	*Aldosterone *Androstenedione *Corticosterone *Cortisol *Cortisone *11-Deoxycortisol *DHEA *DHEAS *Dihydrotestosterone (DHT) * Estradiol *17 α -Hydroxyprogesterone *Progesterone *Testosterone, Total	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days
G193	STEROID PANEL FOR PCOS / CAH DIFFERENTIATION	*Androstenedione *11-Deoxycortisol * DHEA *17- α -Hydroxyprogesterone * Testosterone, Total	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days
G195	STEROID PANEL 1, 5 STEROIDS	*Aldosterone *Corticosterone *Cortisol *Cortisone*11-Deoxycortisol	2 mL (1 mL min) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen. Diurnal variation present; take sample between 8-10 am OR 4-6 pm. Specify time and date on specimen container and test request form. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 4 pm; Report 2 days

C031	STONE ANALYSIS WITH PICTURE		Submit stone (s) in a clean screw capped container. Ship at room temperature or refrigerated.	FTIR, Digital Photography	Use this test code for both Kidney & Gall stones.	Sample Daily by 9 am; Report Next day
T013	STOOL EXAMINATION, HANGING DROP PREPARATION		5 mL (3 mL min.) stool in a leak-proof screw capped container. Specimen must reach LPL within 1 hour of collection at room temperature.	Light microscopy	Also see CULTURE, CHOLERA.	Sample Daily by 5.30 pm; Report 2 hours
T001	STOOL EXAMINATION, ROUTINE; STOOL, R/E		5g (3g min.) stool in a leak-proof screw capped container. Do not send stool specimens in match boxes or absorbant material. Outstation samples not accepted.	Light microscopy		Daily
T002	STOOL EXAMINATION, OCCULT BLOOD		5g (3g min.) stool in a leak-proof screw capped container. Do not use match box or absorbant material. Patient should avoid meat products, pain killers, iron tablets, Vitamin C (Ascorbic Acid) and hard tooth brush for 3 days prior to specimen collection. Outstation specimens not accepted.	Standard Guaiac		Daily
T010	STOOL EXAMINATION, pH & REDUCING SUBSTANCES		5g (3g min.) stool in a leak-proof screw capped container. Do not use matchbox or absorbant material. Outstation samples not accepted.	pH Strip test, Benedict's Test		Daily
S149	STREPTOCOCCUS GROUP B ANTIGEN DETECTION		2 mL (1 mL min.) serum from 1 SST OR 2 mL (1 mL min.) Cerebrospinal fluid OR 10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Latex Agglutination	Also see ANTISTREPTOLYSIN - O, ANTI DNase B, and STREPTOZYME.	Daily

RW175	STREPTOZYME		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Slide Hemagglutination	Also see ANTISTREPTOLYSIN - O, ANTI DNase B, and STREPTOCOCCUS GROUP B ANTIGEN DETECTION.	Sample by 7th of the month; Report after 2–3 weeks
S100	STRIATED / SKELETAL MUSCLE ANTIBODY, IFA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample by Sat / Tue 9 am; Report 4 days
S174	STRIATED / SKELETAL MUSCLE ANTIBODY, IFA IN DILUTIONS		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred.	Immunofluorescence Assay		Sample by Sat / Tue 9 am; Report 4 days
G171	SUCCINYLACETONE, URINE		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Do not use preservatives. Ship refrigerated or frozen. Provide brief clinical history	Gas Chromatography / Mass Spectrometry		Sample Daily by 5 pm; Report 5 days
G163	SUCCINYLACETONE, BLOOD		1 drop of finger /heel prick fresh blood on special filter paper available from LPL. Ship refrigerated. DO NOT FREEZE.	LC-MS/MS		Sample by Mon / Wed / Fri 9 am; Report Wed / Fri / Mon
H226	SUDAN BLACK B STAIN, BLOOD / BONE MARROW		3 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 2 mL (1 mL min.) Bone Marrow from 1 Green Top (Sodium heparin) tube. Ship refrigerated. DO NOT FREEZE. Clinical history mandatory.	Cytochemical staining		Sample Daily by 6 pm; Report 3 days

J058	SUDAN IV STAIN		5g (3g min.) Stool OR 10 mL (5 mL min.) Urine OR 5 mL (2 mL min.) Body Fluids in a screw capped container. Ship refrigerated. DO NOT FREEZE.	Special stain, Microscopy		Sample Daily by 1 pm; Report Next day
S252	SYPHILIS SEROLOGY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Sample by Mon/Thu 9 am; Report Same day
G160	TACROLIMUS; FK506		3 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated . Ideal Sampling Time: Just before the next dose (Trough levels). Provide date and time of sampling, amount & time of dose.	LC-MS/MS		Sample Daily by 5 pm; Report Same day
N111	TACROLIMUS GENOTYPING (CYP3A5)		3 ml (2 ml min.) whole blood in 1 Lavender top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed / Sat
RW178	T4 ANTIBODY		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	Radiobinding Assay		Sample by 7th of the month; Report after 2–3 weeks

RW176	T3 ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoassay	Nonthyroidal chronic diseases, nicotinic acid and nutrition may decrease results. Thyrotoxicosis may not always cause abnormal results. Binding protein variations may cause variations in results. Increases may be associated with prÇ	Sample by 7th of the month; Report after 2–3 weeks
RW179	TBG; THYROXINE BINDING GLOBULIN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescence		Sample by 7th of the month; Report after 2–3 weeks
RW180	TBII; THYROTROPIN - BINDING INHIBITORY IMMUNOGLOB ULIN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radioreceptor Assay		Sample by 7th of the month; Report after 2–3 weeks
OS193	T-CELL GENE REARRANGEME NT MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Sat 11 am; Report 9 days

RW182	TESTOSTERONE FREE, BIOAVAILABLE & TOTAL		3 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Ship refrigerated or frozen.	LC-MS/MS, Spectrophotometry		Sample by 7th of the month; Report after 2–3 weeks
Z384	TESTOSTERONE PANEL, TOTAL & FREE		3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age and sex on test request form.	Radioimmunoassay, Chemiluminescent Immunoassay		Sample by Mon through Sat 9 am; Report Same day
Z083	TESTOSTERONE STIMULATION BY HCG		2 mL (1 mL min.) serum from 1 SST for each timed specimen. Draw fasting baseline specimen on day 1; administer HCG, 5000 IU IM. Draw second specimen for Testosterone after 72 hours. Ship refrigerated or frozen. Both containers must be clearly marked with time and date drawn, and submitted with 1 test request form.	Chemiluminescent Immunoassay	Test performed is Testosterone Total.Refer to Y Codes:Y114 Basal Y115 Post HCG 72 hoursY116 CustomizedInjection HCG is not supplied by LPL.	Daily
Z196	TESTOSTERONE, FREE & SEX HORMONE BINDING GLOBULIN		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age and sex on test request form.	Radioimmunoassay, CLIA		Sample by Mon / Thu 9 am; Report Same Day
R010	TESTOSTERONE, TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age and sex on test request form.	Chemiluminescent Immunoassay		Daily
R009	TESTOSTERONE, FREE		3 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Specify age and sex on test request form.	Radioimmunoassay		Sample by Mon through Sat 9 am; Report Same day

RW183	TETANUS ANTITOXID		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
R090	T3, FREE; FT3		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	This is the recommended test for T3 assay as it detects the biologically active fraction in circulating blood.	Daily
R091	T4, FREE; FT4		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	This is the recommended test for T4 assay as it detects the biologically active fraction in circulating blood.	Daily
Z194	T4, FREE & T3, FREE; FT4 & FT3		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
Z195	T4, FREE & TSH ULTRASENSITIVE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
N093	THALASSEMIA BETA, MUTATION ANALYSIS		4 mL (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	PCR, Reverse Hybridization		Sample by Mon 11 am; Report Saturday

Z049	THALASSEMIA PROFILE	*Iron Studies *Complete Blood Count *Hemoglobin HPLC/Electrophoresis	2 mL (1 mL min.) serum from 1 SST AND 5 mL (3.5 mL min.) whole blood from 2 Lavender Top (EDTA) tubes. Ship refrigerated. DO NOT FREEZE.	Spectrophotometry, Electrical Impedance, VCS & HPLC	HPLC is the method of choice. Hemoglobin electrophoresis (Alkaline pH) on cellulose acetate will be carried out, if necessary, in selected cases, at no additional charge and reported after 3 days.	Sample Daily by 9 am; Report Next day
C043	THALLIUM, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in the kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report next day

C042	THALLIUM, BLOOD		4 mL (2 mL min.) whole blood in 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated. DO NOT FREEZE. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C044	THALLIUM, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

P017	THROAT SWAB, ROUTINE		Submit throat swab obtained by rubbing over affected / localized lesion in aerobic transport medium available from LPL AND 2 air dried smears on sterile glass slides. Ship refrigerated	Light microscopy, Giemsa stain	Does not include Papanicolaou Stain, ZN Stain and Auramine Stain. For AFB stain use test code P030, AFB STAIN, MISCELLANEOUS.	Daily
H038	THROMBIN TIME; TT		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical Clot Detection		Daily

Z497	THROMBOPHILIA COMPREHENSIVE PROFILE	*Protein C, Functional *Protein S, Functional *Antithrombin Activity Functional *Factor V Leiden Mutation *Prothrombin Gene Mutation *MTHFR gene mutation *Lupus anticoagulant	6 mL Whole blood in 2 Blue Top (Sodium Citrate) tubes. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW AND 5 mL (3 mL min.) whole blood in 2 Lavender Top (EDTA) tubes. Ship refrigerated. DO NOT FREEZE. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical & Photo-optical Clot Detection, PCR	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Wed 6 pm; Report on Sat
Z319	THROMBOPHILIA PROFILE	*Protein C, Functional *Protein S, Functional *Antithrombin Activity Functional *Factor V Leiden Mutation *Prothrombin Gene Mutation *MTHFR gene mutation	3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen. DO NOT THAW AND 5 mL (3 mL min.) whole blood in 2 Lavender Top (EDTA) tubes. Ship refrigerated. DO NOT FREEZE. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical & Photo-optical Clot Detection, PCR	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by Wed 6 pm; Report on Sat

H246	THROMBOTIC RISK SCREEN	*Protein C Functional *Protein S Functional *Antithrombin Activity Functional *APCR	6 mL Whole blood in 2 Blue Top (Sodium Citrate) tubes. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Electromechanical & Photo-optical Clot Detection	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by 2nd / 4th Sat 11 am; Report Next Friday
Z045	THYROID PROFILE, FREE	*Free T3 *Free T4 *TSH, Ultrasensitive	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
F012	THYROGLOBULIN; Tg		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. At least 6 weeks should elapse after thyroidectomy or needle biopsy, or I131 therapy before specimen collection.	Chemiluminescent Immunoassay		Daily

Z284	THYROID PROFILE, TOTAL	*T3 Total *T4 Total *TSH	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Not the recommended test for Thyroid disorders. Use Test Code Z045 Thyroid Profile Free for better assessment of Thyroid Disorders.	Daily
Z120	THYROID, COMPREHENSIV E PANEL	*Thyroid Profile, Free * Antithyroid antibodies	5 mL (3 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Daily
S251	TISSUE TRANSGLUTAMI NASE (tTG) ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S228	TISSUE TRANSGLUTAMI NASE TTG-DGP SCREEN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen	Enzyme Immunoassay		Sample by Thu 9 am; Report Same day
S196	TISSUE TRANSGLUTAMI NASE (tTG) ANTIBODY, IgA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen	Enzyme Immunoassay		Sample Daily by 9 am; Report Same day
P083	TLC (TOTAL LEUCOCYTE COUNT), FLUID		5 mL (2 mL min.) Fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Manual count		Daily
H003	TLC (TOTAL LEUCOCYTE COUNT)		3 mL (1 mL min.) whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE.	Electrical Impedence, VCS		Daily

RW253	TOPIRAMATE		3 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Draw sample at Peak (2-4 hrs after dose) or Trough (0.5-1 hr before dose) at steady state. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
Z026	TORCH PANEL, IgM	*Toxoplasma *Rubella *Cytomegalovirus *Herpes Simplex 1+2	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
B155	TORCH PANEL AVIDITY, IgG	*Toxoplasma *Rubella *Cytomegalovirus	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample Daily by 9 am. Report Same day.
Z776	TORCH PANEL, IgG & IgM EXTENDED	*Toxoplasma IgG & IgM *Rubella IgG & IgM *Cytomegalovirus IgG & IgM *Herpes Simplex 1 IgG & IgM *Herpes Simplex 2 IgG & IgM	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA, EIA		Sample by Tue / Fri 9 am; Report Same day

Z778	TORCH PANEL, IgM EXTENDED	*Toxoplasma IgM *Rubella IgM *Cytomegalovirus IgM *Herpes Simplex 1 IgM * Herpes Simplex 2 IgM	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA, EIA		Sample by Tue / Fri 9 am; Report Same day
OS194	TORCH DNA DETECTION	*Toxoplasma* Rubella* Cytomegalovirus * Herpes simplex*	4 mL (2 mL min.) Plasma from 1 Lavender top (EDTA) tube OR 3 mL (2 mL min.) CSF/ Amniotic fluid in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	PCR		Sample by Wed/Sat 11 am; Report Fri/Tue
Z028	TORCH PANEL, IgG & IgM	*Toxoplasma *Rubella *Cytomegalovirus *Herpes Simplex 1+2	5 mL (3 mL min.) serum from 2 SST's. Ship refrigerated or frozen.	CLIA		Daily
Z777	TORCH PANEL, IgG EXTENDED	*Toxoplasma IgG *Rubella IgG *Cytomegalovirus IgG *Herpes Simplex 1 IgG * Herpes Simplex 2 IgG	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA, EIA		Sample by Tue / Fri 9 am; Report Same day
Z027	TORCH PANEL, IgG	*Toxoplasma *Rubella *Cytomegalovirus *Herpes Simplex 1+2	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily

RW186	TOXICOLOGY SCREEN, BLOOD; DRUG SCREEN		5 mL (3 mL min.) serum from 2 Red Top (No Additive) tubes . Do not use SST gel barrier tubes. Ship refrigerated or frozen.	Gas Chromatography, Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
Z197	TOXOPLASMA ANTIBODIES PANEL, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
S015	TOXOPLASMA ANTIBODY, IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
OS195	TOXOPLASMA DNA DETECTION , PCR		4 mL(2 mL min.)Whole blood from 1 Lavender top (EDTA) tube OR 3 mL (2 mL min.) Amniotic fluid /CSF/ Urine /Cord blood in a sterile screw capped container OR Chorionic Villus biopsy in normal saline. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am; Report 8 days
S039	TOXOPLASMA ANTIBODY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA		Daily
B157	TOXOPLASMA AVIDITY, IgG		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample Daily by 9 am; Report Same day.
S001	TPHA; TREPONEMA PALLIDUM HEMAGGLUTIN ATION		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Hemagglutination		Sample by Mon / Wed / Fri 9 am; Report Same day

N071	TPMT (THIOPURINE METHYL TRANSFERASE) GENOTYPING		4 mL (2 mL min.) Whole blood in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N063	t(12;21) (p13;q22); ETV6-RUNX1 (TEL-AML1), PCR QUALITATIVE		3 mL (2 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N066	t(9;11) (p21-22;q23); (MLL-AF9), PCR QUALITATIVE		3 mL (2 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N065	t(11;19) (q23;p13.3); (MLL-ENL), PCR QUALITATIVE		3 mL (2 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N062	t(1;19) (q23;p13.3); TCF3- BX1(E2A-PBX1), PCR QUALITATIVE		3 mL (2 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat
N064	t(4;11) (q21;q23); (MLL-AF4), PCR QUALITATIVE		3 mL (2 mL min.) Whole blood / Bone marrow in 1 Lavender Top (EDTA) tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat

H128	TRANSFERRIN		2 mL (1 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	Immunoturbidimetry		Daily
RW177	T3, REVERSE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	LC/MS/MS		Sample by 7th of the month; Report after 2–3 weeks
RW255	TRICYCLIC ANTIDEPRESSANTS SCREEN, SERUM	Therapeutic drug monitoring for Amitriptyline, Nortriptyline, Imipramine, Desipramine & Doxepin.	3 mL (2 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: Trough levels > 12 hours after dose at a steady state.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
P055	TRIGLYCERIDES, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen. Specify type of fluid on sample container and test request form.	GPO-POD		Daily
B026	TRIGLYCERIDES, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Minimum 12 hours overnight fasting is mandatory.	GPO-POD	Also see LIPID PROFILE TESTS.	Daily
RW188	TRIMIPRAMINE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube OR 2 mL (1 mL min.) Plasma from 1 Green Top (Sodium Heparin) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
B110	TROPONIN - I, HIGH SENSITIVE (hsTnl)		2ml (1 ml min) serum in 1 SST. Ship refrigerated or frozen.	CMIA	Serum is an acceptable sample.	Daily

B085	TROPONIN - T, HIGH SENSITIVE		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Electrochemiluminescence		Daily
RW189	TRYPsin		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoassay	Patients who have undergone in-vivo isotope therapy should refrain from testing for at least 6 weeks.	Sample by 7th of the month; Report after 2–3 weeks
RW191	TRYPTOPHAN, QUANTITATIVE, PLASMA		2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Separate plasma and freeze within 2 hours of collection. Ship frozen. DO NOT THAW. Overnight fasting is mandatory.	Liquid Chromatography, Mass Spectrometry		Sample by 7th of the month; Report after 2–3 weeks
RW190	TRYPTase		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum within 2 hours of collection. Do not use SST gel barrier tubes Ship refrigerated or frozen.	Fluorescence Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
RW192	TSH ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Radioimmunoassay		Sample by 7th of the month; Report after 2–3 weeks
R159	TSH, NEWBORN SCREEN		1 drop of heel prick blood each on 3 spots of filter paper available from LPL. LPL. Clinical details and drug history must accompany sample. Ship refrigerated or frozen.	Fluoroimmunoassay		Sample by Tue / Thu / Sat 5 pm; Report 2 days
S229	TSH RECEPTOR ANTIBODY		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	ECLIA		Sample Daily by 4 pm; Report Same day

Z274	TSH STIMULATION BY TRH (THYROTROPIN RELEASING HORMONE)		2 mL (1 mL min.) serum from 1 SST for each timed specimen. Draw fasting baseline specimen; administer TRH as IV bolus 7 ug/kg body weight (maximum 500 ug) over 30 sec. Draw additional specimens at 30 & 60 min. Ship refrigerated or frozen. All containers must be clearly marked with time drawn. Submit all specimens with 1 test request form.	Chemiluminescent Immunoassay	Refer to Y Codes:Y117 BasalY118 30 min.Y119 60 min.Y120 Customized Injection TRH is not supplied by LPL	Daily
R052	TSH (THYROID STIMULATING HORMONE)		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Sample to be given at the same time in followup cases	Daily
R092	TSH (THYROID STIMULATING HORMONE), ULTRASENSITIV E		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Sample to be given at the same time in followup cases	Daily
RW194	TSI; THYROID STIMULATING IMMUNOGLOB ULIN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunoassay		Sample by 7th of the month; Report after 2–3 weeks

R050	T3, TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Not the recommended test for Thyroid disorders as it is not the active form of the hormone. Use Test Code R090 T3, Free for better assessment of T3 levels.	Daily
R051	T4, TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Not the recommended test for Thyroid disorders as it is not the active form of the hormone. Use Test Code R091 T4, Free for better assessment of Thyroid disorders.	Daily

Z392	T3 & T4, TOTAL		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay	Not the recommended test for Thyroid disorders as it is not the active form of the hormone. Use Test Code Z194 FT3 & FT4 for better assessment of Thyroid disorders.	Daily
B146	TUMOUR NECROSIS FACTOR (TNF), ALPHA		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescent Immunoassay		Sample by Mon 4 pm; Report Same day
H062	TYPHI DOT		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunochromatography		Daily
OS196	TYPHOID DNA DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 AM; REPORT 8 DAYS
RW195	TYPHUS FEVER (RICKETTSIA) ANTIBODIES, IgG & IgM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Immunofluorescence Assay		Sample by 7th of the month; Report after 2–3 weeks
OS197	UGT1A1 GENE POLYMORPHISM (NUCLEOTIDE 'TA' REPEATS) DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am; Report 14 days

G180	UREA CYCLE DISORDER PANEL		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. No preservative required. AND 2 mL (1 mL min.) plasma from 1 Green Top (Sodium Heparin) tube. Transfer plasma into a sterile screw capped vial. Ship refrigerated or frozen. Clinical details and drug history must accompany sample.	LC-MS/MS, GC/MS		Sample Daily by 5 pm; Report 5 days
P087	UREA, FLUID		2 mL (1 mL min.) fluid in a sterile screw capped container. Ship refrigerated or frozen.	Urease UV		Daily
U028	UREA, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine. Do not use preservatives. Ship refrigerated or frozen. Refrigerate during and after collection. Record 24-hour volume on container and test request form.	Urease UV		Daily
U026	UREA NITROGEN, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine. Do not use preservatives. Ship refrigerated or frozen. Refrigerate during and after collection. Record 24 hour volume on test request form and urine container.	Urease UV		Daily
U102	UREA NITROGEN, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Urease UV		Daily
U076	UREA, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Urease UV		Daily
B021	UREA, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Urease UV		Daily
B024	URIC ACID, SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Uricase		Daily

U010	URIC ACID, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine. Do not use preservatives. Ship refrigerated or frozen. Record 24 hour volume on test request form and urine container.	Uricase		Daily
U078	URIC ACID, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine. Ship refrigerated or frozen.	Uricase		Daily
U001	URINE EXAMINATION, ROUTINE; URINE, R/E		10 mL (5 mL min.) aliquot of first morning urine in a sterile screw capped container. Ship refrigerated.	Automated Strip test, Chemical, Light microscopy		Daily
U065	URINE FOR EOSINOPHILS		10 mL (5 mL min.) aliquot of first morning urine in a sterile screw capped container. Ship refrigerated.	Light microscopy		Daily
U082	URINE, GLUCOSE		10 mL (5 mL min.) aliquot of first morning urine in a sterile screw capped container. Ship refrigerated.	Automated strip test / Benedicts Test	For Post Prandial urine, use test code U002	Daily
U066	URINE, DYSMORPHIC RBC		10 mL (5 mL min.) aliquot of random urine in a sterile screw capped container. Ship refrigerated. DO NOT FREEZE.	Phase Contrast Microscopy		Sample Daily by 2 pm; Report Same day
U080	URINE, REDUCING SUBSTANCES		10 mL (5 mL min.) aliquot of first morning urine in a sterile screw capped container. Ship refrigerated.	Dipstick, Benedict's		Daily
S068	U1RNP ANTIBODIES		2 mL (1 mL min.) serum from 1 SST. Overnight fasting is preferred. Ship refrigerated or frozen.	Enzyme Immunoassay	Also See EXTRACTABLE NUCLEAR ANTIGENS.	Sample by Tue / Fri 9 am; Report Same day
U061	UROBILINOGEN, QUALITATIVE, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine without preservative, in a dark coloured container. Ship refrigerated or frozen. Keep specimen refrigerated during collection and protect from light.	Dipstick		Daily

D031	VALPROIC ACID; VALPROATE		2 mL (1 mL min.) serum from 1 Red Top (No Additive) tube. Do not use SST gel barrier tubes. Ship refrigerated or frozen. Recommended specimen: just before the next dose is due at steady state.	Chemiluminescent Immunoassay		Daily
D053	VANCOMYCIN		2 ml (1 ml. min) serum from 1 SST. Ship refrigerated or frozen.	CMIA		Sample Daily by 4 pm; Report Same Day
M220	VANCOMYCIN RESISTANT ENTEROCOCCI (VRE) SCREENING		Rectal / Wound SWAB / Pure ISOLATE	Chromogenic		Sample Daily by 1 pm; Report Next day
Z199	VARICELLA ZOSTER VIRUS (VZV) ANTIBODIES PANEL, IgG & IgM	*Chicken Pox * Herpes Zoster * Shingles	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme immunoassay		Sample by Mon / Thu 9 am; Report Same day
S077	VARICELLA ZOSTER VIRUS (VZV) ANTIBODY, IgG	*Chicken Pox * Herpes Zoster *Shingles	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Thu 9 am; Report Same day
S074	VARICELLA ZOSTER VIRUS (VZV) ANTIBODY, IgM	*Chicken Pox *Herpes Zoster *Shingles	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Enzyme Immunoassay		Sample by Mon / Thu 9 am; Report Same day

N128	VARICELLA ZOSTER VIRUS (VZV), QUALITATIVE PCR		3 ml (2 mL min.) whole blood from 1 Lavender Top (EDTA) tube OR 3 ml (2 mL min.) CSF. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon/Thu 11 am; Report Wed/Sat
Z613	VASCULITIS PANEL	*ANA *ANCA- MPO *ANCA- PR3 *Anti - ds- DNA	3 mL (2 mL min.) serum from 1 SST. Ship refrigerated or frozen. Overnight fasting is preferred	EIA		Sample by Tue / Thu / Sat 9 am, Report Same day
RW196	VASOACTIVE INTESTINAL POLYPEPTIDE; VIP		3 mL (1.5 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma within 2 hours of collection in a plastic transport vial and freeze immediately. Ship frozen. DO NOT THAW.	Extraction, Immunoassay		Sample by 7th of the month; Report after 2–3 weeks
S011	VDRL, CSF		2 mL (1 mL min.) Cerebrospinal fluid in a sterile screw capped container. Ship refrigerated or frozen.	Slide Flocculation	Also see STD Panel	Daily
S008	VDRL (RPR), SERUM		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Slide Flocculation	Also see STD Panel	Daily
N018	VENOUS THROMBOSIS RISK	*Factor V Leiden Mutation *Prothrombin Gene Mutation *MTHFR Gene Mutation	5 mL (3 mL min.) whole blood in 2 Lavender Top (EDTA) tubes. Ship refrigerated. DO NOT FREEZE.	Real Time PCR		Sample by Mon / Thu 11 am; Report Wed / Sat

S250	VGKC (VOLTAGE GATED POTASSIUM CHANNEL) , ANTIBODY		2 mL (1 mL min) serum from 1 SST. Ship refrigerated of frozen. Overnight fasting is preferred.	Immunofluorescence	VGKCs are a family of voltage-gated potassium channels & are membrane proteins - Leucine- rich glioma inactivated protein 1(LGI1) & Contactin- associated protein 2 (CASPR2), responsible for controlling the cell membrane potential. These antigens play a direÇ	Sample by Monday 11 am. Report 5 days
RW197	VISCOSITY, SERUM		3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Kinetic Viscometer		Sample by 7th of the month; Report after 2–3 weeks
RW200	VITAMIN B12 BINDING CAPACITY, UNSATURATED; TRANSCOBALA MIN		2 mL (1 mL min.) serum from 1 SST. Separate serum within 1 hour of collection and transfer to a plastic transport vial. Ship refrigerated or frozen.	Radiobinding Assay		Sample by 7th of the month; Report after 2–3 weeks

R174	VITAMIN B12, ACTIVE; HOLOTRANSCO BALAMIN		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CMIA		Sample by Mon / Thu 4 pm ; Report Same day
R119	VITAMIN D 25 - HYDROXY	This assay measures both D2 and D3 components.	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Chemiluminescence	This assay measures both D2 and D3 components.	Daily
Z617	VITAMIN D COMPREHENSIV E	*Calcium *Phosphorus *Alkaline Phosphatase *Vitamin D 25 Hydroxy	3 mL (2 mL min.) serum from 1 SST. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship refrigerated or frozen.	CLIA, Spectrophotometry		Daily
Z616	VITAMIN D SCREEN	*Calcium *Vitamin D 25 Hydroxy	3 mL (2 mL min.) serum from 1 SST. Avoid use of tourniquet and fist clenching during sampling. Overnight fasting is preferred. Ship refrigerated or frozen.	CLIA, Spectrophotometry		Daily
RW199	VITAMIN B1 (THIAMINE)		3 mL (2 mL min.) Whole blood / Plasma from 1 Lavender Top (EDTA) tube OR 2 mL (1 mL min.) serum from one Red Top (No Additive) tube Transfer whole blood to plastic transport vial and freeze immediately. Separate plasma / serum within 4 hours of collection and freeze. Ship frozen. DO NOT THAW. Wrap in aluminium foil to protect from light. Overnight fasting is preferred.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
RW204	VITAMIN K		4 mL (2 mL min.) plasma from 2 Green Top (Sodium Heparin) tubes. Separate plasma from cells within 1 hour of collection and FREEZE. Ship frozen. DO NOT THAW. Wrap in aluminium foil to protect from light. Overnight fasting is preferred.	High Performance Liquid Chromatography		Sample by 7th of the month; Report after 2–3 weeks

RW201	VITAMIN B2; RIBOFLAVIN		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma from cells within 1 hour of collection and transfer to a plastic transport vial. Freeze immediately. Ship frozen. DO NOT THAW. Wrap in aluminium foil to protect from light.	High Performance Liquid Chromatography, Fluorometric detection		Sample by 7th of the month; Report after 2–3 weeks
R117	VITAMIN E; TOCOPHEROL		3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum from cells within 1 hour of collection. Ship refrigerated or frozen. Protect from light by wrapping in aluminium foil. Overnight fasting is mandatory.	High Performance Liquid Chromatography		Sample by Fri 9 am; Report Tue
RW202	VITAMIN B6; PYRIDOXINE		2 mL (1 mL min.) plasma from 1 Lavender Top (EDTA) tube. Separate plasma and freeze immediately. Ship frozen. DO NOT THAW. Wrap tube in aluminium foil to protect from light. Avoid vitamin supplements & alcohol 24 hours prior to sample collection. Overnight fasting is mandatory.	LC-MS/MS		Sample by 7th of the month; Report after 2–3 weeks
RW203	VITAMIN C; ASCORBIC ACID		2 mL (1 mL min.) serum from 1 SST. Separate serum within 3 hours of collection in a plastic transport vial and freeze. Ship frozen. DO NOT THAW. Wrap in aluminium foil to protect from light. Avoid vitamin supplements & fruits 24 hours prior to sample collection. Overnight fasting is preferred.	Capillary Electrophoresis		Sample by 7th of the month; Report after 2–3 weeks
R118	VITAMIN D 1, 25-DIHYDROXY	This assay measures both D2 and D3 components.	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	CLIA	This assay measures both D2 and D3 components.	Sample by Tue / Thu / Sat 9 am; Report Same day

R096	VITAMIN A (RETINOL)		3 mL (1.5 mL min.) serum from 1 Red Top (No Additive) tube. Separate serum within 1 hour of collection. Ship refrigerated or frozen. Wrap in aluminium foil to protect from light. Overnight fasting is mandatory.	High Performance Liquid Chromatography	Also see CAROTENE.	Sample by Fri 9 am; Report on Tue
R047	VITAMIN B12; CYANOCOBALA MIN		3 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen. Wrap in aluminium foil to protect from light. Overnight fasting is mandatory.	Chemiluminescent Immunoassay	Also see FOLATE, SERUM.	Daily
G183	VITAMIN D, ULTRASENSITIV E	25-Hydroxy Vitamin D2, D3 & 3 epi 25 Hydroxy Vitamin D3	3 mL (1.5 mL min.) serum from 1 SST. Ship refrigerated or frozen.	LC-MS/MS		Sample Daily by 5 pm; Report Same day
U006	VMA (VANILMANDELI C ACID), 24 HOUR URINE		50 mL (10 mL min.) aliquot of 24 hour urine. Collect urine with appropriate volume of 50% HCl to maintain pH between 1-2. Do not use concentrated HCl. Record 24 hour volume on test request form and urine container. Ship refrigerated or frozen. Patient should strictly avoid Theophylline, chocolate, vanilla, banana, alcoholic beverages, tea/coffee, tobacco and strenuous exercise at least 72 hrs before and during specimen collection.	High Performance Liquid Chromatography	Also see PHEOCHROMOCYTO MA PROFILE.	Sample by Wed / Sat 9 am; Report 2 days
U030	VMA (VANILMANDELI C ACID), RANDOM URINE		50 mL (10 mL min.) aliquot of random urine. Add 50% HCl drop by drop to bring the pH between 1–2. Do not use concentrated HCl. Ship refrigerated or frozen. Patient should strictly avoid Theophylline, chocolate, vanilla, banana, alcoholic beverages, tea/coffee, tobacco and strenuous exercise at least 72 hrs before and during specimen collection.	High Performance Liquid Chromatography	Also see PHEOCHROMOCYTO MA PROFILE.	Sample by Wed / Sat 9 am; Report 2 days

H130	von WILLEBRAND FACTOR ANTIGEN; vWF Ag		3 mL Whole blood in 1 Blue Top (Sodium Citrate) tube. Mix thoroughly by inversion. Transport to Lab within 4 hrs. If this is not possible, make PPP within 1 hour of collection as follows: Centrifuge sample at 3600 rpm for 15 min & transfer supernatant to a clean plastic tube. Centrifuge this supernatant again at 3600 rpm for 15 mins & finally transfer the supernatant (PPP) to 1 labelled, clean plastic screw capped vial. FREEZE IMMEDIATELY. Ship frozen.DO NOT THAW. Overnight fasting is preferred. Duly filled Coagulation Requisition Form is mandatory.	Enzyme Immunoassay	It is recommended that patient discontinues Heparin for 1 day and Oral Anticoagulants for 7 days prior to sampling as these drugs may affect test results. Discontinuation should be with prior consent from the treating Physician.	Sample by 1st & 3rd Mon; Report Next day
OS199	WARFARIN SENSITIVITY DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR Sequencing		Sample by Sat 11 am; Report 8 days
S002	WEIL FELIX TEST	*Proteus Antigen OX19, OX2 & OXK	2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Latex Agglutination.	Also see BRUCELLA ANTIBODY, FEBRILE AGGLUTININS PANEL & WIDAL TEST.	Sample Daily by 3.30 pm; Report Next day
S026	WIDAL SLIDE AGGLUTINATION TEST		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Slide Agglutination	Also see BRUCELLA ANTIBODY, FEBRILE AGGLUTININS PANEL & WEIL FELIX TEST.	Daily

S199	WIDAL TUBE AGGLUTINATION TEST		2 mL (1 mL min.) serum from 1 SST. Ship refrigerated or frozen.	Tube Agglutination	Also see BRUCELLA ANTIBODY, FEBRILE AGGLUTININS PANEL & WEIL FELIX TEST.	Sample Daily by 3.30 pm; Report Next day
OS200	WILSON DISEASE ATP7B GENE MUTATION DETECTION		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	PCR		Sample Daily by 11 am ; Report 30 days
Z834	WILSON DISEASE PANEL	*CBC *LFT *Ceruloplasmin* Copper. Serum *Copper. Random Urine	3 mL (1 mL min.) Whole blood in 1 Lavender Top (EDTA) tube AND 2 mL (1.5 mL min.) Serum from 1 SST AND 3 mL (1.0 mL min.) Serum from 1 Metal free White top (Z- No additive) tube available from LPL AND 10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated. Use powderless gloves during specimen collection.	Spectrophotometry, Nephelometry, Electrical Impedence, VCS, ICPMS	Use only specified containers available from LPL. If barium containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Mon through Sat 6 pm ; Report Next day
OS201	WT-1 MUTATION DETECTION		4 mL (2 mL min.) Whole Blood/ Bone marrow in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE. Clinical history is mandatory	PCR Sequencing		Sample by Wed 9 a.m report after 8 days

OS202	Y- CHROMOSOME MICRODELETION DETECTION , PCR		4 mL (2 mL min.) Whole Blood in 1 Lavender top (EDTA) Tube. Ship refrigerated. DO NOT FREEZE.	Multiplex PCR		Sample by Wed 9 a.m report after 9 days
M135	YEAST (CANDIDA SPECIES) IDENTIFICATION		Positive culture growth. Ship refrigerated.	Rapid Automated Identification		Sample Daily by 3.30 pm; Report 3 days
M136	YEAST (CANDIDA SPECIES) SUSCEPTIBILITY TESTING		Positive culture growth. Ship refrigerated.	Rapid Automated Susceptibility testing		Sample Daily by 3.30 pm; Report 3 days
M216	YERSINIA ENTEROCOLITIC A ANTIGEN		5 g fresh STOOL in a sterile screw capped container. Ship refrigerated only. DO NOT FREEZE.	ICT		Sample Daily by 1 pmReport Same day
C033	ZINC, 24 HOUR URINE		10 mL (5 mL min.) aliquot of 24 hour urine collected in an acid washed (metal free) container kit available from LPL. Mix thoroughly. Transfer 10 mL aliquot to vial provided in the kit. Measure 24 hour total volume and record on vial and test request form. Ship refrigerated or frozen. Patients should refrain from sea/ fresh water fish 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

C035	ZINC, RANDOM URINE		10 mL (5 mL min.) aliquot of random urine collected in an acid washed (metal free) container available from LPL. Ship refrigerated or frozen. Patients should refrain from sea/ fresh water fish 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day
C036	ZINC, SERUM / PLASMA		3 mL (1 mL min.) serum from 1 metal free White Top (Z-No Additive) tube OR 3 mL (1 mL min.) plasma from 1 metal free Royal Blue Top (K2EDTA) tube available from LPL. Ship refrigerated or frozen. Patient should avoid sea / fresh water fish 3 days prior to specimen collection. Use powderless gloves during specimen collection.	ICPMS	Use only specified containers available from LPL. If gadolinium or iodine containing contrast media has been administered, the specimen should not be collected for 96 hours.	Sample Daily by 6 pm; Report Next day

M191	Z-N STAIN; ZIEHL NEELSEN STAIN		Miscellaneous specimens (Sputum, Urine, Pus, Body Fluids like pleural, peritoneal, pericardial, gastric, bronchial, semen etc.) in a sterile, screw capped container OR Swabs in sterile normal saline OR Bone marrow in 1 Green Top (Sodium Heparin) tube, OR smeared onto 4 clean glass slides fixed immediately in 95% ethanol for 30 minutes and air dried. For gastric lavage , patient should be fasting for previous 12 hours. Transport the gastric lavage specimens within 4 hrs of collection. Specify type of specimen on container and test request form. Ship refrigerated.	Conventional microscopic examination	For tissue specimens use test code J050 AFB Stain	Sample Daily by 3.30 pm; Report Same day
------	--------------------------------------	--	---	--	---	--